AlaTOP Allergy Screen

DPC®

IMMULITE[®] 2000 AlaTOP Allergy Screen

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2000 Analyzer — for the qualitative detection of IgE antibodies specific to inhalant allergens in serum, as an aid in the differential diagnosis of atopic allergy.

Catalog Number: **L2KAT2** (200 tests)

Test Code: **TOP** Color: **Dark Gray**

CDC Analyte Identifier Code: 2805

CDC Test System Identifier Code: 10418

CLIA Complexity Category: Moderate

Summary and Explanation

In patients with extrinsic asthma, hay fever or atopic eczema, symptoms may appear immediately after exposure to specific allergens. Allergic response of this type (atopic or anaphylactic) is a function of serum reagins, principally immunoglobulin E (IgE). IgE molecules act as points of contact between the allergen and specialized cells that release histamine and other agents upon exposure to the allergen; this initiates the events which we recognize as allergic reactions.^{1,2} When evaluated in the light of other clinical and laboratory findings, *in vitro* allergy screening tests followed by allergen-specific IgE tests can help the physician identify the allergen (or allergens) to which an individual is sensitive.^{5,6,7,8}

Principle of the Procedure

IMMULITE 2000 AlaTOP Allergy Screen is a chemiluminescent enzyme-labeled sequential immunoassay, based on liquid ligand-labeled allergens, monoclonal antibodies, and separation by anti-ligand-coated solid phase. AlaTOP assays use a patented technology (U.S. Patent No. 4,778,751) exploiting liquid-phase kinetics in a bead format.^{3,4} It represents a significant advance over conventional methods relying on allergens attached to a solid-phase support, such as a paper disk.

AlaTOP is unique in its use of allergens in a liquid format. The allergens are covalently bound to a soluble polymer/copolymer matrix, which in turn is labeled with a ligand; anti-ligand is coated on the polystyrene bead to capture the ligand-labeled allergens.

Incubation Cycles: 2 × 30 minutes.

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 AlaTOP Allergy Screen has not been tested with all possible variations of tube types.

Volume Required: 50 µL serum.

Storage: 3 days at 2–8°C,^{9,10} or 6 months at –20°C.¹⁰

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis

B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

AlaTOP Bead Pack (L2AT12)

With barcode. 200 beads, coated with anti-ligand. Stable at 2–8°C until expiration date.

L2KAT2: 1 pack.

AlaTOP Reagent Wedge (L2ATA2)

With barcode. 11.5 mL allergen screen, with preservative. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to monoclonal murine anti-IgE antibodies in buffer, with preservative. Stable at 2–8°C until expiration date.

L2KAT2: 1 wedge.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

AlaTOP Adjustor (L2ATL)

2 mL IgE-free human serum. Stable at 2–8°C for 30 days after opening. *Do not freeze*

L2KAT2: 1 vial.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Alatop Allergy Screen Negative Control (L2ATC1)

2 mL IgE-free human serum. Stable at 2–8°C for 4 months after opening, or until the expiration date marked on the label. *Do not freeze. Do not use if signs of microbiological contamination, such as a*

cloudy appearance, have occurred.

L2KAT2: 1 vial.

Place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

AlaTOP Allergy Screen Positive Control (L2ATC2)

2 mL human serum with IgE specific for the IMMULITE 2000 AlaTOP Allergy Screen assay. Stable at 2–8°C for 4 months after opening, or until the expiration date marked on the label. *Do not freeze. Do not use if signs of microbiological contamination, such as a cloudy appearance, have occurred.*

L2KAT2: 1 vial.

Place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

Also Required

Distilled or deionized water; test tubes.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Operator's Manual.

See the IMMULITE 2000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:

4 weeks.

Quality Control Samples: The controls supplied with the kit should be used to monitor assay performance and the continued applicability of the stored Reactive/Nonreactive Cutoff.

Interpretation of Results

Reactive: A reactive result indicates that antibodies to one or more of the component allergens in the panel are present in the patient sample.

Nonreactive: A nonreactive result indicates nondetectable antibodies to the component allergens.

Limitations

The IMMULITE 2000 AlaTOP procedure provides qualitative (reactive/nonreactive) results only. It does not provide information on the amounts of IgE specific for any component allergens.

As with all diagnostic screening tests, a definitive clinical diagnosis should not be based on the result of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

The following sections contain data *representative* of the IMMULITE 2000 AlaTOP Allergy Screen kit's performance. Results are expressed in kCPS.

Unless otherwise specified, all results were generated on serum samples collected in tubes without anticoagulants, gel separators, or clot-promoting additives.

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Specificity: The antibodies are highly specific for human IgE, with no crossreactivity to other human immunoglobulins (IgA, IgD, IgG, IgM).

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 500 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Method Comparison: The assay was compared to DPC's IMMULITE AlaTOP (LKAT) on 110 patient specimens.

		IML 2000 AlaTOP	
		Nonreactive	Reactive
IML AlaTOP	React	0	64
	Non react	43	3

Agreement: 97.3%
Sensitivity: 100%; Specificity: 93.5%

References

- 1) Halpern GM. Markers of human allergic disease. J Clin Immunoassay 1983;6:131-8.
- 2) Wide L, Bennich H, Johansson SGO. Diagnosis of allergy by an *in vitro* test for allergen antibodies. Lancet 1967;2:1105-7.
- 3) El Shami AS, Alaba O. Liquid-phase *in vitro* allergen-specific IgE assay with *in situ* immobilization. Adv Biosci 1989;74:191-201.
- 4) Alaba O, El Shami AS. Evaluation of non-specific IgE binding: comparison of two *in vitro* allergen assays. Adv Biosci 1989;74:203-14.
- 5) Lichtenstein LM, et al. IgE antibody measurements in ragweed hay fever: relationship to clinical severity and the results of immunotherapy. J Clin Invest 1973;52:472-82.
- 6) Lowenstein H. Cross reactions among pollen antigens. Allergy 1980;35:198-200.
- 7) Weber RW, Nelson HS. Pollen allergens and their interrelationships. Clin Rev Allergy 1985;3:291-318.
- 8) Merrett J, Merrett TG. Screening for IgE-mediated allergy. J Clin Immunoassay 1993;16:164-72.
- 9) Bazaral M, Hamburger RN. Standardization and stability of immunoglobulin E (IgE). J Allergy 1972;49:189-91.
- 10) Tietz NW, editor. Clinical guide to laboratory tests. 3rd ed. Philadelphia: WB Saunders, 1995:358.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782 or 973.927.2828
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

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O Sistema de Qualidade da Diagnostic Products Corporation está registado sob ISO 13485:2003.
Diagnostic Products Corporation
EC REP DPC Biermann GmbH