 IMMULITE[®]
2000

EPO

DPC[®]

IMMULITE[®] 2000 EPO

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2000 Analyzer — for the quantitative measurement of erythropoietin (EPO) in serum or heparinized plasma, as an aid in the diagnosis of anemias and polycythemias.

Catalog Number: **L2KEP2** (200 tests), **L2KEP6** (600 tests)

Test Code: **EPO** Color: **Light Green**

CDC Analyte Identifier Code: 1611

CDC Test System Identifier Code: 10418

CLIA Complexity Category: Moderate

Summary and Explanation

Erythropoietin (EPO) is a glycoprotein hormone consisting of 165 amino acids, with four complex carbohydrate chains attached to the peptide at four linkage sites.¹ It has a molecular weight of 36,000 daltons, 40% of this attributed to the carbohydrate chains. EPO is the primary regulator of erythropoiesis, stimulating the proliferation and differentiation of erythroid precursor cells in bone marrow. In mammals, the fetal liver produces nearly all of the hormone; in adults, hepatic production drops to under 10% and renal secretion accounts for over 90%.^{2,3} The production site is believed to be the proximal renal tubular cells or the peritubular capillary endothelial cells of the renal cortex and outer medulla. The clearance of circulating EPO has not been fully explained, but it is accomplished, in small part, by urinary excretion, and possibly also by hepatic elimination and by uptake into target cells in bone marrow.

EPO adjusts red blood cell production to meet the tissue oxygen demand.¹ It exerts its effect in a complex feedback system, in which renal secretion of the hormone is controlled by an oxygen sensor in the kidney that responds to the partial pressure of oxygen in blood. Under conditions of increased peripheral oxygen, EPO levels diminish. This is seen after correction of hypoxia in healthy subjects (as in descent from a high elevation) and after hypertransfusion.

Anemias may be divided into two categories with respect to EPO levels in blood: those that are primary to EPO levels and those that are secondary. Primary anemias are characterized by an increase of EPO in the blood to attempt to restore red blood cell production levels to normal. Examples of anemias in which EPO levels are elevated include iron deficiency anemia, reduction of blood flow to the kidney (as in blood loss) and hemoglobinopathies with increased affinity of hemoglobin for oxygen.¹ The EPO production rate is seen to increase exponentially with the decrease in available oxygen and with falling hematocrit in nonrenal anemias; in the latter, EPO levels 1,000 times normal have been reported.²

Anemia can be secondary to inflammation, rheumatoid arthritis, neoplasm, and chronic renal disease. The "secondary anemias" may, however, be at least partly attributable to underproduction of EPO.⁴

A failure to produce sufficient EPO accounts for the moderate to severe anemias observed in end-stage renal disease. Decreased EPO production is attributed to destruction of renal production sites; the renal oxygen sensor may also be affected. Levels of the hormone slightly exceed the reference range at most, and are inadequate to counter the blood loss due to dialysis, shortened red blood cell life, iron and folate deficiency, impaired iron transfer to erythroid progenitor cells, and other challenges faced by such patients. Anephric patients demonstrate especially low EPO levels. A few patients with chronic kidney failure, however, exhibit normal hematocrits or less serious anemia, and elevated EPO. Some of these patients have cystic kidneys or viral hepatitis; in the latter, increased EPO may have resulted from enhanced hepatic production.

Overproduction of red blood cells is called polycythemia. Polycythemias may also be divided into two categories depending on whether the condition is primary or secondary to EPO levels. In polycythemia vera, EPO levels are diminished, and erythropoiesis is primary to and

independent of stimulation by EPO. Variation in EPO values can be as much as tenfold for different patients with the same hematocrit.

Certain other conditions may be characterized by the loss of feedback control of oxygen concentration over EPO production, causing an increase in EPO levels.³ These include renal cell carcinomas, in which 2 percent of patients demonstrate erythrocytosis, and some benign renal lesions, such as single or multiple renal cysts, renal artery stenosis and microvascular abnormalities. In addition, approximately 10 percent of renal transplant patients develop erythrocytosis, sometimes from autologous diseased kidney.

Secondary polycythemia is characterized by elevated EPO levels which lead to increased red blood cell mass. This condition may result from a variety of factors, including defective hemoglobin, smoking, pulmonary fibrosis, cardiac disease, tumors, and kidney stones.⁵

When assaying EPO for the differential diagnosis of polycythemias, the possible overlap of values for secondary erythrocytosis or for polycythemia vera with those in the reference range must be considered.³

Principle of the Procedure

IMMULITE 2000 EPO is a solid-phase, chemiluminescent immunometric assay.

Incubation Cycles: 2 × 30 minutes.

Specimen Collection

There are several reports of diurnal variation of erythropoietin in literature references.⁷⁻⁹ It is important to collect samples at a consistent time of day. Morning samples taken between 7:30 am and 12:00 noon have been recommended.

A variation in values may occur if the sample is not clotted at room temperature (15–28°C).

Because EDTA would have a significant effect on results, it should not be used as an anticoagulant.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed or grossly contaminated samples may give erroneous results.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 EPO has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 100 µL serum or heparinized plasma.

Storage: 7 days at 2–8°C or 2 months at –20°C. Avoid repeated freeze/thaw cycles.

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

EPO Bead Pack (L2EP12)

With barcode. 200 beads, coated with anti-ligand derived from streptavidin. Stable at 2–8°C until expiration date. **L2KEP2:** 1 pack. **L2KEP6:** 3 packs.

EPO Reagent Wedge (L2EPA2)

With barcode. 11.5 mL ligand-labeled murine monoclonal anti-EPO antibody, with preservative. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to goat polyclonal anti-EPO antibody, in buffer, with preservative. Stable at 2–8°C until expiration date. **L2KEP2:** 1 wedge. **L2KEP6:** 3 wedges.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

EPO Adjustors (LEPL, LEPH)

Two vials (Low and High) containing lyophilized human recombinant EPO in a nonhuman serum matrix, with preservative. At least 30 minutes before use: Reconstitute each vial with **4.0 mL** distilled or deionized water. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. Stable at 2–8°C for 30 days after reconstitution, or for 6 months (aliquotted) at –20°C.

L2KEP2: 1 set. **L2KEP6:** 2 sets.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately**Multi-Diluent 2 (L2M2Z, L2M2Z4)**

For the on-board dilution of patient samples. One vial, concentrated (ready-to-use), nonhuman protein/buffer matrix, with preservative. Storage: 30 days (after opening) at 2–8°C or 6 months (aliquotted) at –20°C.

L2M2Z: 25 mL **L2M2Z4:** 55 mL

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2M2Z: 3 labels. **L2M2Z4:** 5 labels.

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

LEPCM: Tri-level EPO Control Module

Also Required

Distilled or deionized water; test tubes; controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Operator's Manual

See the IMMULITE 2000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:

2 weeks.

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of EPO.

Expected Values

A study performed with IMMULITE 2000 EPO on 167 apparently healthy laboratory volunteers yielded a median of 10.2 mIU/mL, a 95th percentile of 19.6 mIU/mL and an absolute range of 3.7 to 31.5 mIU/mL.

Caution should be exercised in interpreting the EPO results for patients with diseases other than anemia. In patients with erythrocytosis due to uncompensated hypoxia, serum immunoreactive EPO is elevated; in those with compensated hypoxia, the serum immunoreactive EPO level is usually within the range of normal, and in patients with polycythemia vera, serum immunoreactive EPO is either normal or low. Thus, while an elevated serum EPO level suggests that erythrocytosis is a secondary phenomenon and a low EPO level supports the possibility of autonomous erythropoiesis, a normal serum EPO level excludes neither hypoxia nor autonomous EPO production as the cause of erythrocytosis.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

Variations in EPO results have been seen in the presence of anti-species antibodies.

No drugs have been tested for assay interference.

The results of this assay should be used in conjunction with information available from clinical evaluations and other diagnostic procedures.

Because results obtained with any EPO assay may differ significantly from any other, it is recommended that any serial testing performed on the same patient over time should be performed with the same EPO test.

Lower EPO levels than expected have been seen with anemias associated with the following conditions: rheumatoid arthritis, acquired immunodeficiency syndrome, cancer, ulcerative colitis, sickle cell disease, and in premature neonates.

After allogeneic bone marrow transplant, impaired erythropoietin response may delay erythropoietin recovery.

Patients with hypergammaglobulinemia associated with multiple myeloma or Waldenstrom's disease have impaired production of erythropoietin in relation to hemoglobin concentration, and this has been linked to increased plasma viscosity.

EPO levels of persons living at high altitudes with erythrocytosis may rapidly fall to normal after returning to low altitudes.

EDTA plasma should not be used as a sample type. (See Alternate Sample Type section).

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of

interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in mIU/mL.

Calibration Range: Up to 200 mIU/mL (WHO 2nd IRP 67/343).

Analytical Sensitivity: 0.24 mIU/mL

High-dose Hook Effect:

None up to >100,000 mIU/mL.

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three EPO solutions (25, 50 and 75 mIU/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The assay is highly specific for EPO. (See "Specificity" table).

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 512 mg/dL may cause a slight depression of values. (See "Hemolysis" table).

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To determine if samples other than serum samples are suitable for analysis in the IMMULITE 2000 EPO procedure, blood was collected from 22 volunteers into plain, heparinized, EDTA and Becton Dickinson SST[®] vacutainer tubes. Some of the matched sample types were supplemented with EPO in order to obtain results throughout the calibration range of the assay. All samples were analyzed by

the IMMULITE 2000 EPO procedure, with the following results:

(Heparin) = 1.02 (Serum) – 0.7 mIU/mL
r = 0.996

(EDTA) = 0.61 (Serum) – 1.1 mIU/mL
r = 0.986

(SST®) = 0.97 (Serum) + 0.4 mIU/mL
r = 0.998

Means:

81.5 mIU/mL (Serum)
82.6 mIU/mL (Heparin)
48.7 mIU/mL (EDTA)
79.6 mIU/mL (SST®)

Because EDTA would have a significant effect on results, it should not be used as an anticoagulant.

Method Comparison: The assay was compared to DPC's IMMULITE EPO on 160 samples. (Concentration range: approximately up to 200 mIU/mL. See graph.) By linear regression:

(IML 2000) = 1.02 (IML) + 0.1 mIU/mL
r = 0.992

Means:

34.3 mIU/mL (IMMULITE 2000)
33.4 mIU/mL (IMMULITE)

References

1) Eckardt KU, Bauer C. Erythropoietin in health and disease. *Europ J Clin Invest* 1989;19:117-27. 2) Eschbach J, Adamson J. Recombinant human erythropoietin: implications for nephrology. *Am J Kidney Dis* 1988;11:203-9. 3) Koch KM, Kuhn K, Nonnast-Daniel B, Scigalla P, volume editors. Treatment of renal anemia with recombinant human erythropoietin. In: Berlyne GM, Giovannetti S, series editors. Contributions to Nephrology. New York: Karger, 1988; 66:1-15 and 54-61. 4) Mengel et al. Hematology: principles and practice. Chicago: Year Book Medical Publishers Inc., 1973. 5) Braunwald E, Isselbacher KJ, Petersdorf RG, Wilson JD, Martin JB, Fanci AS, editors. Harrison's principles of internal medicine. 11th ed. New York: McGraw-Hill, 1987:149,1599. 6) Raine AEG. Hypertension, blood viscosity, and cardiovascular morbidity in renal failure: implications of erythropoietin therapy. *Lancet* 1988;1:97-9. 7) Miller ME, Garcia JF, Cohen RA, Cronkite EP, Moccia G, Acevado J. Diurnal levels of immunoreactive erythropoietin in normal subjects and subjects with chronic lung disease. *Br J Haematol* 1981;49:189-200. 8) Wide L, Bengtsson C, Birgegard G. Circadian Rhythm of erythropoietin in human serum. *Br J Haematol* 1989;72:85-90. 9) Cahan C, Decker MJ, Arnold JL, Washington LH, Veldhuis JD, Goldwasser E, Strohl DP. Diurnal Variations in serum erythropoietin levels in healthy subjects and sleep apnea patients. *J Appl Physiol*

1992;72:2112-7. 10) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782 or 973.927.2828
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

Manufactured by EURO/DPC Ltd. under a Quality System registered to ISO 13485:2003.

Tables and Graphs

Precision (mIU/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	4.1	0.23	5.6%	0.24	5.9%
2	12.3	0.71	5.8%	0.78	6.3%
3	16.0	0.75	4.7%	0.79	4.9%
4	31.7	1.5	4.7%	1.8	5.7%
5	64.0	3.2	5.0%	3.6	5.6%
6	126	6.4	5.1%	6.9	5.5%
7	184	9.5	5.2%	10.7	5.8%

Specificity

Compound ¹	µg/mL Added ²	% Cross-reactivity ³
Human Serum Albumin	60	ND
α-1 Antitrypsin	100	0.2%
α-Acid Glycoprotein	100	ND
Human α-Globulin	100	ND
Human Transferrin (iron-saturated)	100	ND
Human Transferrin (nonsaturated)	100	ND

ND: Not detectable.⁴

Linearity (mIU/mL)


	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	39.2	—	—
	4 in 8	17.7	19.6	90%
	2 in 8	8.3	9.8	85%
	1 in 8	4.3	4.9	88%
2	8 in 8	58.7	—	—
	4 in 8	26.5	29.4	90%
	2 in 8	13.1	14.7	89%
	1 in 8	6.3	7.3	86%
3	8 in 8	83.8	—	—
	4 in 8	39.2	41.9	94%
	2 in 8	19.8	21.0	94%
	1 in 8	9.7	10.5	92%
4	8 in 8	104	—	—
	4 in 8	49.8	52	96%
	2 in 8	24.2	26	93%
	1 in 8	12.0	13	92%
5	8 in 8	112	—	—
	4 in 8	53.4	56	95%
	2 in 8	26.5	28	95%
	1 in 8	13.3	14	95%
6	8 in 8	132	—	—
	4 in 8	64.8	66.0	98%
	2 in 8	33.1	33.0	100%
	1 in 8	15.2	16.5	92%

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
EURO/DPC LTD
Diagnostic Products Corporation
EC REP DPC Biermann GmbH