 IMMULITE[®]
2000

Ferritin

DPC[®]

IMMULITE® 2000 Ferritin

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2000 Analyzer — for the quantitative measurement of ferritin in serum, as an aid in the clinical diagnosis of iron deficiency and overload.

Catalog Number: **L2KFE2** (200 tests),
L2KFE6 (600 tests)
Test Code: **FER** Color: **Aqua**

CDC Analyte Identifier Code: 1902
CDC Test System Identifier Code: 10418
CLIA Complexity Category: Moderate

Summary and Explanation

The ferritin molecule consists of a protein shell (MW 450,000) and a core of iron.^{1,2} High concentrations are found in liver cells and in erythrocyte recycling centers (RE cells) of the liver, spleen and bone marrow. In these tissues, ferritin serves as the body's principal storehouse for surplus iron, protecting against the toxic effects of excess and maintaining a readily mobilized reserve for erythropoiesis.^{3,4} Ferritin is also present in human plasma, where its concentration is normally a satisfactory index of body iron stores as measured by quantitative phlebotomy, iron absorption studies, liver biopsy and the microscopic examination of bone marrow aspirates for stainable iron deposits.^{5,6}

This relationship with iron stores can be seen in the pattern of serum ferritin values under a variety of physiological and pathological conditions. For individuals in good health the median level, slightly elevated at birth, reaches a low of about 30 ng/mL at six months, with the increase to adult levels taking place after puberty.⁷ In males, the median level continues to climb — from about 70 ng/mL at age eighteen to almost 200 ng/mL twenty-five years later — whereas in females there is a plateau at 35 or 40 ng/mL throughout the childbearing years and a sharp increase thereafter. For adults in good health, the serum ferritin level has been variously reported to range from 20 (\pm 10) up to 300 (\pm 100) ng/mL for men, and 10 (\pm 5) up to 150 (\pm 50) ng/mL for women.

Concentrations below 10 or 15 ng/mL are typical of uncomplicated iron deficiency anemia. For iron overload, values over 300 or 400 ng/mL are the rule, with levels in the 1 000 – 5 000 ng/mL range common in full-blown cases of hemochromatosis.

Clinical applications of the serum ferritin assay have been extensively reviewed.⁸⁻¹¹ It has important roles to play in the diagnosis of iron deficiency and excess, and in the management of conditions and treatments posing a threat to iron balance. It has proved a valuable aid in discriminating iron deficiency anemia from anemias due to other causes¹² and in exposing the disappearance of iron reserves before the onset of anemia. Serial determinations have been employed to monitor, noninvasively, the progressive erosion of iron stores during pregnancy¹³ and in patients undergoing regular dialysis treatment.¹⁴ Both in company with other routine blood tests¹⁵⁻¹⁷ and on its own, the ferritin assay has been used to screen for iron deficiency in a variety of populations, ranging from blood donors^{18,19} to unselected hospital patients.^{20,21} It is also valuable in screening for precirrhotic hemochromatosis²² and other forms of iron overload²³ and in monitoring patients who are receiving regular blood transfusions²⁴ or iron replacement therapy,²⁵ and thus in danger of accumulating excessive iron stores.

Although iron depletion appears to be the only condition associated with reductions in the serum ferritin level, increases are observed not only in the presence of increased iron stores but also in several other situations, including liver disorders, inflammatory conditions, leukemia, Hodgkin's disease and certain other malignancies. Here, increased levels may reflect the escape of ferritin from damaged liver cells, impaired clearance of ferritin from the plasma, synthesis of ferritin by tumor cells, or an expansion of the iron storage compartment induced by ineffective erythropoiesis. Inflammation tends to raise the ferritin level while lowering the serum iron concentration by stimulating increased ferritin production in RE cells, using iron that would otherwise

be released to plasma transport proteins.²⁶⁻²⁸ In this condition and others, the correlation between iron stores and circulating ferritin continues to hold, but with a shift towards higher values — necessitating an adjustment in the reference range if the ferritin assay is still to be used for distinguishing normal from depleted iron reserves.²⁹

Principle of the Procedure

IMMULITE 2000 Ferritin is a solid-phase, two-site chemiluminescent immunometric assay.

Incubation Cycles: 1 × 30 minutes.

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Ferritin has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 10 µL serum.

Storage: 7 days at 2–8°C, or 2 weeks at –20°C.³⁰

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of

transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

Ferritin Bead Pack (L2FE12)

With barcode. 200 beads, coated with monoclonal murine anti-ferritin. Stable at 2–8°C until expiration date.

L2KFE2: 1 pack. **L2KFE6:** 3 packs.

Ferritin Reagent Wedge (L2FEA2)

With barcode. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to polyclonal goat anti-ferritin in buffer, with preservative. Stable at 2–8°C until expiration date.

L2KFE2: 1 wedge. **L2KFE6:** 3 wedges.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

Ferritin Adjustors (LFEL, LFEH)

Two vials (Low and High), 2.5 mL each, of ferritin in a human protein-based matrix, with preservative. Stable at 2–8°C for 30 days after opening. *Do not freeze.*

L2KFE2: 1 set. **L2KFE6:** 2 sets.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

Multi-Diluent 2 (L2M2Z, L2M2Z4)

For the on-board dilution of patient samples. One vial, concentrated (ready-to-use), nonhuman protein/buffer matrix, with preservative. Storage: 30 days (after opening) at 2–8°C or 6 months (aliquotted) at –20°C.

L2M2Z: 25 mL **L2M2Z4:** 55 mL

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2M2Z: 3 labels **L2M2Z4:** 5 labels

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

Also Required

Distilled or deionized water; test tubes; controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Operator's Manual.

See the IMMULITE 2000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval: 4 weeks.

Quality Control Samples: Use controls or serum pools with at least two levels (low and high) of ferritin.

Expected Values

Reference ranges were generated using IMMULITE 2000 Ferritin in a study involving apparently healthy blood donors who had not donated blood within a year of this sampling.

	Central 95% Range	<i>n</i>
Males	28 – 365 ng/mL	223
Females	5 – 148 ng/mL	193

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitation

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Conversion Factor:
ng/mL × 2.2 → pmol/L

Calibration Range: up to 1,500 ng/mL. Standardized in terms of the WHO 2nd IS 80/578.

Analytical Sensitivity: 0.4 ng/mL.

High-Dose Hook Effect: None up to 150 000 ng/mL.

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three ferritin solutions (1,000, 5,000 and 25,000 ng/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for ferritin.

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of packed red blood cells in concentrations up to 30 μ L/mL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 5,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 15 volunteers into plain, heparinized, EDTA and Becton Dickinson SST[®] vacutainer tubes. Equal volumes of the matched samples were spiked with various concentrations of ferritin, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE 2000 Ferritin procedure.

(Heparin) = 0.95 (Serum) + 1.8 ng/mL
r = 0.998

(EDTA) = 0.93 (Serum) + 1.3 ng/mL
r = 0.998

(SST Tubes) = 0.99 (Plain Tubes) + 1.3 ng/mL
r = 0.997

Means:

61 ng/mL (Serum)
60 ng/mL (Heparin)
58 ng/mL (EDTA)
62 ng/mL (SST Tubes)

Method Comparison: The assay was compared to DPC's IMMULITE Ferritin on 135 samples. (Concentration range: approximately 3 to 1,200 ng/mL. See graph.) By linear regression:

(IML 2000) = 1.15 (IML) – 13.4 ng/mL
r = 0.994

Means:

314 ng/mL (IMMULITE 2000)
285 ng/mL (IMMULITE)

For the 110 data pairs (out of 135) with results by both methods below 500 ng/mL:

(IML 2000) = 1.07 (IML) – 1.4 ng/mL
r = 0.986.

For the 36 data pairs (out of 135) with results by both methods below 100 ng/mL:

(IML 2000) = 0.96 (IML) + 3.7 ng/mL
r = 0.987

References

- 1) Drysdale J, et al. Human isoferritins in normal and disease states. *Semin Hematol* 1977;14:71-88.
- 2) Harrison PM. Ferritin: an iron-storage molecule. *Semin Hematol* 1977;14:55-70.
- 3) Hershko C. Storage iron regulation. *Progress Hematol* 1977;10:105-48.
- 4) Munro H, Linder M. Serum ferritin: diagnostic uses. *Physiol Rev* 1978;58:317-96.
- 5) Lipschitz D, et al. A clinical evaluation of serum ferritin as an index of iron stores. *N Engl J Med* 1974;290:1213-6.
- 6) Bezwoda WR, et al. The relationship between marrow iron stores, plasma ferritin concentrations and iron absorption. *Scand J Haematol* 1979;22:113-20.
- 7) Dallman PR, et al. Iron deficiency in infancy and childhood. *Am J Clin Nutr* 1980;33:86-118.
- 8) Jacobs A, et al. Ferritin in serum. *N Engl J Med* 1975;292:951-6.
- 9) Alfrey CP. Serum ferritin assay. *CRC Crit Rev Clin Lab Sci* 1978;9:179-208.
- 10) Worwood M. Serum ferritin. *CRC Crit Rev Clin Lab Sci* 1979;10:171-204.
- 11) Halliday J, Powell L. *Progress Hematol* 1979;11:229.
- 12) Seiler M, et al. Differentiation of iron deficiency from anemia of chronic disorders: use of serum ferritin assay. *NucCompact* 1978;9:160.
- 13) Ances IG, et al. Serum ferritin as an early determinant of decreased iron stores in pregnant women. *Southern Med J* 1979;72:591-2.
- 14) Eschbach JW, et al. Iron balance in hemodialysis patients. *Ann Intern Med* 1977;87:710-13.
- 15) Hershko C, et al. Serum ferritin and mean corpuscular volume measurement in the diagnosis of β -thalassaemia minor and iron deficiency. *Acta Haematol* 1979;62:236-9.
- 16) Meier F, Beck J. Predictive value theory and the role of serum ferritin in defining iron states. *Ligand Rev* 1980;2:34-8.
- 17) Cook JD, et al. Evaluation of the iron status of a population. *Blood* 1976;48(3):449-55.
- 18) Finch CA, et al. Effect of blood donation on iron stores as evaluated by serum ferritin. *Blood* 1977;50(3):441-7.
- 19) Birgegård G, et al. Serum ferritin levels in male blood donors: relation to number of phlebotomies and iron supplementation. *Vox Sang* 1978;34:65-70.
- 20) Ali MAM, et al. Serum ferritin concentration and bone marrow iron stores: a prospective study. *J Can Med Assoc* 1978;118:945-6.
- 21) Lindstedt G, et al. Serum-ferritin and iron-deficiency anaemia in hospital patients. *Lancet* 1980;205-6.
- 22) Powell LW, et al. Hemochromatosis: 1980 update. *Gastroenterology* 1980;78:374-81.
- 23) Jacobs A. Iron overload – clinical and pathologic aspects. *Semin Hematol* 1977;14:89-108.
- 24) Letsky EA, et al. Serum ferritin in children with thalassaemia regularly transfused. *J Clin Pathol* 1974;27:652-5.
- 25) Kaltwasser JP, et al. Serumferritin als Kontrollparameter bei oraler Eisentherapie. *Deut Med Wochenschr* 1977;102:1150-5.
- 26) Cartwright GE, Lee GR. The anaemia of chronic disorders. *Brit J Haematol* 1971;21:147-52.
- 27) Elin RJ, et al. Effect of induced fever on serum iron and ferritin concentrations in man. *Blood* 1977;49:147-53.

28) Birgegård G, et al. Serum ferritin during infection. Scand J Haematol 1978;21:333-40.
 29) Smith RJ, et al. Serum ferritin levels in the anemia of rheumatoid arthritis. J Rheumatol 1977;4:389-92. 30) Tietz NW, editor. Clinical guide to laboratory tests. 3rd ed. Philadelphia: W.B. Saunders, 1995:234.

Additional References Arosio P, et al. Ferritin: biochemistry and methods of determination. Ligand Quarterly 1981 Fall; 4(3):45-51. • Cook JD. Clinical evaluation of iron deficiency. Semin Hematol 1982;19:6-17. • Ellis D. Serum ferritin compared with other indices of iron status in children and teenagers undergoing maintenance hemodialysis. Clin Chem 1979;25:741-4. • Finch CA, Huebers H. Perspectives in iron metabolism. N Engl J Med 1982;306:1520-8. • Grail A, Hancock BW, Harrison P. Serum ferritin in normal individuals and in patients with malignant lymphoma and chronic renal failure measured with seven different commercial immunoassay techniques. J Clin Pathol 1982;35:1204-12. • Halliday J, Powell L. Iron overload. Semin Hematol 1982;19:42-53. • Kaltwasser JP. Indikation zur Serumferritin-Bestimmung. Deut Med Wochenschr 1980;105:319-21. • Lundberg PE, Lindstedt G, et al. Increase in serum ferritin concentration induced by fasting. Clin Chem 1984;30:161-3. • Valberg L. Plasma ferritin concentrations: their clinical significance and relevance to patient care. J Can Med Assoc 1980;122:1240-7.

Technical Assistance

In the United States, contact DPC's Technical Services department.
 Tel: 800.372.1782 or 973.927.2828
 Fax: 973.927.4101. Outside the United States, contact your National Distributor.

Manufactured by EURO/DPC Ltd. under a Quality System registered to ISO 13485:2003.

Tables and Graphs

Precision (ng/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	13	0.44	3.4%	0.76	5.9%
2	158	4.7	3.0%	6.3	4.0%
3	323	11	3.4%	18	5.6%
4	531	21	4.0%	28	5.3%
5	1,081	57	5.3%	78	7.2%


Linearity (ng/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	16 in 16 ⁵	27	—	—
	8 in 16	13	14	93%
	4 in 16	7.9	6.8	116%
	2 in 16	2.8	3.4	82%
	1 in 16	1.9	1.7	112%
2	16 in 16	58	—	—
	8 in 16	32	29	110%
	4 in 16	16	15	107%
	2 in 16	9.0	7.3	123%
	1 in 16	3.9	3.6	108%
3	16 in 16	110	—	—
	8 in 16	59	55	107%
	4 in 16	28	28	100%
	2 in 16	14	14	100%
	1 in 16	8.1	6.9	117%
4	16 in 16	963	—	—
	8 in 16	479	482	99%
	4 in 16	237	241	98%
	2 in 16	122	120	102%
	1 in 16	61	60	102%
5	16 in 16	1,121	—	—
	8 in 16	571	561	102%
	4 in 16	289	280	103%
	2 in 16	148	140	106%
	1 in 16	73	70	104%

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
EURO/DPC LTD
Diagnostic Products Corporation
EC REP DPC Biermann GmbH