 IMMULITE[®]
2000

Homocysteine

DPC[®]

IMMULITE[®] 2000 Homocysteine

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2000 Analyzer — for the quantitative determination of L-homocysteine in human plasma or serum. This device can assist in the diagnosis and treatment of patients suspected of having hyperhomocysteinemia or homocystinuria.

Catalog Number: **L2KHO2** (200 tests)

Test Code: **HCY** Color: **Dark Gray**

Warning: Specimens from patients who are on drug therapy involving S-adenosyl-methionine may show falsely elevated levels of homocysteine.

Results on specimens obtained from patients taking methotrexate, carbamazepine, phenytoin, nitrous oxide, anti-convulsants, or 6-azauridine triacetate should be interpreted with caution as these substances interfere with homocysteine determination.

Summary and Explanation

Total homocysteine (tHcy) has emerged as an important risk factor in the assessment of cardiovascular disease.^{3-7,12} Hcy, a thiol-containing amino acid, is produced by the intracellular demethylation of methionine. Hcy thereby serves as a pool that can be later scavenged for use in the remanufacture of either methionine through the action of the folate-dependent enzyme methionine synthase or cysteine using the B6 dependent transsulphuration pathway.^{1,2} Hcy in plasma is found primarily in a protein bound form but free, oxidized and disulfide forms are also present.

Highly elevated levels of tHcy are found in patients with homocystinuria, a rare genetic disorder of the enzymes involved in Hcy metabolism.^{1,2,3} Homocystinuria patients exhibit arterial thromboembolism, mental retardation and early arteriosclerosis.¹ Less severe genetic

defects are also associated with moderate levels of Hcy.^{3,4,5}

Homocysteine has been identified as a indicator of cardiovascular disease. A meta-analysis of 27 epidemiological studies has suggested that a 5 µmol/L increase in tHcy could be associated with an odds ratio for coronary artery disease (CAD) of 1.6 for men and 1.8 for women, the same increase in risk as a 0.5 mmol/L increase in cholesterol.⁷ In addition, patients with chronic renal disease complicated by arteriosclerotic cardiovascular disease show elevated tHcy due to the inability of the kidney to remove Hcy from the blood.^{1,2,3}

Principle of the Procedure

Competitive Immunoassay.

The IMMULITE 2000 performs an on-line one-cycle sample pretreatment of patient plasma or serum with S-adenosyl-L-homocysteine (SAH) hydrolase and dithiothreitol (DTT) solution in a reaction tube containing no bead. After a 30-minute incubation, the treated sample is transferred to a second reaction tube containing a SAH-coated polystyrene bead and an alkaline phosphatase-labeled antibody specific for SAH. During a 30-minute incubation, the converted SAH from the sample pretreatment competes with immobilized SAH for binding alkaline phosphatase-labeled anti-SAH antibody. Unbound enzyme conjugate is removed by centrifugal wash. Substrate is added and the procedure continues as described for typical immunoassays in the Operator's Manual.

Incubation Cycles: 2 × 30 minutes. 2 test positions per assay: 1 sample treatment cup; 1 immunoreaction cup.

Cycle 1: Release of bound homocysteine and conversion to SAH.

Cycle 2: Immunoreaction.

Specimen Collection

Heparinized and EDTA plasma are the samples of choice, but serum is also suitable for use. It is important to separate the plasma or serum from the cells as soon as possible after collection, as

synthesis of HCY will take place in red blood cells after sampling. **Samples should be stored on ice between the time of sampling and centrifugation.** Note that storage on ice makes the use of serum samples particularly difficult.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Homocysteine has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 15 μ L plasma or serum.

Storage: 14 days at 2–8°C, or 6 months at –20°C.

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of

potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

Homocysteine Bead Pack (L2HO12)

With barcode. 200 beads, coated with S-adenosyl-L-homocysteine (SAH). Stable at 2–8°C until the expiration date.

L2KHO2: 1 pack.

Homocysteine Reagent Wedge (L2HOA2)

With barcode. One wedge, containing 3 reagents: 15.5 mL bovine S-adenosyl-L-homocysteine hydrolase in buffer, with preservative. 4.5 mL dithiothreitol (DTT) in buffer. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to murine monoclonal anti-SAH in buffer. Stable at 2–8°C until expiration date.

L2KHO2: 1 wedge.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

Homocysteine Adjustors (LHOL, LHOH)

Two amber vials (Low and High), 2.0 mL each, of synthetically-derived S-adenosyl-L-homocysteine (SAH) in a protein/buffer matrix. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

L2KHO2: 1 set.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

Homocysteine Sample Diluent (L2HOZ)

For the on-board dilution of high samples. 25 mL of concentrated (ready-to-use) homocysteine-free protein/buffer matrix.

Storage: 30 days (after opening) at 2–8°C or 6 months (aliquotted) at –20°C.

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2HOZ: 3 labels

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

CCCM: A bi-level, nonhuman serum-based Cardiac Marker Control Module, containing homocysteine as one of four different constituents.

Also required

Distilled or deionized water; test tubes; controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Operator's Manual.

See the IMMULITE 2000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval: 4 weeks.

Quality Control Samples: Use controls or plasma or serum pools with at least two levels (low and high) of homocysteine.

Expected Values

Homocysteine levels can vary with age, gender, geographical area and genetic factors, therefore it is important for laboratories to establish their own reference ranges based on their local populations. Literature suggests a reference range of 5 – 15 µmol/L for adult males and females,¹¹ but also notes that men tend to have higher levels than women, and that postmenopausal women tend to have higher levels than premenopausal.¹²

One hundred and twenty samples from apparently healthy adult male and female volunteers age 22–66 were analyzed

using the IMMULITE 2000 Homocysteine procedure. The samples were drawn in heparinized plasma tubes and kept on ice prior to separation of the plasma from the cells. The median value was 7.7 µmol/L, with a central 95% range of

5.0 – 12 µmol/L.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

Specimens from patients who are on drug therapy involving S-adenosyl-methionine may show falsely elevated levels of homocysteine.

Although analysis of the parent compounds of carbamazepine, phenytoin, 6-azauridine, and anthopterin indicate no cross-reactivity, specimens obtained from patients treated with these drugs, as well as with methotrexate, nitrous oxide, and other anticonvulsants, should be interpreted with caution as these substances have been shown to interfere in some homocysteine assays.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in µmol/L. (Unless otherwise noted, all were generated on heparinized plasma samples.)

Calibration Range: 2 to 50 µmol/L.

Analytical Sensitivity: 0.5 µmol/L

Precision: Samples were repeatedly assayed in quadruplicate over the course of several days, for a total of 20 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three homocysteine solutions (125, 250, 500 µmol/L), were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for homocysteine. (See "Specificity" table.)

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 512 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To determine whether serum can be used in the IMMULITE 2000 Homocysteine procedure, blood was collected on ice from 34 laboratory volunteers into heparinized, EDTA and plain vacutainer tubes. Samples were separated from the cells, and matched samples were spiked with homocysteine and then assayed by the IMMULITE 2000 Homocysteine procedure, with the following results.

(EDTA) = 0.98 (Heparin) + 0.85 µmol/L
r = 0.954

(Serum) = 1.02 (Heparin) – 0.53 µmol/L
r = 0.975

Means:

20.2 µmol/L (Heparin)
20.7 µmol/L (EDTA)
20.1 µmol/L (Serum)

To assess the effect of storage temperature, heparinized, EDTA and plain vacutainer tubes were collected from five volunteers for each tube type. Some tubes were stored at room temperature, while other tubes were kept on ice for various time periods *prior to separation*. The graphs below show the effect of storage

time and temperature for heparin, EDTA and serum. (See graphs 1-3).

Method Comparison 1: The assay was compared to a commercially available manual enzyme immunoassay (Kit A) on 168 plasma samples (Concentration range: approximately 4 to 43 µmol/L. See "Method Comparison 1" graph.) By linear regression:

(IML 2000) = 1.0 (Kit A) + 0.24 µmol/L
r = 0.966

Means:

12.9 µmol/L (IMMULITE 2000)
12.7 µmol/L (Kit A)

Method Comparison 2: The assay was also compared to an in-house HPLC method in use at a reference laboratory in the United States on 95 heparinized plasma samples (Concentration range: approximately 4 to 44 µmol/L. See "Method Comparison 2" graph.) By linear regression:

(IML 2000) = 0.97 (HPLC) + 0.71 µmol/L
r = 0.974

Means:

13.4 µmol/L (IMMULITE 2000)
13.2 µmol/L (HPLC)

Method Comparison 3: The assay was compared to another commercially available immunoassay (Kit B) on 113 plasma samples (Concentration range: approximately 4 to 23 µmol/L. See "Method Comparison 3" graph.) By linear regression:

(IML 2000) = 0.90 (Kit B) – 0.02 µmol/L
r = 0.925

Means:

8.7 µmol/L (IMMULITE 2000)
9.6 µmol/L (Kit B)

References

- 1) Selhub J, Miller JW. The pathogenesis of homocysteinemia: interruption of the coordinate regulation by s-adenosylmethionine of the remethylation and transsulfuration of homocysteine. *Am J Clin Nutr* 1992;55:131-8. 2) Jacques PF, Bostom AG, Williams RR, et al. Relation between folate status, a common mutation in methylenetetrahydrofolate reductase and plasma homocysteine concentrations. *Circulation* 1996;93:7-9. 3) Malinow MR. Plasma homocysteine and arterial occlusive diseases: A mini review. *Clin Chem*. 1995;40:173-76. 4) Clarke R, Daly L, Robinson K, Naughten E, et al. Hyperhomocysteinemia: An independent risk factor for vascular disease. *New Eng J Med* 1991;324:1149-55. 5) Deloughery TG, Evans A,

Sadeghi A, et al. Common mutation in methylenetetrahydrofolate reductase: Correlation with homocysteine metabolism and late-onset vascular disease. *Circulation* 1996;94:3074-78. 6) Bostom AG, Selhub J. Homocysteine and arteriosclerosis: Subclinical and clinical disease associations. *Circulation* 1999;99:2361-63. 7) Boushey CJ, Beresford SAA, et al. A quantitative assessment of plasma homocysteine as a risk factor for vascular disease. *JAMA* 1995;274(13):1049-57. 8) Guttormsen AB, Svarstad E, et al. Elimination of homocysteine in subjects with end-stage renal failure. *Irish J Med Sci* 1995;164:8. 9) Bostom AG, Lathrop, L. Hyper-homocysteinemia in end-stage renal disease (ESDR): Prevalence, etiology and potential relationship to arteriosclerotic outcomes. *Kidney Int* 1997;52:10-20. 10) Refsum H, Ueland PM. Clinical significance of pharmacological modulation of homocysteine metabolism. *TIPS* 1990;11:411-16. 11) Ueland PM, Refsum H, Stabler SP, et al. Total homocysteine in plasma or serum: Methods and clinical applications. *Clin Chem* 1993;39:1764-79. 12) Nehler MR, Taylor LM, Porter JM. Homocysteine as a risk factor for atherosclerosis: A Review. *Cardiovasc Pathol* 1997;6:1-9.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782 or 973.927.2828
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

The Quality System of Diagnostic Products Corporation is registered to ISO 13485:2003.

Tables and Graphs

Precision (µmol/L)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	3.94	0.29	7.4%	0.41	10.4%
2	9.88	0.47	4.8%	0.75	7.6%
3	11.1	0.38	3.4%	0.46	4.1%
4	25.9	1.06	4.1%	1.31	5.1%

Specificity (µmol/L)

Compound ¹	Amount Added ² µmol/L	% Cross reactivity ³
Adenosine	5000	ND
S-adenosyl-L-methionine	500	0.6%
Cystathionine	500	6.1%
L-Cysteine	100000	ND
Gluthathione	100000	ND

ND: Not detectable.⁴

Linearity (µmol/L)


	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	19.2	—	—
	4 in 8	9.47	9.60	99%
	2 in 8	4.94	4.80	103%
	1 in 8	2.11	2.40	88%
2	8 in 8	25.2	—	—
	4 in 8	13.8	12.6	110%
	2 in 8	6.15	6.30	98%
	1 in 8	3.34	3.15	106%
3	8 in 8	31.6	—	—
	4 in 8	15.7	15.8	99%
	2 in 8	7.44	7.90	94%
	1 in 8	3.37	3.95	85%
4	8 in 8	47.5	—	—
	4 in 8	24.2	23.7	102%
	2 in 8	11.5	11.9	97%
	1 in 8	4.91	5.94	83%

Recovery (µmol/L)

	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	7.41	—	—
	A	13.3	13.3	100%
	B	20.2	19.5	104%
	C	34.4	32.0	108%
2	—	9.96	—	—
	A	14.7	15.7	94%
	B	22.6	22.0	103%
	C	36.2	34.5	105%
3	—	13.3	—	—
	A	16.9	18.9	89%
	B	26.0	25.1	104%
	C	32.0	37.6	85%

	Solution ¹	Observed ²	Expected ³	%O/E ⁴
4	—	20.2	—	—
	A	27.7	25.4	109%
	B	33.0	31.7	104%
	C	44.5	44.2	101%


Method Comparison 1

$$(IML\ 2000) = 1.0\ (Kit\ A) + 0.024\ \mu\text{mol/L}$$

$$r = 0.966$$


Method Comparison 2

$$(IML\ 2000) = 0.97\ (HPLC) + 0.71\ \mu\text{mol/L}$$

$$r = 0.974$$

Method Comparison 3


$$(IML\ 2000) = 0.90\ (Kit\ B) - 0.02\ \mu\text{mol/L}$$


$$r = 0.925$$


Effect of Storage Temperature on Different Sample Types

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O Sistema de Qualidade da Diagnostic Products Corporation está registado sob ISO 13485:2003.
EC REP DPC Biermann GmbH