 IMMULITE[®]
2000

Albumin

DPC[®]

IMMULITE[®] 2000 Albumin

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2000 Analyzer — for the detection and quantitation of microalbuminuria, that is, concentrations of albumin in human urine not detectable by conventional dipstick methods. It is intended for use as an aid in monitoring the development of incipient diabetic nephropathy.

Catalog Number: L2KHA2 (200 tests), **L2KHA6** (600 tests)

Test Code: **ALB** Color: **Brown**

CDC Analyte Identifier Code: 4019
CDC Test System Identifier Code: 10418
CLIA Complexity Category: Moderate

Summary and Explanation

Microalbuminuria is a condition characterized by increased urinary excretion of albumin in the absence of overt nephropathy.^{9,10,21} Conventional dipstick and acid precipitation tests for detecting protein in urine lack the sensitivity required to delineate this condition. Dip-sticks may yield negative or trace results even when the albumin excretion rate is 10 or 20 times normal; and the rate must increase to 200 or 300 micrograms per minute ($\mu\text{g}/\text{min}$) before nephropathy becomes clinically apparent as persistent proteinuria.^{14,21} Published studies based on highly sensitive immunoassays for albumin have established the upper limit of normal for adults as approximately 12 or 15 $\mu\text{g}/\text{min}$; accordingly, investigators have typically fixed upon a range extending from 15 or 30 $\mu\text{g}/\text{min}$ to about 150 $\mu\text{g}/\text{min}$ as an operational definition of microalbuminuria.^{11,12,20}

Measurement of urinary albumin has been used in conjunction with other assays to help localize the site of renal impairment. Thus, elevated urinary albumin levels are to be expected when the glomerular filtration process is overactive; impaired tubular resorption, on the other hand, leads to increased urinary levels of smaller molecules like β_2 -microglobulin.¹⁶

A variety of factors other than glomerular integrity — notably exercise — may have an influence on reference ranges for urinary albumin. Some reports indicate that ranges based on overnight specimens are lower than ranges based on daytime or 24-hour specimens.^{2,6,17} It is tempting to explain this in terms of differences in posture and activity which the collection times entail. But other studies have demonstrated that exercise, within limits, has no significant effect on the albumin excretion rate in healthy adults.^{1,7,15}

This does not hold true for individuals with diabetes or essential hypertension.^{1,8,13,15,18,24} Although newly diagnosed diabetics may have baseline levels within the reference range, their urinary albumin levels are likely to show an increased responsiveness to exercise, a responsiveness which can be eliminated in the short term by careful metabolic control.^{4,5,19,25}

Interest in measuring subclinical elevations in the albumin excretion rate has focused on individuals with an already established diagnosis of diabetes or essential hypertension.³ Providing proper care is taken to minimize the influence of exercise and poor metabolic control on the albumin excretion rate, the urinary albumin level has proved to be an excellent predictor of progression to overt nephropathy in both insulin-dependent and non-insulin-dependent diabetes.^{11,12,14,20,23} Thus, according to the literature, it is very unlikely that a diabetic with a resting level of less than 15 $\mu\text{g}/\text{min}$ will progress within a decade to overt clinical nephropathy; whereas diabetics who have resting albumin levels greater than 15 $\mu\text{g}/\text{min}$ even while in good metabolic control can be expected to show yearly increases of (on average) 25 $\mu\text{g}/\text{min}$.¹²

Principle of the Procedure

IMMULITE 2000 Albumin is a solid-phase, competitive chemiluminescent enzyme immunoassay.

Incubation Cycles: 1 × 30 minutes.

Specimen Collection

Collect a timed 3-hour, overnight or 24-hour urine sample — without preservative — recording the time, duration and total volume of the collection. Mix well and retain an aliquot for analysis. Samples should not be acidified.

The urine should be cleared by filtration or centrifugation before use. (Bloody specimens are unsuitable for use.) Remove the supernatant urine from the sediment and mix by *gentle* swirling or inversion. Use disposable-tip pipets.

Volume Required: 25 μ L urine.

Storage: 14 days at 2–8°C²⁸
Do *not* freeze.^{27,28}

Before assay, allow the samples to come to room temperature (15–28°C). Samples which are cloudy or have particulate material should be clarified by low-speed centrifugation.

Dilution of High Samples: All samples expected to have levels above the assay's calibration range should be diluted on-board.

Consult the Summary and Explanation section and the References for information regarding the effect of exercise and other factors on urinary albumin determinations.

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

Albumin Bead Pack (L2HA12)

With barcode. 200 beads, coated with monoclonal murine anti-HSA (human serum albumin). Stable at 2–8°C until expiration date.

L2KHA2: 1 pack. **L2KHA6:** 3 packs.

Albumin Reagent Wedge (L2HAA2)

With barcode. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to HSA (human serum albumin) in buffer. Stable at 2–8°C until expiration date.

L2KHA2: 1 wedge. **L2KHA6:** 3 wedges.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

Albumin Adjustors (LHAL, LHAH)

Two vials (Low and High) of lyophilized albumin in a protein/buffer matrix, with preservative. At least 30 minutes before use: Reconstitute each vial with **2.0 mL** distilled or deionized water. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. Stable at 2–8°C for 30 days after reconstitution, or for 6 months (aliquotted) at –20°C.

L2KHA2: 1 set. **L2KHA6:** 2 sets.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

Multi-Diluent 2 (L2M2Z, L2M2Z4)

For on-board dilution of high samples.

One vial, concentrated (ready-to-use), nonhuman protein/buffer matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

L2M2Z: 25 mL **L2M2Z4:** 55 mL

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by

the on-board reader.
L2M2Z: 3 labels **L2M2Z4:** 5 labels

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes
(16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

LHACM: Bi-level, protein-based albumin control.

Also Required

Distilled or deionized water; test tubes; controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Operator's Manual.

See the IMMULITE 2000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:
2 weeks.

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of albumin.

Expected Values

Results for albumin are usually reported as excretion rates, in micrograms of albumin per minute ($\mu\text{g}/\text{min}$). This is obtained by multiplying the concentration in $\mu\text{g}/\text{mL}$ by the total volume of the urine collection (in milliliters) and then dividing by the duration of the collection (in minutes).

Based on its relationship to DPC's IMMULITE Albumin (and DPC's Double Antibody Albumin), IMMULITE 2000 Albumin can be expected to have essentially the same reference range limits for urine samples, expressed as micrograms of albumin excreted per minute ($\mu\text{g}/\text{min}$):

nondetectable to 18 $\mu\text{g}/\text{min}$,

with a median of approximately 6 $\mu\text{g}/\text{min}$.

A total of 55 timed, 3-hour morning urine samples were collected from healthy laboratory volunteers at two separate laboratories and assayed by the DPC's

Double Antibody Albumin kit. 95% had excretion rates of less than 13 $\mu\text{g}/\text{min}$. 24-hour as well as 3-hour specimens were collected from the 20 volunteers at one of the two laboratories. 90% of the results fell below 11.5 and 12.5 $\mu\text{g}/\text{min}$ for the 24-hour and 3-hour specimens, respectively.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

There are other potential causes of microalbuminuria besides incipient diabetic nephropathy. Consideration must be given to exercise, poor diabetic control and nondiabetic renal or systemic diseases, including hypertension. Urinary tract infections and congestive heart disease are also possible causes of subclinical elevations in the albumin excretion rate.

Samples with a pH less than 4 or greater than 8 may yield results which are respectively too high or too low. Acidified samples are accordingly unsuitable for use.

The test should not be performed if the sample exhibits significant bacterial growth or if the patient shows signs of a urinary tract infection.

Bloody specimens are unsuitable for use, even if cleared by centrifugation, since blood is a likely sign of contamination — e.g. from menstrual flow — and albumin circulates in plasma at levels approximately 4,000 times those normally encountered in urine.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in $\mu\text{g}/\text{mL}$. (Unless otherwise noted, all results were generated on urine samples.)

Calibration Range: 2.5 – 60 $\mu\text{g}/\text{mL}$

Analytical Sensitivity: 1.0 $\mu\text{g}/\text{mL}$

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three albumin solutions (100, 200 and 400 µg/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for human albumin. (See "Specificity" table.)

Method Comparison: The assay was compared to DPC's IMMULITE Albumin on 117 urine samples (Concentration range: approximately 2.5 to 60 µg/mL. See graph.) By linear regression:

$$(IML\ 2000) = 1.08 (IML) - 1.4\ \mu\text{g/mL}$$
$$r = 0.988$$

Means:

19.6 µg/mL (IMMULITE 2000)

19.4 µg/mL (IMMULITE)

References

1) Dahlquist G, et al. Effect of metabolic control and duration on exercise-induced albuminuria in diabetic teenagers. *Acta Paediatr Scand* 1983;72:895-902. 2) Davies AG, et al. Urinary albumin excretion in school children. *Arch Dis Child* 1984;59:625-30. 3) Hostetter TH. Diabetic nephropathy. *New Engl J Med* 1985;312:642-4. 4) Kroc Collaborative Study Group. Blood glucose control and the evolution of diabetic retinopathy and albuminuria; a preliminary multicenter trial. *New Engl J Med* 1984;311:365-72. 5) Microalbuminuria and renal hyperfiltration. In: Gambino R, editor. *Lab Report for Physicians* 1984; 6(8): 57-61. 6) Mogensen CE. Urinary albumin excretion in early and long-term juvenile diabetes. *Scand J Clin Lab Invest* 1971;28:183-93. 7) Mogensen CE, Vittinghus E, Sølling K. Abnormal albumin excretion after two provocative renal tests in diabetes: physical exercise and lysine injection. *Kidney Int* 1979;16:385-93. 8) Mogensen CE. Antihypertensive treatment inhibiting the progression of diabetic nephropathy. *Acta Endocrinol* 1980;94 (suppl. 238):103-11. 9) Mogensen CE, et al. Early changes in kidney function, blood pressure and the stages in diabetic nephropathy. In: Keen H, Legrain M, editors. *Prevention and treatment of diabetic nephropathy*. Boston: MTP Press, 1983: 57-83, 99-105. 10) Mogensen CE, Christensen CK, Vittinghus E. The stages in diabetic renal disease, with emphasis on the stage of incipient diabetic nephropathy. *Diabetes* 1983;32 (Suppl. 2):64-78. 11) Mogensen CE. Microalbuminuria predicts clinical proteinuria and early mortality in maturity-onset diabetes. *New Engl J Med* 1984;310:356-60. 12) Mogensen CE, Christensen CK. Predicting diabetic nephropathy in insulin-dependent

patients. *New Engl J Med* 1984;311:89-93. 13) Mohamed A, et al. Response of urinary albumin to submaximal exercise in newly diagnosed non-insulin dependent diabetes. *Br Med J* 1984;288:1342-3. 14) Parving H-H, et al. Early detection of patients at risk of developing diabetic nephropathy; a longitudinal study of urinary albumin excretion. *Acta Endocrinol* 1982;100:550-5. 15) Pedersen EB, Mogensen CE, Larsen JS. Effects of exercise on urinary excretion of albumin and β₂-microglobulin in young patients with mild essential hypertension without treatment and during long-term propranolol treatment. *Scand J Clin Lab Invest* 1981;41:493-8. 16) Ratge D, Wissner H. Urinary protein profiling by high performance gel permeation chromatography. *J Chromatogr* 1982;230:47-56. 17) Rowe DJF, et al. Effect of glycaemic control and duration of disease on overnight albumin excretion in diabetic children. *Br Med J* 1984;289:957-9. 18) Viberti GC, et al. Increased glomerular permeability to albumin induced by exercise in diabetic subjects. *Diabetologia* 1978;14:293-300. 19) Viberti GC, et al. Effect of control of blood glucose on urinary excretion of albumin and β₂-microglobulin in insulin-dependent diabetes. *New Engl J Med* 1979;300:638-41. 20) Viberti GC, et al. Microalbuminuria as a predictor of clinical nephropathy in insulin-dependent diabetes mellitus. *Lancet* 1982;(26 Jun):1430-2. See also *ibid* 1982;ii(7 Aug):330, 1982;ii(30 Oct):985, 1982;ii(18 Dec):1407 and 1983;ii(12 Feb):352. 21) Viberti GC, et al. Proteinuria in diabetes mellitus: role of spontaneous and experimental variation of glycemia. *Kidney Int* 1982;21:714-20. 22) Viberti GC, Mogensen CE, et al. Urinary excretion of albumin in normal man: the effect on water loading. *Scand J Clin Lab Invest* 1982;42:147-51. 23) Viberti GC, et al. Long term correction of hyperglycemia and progression of renal failure in insulin dependent diabetes. *Br Med J* 1983;286:598-602. 24) Vittinghus E, Mogensen CE. Graded exercise and protein excretion in diabetic man and the effect of insulin treatment. *Kidney Int* 1982;21:725-9. 25) Wiseman MJ, et al. Effect of blood glucose control on increased glomerular filtration rate and kidney size in insulin-dependent diabetes. *New Engl J Med* 1985;312:617-21. 26) Rosenstock J, Raskin P. Early diabetic nephropathy: assessment and potential therapeutic interventions. *Diabetes Care* 1986;9:529-45. 27) Erving LD, et al. Screening of microalbuminuria in patients with diabetes mellitus: frozen storage of urine samples decreases their albumin content. *Clin Chem* 1989;35(2):308-10. 28) Osberg I, et al. Effects of storage time and temperature on measurement of small concentrations of albumin in urine. *Clin Chem* 1990;36(8):1428-30.

Technical Assistance

In the United States, contact DPC's Technical Services department.

Tel: 800.372.1782 or 973.927.2828
 Fax: 973.927.4101. Outside the United States, contact your National Distributor.

Manufactured by EURO/DPC Ltd. under a Quality System registered to ISO 13485:2003.

Tables and Graphs

Precision (µg/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	7.6	0.52	6.8%	0.55	7.2%
2	15.6	1.01	6.5%	1.16	7.4%
3	18.2	1.05	5.8%	1.81	9.9%
4	35.1	2.61	7.4%	3.32	9.5%
5	43.5	3.58	8.2%	4.08	9.4%

Linearity (µg/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8	8.21	—	—
	4 in 8	4.22	4.11	103%
	2 in 8	2.25	2.05	110%
	1 in 8	1.08	1.03	105%
2	8 in 8	15.0	—	—
	4 in 8	7.38	7.48	99%
	2 in 8	3.53	3.74	94%
	1 in 8	2.12	1.87	113%
3	8 in 8	21.5	—	—
	4 in 8	11.9	10.8	110%
	2 in 8	5.87	5.38	109%
	1 in 8	2.23	2.69	83%
4	8 in 8	23.5	—	—
	4 in 8	11.9	11.7	102%
	2 in 8	5.95	5.88	101%
	1 in 8	2.51	2.94	85%
5	8 in 8	33.8	—	—
	4 in 8	19.1	16.9	113%
	2 in 8	9.90	8.45	117%
	1 in 8	4.45	4.23	105%
6	8 in 8	38.9	—	—
	4 in 8	20.7	19.4	107%
	2 in 8	10.4	9.72	107%
	1 in 8	5.09	4.86	105%

Recovery (µg/mL)


	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	11.1	—	—
	A	16.8	15.5	108%
	B	23.0	20.5	112%
	C	28.3	30.5	93%
2	—	12.2	—	—
	A	17.5	16.6	105%
	B	23.1	21.6	107%
	C	32.3	31.6	102%
3	—	32.4	—	—
	A	39.1	35.8	109%
	B	39.9	40.8	98%
	C	47.5	50.8	94%
4	—	34.0	—	—
	A	35.8	37.3	96%
	B	39.5	42.3	93%
	C	54.3	52.3	104%

Specificity

Compound ¹	µg/mL Added ²	% Cross-reactivity ³
Alpha-fetoprotein	20	0.9%
Apotransferrin	100	ND
Urea	625	ND
Human IgG	500	ND
Creatinine	2,500	ND
Bovine serum albumin	10,000	ND
Transferrin	100	ND
Conjugated Bilirubin	50	ND
	20	ND
	10	ND
Gentamycin sulfate	120	ND
Glucose	5,000	ND

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
EURO/DPC LTD
Diagnostic Products Corporation
EC REP DPC Biermann GmbH