 IMMULITE<sup>®</sup>  
2000

## **Insulin**

**DPC<sup>®</sup>**

# IMMULITE<sup>®</sup> 2000 Insulin

---

## English

---

**Intended Use:** For *in vitro* diagnostic use with the IMMULITE 2000 Analyzer — for the quantitative measurement of insulin in serum or heparinized plasma, for the management of diabetes.

Catalog Number: **L2KIN2** (200 tests),  
**L2KIN6** (600 tests)  
Test Code: **INS** Color: **Orange**

### Summary and Explanation

Human insulin is a polypeptide hormone originating in the beta cells of the pancreas and serving as a principal regulator for the storage and production of carbohydrates. Its secretion is normally stimulated by increases in the amount of glucose in circulation. This leads to higher insulin levels and more rapid tissue-assimilation of glucose — followed by a decline in the insulin level as the glucose level subsides.

In a number of conditions, notably insulinoma and diabetes, this relationship is impaired. Insulin tends to circulate at inappropriately high levels in patients with insulin-secreting pancreatic tumors; such tumors can thus be a cause of hypoglycemia. Accordingly, insulin immunoassays — used sometimes in connection with provocative doses of tolbutamide or calcium — play an essential role in the identification (and localization) of insulinomas. The finding of fasting hypoglycemia in association with an *inappropriately high* serum insulin concentration is considered diagnostic.

Insulin levels do not figure in the subclassification of diabetes worked out by the National Diabetes Data Group. Nevertheless, when obtained in the course of a glucose tolerance test, they appear to be of some prognostic value in predicting the benefits of insulin therapy and the likelihood of progression to insulin-dependence and the complications (such as retinopathy) characteristic of diabetes.

The application of insulin immunoassay to patients already undergoing insulin

therapy is complicated by the fact that such therapy typically leads to the formation of anti-insulin antibodies capable of interfering with the assay. Some investigators have sought therefore to measure insulin in urine, or in serum samples subjected to column chromatography or PEG precipitation. But the measurement of "free" insulin remains of limited interest as a technique for monitoring insulin therapy in the absence of statistics establishing therapeutic or toxic ranges. So far it appears that glucose control in diabetics cannot in general be achieved by normalizing the insulin profile. Nor is it known at what point abnormally high insulin levels become dangerous.

### Principle of the Procedure

IMMULITE 2000 Insulin is a solid-phase, two-site chemiluminescent immunometric assay.

**Incubation Cycles:** 1 × 60 minutes.

### Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

**EDTA tubes** should *not* be used in the IMMULITE 2000 Insulin procedure.

In collecting samples for insulin determinations, it is important to avoid hemolysis, which can lead to spuriously low results.<sup>21,22</sup>

Icteric, lipemic, or grossly contaminated samples may give erroneous results.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Insulin has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

**Volume Required:** 100 µL serum or heparinized plasma.

**Storage:** 7 days at 2–8°C, or 3 months at –20°C.<sup>20</sup>

## Warnings and Precautions

For *in vitro* diagnostic use.

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight. (See insert.)

**Water:** Use distilled or deionized water.

## Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

### Insulin Bead Pack (L2IN12)

With barcode. 200 beads, coated with monoclonal murine anti-insulin antibody. Stable at 2–8°C until expiration date. **L2KIN2:** 1 pack. **L2KIN6:** 3 packs.

### Insulin Reagent Wedge (L2INA2)

With barcode. 18.5 mL alkaline phosphatase (bovine calf intestine) conjugated to polyclonal sheep anti-insulin antibody and alkaline phosphatase (bovine calf intestine) conjugated to

monoclonal murine anti-insulin antibody in buffer, with preservative. Stable at 2–8°C until expiration date.

**L2KIN2:** 1 wedge. **L2KIN6:** 3 wedges.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

### Insulin Adjustors (LINL, LINH)

Two vials (Low and High), lyophilized insulin in a nonhuman serum matrix, with preservative. At least 30 minutes before use, reconstitute each vial by adding **4.0 mL** distilled or deionized water. Mix by *gentle* swirling or inversion. Aliquot and freeze: stable at –20°C for 60 days after reconstitution.

**L2KIN2:** 1 set. **L2KIN6:** 2 sets.

### Insulin Controls (LINC1, LINC2)

Two vials of lyophilized insulin in a nonhuman serum matrix, with preservative. At least 30 minutes before use, reconstitute each vial by adding **4.0 mL** distilled or deionized water. Mix by *gentle* swirling or inversion. Aliquot and freeze: stable at –20°C for 60 days after reconstitution.

**L2KIN2:** 1 set. **L2KIN6:** 2 sets.

Refer to the control insert for concentration levels.

Before running adjustors or controls, place the appropriate Aliquot Labels (supplied with the kit) on test tubes, so that the barcodes can be read by the on-board reader.

## Kit Components Supplied Separately

### Insulin Sample Diluent (L2INZ)

For on-board dilution of high samples. 25 mL of concentrated (ready-to-use) processed, insulin-free nonhuman serum matrix, with preservative. Storage: 30 days (after opening) at 2–8°C or 6 months (aliquotted) at –20°C.

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 x 100 mm test tube, so that the barcodes can be read by the on-board reader.

**L2INZ:** 3 labels

**L2SUBM:** Chemiluminescent Substrate  
**L2PWSM:** Probe Wash  
**L2KPM:** Probe Cleaning Kit  
**LRXT:** Reaction Tubes (disposable)  
**L2ZT:** 250 Sample Diluent Test Tubes  
(16 × 100 mm)  
**L2ZC:** 250 Sample Diluent Tube Caps

Also required  
Distilled or deionized water; test tubes.

### Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Operator's Manual.

See the IMMULITE 2000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

**Recommended Adjustment Interval:**  
4 weeks.

**Quality Control Samples:** Use Insulin Controls supplied with the kit.

### Expected Values

A study performed on 83 fasting, apparently healthy laboratory volunteers yielded a median of 9.3  $\mu\text{IU/mL}$  and a lower 95% range of up to 29.1  $\mu\text{IU/mL}$ .

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

### Limitations

Circulating anti-insulin antibodies are often found in patients who have been treated with nonhuman forms of insulin. If present, these antibodies may interfere with the assay.

For individuals who are significantly overweight, fasting insulin levels are typically somewhat higher than for adults of normal weight.

EDTA plasma has an effect on the measurement of insulin in the IMMULITE 2000 Insulin procedure.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-

33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

### Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in  $\mu\text{IU/mL}$ . (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

**Conversion Factor:**  
 $\mu\text{IU/mL} \times 1 \rightarrow \text{mIU/L}$   
 $\mu\text{IU/mL} \times 7.217 \rightarrow \text{pmol/L}$

**Calibration Range:** up to 300  $\mu\text{IU/mL}$   
(WHO NIBSC 1st IRP 66/304).

**Analytical Sensitivity:** 2  $\mu\text{IU/mL}$

**High-dose Hook Effect:** None up to 70,000  $\mu\text{IU/mL}$ .

**Precision:** Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

**Linearity:** Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

**Recovery:** Samples spiked 1 to 19 with three insulin solutions (280, 960 and 2,000  $\mu\text{IU/mL}$ ) were assayed. (See "Recovery" table for representative data.)

**Specificity:** The antibody is specific for insulin. (See "Specificity" table.)

**Alternate Sample Type:** To assess the effect of alternate sample types, blood was collected from 20 volunteers into plain glass and plastic serum tubes, plastic heparinized and EDTA tubes, sodium fluoride and plastic gel barrier tubes (SST<sup>®</sup>). All tubes were from Becton Dickinson. Eight samples were spiked with insulin to provide values throughout the calibration range of the assay. All samples were assayed by the IMMULITE 2000

Insulin procedure, with the following results. By linear regression:

(Serum Plastic) = 0.99 (Serum Glass) – 0.64 µIU/mL  
r = 0.999

(Heparin Plastic) = 1.03 (Serum Glass) + 1.23 µIU/mL  
r = 0.999

(SST Plastic) = 1.01 (Serum Glass) + 0.51 µIU/mL  
r = 0.999

(NaFI Plastic) = 0.88 (Serum Glass) – 7.1 µIU/mL  
r = 0.995

(EDTA Plastic) = 0.26 (Serum Glass) + 0.48 µIU/mL  
r = 0.977

Means:

106 µIU/mL (Serum Glass)  
104 µIU/mL (Serum Plastic)  
110 µIU/mL (Heparin Plastic)  
108 µIU/mL (SST Plastic)  
86.4 µIU/mL (NaFI Plastic)  
27.3 µIU/mL (EDTA Plastic)

EDTA plasma is not suitable for use.  
Sodium fluoride causes a depression in values.

**Bilirubin:** Severe icterus (bilirubin up to 200 mg/L) may cause a depression of values. (See "Bilirubin" table.)

**Lipemia:** Presence of triglycerides may cause a depression of values. (See "Lipemia" table.)

**Method Comparison:** The assay was compared to DPC's Coat-A-Count Insulin on 112 patient samples. (Concentration range: approximately 2 to 77 µIU/mL. See graph.) By linear regression:

(IML 2000) = 0.87 (CAC) + 1.7 µIU/mL  
r = 0.95

Means:

22.4 µIU/mL (IMMULITE 2000)  
23.6 µIU/mL (Coat-A-Count)

## References

1) Bates HM. Insulinoma and pheochromocytoma. Lab Management 1983 (May);21:11–12, 15. 2) Bennett PH. The diagnosis of diabetes; new international classification and diagnostic criteria. Ann Rev Med 1983;34:295–309. 3) Berntorp K et al. Relation between plasma insulin and blood glucose in a cross-sectional population study of the oral glucose tolerance test. Acta Endocrinol 1983;102:549–56.4) Boehm TM, Lebovitz HE. Statistical analysis of glucose and insulin responses to intravenous tolbutamide: evaluation of hypoglycemic and hyperinsulinemic states. Diabetes Care 1979;2:479–90. 5) Fajans SS, Floyd JC. Diagnosis and medical management of insulinomas. Ann Rev Med 1979;30:313–29.

6) Hirata Y et al. Autoimmunity in diabetes mellitus. Gunma Symposia Endocrinol 1980;17:81–91. 7) Kahn CR, Rosenthal AS. Immunologic reactions to insulin: insulin allergy, insulin resistance and the autoimmune insulin syndrome. Diabetes Care 1979;2:283–95. 8) Kaplan EL, Rubinstein AH et al. Calcium infusion: a new provocative test for insulinomas. Annals Surgery 1979;190:501–7. 9) Malone JI, Root AW. Plasma free insulin concentrations: keystone to effective management of diabetes mellitus in children. J Pediatr 1981;99:862–7. See also *ibid* 1983;102:800–2. 10) Marschner I et al. Group experiments on the radioimmunological insulin determination. Horm Metabol Res 1974;6:293–6. 11) National Diabetes Data Group. Classification and diagnosis of diabetes mellitus and other categories of glucose intolerance. Diabetes 1979;28:1039–57. 12) Pfeifer MA et al. Insulin secretion in diabetes mellitus. Am J Med 1981;70:579–88. 13) Reeves WG. Insulin antibody determination: theoretical and practical considerations. Diabetologia 1983;24:399–403. 14) Sachse G, Willms B. Klinische Bedeutung der Seruminsulinbestimmung: Diagnose und Prognose bei drohendem Sekundrversagen der Sulfonylharnstofftherapie. Medizin Klinik 1979;74:1635–40. 15) Schade DS, Eaton RP, Spencer W. Normalization of plasma insulin profiles in diabetic subjects with programmed insulin delivery. Diabetes Care 1980;3:9–14. 16) Service FJ et al. Insulinoma: clinical and diagnostic features of 60 consecutive cases. Mayo Clinic Proceed 1976;51:417–29. 17) Turkington RW, Estkowski A, Link M. Secretion of insulin or connecting peptide; a predictor of insulin dependence of obese diabetics. Arch Intern Med 1982;142:1102–5. 18) Turkington RW, Weindling HK. Insulin secretion in the diagnosis of adult-onset diabetes mellitus. JAMA 1978;240:833–36. See also *ibid* 1979;241:1462–3. 19) Yalow R, Bauman WA. Plasma insulin in health and disease. In: Ellenberg M, Rifkin H, editors. Diabetes mellitus; theory and practice. New York: Excerpta Medica, 1983: 119–50. 20) Burtis CA, Ashwood ER, editors. Tietz textbook of clinical chemistry. 2nd ed. Philadelphia: WB Saunders, 1994:943–4. 21) Sapin R, Ongagna JC, Gasser F, Grucker D. Insulin measurements in haemolysed serum; influence of insulinase inhibitors. Clin Chim Acta 1998;274:111–7. 22) Chevenne D, Letailleur A, Trivin F, Porquet D. Effect of hemolysis on the concentration of insulin in serum determined by RIA and IRMA. Clin Chem 1998;44:354–6.

## Technical Assistance

In the United States, contact DPC's Technical Services department.  
Tel: 800.372.1782 or 973.927.2828  
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

## Tables and Graphs

### Precision ( $\mu\text{U/mL}$ )

	Mean <sup>3</sup>	Within-Run <sup>1</sup>		Total <sup>2</sup>	
		SD <sup>4</sup>	CV <sup>5</sup>	SD	CV
1	7.67	0.42	5.5%	0.56	7.3%
2	12.5	0.50	4.0%	0.62	4.9%
3	17.2	0.57	3.3%	0.71	4.1%
4	26.4	1.04	3.9%	1.32	5.0%
5	100	3.75	3.8%	4.20	4.2%
6	291	10.8	3.7%	15.5	5.3%

### Specificity

Compound <sup>1</sup>	Amount Added <sup>2</sup> ng/mL	% Cross reactivity <sup>3</sup>
C-peptide	1,000	ND
Glucagon	1,000	ND
Pro-insulin	10	8%

ND: Not detectable<sup>4</sup>.

### Linearity ( $\mu\text{U/mL}$ )

	Dilution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	8 in 8 <sup>5</sup>	149	—	—
	4 in 8	70.4	74.5	94%
	2 in 8	38.5	37.3	103%
	1 in 8	20.0	18.6	108%
2	8 in 8	157	—	—
	4 in 8	69.0	78.5	88%
	2 in 8	37.2	39.2	95%
	1 in 8	19.9	19.6	101%
3	8 in 8	200	—	—
	4 in 8	83.4	100	83%
	2 in 8	45.6	50.0	91%
	1 in 8	25.6	25.0	102%
4	8 in 8	239	—	—
	4 in 8	102	120	85%
	2 in 8	52.2	59.8	87%
	1 in 8	30.2	29.9	101%
5	8 in 8	269	—	—
	4 in 8	124	135	92%
	2 in 8	64.0	67.3	95%
	1 in 8	37.2	33.6	110%

### Recovery (µIU/mL)

	Solution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	—	11.8	—	—
	A	25.2	25.2	100%
	B	54.5	59.2	92%
	C	99.1	111	89%
2	—	15.6	—	—
	A	28.2	28.8	98%
	B	56.1	62.8	89%
	C	94.2	115	82%
3	—	16.4	—	—
	A	28.6	29.6	97%
	B	59.3	63.6	93%
	C	101	116	87%
4	—	24.5	—	—
	A	36.8	37.3	99%
	B	64.3	71.3	90%
	C	108	123	88%
5	—	27.9	—	—
	A	39.8	40.5	98%
	B	70.5	74.5	95%
	C	112	127	89%
6	—	39.8	—	—
	A	47.5	51.6	92%
	B	75.2	85.6	88%
	C	115	138	83%


### Bilirubin

Unspiked <sup>1</sup>	Conjugated <sup>2</sup>		Unconjugated <sup>3</sup>	
	100 mg/L	200 mg/L	100 mg/L	200 mg/L
1	7.6	5.8	5.7	5.5
2	15.1	12.5	12.8	12.5
3	20.6	16.4	16.9	16.8
4	69.1	58.6	55.6	56.8
5	266	217	209	213

### Lipemia

	Triglycerides Added mg/dL	Observed	Expected	%O/E
1	—	6.22	—	—
	1,000	5.32	5.91	90%
	2,000	4.96	5.60	89%
	3,000	4.28	5.29	81%
2	—	16.5	—	—
	1,000	13.8	15.7	88%
	2,000	12.2	14.9	82%
	3,000	10.6	14.0	76%
3	—	24.7	—	—
	1,000	22.9	23.5	98%
	2,000	21.0	22.2	94%
	3,000	19.1	21.0	91%
4	—	91.0	—	—
	1,000	86.9	86.5	101%
	2,000	84.7	81.9	103%
	3,000	69.8	77.4	90%
5	—	273	—	—
	1,000	251	259	97%
	2,000	220	246	90%
	3,000	198	232	85%

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
**EURO/DPC** LTD
Diagnostic Products Corporation
EC REP DPC Biermann GmbH