 IMMULITE[®]
2000

C-Peptide

DPC[®]

IMMULITE[®] 2000 C-Peptide

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2000 Analyzer — for the quantitative measurement of C-peptide in serum, heparinized plasma, or urine, as an aid in the diagnosis and treatment of patients with abnormal insulin secretion.

Catalog Number: **L2KPE2** (200 tests),
L2KPE6 (600 tests)

Test Code: **CPE** Color: **Dark Gray**

CDC Analyte Identifier Code: 1040
CDC Test System Identifier Code: 10418
CLIA Complexity Category: Moderate

Summary and Explanation

Human C-peptide is a 31 amino acid chain with a molecular mass of approximately 3,020 daltons. Metabolically inert, it originates in the pancreatic β -cells as a by-product of the enzymatic cleavage of proinsulin to insulin.^{1,2,5} In this process, insulin and C-peptide are split from the prohormone and secreted into the portal circulation in equimolar concentrations.^{4,5,7} It is this fact which underlies the clinical interest in plasma determinations of C-peptide.

Within limits, C-peptide levels can serve as a valuable index to insulin secretion. Thus, low C-peptide levels are to be expected where insulin secretion is diminished, as in insulin-dependent diabetes, or suppressed, as a normal response to exogenous insulin; whereas elevated C-peptide levels may result from the increased β -cell activity observed in insulinomas.^{3,4,6,9}

Accordingly, in the differential diagnosis of hypoglycemia, C-peptide determinations can be used to supplement insulin measurements as an index to pancreatic activity in the classic 72-hour fasting test, and as the sole indicator of pancreatic activity where insulin itself is administered to check for suppressibility.^{1,8} In addition, covert self-administration of insulin can be virtually ruled out as the cause of hyperinsulinemia by the finding of an elevated C-peptide level.^{2,3,8,9}

Circulating anti-insulin antibodies are commonly encountered in patients who

have undergone insulin therapy. These would typically interfere with immunoassays for insulin, making it impossible to use insulin measurements in this context to check on residual β -cell activity, even if treatment were temporarily suspended. C-peptide measurements have therefore been used as an alternative in this context, to yield information on the natural history of insulin-dependent diabetes, to indirectly monitor insulin secretion in the presence of anti-insulin antibodies, and to help settle on an appropriate course of treatment.^{3,6,7,10}

C-peptide has also been measured as an additional means for evaluating glucose tolerance and glibenclamide-glucose tests.^{2,3,10}

Principle of the Procedure

IMMULITE 2000 C-Peptide is a solid-phase, competitive chemiluminescent enzyme immunoassay.

Incubation Cycles: 1 × 60 minutes.

Specimen Collection

Serum and Heparinized Plasma:

The patient should be fasting. Collect blood by venipuncture,¹² avoiding hemolysis, into plain tubes (without anticoagulant) or heparinized tubes, noting the time of collection, and separate the serum or plasma from the cells.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

EDTA plasma and sodium fluoride plasma are unsuitable for use.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 C-Peptide has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Storage: Assay within 2-3 hours, or store frozen at -20°C for 1 week.¹¹

Urine – Collection and Storage:

Collect a 24-hour urine, without preservative, keeping the specimen refrigerated at $2-8^{\circ}\text{C}$ during collection. Record the total volume of the collection and retain a well-mixed aliquot for analysis. Before assay, clear the sample by centrifugation or filtration.

Storage: For longer storage, aliquot and freeze: stable at -20°C for 30 days.

Dilution Factor: At least 5. Using a dilution factor of 20 will bring *normal* urine samples within the assay's calibration range. (For urine samples, select 5 or 20 in the Dilution Factor window.)

Volume Required:

75 μL serum, plasma, or urine sample.

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at $2-8^{\circ}\text{C}$. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

C-Peptide Bead Pack (L2PE12)

With barcode. 200 beads, coated with polyclonal rabbit anti-C-peptide. Stable at $2-8^{\circ}\text{C}$ until expiration date.

L2KPE2: 1 pack **L2KPE6:** 3 packs.

C-Peptide Reagent Wedge (L2PEA2)

With barcode. 17.5 mL alkaline phosphatase (bovine calf intestine) conjugated to C-peptide in buffer. Stable at $2-8^{\circ}\text{C}$ until expiration date.

L2KPE2: 1 wedge. **L2KPE6:** 3 wedges

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

C-Peptide Adjustors (LPEL, LPEH)

Two vials (Low and High) of lyophilized C-peptide in buffered human albumin. Reconstitute each vial with **4.0 mL** distilled or deionized water. Let stand for 30 minutes. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. After reconstitution, aliquot and freeze. Stable at -20°C for 6 months. Discard aliquots after use.

L2KPE2: 1 set. **L2KPE6:** 2 sets.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

C-Peptide Sample Diluent (L2PEZ)

For the on-board dilution of urine samples and high serum/plasma samples. 25 mL of concentrated (ready-to-use), processed, C-peptide-free buffered human albumin, with preservative. Stable at $2-8^{\circ}\text{C}$ for 30 days after opening, or for 6 months (aliquotted) at -20°C .

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16×100 mm test tube, so that the barcode can be read by the on-board reader.

L2PEZ: 3 labels.

L2SUBM: Chemiluminescent Substrate
L2PWSM: Probe Wash
L2KPM: Probe Cleaning Kit
LRXT: Reaction Tubes (disposable)
L2ZT: 250 Sample Diluent Test Tubes
(16 × 100 mm)
L2ZC: 250 Sample Diluent Tube Caps
PECM: Tri-level C-Peptide Control Module
Also Required
Distilled or deionized water; test tubes;
controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Operator's Manual.

See the IMMULITE 2000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:
2 weeks.

Quality Control Samples: Use controls or serum pools with at least two levels (low and high) of C-peptide.

Expected Values

Serum and Heparinized Plasma:
Serum samples were collected from 71 fasting laboratory volunteers and analyzed by the IMMULITE 2000 C-Peptide procedure, with the median of 2.1 ng/mL (0.7 nmol/L; 695 pmol/L) and a nonparametric central 95% reference range of

1.1 – 5.0 ng/mL
(0.4 – 1.7 nmol/L; 364 – 1,655 pmol/L)

Urine: Mean ± SD of 79 ± 56 µg/day, with a range of 2 to 260 µg/day, representing the central 95% of the observations.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

Because metabolism of C-peptide differs from that of insulin, C-peptide levels are at best a semi-quantitative index of insulin secretion. The half-life of C-peptide in plasma has been estimated as approximately 30 minutes, compared to

approximately 5 minutes for insulin. Because of the difference in half-life, C-peptide circulates in plasma at a level roughly five times that of insulin, even though the two molecules are secreted in an equimolar ratio. Again, the liver plays a major role in clearing insulin, whereas C-peptide is removed by degradation and elimination mainly through the kidneys. Hepatic and renal complications will therefore affect the circulating C-peptide/insulin ratio.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Conversion Factors:
ng/mL × 0.331 → nmol/L
ng/mL × 331 → pmol/L

Calibration Range: 0.5 – 7 ng/mL
(0.17 – 2.3 nmol/L; 166 – 2,317 pmol/L).
Standardized to WHO 1st IRP 84/510.

Analytical Sensitivity: 0.3 ng/mL
(0.10 nmol/L; 99 pmol/L).

Precision: Serum samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Serum samples spiked 1 to 19 with three C-peptide solutions (15, 54 and 112 ng/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for C-peptide. (See "Specificity" table.)

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 381 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of lipemia in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: No significant effect for heparinized plasma. EDTA plasma and sodium fluoride plasma are unsuitable for use. By linear regression:

(Heparin) = 0.97 (Serum) + 0.51 ng/mL
r = 0.969

Means:
4.7 ng/mL (Serum)
5.1 ng/mL (Heparin)

Method Comparison: The assay was compared to DPC's IMMULITE C-Peptide procedure on 66 samples. (Concentration range: approximately 0.5 – 6.5 ng/mL. See graph.) By linear regression:

(IML 2000) = 0.97 (IML) + 0.16 ng/mL
r = 0.990

Means:
2.84 ng/mL (IMMULITE 2000)
2.76 ng/mL (IMMULITE)

References

1) Beischer W. Proinsulin and C-peptide in humans. In: Fotherby K, Pal S, editors. Hormones in Normal and Abnormal Human Tissues. Vol 3. Berlin: Walter DeGruyter, 1983: 1-43. 2) Beyer J, Krause U, Cordes U. C-peptide: its biogenesis, structure, determination and clinical significance. *Giornale Ital Chem Clin* 1979;4 Suppl 1:9-22. 3) Bonser A, Garcia-Webb P. C-Peptide measurement: methods and clinical utility. *CRC Crit Rev Clin Lab Sci* 1984;19:297-352. 4) Gonen G, Rubenstein AH, Horwitz DL, Blix PM. Clinical significance of C-peptide. In: Baba S, Kaneko T, Yanaihara N, editors. Proinsulin, Insulin, C-Peptide. Amsterdam: Excerpta Medica, 1979: 246-53. 5) Horwitz D, et al. Proinsulin, insulin and C-peptide concentrations in human portal and

peripheral blood. *J Clin Invest* 1975;55:1278-83. 6) Horwitz D, Kuzuya H, Rubenstein AH. Circulating serum C-peptide. *N Engl J Med* 1976;295:207-9. 7) Polonsky KS, Rubenstein AH. C-peptide as a measure of the secretion and hepatic extraction of insulin; pitfalls and limitations. *Diabetes* 1984;33:486-93. 8) Rendell M. Expanding clinical use of C-peptide radioimmunoassay. *Acta Diabetol Lat* 1983;20:105-13. 9) Rubenstein AH, Kuzuya H, Horwitz DL. Clinical significance of circulating C-peptide in diabetes mellitus and hypoglycemic disorders. *Arch Intern Med* 1977;137:625-32. 10) Turkington RW, Estkowski A, Link M. Secretion of insulin or connecting peptide; a predictor of insulin dependence of obese diabetics. *Arch Intern Med* 1982;142:1102-5. 11) Data on file. 12) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782 or 973.927.2828
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

Manufactured by EURO/DPC Ltd. under a Quality System registered to ISO 13485:2003.

Tables and Graphs

Precision (ng/mL)

	Within-Run ¹			Total ²	
	Mean ³	SD ⁴	CV ⁵	SD	CV
1	0.74	0.105	14.1%	0.137	18.6%
2	1.49	0.093	6.3%	0.131	8.8%
3	3.02	0.210	7.0%	0.252	8.4%
4	6.45	0.319	4.9%	0.405	6.3%

Specificity

Compound ¹	Amount Added ²	% Cross-reactivity ³
Insulin	200 µIU/mL	ND
Glucagon	15,000 ng/mL	ND
Secretin	15,000 ng/mL	ND
Proinsulin	10 ng/mL	17%

ND: not detectable.⁴


Linearity (ng/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	4 in 4 ⁵	1.8	—	—
	2 in 4	0.82	0.9	91%
	1 in 4	0.41	0.45	91%
2	8 in 8	3.7	—	—
	4 in 8	1.8	1.9	95%
	2 in 8	0.87	0.93	94%
	1 in 8	0.39	0.46	85%
3	8 in 8	4.9	—	—
	4 in 8	2.6	2.5	104%
	2 in 8	1.2	1.2	100%
	1 in 8	0.65	0.61	107%
4	8 in 8	6.2	—	—
	4 in 8	3.2	3.1	103%
	2 in 8	1.6	1.6	100%
	1 in 8	0.69	0.78	88%

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
EURO/DPC LTD
Diagnostic Products Corporation
EC REP DPC Biemann GmbH