 IMMULITE[®]
2000

High Sensitivity CRP

DPC[®]

IMMULITE® 2000 High Sensitivity CRP

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2000 Analyzer — for the quantitative measurement of C-reactive protein (CRP) in serum or plasma, as an aid in the detection and evaluation of the infection, tissue injury, inflammatory disorders and associated diseases.

Catalog Number: **L2KCR2** (200 tests), **L2KCR6** (600 tests)

Test Code: **CRP** Color: **Dark Blue**

CDC Analyte Identifier Code: 1001

CDC Test System Identifier Code: 10418

CLIA Complexity Category: Moderate

Summary and Explanation

CRP is an acute phase β -globulin with a molecular mass of approximately 118,000 daltons.¹ CRP is highly conserved, composed of five identical cyclic globular subunits, and is classified as a member of the pentraxin superfamily of proteins.^{1,2}

CRP levels in serum or plasma may rise during general, nonspecific response to infectious and non-infectious inflammatory conditions such as rheumatoid arthritis, cardiovascular disease and peripheral vascular disease.¹⁷ CRP is synthesized in the liver and is normally present as a trace constituent of serum or plasma.

While the precise *in vivo* functions of CRP during the inflammation state are not known, there is considerable evidence indicating a role in the recognition and elimination of foreign pathogens, as well as potentially toxic endogenous substances released from damaged tissue, by assisting humoral and cellular immunity.^{1,3} In normal individuals, hepatic cells produce CRP at constitutive levels of less than 1 mg/dL.³

Within 24 to 48 hours following acute tissue damage, production dramatically rises to approximately 1,000 times the constitutive level. CRP levels may remain elevated for several days before returning to normal.^{1,3,4}

As elevated CRP values are always associated with pathological changes, the CRP assay provides useful information for the diagnosis, therapy, and monitoring of inflammatory conditions and associated diseases.⁷⁻⁹ Additionally, measurement of CRP by high-sensitivity CRP assays may add to the predictive value of other markers used to assess the risk of cardiovascular and peripheral vascular diseases.¹¹⁻¹⁶

Principle of the Procedure

IMMULITE 2000 High Sensitivity CRP is a solid-phase, chemiluminescent immunometric assay.

Incubation Cycles: 1 × 30 minutes.

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 High Sensitivity CRP has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 5 μ L serum or plasma (heparinized or EDTA).

Automatic Predilution Factor: 100

Storage: 3 days at 2–8°C, or 2 months at –20°C.

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

CRP Bead Pack (L2CR12)

With barcode. 200 beads, coated with anti-ligand. Stable at 2–8°C until expiration date.

L2KCR2: 1 pack. **L2KCR6:** 3 packs.

CRP Reagent Wedge (L2CRA2)

With barcode. 11.5 mL ligand-labeled anti-CRP murine monoclonal antibody and alkaline phosphatase (bovine calf intestine) conjugated to rabbit polyclonal anti-CRP antibody in buffer. Stable at 2–8°C until expiration date.

L2KCR2: 1 wedge. **L2KCR6:** 3 wedges.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

CRP Adjustors (LCRL, LCRH)

Two vials (Low and High), 2.0 mL each, of CRP in a protein/buffer matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

L2KCR2: 1 set. **L2KCR6:** 2 sets.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

CRP Sample Diluent (L2CRZ)

For the on-board dilution of patient samples. 50 mL of a concentrated (ready-to-use) CRP-free protein/buffer matrix. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

L2KCR2: 1 vial. **L2KCR6:** 3 vials.

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2KCR2: 3 labels. **L2KCR6:** 5 labels.

Kit Components Supplied Separately

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

LCRCM: Tri-level CRP Control Module

Also Required

Distilled or deionized water; test tubes; controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Operator's Manual.

See the IMMULITE 2000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval: 2 weeks.

Quality Control Samples: Use controls or pools with at least two levels (low and high) of CRP.

Expected Values

Based on its relationship to DPC's IMMULITE High Sensitivity CRP (see Method Comparison), the assay can be expected to have essentially the same reference ranges.

A study performed on 100 apparently healthy laboratory volunteers yielded a median of 0.14 mg/dL and an upper 97.5th percentile of 1.1 mg/dL.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

For diagnostic purposes, CRP results should be used in conjunction with other test results, the overall clinical presentation to the physician, and all other appropriate information.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in mg/dL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Conversion Factor: mg/dL \times 10 \rightarrow mg/L

Calibration Range: Up to 15 mg/dL (150 mg/L) WHO 1st IS 85/506.

Analytical Sensitivity: 0.01 mg/dL (0.1 mg/L).

Functional Sensitivity: <0.02 mg/dL (<0.2 mg/L).

High Dose Hook Effect: None up to 400 mg/dL.

Precision: Samples were repeatedly assayed in quadruplicate over the course

of several days, for a total of 20 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with 3 CRP solutions (60, 120, and 225 mg/dL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for CRP. (See "Specificity" table.)

Bilirubin: Presence of conjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 512 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 19 volunteers into plain, heparinized, EDTA and Becton Dickinson SST[®] vacutainer tubes. Equal volumes of the matched samples were spiked with various concentrations of CRP, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE 2000 High Sensitivity CRP procedure.

(Heparin) = 1.12 (Serum) – 0.09 mg/dL
r = 0.995

(EDTA) = 1.13 (Serum) – 0.14 mg/dL
r = 0.996

(SST) = 1.11 (Plain Tubes) – 0.07 mg/dL
r = 0.995

Means:

6.41 mg/dL (Serum)
7.05 mg/dL (Heparin)
7.09 mg/dL (EDTA)
6.95 mg/dL (SST)

Method Comparison 1: The assay was compared to DPC's IMMULITE High Sensitivity CRP on 435 samples. (Concentration range: up to approximately 15 mg/dL. See graph.) By linear regression:

(IML 2000) = 1.00 (IML) – 0.003 mg/dL
r = 0.993

Means:
0.99 mg/dL (IMMULITE 2000)
0.99 mg/dL (IMMULITE)

Method Comparison 2: A second method comparison study was performed using another commercially available method (Kit A) on 199 samples. (Concentration range: up to approximately 2.5 mg/dL. See graph.) By linear regression:

$$(IML\ 2000) = 1.01 (\text{Kit A}) - 0.01\ \text{mg/dL}$$

$$r = 0.994$$

Means:
0.23 mg/dL (IMMULITE 2000)
0.24 mg/dL (Kit A)

References

1) Schultz DR, Arnold PI. Properties of four acute phase proteins: CRP, serum amyloid A protein, β -acid glycoprotein, and fibrinogen. *Semin Arthritis Rheum* 1990;20:129-47. 2) Osmand AP, Friedenson B, Gewurz H, et al. Characterization of CRP and the complement subcomponent C1t as homologous proteins displaying cyclic pentameric symmetry (pentraxins). *Proc Nat Acad Sci USA* 1977;74:739-43. 3) Tietz NW, editor. *Clinical guide to laboratory tests*. 3rd ed. Philadelphia: WB Saunders, 1995:358. 4) Kushner I. The phenomenon of acute phase response. *Ann N Y Acad Sci* 1982;389:39-48. 5) National Committee for Clinical Laboratory Standards. *Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard*. 4th ed. NCCLS Document H3-A4. Wayne, PA: NCCLS, 1998. 6) Kushner I, Rzewnicki DL. The acute phase response: general aspects. *Baillieres Clin Rheumatol* 1994;8:513-30. 7) Dixon JS, Bird HA, Sitton NG, et al. CRP in the serial assessment of disease activity in rheumatoid arthritis. *Scand J Rheumatol* 1984;13:39-44. 8) Hind CR, Pepys MB. The role of CRP measurement in clinical practice. *Int Med* 1984;5:112-51. 9) van Leeuwen MA, van Rijswijk MH. Acute phase proteins in the monitoring of inflammatory disorders. *Baillieres Clin Rheumatol* 1994;8:531-52. 10) Ridker PM, Cushman M, Stampfer MJ, et al. Plasma concentration of CRP and risk of developing peripheral vascular disease. *Circulation* 1998;97:425-8. 11) Ridker PM, Buring JE, Shih J, et al. Prospective study of CRP and the risk of future cardiovascular events among apparently healthy women. *Circulation* 1998;98:731-3. 12) Ridker PM, Glynn RJ, Hennekens CH. CRP adds to the predictive value of total and HDL cholesterol in determining risk of first myocardial infarction. *Circulation* 1998;97:2007-11. 13) Ridker PM, Rifai N, Pfeffer MA, et al. Inflammation, pravastatin, and the risk of coronary events after myocardial infarction in patients with average cholesterol levels. *Cholesterol and Recurrent Events (CARE) Investigators. Circulation* 1998;98:839-44. 14)

Tracy RP, Lemaitre RN, Psaty BM, et al. Relationship of CRP to risk of cardiovascular disease in the elderly. Results from the Cardiovascular Health Study and the Rural Health Promotion Project. *Arterioscler Thromb Vasc Biol* 1997;17:1121-7. 15) Macy EM, Hayes TE, Tracy RP. Variability in the measurement of CRP in healthy subjects: implications for reference intervals and epidemiological applications. *Clin Chem* 1997;43:52-8. 16) Ridker PM, Cushman M, Stampfer MJ, et al. Inflammation, aspirin, and the risk of cardiovascular disease in apparently healthy men. *N Engl J Med* 1997;336:973-9. Erratum: *ibid* 1997;337:356. 17) Du Clos TW. Function of CRP. *Ann Med* 2000; 32:274-8.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782 or 973.927.2828
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

The Quality System of Diagnostic Products Corporation is registered to ISO 13485:2003.

Tables and Graphs

Precision (mg/dL)

	Within-Run ¹			Total ²	
	Mean ³	SD ⁴	CV ⁵	SD	CV
1	0.023	0.002	8.7%	0.002	8.7%
2	0.085	0.004	4.7%	0.006	7.1%
3	0.32	0.009	2.8%	0.010	3.1%
4	1.23	0.035	2.8%	0.040	3.3%
5	2.18	0.074	3.4%	0.082	3.8%
6	9.37	0.48	5.1%	0.49	5.2%

Specificity

Compound ¹	mg/dL Added ²	% Cross reactivity ³
Human Serum Albumin	1,000	ND
Human IgG	1,000	ND
Transferrin	1,000	ND

ND: Not detectable.⁴


Linearity (mg/dL)


	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	0.87	—	—
	4 in 8	0.41	0.44	93%
	2 in 8	0.20	0.22	91%
	1 in 8	0.10	0.11	91%
2	8 in 8	2.31	—	—
	4 in 8	1.15	1.16	99%
	2 in 8	0.59	0.58	102%
	1 in 8	0.30	0.29	103%
3	8 in 8	7.22	—	—
	4 in 8	3.67	3.61	102%
	2 in 8	1.76	1.81	97%
	1 in 8	0.90	0.90	100%
4	8 in 8	7.86	—	—
	4 in 8	3.95	3.93	101%
	2 in 8	1.82	1.97	92%
	1 in 8	1.00	0.98	102%
5	8 in 8	13.2	—	—
	4 in 8	6.56	6.60	99%
	2 in 8	3.24	3.30	98%
	1 in 8	1.56	1.65	95%

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O Sistema de Qualidade da Diagnostic Products Corporation está registado sob ISO 13485:2003.
EC REP DPC Biermann GmbH