 IMMULITE[®]
2000

Carbamazepine

DPC[®]

IMMULITE[®] 2000 Carbamazepine

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2000 Analyzer — for the quantitative measurement of carbamazepine in serum or heparinized plasma, as an aid in monitoring the therapeutic administration of this drug.

Catalog Number: **L2KCB2** (200 tests),
L2KCB6 (600 tests)
Test Code: **CMP** Color: **Light Gray**

Summary and Explanation

Carbamazepine is an iminostilbene derivative used for the treatment of epilepsy; trigeminal neuralgia; and simple and complex, partial and generalized convulsive seizures.¹⁻³ It is administered either alone or in combination with other antiepileptic drugs.

Carbamazepine is rapidly absorbed into the bloodstream, where it is highly protein bound (60 to 80 percent).³ The drug is metabolized by hepatic oxidation to the 10,11-epoxide form, which is as active as the parent drug. The epoxide form is further metabolized to the 10,11-dihydroxide, which is eliminated in the urine.^{1,4} The presence of other drugs can significantly alter the rate of epoxide formation: phenytoin, phenobarbital and felbamate induce the metabolic enzyme activity, whereas erythromycin and propoxyphene inhibit this activity.¹⁻⁴

The therapeutic range for carbamazepine is between 4 and 12 µg/mL (between 17 and 51 µmol/L), with toxic levels above 15 µg/mL.³ However, individual differences in absorption and metabolism can result in highly variable responses to carbamazepine. Monitoring levels assists the physician in optimizing dosage and minimizing toxic side effects for each patient.

Principle of the Procedure

IMMULITE 2000 Carbamazepine is a solid-phase, competitive chemiluminescent enzyme immunoassay.

Incubation Cycles: 1 × 30 minutes.

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Carbamazepine has not been tested with all possible variations of tube types.

Volume Required: 20 µL serum or plasma.

Storage: 2 days at 2–8°C, or 2 months at –20°C.

Note: The time of collection will depend on whether the desired measurement is for peak or trough values.⁵

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of

potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

Carbamazepine Bead Pack (L2CB12)

With barcode. 200 beads, coated with monoclonal murine anti-carbamazepine antibody. Stable at 2–8°C until expiration date.

L2KCB2: 1 pack. **L2KCB6:** 3 packs.

Carbamazepine Reagent Wedge (L2CBA2)

With barcode. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to carbamazepine in buffer, with preservative. Stable at 2–8°C until expiration date.

L2KCB2: 1 wedge. **L2KCB6:** 3 wedges.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

Carbamazepine Adjustors (LCBL, LCBH)

Two vials (Low and High), 2.0 mL each, of carbamazepine in processed human serum, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

L2KCB2: 1 set. **L2KCB6:** 2 sets.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

SDCM: Bi-level, human serum-based control containing carbamazepine

Also Required

Distilled or deionized water; test tubes; controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Operator's Manual.

See the IMMULITE 2000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval: 2 weeks.

Quality Control Samples: Use controls or serum pools with at least two levels (low and high) of carbamazepine.

Therapeutic Range

The following results are reported in the literature, with the therapeutic range for carbamazepine between 4 and 12 µg/mL (17 and 51 µmol/L) and toxic levels above 15 µg/mL (63 µmol/L) for trough level in adults and children.³

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

A complete understanding of the information contained in this package insert is necessary for proper use of this product.

The specificity table should be consulted for a list of substances which have measurable crossreactivity in this assay. Other medications and/or food not found on this list may nevertheless crossreact in the assay.

The following should be considered in the interpretation of the carbamazepine concentration:

- All concurrent medications taken by the patient
- The time of sampling in relation to the dosing history
- Existing diseases or conditions

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing

interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in µg/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Conversion Factor:
µg/mL × 4.23 → µmol/L

Calibration Range: 1.25 – 20 µg/mL
(5.29 – 85 µmol/L).

Analytical Sensitivity: 0.2 µg/mL
(0.9 µmol/L).

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three carbamazepine solutions (24, 85 and 170 µg/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for carbamazepine. (See "Specificity" table.)

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 512 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: Samples (*n* = 20) were collected into plain and heparinized vacutainer tubes. All samples were spiked with carbamazepine and then assayed by the IMMULITE 2000 Carbamazepine procedure, with the following results.

(Heparin) = 0.99 (Serum) – 0.08 µg/mL
r = 0.999

Means:
7.40 µg/mL (Serum)
7.24 µg/mL (Heparin)

Method Comparison 1: The assay was compared to DPC's IMMULITE Carbamazepine on 97 samples. (Concentration range: approximately 1.3 to 15 µg/mL. See graph 1.) By linear regression:

(IML 2000) = 0.99 (IML) – 0.07 µg/mL
r = 0.990

Means:
7.6 µg/mL (IMMULITE 2000)
7.8 µg/mL (IMMULITE)

Method Comparison 2: The assay was compared to Kit A on 45 samples. (Concentration range: approximately 2 to 12 µg/mL. See graph 2.) By linear regression:

(IML 2000) = 1.08 (Kit A) – 0.05 µg/mL
r = 0.976

Means:
8.2 µg/mL (IMMULITE 2000)
7.6 µg/mL (Kit A)

Method Comparison 3: The assay was compared to Kit B on 71 samples. (Concentration range: approximately 2 to 14 µg/mL. See graph 3.) By linear regression:

(IML 2000) = 1.01 (Kit B) – 0.58 µg/mL
r = 0.970

Means:
8.0 µg/mL (IMMULITE 2000)
8.5 µg/mL (Kit B)

References

- 1) Drug monitoring data pocket guide II. 2nd ed. Washington, DC: AACC Press, 1994: 49-51.
- 2) McNamara JO. Drugs effective in the treatment of the epilepsies. In: GL Hardman, LL Limbird, PB Molinoff, RW Ruddon, editors. Goodman & Gilman's the pharmacological basis of

therapeutics. 9th ed. New York: McGraw-Hill, 1996: 461-86. 3) Moyer TP, Pippenger CE. Therapeutic drug monitoring. In: Burtis CA, Ashwood ER, editors. Tietz textbook of clinical chemistry. 2nd ed. Philadelphia: WB Saunders, 1994: 1094-154. 4) Kutt H. Carbamazepine: chemistry and methods of determination. Complex partial seizures and their treatment. Adv Neurol 1975;11:249-61. 5) Jacobs DS, Kaster BL Jr, Demott WR, Wolfson WL. Laboratory test handbook. Stowe, OH: Lexi-Compl. Mosby, 1990: 782.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782 or 973.927.2828
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

Manufactured by EURO/DPC Ltd. under a Quality System registered to ISO 13485:2003.

Tables and Graphs

Precision (µg/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV% ⁵	SD	CV
1	4.3	0.15	3.5%	0.26	6.0%
2	8.8	0.29	3.3%	0.54	6.1%
3	10	0.32	3.2%	0.59	5.9%
4	16	0.56	3.5%	1.04	6.5%
5	20	0.89	4.5%	1.16	5.8%

Linearity (µg/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	8.5	—	—
	4 in 8	4.0	4.3	93%
	2 in 8	2.0	2.1	95%
	1 in 8	<1.25	1.1	—
2	8 in 8	11	—	—
	4 in 8	5.3	5.5	96%
	2 in 8	2.5	2.8	89%
	1 in 8	1.3	1.4	93%
3	8 in 8	13	—	—
	4 in 8	6.0	6.5	92%
	2 in 8	3.2	3.3	97%
	1 in 8	1.6	1.6	100%
4	8 in 8	14	—	—
	4 in 8	6.8	7.0	97%
	2 in 8	3.3	3.5	94%
	1 in 8	1.6	1.8	89%
5	8 in 8	16	—	—
	4 in 8	7.9	8.0	99%
	2 in 8	3.8	4.0	95%
	1 in 8	1.9	2.0	95%

Recovery (µg/mL)


	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	ND	—	—
	A	1.2	1.2	100%
	B	4.5	4.3	105%
	C	8.8	8.5	104%
2	—	1.4	—	—
	A	2.4	2.5	96%
	B	5.7	5.6	102%
	C	10.4	9.8	106%
3	—	3.3	—	—
	A	4.2	4.4	95%
	B	7.6	7.4	103%
	C	11.9	11.7	102%
4	—	5.4	—	—
	A	6.3	6.3	100%
	B	9.4	9.4	100%
	C	14.3	13.6	105%
5	—	7.2	—	—
	A	7.9	8.1	98%
	B	11.0	11.1	99%
	C	15.3	15.4	99%

Specificity



Compound ¹	µg/mL Added ²	% Cross reactivity ³
Carbamazepine	20	99.5%
Carbamazepine-10,11-epoxide	20	8.8%
Clonazepam	0.12	ND
Diazepam	2.5	ND
Ethosuximide	100	ND
Imipramine	1.0	ND
Phenobarbital	40	0.1%
Phenylbutazone	45	ND
Phenytoin	20	ND
5-(4-Hydroxyphenyl)-5-phenyl-hydantoin	20	ND
Secobarbital	15	ND
Theophylline	20	ND
Valproic Acid	100	ND

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
EURO/DPC LTD
Diagnostic Products Corporation
EC REP DPC Biermann GmbH