 IMMULITE[®]
2000

Theophylline

DPC[®]

IMMULITE[®] 2000 Theophylline

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2000 Analyzer — for the quantitative measurement of theophylline in serum, as an aid in monitoring the therapeutic administration of this bronchodilator.

Catalog Numbers: **L2KTN2** (200 tests), **L2KTN6** (600 tests)

Test Code: **THE** Color: **Light Green**

CDC Analyte Identifier Code: 6104

CDC Test System Identifier Code: 10418

CLIA Complexity Category: Moderate

Summary and Explanation

Theophylline is a methylxanthine derivative which is widely used clinically for its bronchodilator activity in the treatment of respiratory conditions such as asthma. As a bronchodilator, theophylline appears to exert its effect on the smooth muscle of the bronchi. Theophylline is a competitive inhibitor of the enzyme cyclic nucleotide phosphodiesterase, an enzyme that catalyzes the conversion of cyclic adenosine monophosphate (cAMP) to 5'-adenosine monophosphate (5'-AMP).¹ Mechanisms other than phosphodiesterase inhibition have been suggested to mediate the bronchodilator effect of theophylline.² Other methylxanthines which have measurable crossreactivity — such as caffeine, theobromine, and 8-chlorotheophylline — occur in coffee, tea, soft drinks, chocolate, and certain over-the-counter medications such as dimenhydrinate.

The importance of measuring theophylline concentrations is related to the relatively narrow therapeutic index for theophylline as well as the inter- and intra-patient variability in the disposition of theophylline. There exists a window of concentrations below which theophylline levels are subtherapeutic, and above which there is an increased incidence of toxic side effects, such as gastrointestinal disturbances, tachycardia, cardiac arrhythmias, convulsions and coma.^{3,4,5}

Patients with the same theophylline levels will not necessarily experience the same therapeutic or toxic effects. Variability in the disposition of theophylline results in widely varying serum concentrations in patients receiving the same dose.^{6,7} Theophylline is excreted primarily via metabolism in the liver as 1-methyluric acid, 1,3-dimethyluric acid, and methylxanthines.⁸ Factors which affect liver metabolism — such as smoking, extremes of age, concomitant drug therapy, alcoholism, and diet — can lead to either subtherapeutic or toxic theophylline concentrations in a patient receiving a dosage which would be effective in the absence of these conditions.

Another source of variation results from the specific dosage formulation employed and the route of administration.⁹ This variation exists in part due to the low solubility characteristics of theophylline as well as to the wide variety of sustained-release dosage formulations which are used. Intravenously, theophylline is administered as aminophylline, which is a complex of two theophylline molecules with one mole of ethylenediamine to enhance solubility.

Therapeutic levels are generally regarded to be between 10 and 20 µg/mL,¹⁰ although levels between 5 and 10 µg/mL may still be effective in certain individuals.¹¹ Toxic levels are generally considered to be greater than 20 µg/mL.

Circulating theophylline levels generally reach a plateau or steady state upon administration of a fixed daily dosage. A serum theophylline level should not be assumed to reflect the steady-state concentration unless a continuous fixed dosage has been administered for at least 3 to 5 half-lives for theophylline.¹⁵ Although the average *in vivo* half-life for theophylline in adults is 4.5 hours, it varies widely in this population (3 to 9.5 hours) and is significantly shorter in children (mean 3.6 hours, range 1.5 to 9.5 hours).¹⁰ The following diseases and conditions have been implicated in altering the pharmacokinetics of theophylline: congestive heart failure, respiratory viral infections, liver disease, and renal

disease.¹⁶ Concomitant use or discontinuation of the following drugs may also affect theophylline pharmacokinetics: cimetidine, barbiturates, erythromycin, troleandomycin, phenytoin and allopurinol.^{16,17} Consult the manufacturer's package insert for the theophylline drug product and the potentially interfering drug product for further information.

Principle of the Procedure

IMMULITE 2000 Theophylline is a solid-phase, enzyme-labeled, competitive chemiluminescent immunoassay.

Incubation Cycles: 1 × 30 minutes.

Specimen Collection

The time of collection is important. In order to obtain theophylline samples which accurately reflect pharmacological levels in the patient, guidelines should be followed with respect to sampling times and the interpretation of results. Samples should not be taken immediately following dosage, during absorption and distribution of the drug. Following oral dosage with solid formulation, samples should be drawn 4 hours post-dose to allow for complete absorption and distribution.¹² Following intravenous administration or oral administration with a liquid formulation, samples may be drawn 1 hour post-dose.¹³ In premature infants given rectal administration of aminophylline, samples may be drawn two hours post-dose.¹⁴ Consistent sampling times with respect to dosing appropriate for each clinical setting should be employed.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Theophylline has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 50 µL serum.

Storage: 7 days at 2–8°C or 2 months at –20°C.

SST Tubes: Manufacturer's studies have shown that Becton Dickinson SST[®] vacutainer tubes may cause elevations in values for therapeutic drug assays.

Dilution of High Samples: All samples expected to have levels above the assay's calibration range should be diluted before assay with an untreated patient sample.

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

Theophylline Bead Pack (L2TN12)

With barcode. 200 beads, coated with polyclonal rabbit anti-theophylline antibody. Stable at 2–8°C until expiration.

date.

L2KTN2: 1 pack. **L2KTN6:** 3 packs.

Theophylline Reagent Wedge (L2TNA2)

With barcode. 11.5 mL of alkaline phosphatase (bovine calf intestine) conjugated to theophylline, in buffer.

Stable at 2–8°C until expiration date.

L2KTN2: 1 wedge. **L2KTN6:** 3 wedges.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

Theophylline Adjustors (LTNL, LTNH)

Two vials (Low and High), 4 mL each, of theophylline in processed human serum, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

L2KTN2: 1 set. **L2KTN6:** 2 sets.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components

Supplied Separately

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

SDCM: Bi-level, multi-constituent control

Also Required

Distilled or deionized water; test tubes; controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Operator's Manual.

See the IMMULITE 2000 Operator's Manual for: preparation, setup, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:
2 weeks.

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of theophylline.

Therapeutic Range

As is indicated in the Summary and Explanation section, therapeutic levels are generally regarded to be between 10 and 20 µg/mL,¹⁰ although levels between 5 and 10 µg/mL may still be effective in certain individuals.¹¹ Toxic levels are generally considered to be greater than 20 µg/mL.

Limitations

A complete understanding of the information contained in this package insert is necessary for proper use of this product.

The specificity table should be consulted for a list of substances which have measurable crossreactivity in this assay. Other medications and/or food not found on this list may nevertheless crossreact in the assay.

The following should be considered in the interpretation of the theophylline concentration:

- All concurrent medications taken by the patient
- The time of sampling in relation to the dosing history
- Existing diseases or conditions

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in µg/mL. (Unless

otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Calibration Range: 1 – 60 µg/mL

Analytical Sensitivity: 0.3 µg/mL

Precision: Samples were processed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three theophylline solutions (125, 250 and 500 µg/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for theophylline. (See "Specificity" table.)

Method Comparison: The assay was compared to DPC's IMMULITE Theophylline assay on 155 patient samples. (Concentration range: approximately 1 to 60 µg/mL. See graph.)
By linear regression:

(IML 2000) = 0.98 (iML) + 0.73 µg/mL
r = 0.992

Means:
27.0 µg/mL (IMMULITE 2000)
26.7 µg/mL (IMMULITE)

Alternate Sample Types: To assess the effect of alternate sample types, blood was collected from 40 volunteers into plain, heparinized, EDTA and Becton Dickinson SST[®] vacutainer tubes. Equal volumes of the matched samples were spiked with various concentrations of theophylline, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE 2000 Theophylline procedure. By linear regression:

(Heparin) = 0.86 (Serum) + 1.8 µg/mL
r = 0.970

(EDTA) = 0.90 (Serum) + 2.5 µg/mL
r = 0.964

(SST) = 1.09 (Serum) + 0.5 µg/mL
r = 0.990

Means:
28.8 µg/mL (Serum)
26.5 µg/mL (Heparin)
28.4 µg/mL (EDTA)
31.8 µg/mL (SST)

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations

up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 512 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of tryglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

References

Butcher RW, Sutherland EW. Adenosine 3',5'-phosphate in biological materials. *J Biol Chem* 1962;237:1244-50. 2) Isles AF, MacLeod SM, Levison H. Theophylline. New thoughts about an old drug. *Chest* 1982;82:49S-54S. 3) Hendeles L, et al. Frequent toxicity from IV aminophylline infusions in critically ill patients. *Drug Intel Clin Pharmacol* 1977;11:12-8. 4) Jacobs MH, Senior RM, Kessler G. Clinical experience with theophylline. *JAMA* 1976;235:1983-6. 5) Zwillich CW, Sutton FD, Neff TA, et al. Theophylline-induced seizures in adults; correlation with serum concentrations. *Ann Intern Med* 1975;82:784-7. 6) Weinberger MW, et al. Intravenous aminophylline dosage; use of serum theophylline measurement for guidance. *JAMA* 1976;235:2110-3. 7) Kordash TR, Van Dellen RG, McCall JT. Theophylline concentrations in asthmatic patients after administration of aminophylline. *JAMA* 1977;238:139-41. 8) Piafsky KM, Ogilvie RI. Dosage of theophylline in bronchial asthma. *New Engl J Med* 1975;292:1218-22. 9) Weinberger M, Hendeles L. Slow-release theophylline; rationale and basis for product selection. *New Engl J Med* 1983;308:760-4. 10) Webb-Johnson DC, Andrews JL. Bronchodilator therapy. *New Engl J Med* 1977;297:476-82 and 758-64. 11) Lesko LJ. In: Benet LZ, et al, editors. *Pharmacokinetic basis for drug treatment*. New York: Raven Press, 1984: 321-424. 12) Pribor HC, Morrell G. *Lab Management* 1977(Feb):14. 13) Hendeles L, Weinberger M, Johnson G. Monitoring serum theophylline levels. *Clin Pharmacokinetics* 1978;3:294-312. 14) Neese AL, Soyka LF. Development of radioimmunoassay for theophylline. *Clin Pharmacol Therapeut* 1977;21:633-42. 15) Rowe DJF, et al. The clinical use and measurement of theophylline. *Ann Clin Biochem* 1988;25:4-26. 16) Gibaldi M, Perrier D. *Pharmacokinetics*. In: *Drugs and the pharmaceutical sciences*. New York: Marcel Dekker, 1975: 97. 17) Reitberg DP, Bernhard H, Schentag JJ. Alteration of theophylline clearance and half-life by cimetidine in normal volunteers. *Ann Intern Med* 1981;95:582-5.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782 or 973.927.2828
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

Manufactured by EURO/DPC Ltd. under a Quality System registered to ISO 13485:2003.

Tables and Graphs

Precision (µg/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	0.99	0.12	12.1%	0.13	13.4%
2	5.32	0.47	8.9%	0.42	7.9%
3	12.2	0.73	6.0%	0.72	6.0%
4	22.0	1.86	8.5%	1.88	8.5%
5	47.4	1.95	4.1%	2.81	5.9%

Linearity (µg/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	8.53	—	—
	4 in 8	4.72	4.27	111%
	2 in 8	2.07	2.13	97%
	1 in 8	1.28	1.07	120%
2	8 in 8	12.0	—	—
	4 in 8	6.20	6.0	103%
	2 in 8	3.05	3.0	102%
	1 in 8	1.40	1.5	94%
3	8 in 8	25.3	—	—
	4 in 8	12.3	12.7	97%
	2 in 8	6.69	6.33	106%
	1 in 8	3.12	3.16	99%
4	8 in 8	39.7	—	—
	4 in 8	19.9	19.9	100%
	2 in 8	9.58	9.93	96%
	1 in 8	5.70	4.96	115%
5	8 in 8	38.5	—	—
	4 in 8	19.3	19.3	100%
	2 in 8	9.40	9.63	98%
	1 in 8	5.52	4.81	115%

Recovery (µg/mL)


	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	1.60	—	—
	A	7.94	7.77	102%
	B	13.9	14.0	99%
	C	24.4	26.5	92%
2	—	2.62	—	—
	A	9.10	8.74	104%
	B	14.9	15.0	99%
	C	26.1	27.5	95%
3	—	5.71	—	—
	A	11.5	11.7	98%
	B	17.1	17.9	96%
	C	27.9	30.4	92%
4	—	10.3	—	—
	A	15.1	16.0	94%
	B	23.4	22.3	105%
	C	33.2	34.8	95%
5	—	17.9	—	—
	A	22.6	23.3	97%
	B	28.9	29.5	98%
	C	38.2	42.0	91%

Specificity

Compound ¹	µg/mL Added ²	Apparent µg/mL ³	% Cross-reactivity ⁴
Caffeine	1,000	3.62	0.4%
Acetamidophenol	1	ND	ND
3,7 Dimethyluric	1	ND	ND
1-Methylxanthine	1	ND	ND
1-Methyluric Acid	1	ND	ND
3-Methyluric Acid	1	ND	ND
1,3,7-Trimethyluric Acid	1	ND	ND
Phenobarbital	1	ND	ND
Theobromine	1	ND	ND

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
EURO/DPC LTD
Diagnostic Products Corporation
EC REP DPC Biermann GmbH