 IMMULITE<sup>®</sup>  
2000

## **Valproic Acid**

# IMMULITE<sup>®</sup> 2000 Valproic Acid

---

## English

---

**Intended Use:** For *in vitro* diagnostic use with the IMMULITE 2000 Analyzer — for the quantitative measurement of valproic acid in serum or heparinized plasma, as an aid in monitoring drug therapy.

Catalog Number: **L2KVA2** (200 tests), **L2KVA6** (600 tests)

Test Code: **VAL** Color: **Aqua**

### Summary and Explanation

Valproic acid is used primarily in the treatment of simple and complex absence seizures.<sup>1</sup> It can be administered alone or in combination with other drugs such as phenobarbital or phenytoin. Valproic acid increases the concentration of  $\gamma$ -aminobutyric acid (GABA) by inhibiting GABA transaminase.<sup>2</sup> GABA is a potent inhibitor of pre- and postsynaptic discharges in the central nervous system.

Valproic acid circulates mainly bound to protein. It competes with phenytoin for protein-binding sites. Valproic acid also inhibits the renal clearance of phenobarbital. When these two drugs are used in combination, the concentration of phenobarbital can be increased by up to 40 percent.<sup>3,4</sup>

Valproic acid is rapidly and almost completely absorbed after oral administration, reaching its peak concentration after 1 to 2 hours. In adults, its half-life is 16 hours.<sup>4</sup> The therapeutic range for valproic acid is 50 to 100  $\mu\text{g/mL}$ . Levels above 100  $\mu\text{g/mL}$  are considered toxic. Plasma concentrations should be monitored to help maintain an effective dosage and prevent side effects caused by an excessive concentration.<sup>4</sup>

### Principle of the Procedure

IMMULITE 2000 Valproic Acid is a solid-phase, competitive chemiluminescent enzyme immunoassay.

**Incubation Cycles:** 1  $\times$  30 minutes.

### Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Valproic Acid has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

**Volume Required:** 10  $\mu\text{L}$  serum or heparinized plasma.

**Storage:** 2 days at 2–8°C, or 1 month at –20°C.<sup>5</sup>

### Warnings and Precautions

For *in vitro* diagnostic use.

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight. (See insert.)

**Water:** Use distilled or deionized water.

## Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

### Valproic Acid Bead Pack (L2VA12)

With barcode. 200 beads, coated with monoclonal murine anti-valproic acid. Stable at 2–8°C until expiration date. **L2KVA2:** 1 pack. **L2KVA6:** 3 packs.

### Valproic Acid Reagent Wedge (L2VAA2)

With barcode. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to 6-amino-2-propylhexanoic acid in a human protein buffer/matrix. Stable at 2–8°C until expiration date. **L2KVA2:** 1 wedge. **L2KVA6:** 3 wedges.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

### Valproic Acid Adjustors (LVAL, LVAH)

Two vials (Low and High), 2.0 mL each, of valproic acid in processed human serum, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

**L2KVA2:** 1 set. **L2KVA6:** 2 sets.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

## Kit Components Supplied Separately

**L2SUBM:** Chemiluminescent Substrate

**L2PWSM:** Probe Wash

**L2KPM:** Probe Cleaning Kit

**LRXT:** Reaction Tubes (disposable)

**SDCM:** Bi-level, human serum-based multi-constituent control module

Also Required

Distilled or deionized water; test tubes; controls.

## Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Operator's Manual.

See the IMMULITE 2000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

### Recommended Adjustment Interval:

2 weeks.

**Quality Control Samples:** Use controls or serum pools with at least two levels (low and high) of valproic acid.

## Therapeutic Range

The commonly accepted therapeutic range for valproic acid is 50 – 100 µg/mL (347 – 693 µmol/L). Toxicity has been reported at concentrations greater than 100 µg/mL.<sup>4</sup>

Consider these limits as *guidelines* only.

## Limitations

A complete understanding of the information contained in this package insert is necessary for proper use of this product.

The specificity table should be consulted for a list of substances which have measurable crossreactivity in this assay. Other medications and/or food not found on this list may nevertheless crossreact in the assay.

The following should be considered in the interpretation of the valproic acid concentration:

- All concurrent medications taken by the patient
- The time of sampling in relation to the dosing history
- Existing diseases or conditions

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic

purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

### Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in µg/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

**Conversion Factor:**  
µg/mL × 6.93 → µmol/L

**Calibration Range:** 10 – 200 µg/mL  
(69 – 1,386 µmol/L).

**Analytical Sensitivity:** 1.5 µg/mL  
(10.4 µmol/L).

**Precision:** Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

**Linearity:** Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

**Recovery:** Samples spiked 1 to 19 with three valproic acid solutions (525, 1,050 and 2,300 µg/mL) were assayed. (See "Recovery" table for representative data.)

**Specificity:** The antibody is highly specific for Valproic Acid. (See "Specificity" table.)

**Bilirubin:** Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

**Hemolysis:** Presence of hemoglobin in concentrations up to 512 mg/dL has no effect on results, within the precision of the assay.

**Lipemia:** Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

**Alternate Sample Type:** Samples ( $n = 15$ ) were collected into plain and heparinized vacutainer tubes. All samples were spiked with valproic acid and then assayed by the IMMULITE 2000 Valproic Acid procedure.

(Heparin) = 1.04 (Serum) – 0.13 µg/mL  
 $r = 0.970$

Means:  
82 µg/mL (Serum)  
85 µg/mL (Heparin)

**Method Comparison:** The assay was compared to DPC's IMMULITE Valproic Acid assay on 128 samples. (Concentration range: approximately 13 to 200 µg/mL. See graph.) By linear regression:

(IML 2000) = 1.0 (IML) – 1.8 µg/mL  
 $r = 0.997$

Means:  
79 µg/mL (IMMULITE 2000)  
80 µg/mL (IMMULITE)

### References

- 1) Gerson B, editor. Essentials of therapeutic drug monitoring. New York: IGAKU-SHOIN, 1983: 215-26.
- 2) Boeckx RL. In: Moyer TP, Boeckx RL, editors. Applied therapeutic drug monitoring. Washington DC: The American Association for Clinical Chemistry, 1984: 79-81.
- 3) Drug Monitoring Data Pocket Guide II, 2nd ed. Washington DC: AACC Press, 1994: 133-5.
- 4) Moyer TP. Therapeutic drug monitoring. In: Burtic AC, Ashwood ER, editors. Tietz textbook of clinical chemistry, 3rd ed. Philadelphia: WB Saunders, 1999; 873,876-7.
- 5) Quest Diagnostics. The 1997/1998 Quest Diagnostics Nichols Institute reference manual. San Juan Capistrano: Quest Diagnostics, 1997: 271.

### Technical Assistance

In the United States, contact DPC's Technical Services department.  
Tel: 800.372.1782 or 973.927.2828  
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

The Quality System of Diagnostic Products Corporation is registered to ISO 13485:2003.

---

## Tables and Graphs

---

### Precision (µg/mL)

	Mean <sup>3</sup>	Within-Run <sup>1</sup>		Total <sup>2</sup>	
		SD <sup>4</sup>	CV <sup>5</sup>	SD	CV
1	15.6	0.6	3.8%	1.1	7.1%
2	34.3	1.4	4.1%	2.3	6.7%
3	55.5	1.5	2.7%	3.0	5.4%
4	110	3.0	2.7%	5.1	4.6%

**Linearity (µg/mL)**

	Dilution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	8 in 8	103	—	—
	4 in 8	49.2	51.5	96%
	2 in 8	21.5	25.8	83%
	1 in 8	10.8	12.9	84%
2	8 in 8	128	—	—
	4 in 8	63.0	64.0	98%
	2 in 8	29.8	32.0	93%
	1 in 8	14.1	16.0	88%
3	8 in 8	150	—	—
	4 in 8	76.8	75.0	102%
	2 in 8	37.0	37.5	99%
	1 in 8	15.9	18.8	85%
4	8 in 8	167	—	—
	4 in 8	86.0	83.5	103%
	2 in 8	42.5	41.8	102%
	1 in 8	17.5	20.9	84%
5	8 in 8	185	—	—
	4 in 8	99.4	92.5	107%
	2 in 8	45.5	46.3	98%
	1 in 8	21.5	23.1	93%

**Recovery (µg/mL)**


	Solution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	—	ND	—	—
	A	25.6	26.3	98%
	B	52.5	52.5	100%
	C	112	115	97%
2	—	10.1	—	—
	A	35.9	35.8	100%
	B	65.4	62.1	105%
	C	125	125	100%
3	—	14.8	—	—
	A	40.9	40.3	101%
	B	69.1	66.6	104%
	C	129	129	100%
4	—	25.2	—	—
	A	52.1	50.2	104%
	B	80.2	76.5	105%
	C	143	139	103%
5	—	50.2	—	—
	A	77.5	73.9	105%
	B	110	100	110%
	C	165	163	102%

**Specificity**

Compound <sup>1</sup>	µg/mL Added <sup>2</sup>	% Cross reactivity <sup>3</sup>
Carbamazepine	140	ND
Clonazepam	1.2	16%
Diazepam	25	2%
Ethosuximide	120	ND
2-propyl Glutaric Acid	250	11%
2- <i>n</i> -propyl-4 Pentenoic Acid	100	20%
Phenobarbital	200	ND
Phenytoin	160	ND
Sodium Salicylate	100	ND

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O Sistema de Qualidade da Diagnostic Products Corporation está registado sob ISO 13485:2003.
EC REP DPC Biermann GmbH