# **hGH**

**DOUBLE ANTIBODY**

**DPC®**

# Double Antibody hGH

## English

**Intended Use:** Double Antibody Human Growth Hormone (hGH) is an  $^{125}\text{I}$  radioimmunoassay designed for the quantitative measurement of human growth hormone in serum or plasma. It is intended strictly for *in vitro* diagnostic use.

Catalog numbers: **KGHD1** (70 tubes), **KGHD2** (140 tubes).

The 70-tube kit contains less than 3 microcuries (111 kilobecquerels) of radioactive  $^{125}\text{I}$  growth hormone, and the 140-tube kit contains less than 6 microcuries (222 kilobecquerels).

### Summary and Explanation of the Test

Human growth hormone (hGH, somatotropin) is a polypeptide originating in the anterior pituitary. It is 191 amino acids in length and has a molecular mass of approximately 22,000 daltons. Its metabolic effects are primarily anabolic. It promotes protein conservation and engages a wide range of mechanisms for protein synthesis. It also enhances glucose transport and facilitates the buildup of glycogen stores. Another family of peptide hormones, the somatomedins, mediates its cascade of growth-promoting actions.

Measurement of hGH is primarily of interest in the diagnosis and treatment of various forms of inappropriate growth hormone secretion. Clinical disorders of hyposecretion include dwarfism and unattained growth potential. Hypersecretion is associated with gigantism and acromegaly.

Caution must be exercised in the clinical interpretation of growth hormone levels. These vary throughout the day, making it difficult to define a reference range or to judge an individual's status based on single determinations. Many factors are known to influence the rate of growth hormone secretion, including periods of sleep and wakefulness, exercise, stress, hypoglycemia, estrogens, corticosteroids, L-dopa, and others. Because of its

similarity to prolactin and placental lactogen, earlier growth hormone immunoassays were often plagued with falsely high values in pregnant and lactating women.

Because not all acromegalic individuals have elevated baseline levels, suppression tests based on glucose loading are of value in this context. In spite of the induced hyperglycemia, there is rarely a decrease from baseline levels in acromegaly.

Growth hormone-deficient individuals have fasting/resting levels similar to those found in normal individuals. Various challenge tests have therefore been devised to differentiate these groups. Thus, with the onset of deep sleep or after 15 to 20 minutes of vigorous exercise, growth hormone levels normally show a rise. Other tests of growth hormone responsiveness are based on the administration of L-dopa, arginine and insulin. Propanolol or estrogen is sometimes given in conjunction with the primary stimulus to accentuate the response.

A small number of cases of dwarfism have been documented in which both the basal level and the response to challenge testing were normal. Such cases may involve tissue insensitivity to either growth hormone or the somatomedins, or the presence of antibodies or immunoreactive but biologically inactive growth hormone.

### Principle of the Procedure

In the Double Antibody hGH procedure,  $^{125}\text{I}$ -labeled growth hormone competes with growth hormone in the patient sample for sites on growth hormone-specific antibody. After incubation for a fixed time, separation of bound from free is achieved by the PEG-accelerated double-antibody method. The tube is then counted in a gamma counter, the counts being inversely related to the amount of growth hormone present in the patient sample.

**Reagents to Pipet:** 3

**Total Incubation Time:** 2 Hours

**Separation:** The ready-to-use Precipitating Solution combines second antibody and dilute PEG.

**Total Counts at Iodination:**  
approximately 60,000 cpm

## Warnings and Precautions

For *in vitro* diagnostic use.

**Reagents:** Store at 2–8°C in a refrigerator designated for incoming radioactive materials. Dispose of in accordance with applicable laws.

Do not use reagents beyond their expiration dates.

Some components supplied in this kit may contain human source material and/or other potentially hazardous ingredients which necessitate certain precautions.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

**Water:** Use distilled or deionized water.

## Radioactivity

A copy of any radioisotope license certificate (Specific or General) issued to a US customer must be on file with Diagnostic Products Corporation before kits or components containing radioactive material can be shipped. These radioactive materials may be acquired by any customer with the appropriate Specific license. Under a General license these radioactive materials may be acquired only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories and hospitals — and strictly for *in vitro* clinical or laboratory tests not involving external or internal administration of the radioactive material or its radiation to human beings or other animals. Its acquisition, receipt, storage, use, transfer and disposal are all subject to the regulations and a (General or

Specific) license of the U.S. Nuclear Regulatory Commission or a State with which the NRC has entered into an agreement for the exercise of regulatory control.

Handle radioactive materials according to the requirements of your General or Specific license. To minimize exposure to radiation, the user should adhere to guidelines set forth in the National Bureau of Standards publication on the *Safe Handling of Radioactive Materials* (Handbook No. 92, issued March 9, 1964) and in subsequent publications issued by State and Federal authorities.

Wipe up spills promptly and decontaminate affected surfaces. Avoid generation of aerosols. Dispose of solid radioactive waste according to license requirements. General licensees (holders of NRC Form 483) may dispose of solid radioactive waste as nonradioactive waste, after removing labeling. Specific licensees (NRC Form 313) should refer to Title 10, Code of Federal Regulations, Part 20. Licensees in Agreement States should refer to the appropriate regulations of their own state. General licensees may dispose of liquid radioactive waste of the type contained in this product through a laboratory sink drain. Licensees must remove or deface labels from empty containers of radioactive materials before disposal of solid waste. Specific licensees may dispose of small quantities of liquid radioactive waste of the type used in this product through a laboratory sink drain. Refer to the appropriate regulations applicable to your laboratory.

## Materials Supplied: Initial Preparation

Prepare all components at least 10 minutes before use.

### hGH Antiserum (GHD1)

Lyophilized growth hormone antiserum, with preservative. Reconstitute each vial by adding **7.7 mL** distilled or deionized water. Mix by *gentle* inversion. Stable at 2–8°C for 30 days after reconstitution. Color: blue.

**KGHD1:** 1 vial. **KGHD2:** 2 vials.

#### <sup>125</sup>I hGH (GHD2)

Lyophilized iodinated growth hormone, with preservative. Reconstitute each vial by adding **7.7 mL** distilled or deionized water. Mix by *gentle* inversion. Stable at 2–8°C for 30 days after reconstitution, or until the expiration date marked on the vial. Color: yellow.

**KGHD1:** 1 vial. **KGHD2:** 2 vials.

#### hGH Calibrators (GHD3–9)

Seven vials, labeled A through G, of lyophilized calibrators. At least 30 minutes before use, reconstitute the zero calibrator **A with 5.0 mL** distilled or deionized water, and each of the remaining calibrators **B through G with 2.0 mL** distilled or deionized water. Use volumetric pipets and mix by *gentle* inversion. Stable at 2–8°C for 30 days after reconstitution. For longer storage, aliquot and freeze: stable at –20°C for 6 months.

**KGHD1:** 1 set. **KGHD2:** 1 set.

The calibrators have *lot-specific* values of *approximately* 0, 1, 2, 5, 10, 20 and 30 ng/mL *in a human serum matrix*. The assay is standardized in terms of the World Health Organization's International Reference Preparations of human Growth Hormone for immunoassay, numbers 66/217 and 80/505, and the National Institute of Health's reference preparation for hGH NIAMDD-hGH-RP-1. *Refer to the vial labels for exact concentrations in ng/mL*. Intermediate calibration points can be obtained by mixing calibrators in suitable proportions.

#### Precipitating Solution (ACGG)

Goat anti-rabbit gamma globulin (GARGG) and dilute polyethylene glycol (PEG) in saline. The Precipitating Solution is supplied in liquid form, ready to use. Stable at 2–8°C for 30 days after opening. The Precipitating Solution should be thoroughly mixed before use, since a fine precipitate may form after storage at 2–8°C. Color: blue.

**KGHD1:** 1 vial. **KGHD2:** 2 vials.

#### Materials Required But Not Provided

Gamma counter

Centrifuge — preferably refrigerated and capable of at least 3000×g

Vortex mixer

#### Reagent Preparation

Distilled or deionized water

Volumetric pipets: 2.0, 5.0 and 7.7 mL

#### Radioimmunoassay

Plain 12x75 mm polypropylene tubes — available from DPC

Micropipets: 100 µL and 1.0 mL.

Foam decanting rack — available from DPC

Logit-log graph paper — available from DPC (catalog number: ZP797)

A tri-level human serum-based immunoassay control, containing hGH as one of over 25 assayed constituents, is available from DPC (catalog number: CON6).

#### Specimen Collection

Patient must be fasting and at complete rest 30 minutes before blood collection.<sup>6,7</sup>

Collect blood by venipuncture<sup>5</sup> into plain, heparinized or EDTA tubes, avoiding hemolysis, and separate the serum or plasma from the cells. The time of collection should be noted.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Lipemic, hemolyzed or grossly contaminated samples may give erroneous results.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. Double Antibody hGH has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

**Volume Required:** 100 µL of serum or plasma per tube.

**Storage:** Refrigerate immediately. Stable at 2–8°C for 8 hours,<sup>6,7</sup> or for 2 months (aliquotted) at –20°C. (Avoid repeated thawing and freezing.)

Before assay, allow the samples to come to room temperature and mix by *gentle* swirling or inversion. *Do not* attempt to thaw frozen specimens by heating them in a waterbath.

## Radioimmunoassay Procedure

All components except the Precipitating Solution must be at room temperature (15–28°C) before use.

- 1 Label eighteen tubes in duplicate: T (total counts), NSB (nonspecific binding), A (maximum binding) and B through G. Label additional tubes, also in duplicate, for patient samples and controls.

| Calibrator | Approximate ng/mL |
|------------|-------------------|
| A (MB)     | 0                 |
| B          | 1                 |
| C          | 2                 |
| D          | 5                 |
| E          | 10                |
| F          | 20                |
| G          | 30                |

**Note:** The calibrators are *lot-specific*. Refer to the vial labels for values in ng/mL.

- 2 Pipet **100 µL** of the zero calibrator A into the NSB and A tubes, and **100 µL** of each of the remaining calibrators B through G into correspondingly labeled tubes. Pipet **100 µL** of each patient sample (serum or plasma) and **100 µL** of each control into the tubes prepared.

Patient samples expected to contain hGH concentrations greater than the highest calibrator (30 ng/mL) should be diluted in the zero calibrator before assay.

It is good practice to use a *disposable-tip* micropipet, changing the tip between samples, in order to avoid carryover contamination.

- 3 Add **100 µL** of hGH Antiserum (Blue) to all tubes except the NSB and T tubes. Vortex.

All tubes except the NSB and T tubes should be blue-green. A repeating dispenser is recommended for the 100 µL additions of antiserum and tracer.

- 4 Incubate for **1 hour** at room temperature (15–28°C).  
The incubation may be extended up to 90 minutes at room temperature providing the second incubation (step 6) is similarly modified.

- 5 Add **100 µL** of <sup>125</sup>I hGH (Yellow) to all tubes. Vortex.

All tubes except the NSB and T tubes should be green. Remove the T tubes for counting at step 10; they require no further processing.

- 6 Incubate for **1 hour** at room temperature (15–28°C).

This incubation step should be handled in the same fashion as step 4 above.

- 7 Add **1.0 mL** of *cold* Precipitating Solution (Blue) to all tubes. Vortex.

All tubes should be blue, except for the T tubes which should be yellow. Before dispensing, the Precipitating Solution should be thoroughly mixed (without foaming) until all visible precipitate has disappeared.

- 8 Centrifuge for **15 minutes** at 3000xg.

- 9 Decant (or aspirate) the supernatant, retaining the precipitate for counting.

Using a foam decanting rack, gently tap the tubes on absorbant paper and blot the rims, to remove residual droplets.

- 10 Count each tube for **1 minute**.

## Calculation of Results

To calculate growth hormone concentrations from a logit-log representation of the calibration curve, first calculate for each pair of tubes the average NSB-corrected counts per minute:

Net Counts = (Average CPM) minus (Average **NSB** CPM)

Then determine the binding of each pair of tubes as a percent of maximum binding (MB), with the NSB-corrected counts of the A tubes taken as 100%:

Percent Bound = (Net Counts / Net **MB** Counts) × 100

Using logit-log graph paper, plot Percent Bound on the vertical (probability) axis against Concentration on the horizontal (logarithmic) axis for each of the nonzero calibrators, and draw a straight line approximating the path of these points. Results for the unknowns may then be read from the line by interpolation.

It is good practice to inspect results for agreement within replicates, and to construct a graph of the calibration curve (even if the calculations are handled by computer) as a visual check on the appropriateness of the transformation used and as a way to detect deviant calibration points. We also recommend keeping track of these data reduction parameters:

$T$  = Total Counts (as counts per minute)

$\%NSB = 100 \times (\text{Average NSB Counts} / \text{Total Counts})$

$\%MB = 100 \times (\text{Net Counts} / \text{Total Counts})$

And the 20, 50 and 80 percent "intercepts," where

20% = Concentration at 20 Percent Bound, etc.

Note that other approaches, e.g. a sound implementation of the 4-parameter logistic, are also acceptable. See Dudley RA, et al. Guidelines for immunoassay data reduction. Clin Chem 1985;31:1264-71.

**Example Run:** For illustration only, not for calculating results from another run. Because the calibrator values are *lot-specific*, concentrations listed in the right-most column may not match the values of the calibrators supplied in your shipment. (See "Example Run" table.)

## Quality Control

Controls or serum pools with at least two hGH concentration levels — low and high — should routinely be assayed as unknowns.

Report patient results only if the control results for that assay meet your laboratory's established criteria for acceptability.

It is good laboratory practice to record for each assay the lot numbers of the components used, as well as the dates when they were first reconstituted or opened. We also recommend charting control results from day to day — as described, for example, in Westgard JO, et al. A multi-rule chart for quality control. Clin Chem 1981;27:493-501. Note that repeat samples can serve as a valuable additional tool for monitoring interassay precision, and that pairs of control tubes can be spaced throughout the assay to help verify the absence of significant drift.

**Centrifugation:** The procedure calls for centrifuging at 3000xg for 15 minutes. Lower accelerations are satisfactory only if the centrifugation time is increased appropriately — for example, to 30 minutes at 1500xg. A high-speed, refrigerated centrifuge is desirable but not essential. Use the formulas below to calculate the acceleration of your centrifuge at a given speed, or the speed (in revolutions per minute) required to achieve a desired g force.

$g \text{ force} = 28.38 \times (\text{rpm} / 1,000)^2 \times \text{radius [inches]}$

$\text{rpm} = 187.7 \times \sqrt{(g \text{ force} / \text{radius [inches]})}$

or

$g \text{ force} = 11.18 \times (\text{rpm} / 1,000)^2 \times \text{radius [cm]}$

$\text{rpm} = 299.07 \times \sqrt{(g \text{ force} / \text{radius [cm]})}$

## Expected Values

A reference range study was performed at a large hospital laboratory in which fasting samples from 72 adult volunteers were assayed by the Double Antibody hGH procedure. The results of the study suggest a reference range for adults of:

Nondetectable to 7 ng/mL.

Laboratories should consider these results *as guidelines only*. Each laboratory should establish its own reference ranges.

## Performance Data

See Tables and Graphs for data *representative* of the Double Antibody hGH kit's performance. Results are expressed in ng/mL.

### Conversion Factor:

$\text{ng/mL} \times 2 \rightarrow \text{mIU/L [WHO IRP 66/217]}$

$\text{ng/mL} \times 2.6 \rightarrow \text{mIU/L [WHO 1<sup>st</sup> IRP 80/505]}$

### Calibration Range:

Approximately 1–30 ng/mL.

The calibrator values are *lot-specific*.

**Analytical Sensitivity:** 0.9 ng/mL

### Intraassay Precision (Within-Run):

Statistics were calculated for each of six samples from the results of 20 pairs of tubes in a single assay. (See "Intraassay Precision" table.)

### Interassay Precision (Run-to-Run):

Statistics were calculated for each of six samples from the results of pairs of tubes

in 20 different assays. (See "Interassay Precision" table.)

**End-of-Run Effect:** None up to approximately 150 tubes. (See "End-of-Run Effect" table.)

**Specificity:** The antibody is highly specific for human growth hormone. (See "Specificity" table.)

**Linearity:** Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

**Recovery:** Samples spiked 1 to 19 with four hGH solutions (50, 100, 200 and 400 ng/mL) were assayed. (See "Recovery" table for representative data.)

**Bilirubin:** Presence of bilirubin in concentrations up to 100 mg/L has no effect on results, within the precision of the assay.

**Hemolysis:** Results obtained on hemolyzed samples should be interpreted with caution, since there is a trend — though not a uniform one — towards lower values as a function of increasing hemoglobin concentration. (See "Hemolysis" table.)

**Free Fatty Acid:** Presence of oleic acid in concentrations up to 10 mmol/L has no clinically significant *in vitro* effect on the results.

**Alternate Sample Type:** To determine whether anticoagulants interfere with the assay, blood was collected from 20 healthy volunteers into plain, heparinized and EDTA vacutainer tubes. All samples were assayed by the Double Antibody hGH procedure. By linear regression:

(Heparin) =  $1.02 (\text{Serum}) + 0.04 \text{ ng/mL}$   
 $r = 0.998$

(EDTA) =  $0.97 (\text{Serum}) + 0.09 \text{ ng/mL}$   
 $r = 0.995$

Means:

3.8 ng/mL (Heparin)  
3.7 ng/mL (Serum)  
3.7 ng/mL (EDTA)

**Method Comparison:** The assay was compared to a commercially available radioimmunoassay (Kit A) on 145 patient serum samples. (Concentration range: nondetectable to approximately 7.4 ng/mL by the Double Antibody procedure.) By linear regression:

(Kit A) =  $1.13 (\text{DAb}) + 0.13 \text{ ng/mL}$   
 $r = 0.985$

Means:

1.72 ng/mL (Double Antibody)  
2.07 ng/mL (Kit A)

The slightly lower results obtained with the Double Antibody kit are probably due to the fact that it employs calibrators in a human serum matrix for maximal compatibility with typical patient samples, whereas Kit A makes use of a bovine-buffer matrix.

## References

1) Daughaday WM, Cryer PE. Growth hormone hypersecretion and acromegaly. Hospital Practice. no. 8 August 1978;13:75-80. 2) Raite S. The standards for human growth hormone assays. In: Laron Z, Butenandt O, editors. Evaluation of Growth Hormone Secretion. Basel: Karger, 1983: 162-9. 3) Van Wyk JJ. Diagnosis of growth hormone deficiency and other causes of growth failure. The Endocrine Society: Syllabus for the 35th Annual Postgraduate Assembly October 1983: 90-106. 4) Van Wyk JJ, Underwood L. Growth hormone, somatomedins and growth failure. Hospital Practice no. 8 August 1978;13:57-67. 5) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 6) Tietz NW, editor. Clinical guide to laboratory tests. 3rd ed. Philadelphia: W.B. Saunders, 1995: 300. 7) Burtis CA, Ashwood ER, editors. Tietz textbook of clinical chemistry. 2nd ed. Philadelphia: W.B. Saunders, 1994: 1670.

## Technical Assistance

In the United States, contact DPC's Technical Services department.  
Tel: 800.372.1782, 310.645.8200  
Fax: 310.645.9999. To place an order:  
Tel: 800.372.1782. Fax: 800.234.4372.  
Outside the United States, contact your National Distributor.

The Quality System of Diagnostic Products Corporation is registered to ISO 13485:2003.

## Tables and Graphs

### Example Run

| Tube <sup>1</sup>      | Duplicate CPM <sup>2</sup> | Average CPM <sup>3</sup> | Net CPM <sup>4</sup> | Percent Bound <sup>5</sup> | Approx. hGH ng/mL <sup>6</sup> |
|------------------------|----------------------------|--------------------------|----------------------|----------------------------|--------------------------------|
| T <sup>7</sup>         | 51,017<br>51,434           | 51,226                   |                      |                            |                                |
| NSB <sup>8</sup>       | 829<br>986                 | 908                      | 0                    |                            |                                |
| A<br>(MB) <sup>9</sup> | 23,511<br>23,725           | 23,618                   | 22,710               | 100%                       | 0                              |
| B                      | 21,961<br>22,020           | 21,991                   | 21,083               | 92.8%                      | 1                              |
| C                      | 20,094<br>20,325           | 20,210                   | 19,302               | 85.0%                      | 2                              |
| D                      | 15,001<br>15,108           | 15,055                   | 14,147               | 62.3%                      | 5                              |
| E                      | 8,910<br>9,030             | 8,971                    | 8,063                | 35.5%                      | 10                             |
| F                      | 4,630<br>4,724             | 4,678                    | 3,770                | 16.6%                      | 20                             |
| G                      | 3,400<br>3,490             | 3,446                    | 2,538                | 11.2%                      | 30                             |
| Unknowns <sup>10</sup> |                            |                          |                      |                            |                                |
| X1                     | 21,826<br>22,009           | 21,918                   | 21,010               | 92.5%                      | 1.2                            |
| X2                     | 13,733<br>13,785           | 13,760                   | 12,852               | 56.6%                      | 5.7                            |
| X3                     | 8,717<br>8,794             | 8,756                    | 7,848                | 34.6%                      | 10.6                           |

Quality Control Parameters:<sup>11</sup>

T<sup>7</sup> = 51,226 cpm

%NSB<sup>8</sup> = 1.8%

%MB<sup>9</sup> = 44%

20% Intercept<sup>12</sup> = 17.7 ng/mL

50% Intercept = 6.8 ng/mL

80% Intercept = 2.6 ng/mL

### Intraassay Precision (ng/mL)

|   | Mean <sup>1</sup> | SD <sup>2</sup> | CV <sup>3</sup> |
|---|-------------------|-----------------|-----------------|
| 1 | 1.1               | 0.065           | 5.9%            |
| 2 | 1.7               | 0.048           | 2.8%            |
| 3 | 2.3               | 0.070           | 3.0%            |
| 4 | 5.3               | 0.086           | 1.6%            |
| 5 | 6.0               | 0.092           | 1.5%            |
| 6 | 10.9              | 0.167           | 1.5%            |

### Interassay Precision (ng/mL)

|   | Mean <sup>1</sup> | SD <sup>2</sup> | CV <sup>3</sup> |
|---|-------------------|-----------------|-----------------|
| 1 | 1.1               | 0.091           | 8.3%            |
| 2 | 1.6               | 0.084           | 5.3%            |
| 3 | 2.3               | 0.101           | 4.4%            |
| 4 | 5.3               | 0.122           | 2.3%            |
| 5 | 5.7               | 0.102           | 1.8%            |
| 6 | 10.7              | 0.365           | 3.4%            |

### End-of-Run Effect (ng/mL)

|   | Tube's<br>31–42 | Tubes<br>83–94 | Tubes<br>135–146 | Tubes<br>187–198 |
|---|-----------------|----------------|------------------|------------------|
| 1 | 1.6             | 1.6            | 1.7              | 1.8              |
| 2 | 1.6             | 1.6            | 1.8              | 1.8              |
| 3 | 2.5             | 2.5            | 2.5              | 2.7              |
| 4 | 3.9             | 3.9            | 4.2              | 4.1              |
| 5 | 8.7             | 8.9            | 8.7              | 8.8              |
| 6 | 12.6            | 12.7           | 13.1             | 12.8             |

### Specificity

| Compound <sup>1</sup> | ng/mL<br>Added <sup>2</sup> | % Cross-<br>reactivity <sup>3</sup> |
|-----------------------|-----------------------------|-------------------------------------|
| Prolactin             | 40,000                      | ND                                  |
|                       | 10,000                      | ND                                  |
|                       | 1,000                       | ND                                  |
|                       | 200                         | ND                                  |
| HPL                   | 100,000                     | ND                                  |
| TSH                   | 80,000                      | 0.006%                              |
| FSH                   | 5,340                       | ND                                  |
| LH                    | 5,340                       | ND                                  |
| HCG, $\beta$ -subunit | 10,000                      | ND                                  |

ND: not detectable<sup>4</sup> (crossreactivity of less than 0.006%<sup>5</sup>).

### Hemolysis

| Packed Red Blood Cells <sup>1</sup> |                       |               |               |               |
|-------------------------------------|-----------------------|---------------|---------------|---------------|
|                                     | Unspiked <sup>2</sup> | 10 $\mu$ L/mL | 20 $\mu$ L/mL | 30 $\mu$ L/mL |
| 1                                   | 2.2                   | 1.8           | 1.7           | 1.5           |
| 2                                   | 5.1                   | 4.6           | 4.5           | 4.2           |
| 3                                   | 4.9                   | 5.1           | 4.9           | 4.5           |

## Linearity (ng/mL)

|   | Dilution <sup>1</sup> | Observed <sup>2</sup> | Expected <sup>3</sup> | %O/E <sup>4</sup> |
|---|-----------------------|-----------------------|-----------------------|-------------------|
| 1 | 8 in 8 <sup>5</sup>   | 3.2                   | —                     | —                 |
|   | 4 in 8                | 1.7                   | 1.6                   | 106%              |
|   | 2 in 8                | 0.9                   | 0.8                   | 113%              |
|   | 1 in 8                | 0.4                   | 0.4                   | 100%              |
| 2 | 8 in 8                | 8.5                   | —                     | —                 |
|   | 4 in 8                | 4.2                   | 4.3                   | 98%               |
|   | 2 in 8                | 2.2                   | 2.1                   | 105%              |
|   | 1 in 8                | 1.1                   | 1.1                   | 100%              |
| 3 | 8 in 8                | 8.8                   | —                     | —                 |
|   | 4 in 8                | 4.5                   | 4.4                   | 102%              |
|   | 2 in 8                | 2.3                   | 2.2                   | 105%              |
|   | 1 in 8                | 1.2                   | 1.1                   | 109%              |

## Recovery (ng/mL)

|   | Solution <sup>1</sup> | Observed <sup>2</sup> | Expected <sup>3</sup> | %O/E <sup>4</sup> |
|---|-----------------------|-----------------------|-----------------------|-------------------|
| 1 | —                     | 0.0                   | —                     | —                 |
|   | A                     | 2.5                   | 2.5                   | 100%              |
|   | B                     | 4.9                   | 5.0                   | 98%               |
|   | C                     | 9.8                   | 10.0                  | 98%               |
|   | D                     | 22.6                  | 20.0                  | 113%              |
| 2 | —                     | 0.0                   | —                     | —                 |
|   | A                     | 2.5                   | 2.5                   | 100%              |
|   | B                     | 4.9                   | 5.0                   | 98%               |
|   | C                     | 10.2                  | 10.0                  | 102%              |
|   | D                     | 20.2                  | 20.0                  | 101%              |
| 3 | —                     | 2.9                   | —                     | —                 |
|   | A                     | 5.3                   | 5.3                   | 100%              |
|   | B                     | 7.8                   | 7.8                   | 100%              |
|   | C                     | 12.0                  | 12.8                  | 94%               |
|   | D                     | 22.3                  | 22.8                  | 98%               |

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O Sistema de Qualidade da Diagnostic Products Corporation está registado sob ISO 13485:2003.
EC REP DPC Biermann GmbH