Aldosterone

DPC[®]

Coat-A-Count[®] Aldosterone

English

Intended Use: Coat-A-Count Aldosterone is a solid-phase ¹²⁵I radioimmunoassay designed for the quantitative measurement of aldosterone levels in unextracted serum. Coupled with an ethyl acetate extraction procedure, it is also suitable for assaying urine samples. The kit is intended strictly for *in vitro* diagnostic use as an aid in clinical diagnosis.

Catalog numbers: **TKAL1** (100 tubes), **TKAL2** (200 tubes)

The 100-tube kit contains less than 3.5 microcuries (130 kilobecquerels) of radioactive ¹²⁵I aldosterone, and the 200-tube kit contains less than 7 microcuries (259 kilobecquerels).

Summary and Explanation

Normal serum or plasma levels of aldosterone are dependent on the sodium intake and whether the patient is standing or recumbent. High sodium intake will tend to suppress serum aldosterone, whereas low sodium intake will elevate serum aldosterone.

Principle of the Procedure

The Coat-A-Count procedure is a solid-phase radioimmunoassay, based on aldosterone-specific antibody immobilized to the wall of a polypropylene tube. ¹²⁵I-labeled aldosterone competes for a fixed time with aldosterone in the patient sample for antibody sites. The tube is then decanted, to separate bound from free, and counted in a gamma counter. The amount of aldosterone present in the patient sample is determined from a calibration curve.

Reagents to Pipet: 1

Total Incubation Time: 18 Hours at room temperature (15–28°C).

Total Counts at Iodination: approximately 40,000 cpm

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C in a refrigerator designated for incoming radioactive materials. Dispose of in accordance with applicable laws.

Do not use reagents beyond their expiration dates.

Some components supplied in this kit may contain human source material and/or other potentially hazardous ingredients which necessitate certain precautions.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Water: Use distilled or deionized water.

Radioactivity

A copy of any radioisotope license certificate (Specific or General) issued to a US customer must be on file with Diagnostic Products Corporation before kits or components containing radioactive material can be shipped. These radioactive materials may be acquired by any customer with the appropriate Specific license. Under a General license these radioactive materials may be acquired only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories and hospitals — and strictly for *in vitro* clinical or laboratory tests not involving external or internal administration of the radioactive material or its radiation to human beings or other animals. Its acquisition, receipt, storage, use, transfer and disposal are all subject to the regulations and a (General or Specific) license of the U.S. Nuclear Regulatory Commission or a State with which the NRC has entered into an

agreement for the exercise of regulatory control.

Handle radioactive materials according to the requirements of your General or Specific license. To minimize exposure to radiation, the user should adhere to guidelines set forth in the National Bureau of Standards publication on the *Safe Handling of Radioactive Materials* (Handbook No. 92, issued March 9, 1964) and in subsequent publications issued by State and Federal authorities.

Wipe up spills promptly and decontaminate affected surfaces. Avoid generation of aerosols. Dispose of solid radioactive waste according to license requirements. General licensees (holders of NRC Form 483) may dispose of solid radioactive waste as nonradioactive waste, after removing labeling. Specific licensees (NRC Form 313) should refer to Title 10, Code of Federal Regulations, Part 20. Licensees in Agreement States should refer to the appropriate regulations of their own state. General licensees may dispose of liquid radioactive waste of the type contained in this product through a laboratory sink drain. Licensees must remove or deface labels from empty containers of radioactive materials before disposal of solid waste. Specific licensees may dispose of small quantities of liquid radioactive waste of the type used in this product through a laboratory sink drain. Refer to the appropriate regulations applicable to your laboratory.

Materials Supplied: Initial Preparation

Aldosterone Ab-Coated Tubes (TAL1)

Polypropylene tubes coated with antibodies to aldosterone and packaged in zip-lock bags. Store refrigerated and protected from moisture, carefully resealing the bags after opening. Stable at 2–8°C for one year from the date of manufacture. Color: light brown.

TKAL1: 100 tubes. **TKAL2:** 200 tubes.

¹²⁵I Aldosterone (TAL2)

Lyophilized iodinated aldosterone. Reconstitute each vial by adding a measured **110 mL** of distilled or deionized water. Let stand for 10 minutes, then mix by *gentle* inversion. Stable at 2–8°C for 30 days after reconstitution, or until the

expiration date marked on the label.

TKAL1: 1 vial. **TKAL2:** 2 vials.

Aldosterone Calibrators (ALC3–9)

Seven vials, labeled A through G, of lyophilized processed human serum. *At least 30 minutes prior to use*, reconstitute the zero calibrator **A with 6.0 mL** of distilled or deionized water, and each of the remaining calibrators **B through G with 3.0 mL** of distilled or deionized water. Use volumetric pipets and mix by *gentle* swirling. Store frozen: stable at –20°C for 30 days after reconstitution. For maximum integrity, aliquot and freeze.

TKAL1: 1 set. **TKAL2:** 1 set.

The reconstituted calibrators contain, respectively, 0, 25, 50, 100, 200, 600 and 1,200 picograms of aldosterone per milliliter (pg/mL) in processed human serum; equivalently: 0, 69, 139, 278, 555, 1,665 and 3,300 picomoles per liter (pmol/L).

Materials Required But Not Provided

Gamma counter — compatible with standard 12x75 mm tubes

Vortex mixer

Reagent Preparation

Distilled or deionized water

Graduated cylinder: 110 mL

Volumetric pipets: 3.0 mL and 6.0 mL

Urine Extraction

See Sample Preparation: Urine, below.

Radioimmunoassay

Plain 12x75 mm polypropylene tubes — for use as NSB tubes, available from DPC
Micropipets: 200 µL and 1.0 mL.

Foam decanting rack — available from DPC.

Logit-log graph paper — available from DPC (catalog number: ZP797)

A tri-level, human serum-based immunoassay control, containing aldosterone as one of over 25 constituents, is available from DPC (catalog number: CON6).

Specimen Collection

Serum:

The patient need not be fasting, and no special preparations are necessary. Collect blood by venipuncture⁶ into plain or heparinized tubes, noting the time of collection, and separate the serum from the cells by centrifugation.

Heparinized plasma yields results which are comparable to serum, but EDTA plasma yields results which are on average about 15% higher.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. Coat-A-Count Aldosterone has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 200 μ L per tube.

Storage: 2–8°C for 7 days, or 2 months frozen at –20°C.

Before assay, allow the samples to come to room temperature and mix by *gentle* swirling or inversion. Aliquot, if necessary, to avoid repeated thawing and freezing. Do not attempt to thaw frozen specimens by heating them in a waterbath.

Urine:

Obtain a 24-hour sample, keeping the specimen refrigerated during collection. Record the total volume. Unless the analysis is to be carried out immediately, reserve a well-mixed aliquot after adding 1 gram of boric acid per 100 mL urine.

Volume Required: 0.25 mL urine (for extraction procedure).

Storage: 7 days at 2–8°C or one month frozen at –20°C.

Sample Preparation: Urine

Precautions: Ethyl acetate is a volatile, flammable organic solvent. Conduct the

evaporation step under a fume hood equipped with an explosion-proof exhaust fan. Avoid open flames, and do not pipet by mouth. The ethyl acetate must be of at least spectrophotometric grade.

- 1 Label one glass or polypropylene tube (uncoated) for each urine sample.

The tubes should have tight-fitting caps and be able to withstand centrifuging at 1500xg

- 2 Pipet **250 μ L** of each urine sample into the appropriate tube.

If the sample is cloudy or if a precipitate has formed, first centrifuge the urine and work with the supernatant.

- 3 **Hydrolysis:** Add **25 μ L** of 3.2 N HCl (not supplied) to every tube. Cap securely and incubate for 24 hours at room temperature *in the dark*.

3.2 N HCl can be made by adding 1.0 mL concentrated HCl (12N) to 2.75 mL distilled water. Do not add water to concentrated acid, since this may cause splattering.

- 4 **Extraction:** Add **2.5 mL** ethyl acetate (not supplied) to every tube. Cap securely.

- 5 Mix by gentle inversion for **60 minutes**.

Use a mechanical rotator set at 15–20 revolutions-per-minute.

- 6 Centrifuge for **5 minutes** at about 1500xg, to separate the two layers.

Any sample partially emulsified should be shaken vigorously and centrifuged again.

- 7 **Evaporation:** Transfer exactly **50 μ L** of the upper (ethyl acetate) phase cleanly into one plain (uncoated) 12x75 mm polypropylene tube.

Suitable tubes are available from DPC. (Do not use polystyrene or antibody-coated tubes.) Pipet directly to the bottom of the tube using a positive-displacement micropipet. The remainder of the ethyl acetate phase may be retained for future use simply by freezing the extraction tube at –20°C; it is not necessary to separate the ethyl acetate phase from the aqueous phase.

- 8 Evaporate to complete dryness under a gentle stream of nitrogen at 37°C.
- 9 Add **0.5 mL** of the zero calibrator A. Thoroughly resuspend the extract by vortexing.

Proceed with the radioimmunoassay procedure, as described below. If not assayed immediately, the extracts should be stored tightly capped at –20°C, either dry (following step 8) or resuspended. Volumes may be scaled up, if desired, providing the proportions specified in the extraction procedure remain unchanged.

Radioimmunoassay Procedure

A single calibration curve provides the basis for determining aldosterone concentrations in serum and in urine. Accordingly, both urinary extracts and unextracted serum samples may be processed in the same assay. The calibrators, once reconstituted, are ready to use; only urine samples must be subjected to extraction, as described in the Sample Preparation: Urine section. All components must be at room temperature (15–28°C) before use.

- 1 **Plain Tubes:** Label four plain (uncoated) 12x75 mm polypropylene tubes T (total counts) and NSB (nonspecific binding) in duplicate.

Because nonspecific binding in the Coat-A-Count procedure is characteristically low, the NSB tubes may be safely omitted without compromising accuracy or quality control.

Coated Tubes: Label fourteen Aldosterone Ab-Coated tubes A (maximum binding) and B through G in duplicate. Label additional antibody-coated tubes, also in duplicate, for controls and patient samples.

Calibrators	pg/mL	ng/dL	pmol/L
A (MB)	0	0	0
B	25	2.5	69
C	50	5	139
D	100	10	278
E	200	20	555
F	600	60	1,665
G	1,200	120	3,300

- 2 Pipet **200 µL** of the zero calibrator A into the NSB and A tubes, and **200 µL** of each remaining calibrator, control and patient sample (serum) into the tubes prepared. Pipet **200 µL** of each urinary *extract* into appropriately labeled tubes.

Patient samples expected to contain aldosterone concentrations greater than the highest calibrator (1,200 pg/mL) should be diluted in the zero calibrator before assay.

It is good practice to use a disposable-tip micropipet, changing the tip between samples, in order to avoid carryover contamination.

- 3 Add **1.0 mL** of ¹²⁵I Aldosterone to every tube. Vortex briefly and gently.

Laboratories equipped with a reliable pipettor-diluter may handle steps 2 and 3 simultaneously, but care should be taken to avoid carryover from sample to sample. No more than 10 minutes should elapse during the dispensing of the tracer. Set the T tubes aside for counting (at step 6); they require no further processing.

- 4 Incubate for **18 hours** at room temperature (15–28°C).

- 5 Decant thoroughly.

Removing all visible moisture will greatly enhance precision. Using a foam decanting rack, decant the contents of all tubes (except the T tubes) and allow them to drain for 2 or 3 minutes. Then strike the tubes sharply on absorbant paper to shake off all residual droplets.

- 6 Count for **1 minute** in a gamma counter.

Calculation of Results

To obtain results in terms of concentration from a logit-log representation of the calibration curve, first calculate for each pair of tubes the average NSB-corrected counts per minute:

Net Counts = (Average CPM) minus (Average **NSB** CPM)

Then determine the binding of each pair of tubes as a percent of maximum binding (MB), with the **NSB**-corrected counts of the A tubes taken as 100%:

Percent Bound = (Net Counts / Net **MB** Counts) × 100

(The calculation can be simplified by omitting the correction for nonspecific binding; samples within range of the calibrators yield virtually the same results when Percent Bound is calculated directly from Average CPM.)

Using logit-log graph paper, plot Percent Bound on the vertical (probability) axis against Concentration on the horizontal (logarithmic) axis for each of the nonzero calibrators, and draw a straight line approximating the path of these points. Results for the unknowns may then be read from the line by interpolation.

Urine Samples: The result in "pg/mL" as read from the calibration curve must be multiplied by 100 to obtain the aldosterone concentration, in picograms per milliliter, of the original, unextracted urine sample. Divide this figure by 1,000, then multiply by the total volume *in liters*; to report the 24-hour aldosterone output in micrograms per day. (A correction factor of 100 is used because the urine samples are twice diluted 1-in-10: first by extracting 0.25 mL urine into 2.5 mL ethyl acetate, then by reconstituting the residue of 0.05 mL in 0.5 mL of zero calibrator. The addition of hydrochloric acid during the hydrolysis step has no effect on the dilution.)

Data Reduction Parameters

It is good practice to inspect results for agreement within replicates, and to construct a graph of the calibration curve (even if the calculations are handled by computer) as a visual check on the appropriateness of the transformation used and as a way to detect deviant calibration points. We also recommend keeping track of these data reduction parameters:

T = Total Counts (as counts per minute)

%NSB = $100 \times (\text{Average NSB Counts} / \text{Total Counts})$

%MB = $100 \times (\text{Net Counts} / \text{Total Counts})$

And the 20, 50 and 80 percent "intercepts," where

20% = Concentration at 20 Percent Bound, etc.

Note that other approaches, e.g. a sound implementation of the 4-parameter logistic, are also acceptable. See Dudley RA, et al.

Guidelines for immunoassay data reduction. Clin Chem 1985;31:1264-71.

Example Run: For illustration only, not for calculating results from another run. (See "Example Run" table.)

Quality Control

Controls or serum pools with at least two aldosterone concentration levels — low and high — should routinely be assayed as unknowns.

Report patient results only if the control results for the assay meet your laboratory's established criteria for acceptability.

It is good laboratory practice to record for each assay the lot numbers of the components used, as well as the dates when they were first reconstituted or opened. We also recommend charting control results from day to day — as described, for example, in Westgard JO, et al. A multi-rule chart for quality control. Clin Chem 1981;27:493-501. Note that repeat samples can serve as a valuable additional tool for monitoring interassay precision, and that pairs of control tubes can be spaced throughout the assay to help verify the absence of significant drift.

Solvent Blanks: For reliable results, the ethyl acetate must be of at least spectrophotometric grade and the nitrogen source, if used, must not contaminate the evaporation tubes with oil or moisture. (The use of a trap containing calcium sulfate will help prevent such contamination.) As a check on the purity of these reagents, which are not supplied with the kit, we recommend processing a "solvent blank" periodically. This involves extracting and assaying distilled or sterile water, in the manner prescribed for urine samples, to verify that it yields an apparent aldosterone concentration of essentially zero, that is, below the detection limit of the assay.

Expected Values

The following aldosterone levels for serum (or plasma) are based on normal sodium intake. Note that normal serum or plasma levels of aldosterone are dependent on the sodium intake and whether the patient is standing or recumbent. High sodium intake will tend to suppress serum

aldosterone, whereas low sodium intake will elevate serum aldosterone.

Standing	4–31 ng/dL
Recumbent	1–16 ng/dL

The following urinary aldosterone levels have been reported.

Normal Diet	6–25 µg/day
Low Salt	17–44 µg/day
High Salt	0–6 µg/day

A complete investigation of aldosterone pathophysiology must also include assessment of plasma renin/angiotensin activity.

Laboratories should consider these results as *guidelines only*. Each laboratory should establish its own reference ranges.

Performance Data

See Tables and Graphs for data *representative* of the Coat-A-Count Aldosterone kit's performance. Aldosterone results in the sections below are expressed as pg/mL. Unless otherwise specified, all results were generated on serum samples.

Conversion Factors:

pg/mL $\times$ 2.775 $\rightarrow$ pmol/L

pg/mL $\div$ 10 $\rightarrow$ ng/dL

Calibration Range: 25 – 1,200 pg/mL

Analytical Sensitivity: 11 pg/mL.

Intraassay Precision (Within-Run):

Statistics were calculated for each of six samples from the results of twenty pairs of tubes in a single assay. (See "Intraassay Precision" table.)

Interassay Precision (Run-to-Run):

Statistics were calculated for each of six samples from the results of pairs of tubes in 20 different assays. (See "Interassay Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

End-of-Run Effect: None up to approximately 400 tubes. (See "End-of-Run Effect" table.)

Recovery: Samples spiked 1 to 19 with three aldosterone solutions (4,576, 8,347

and 15,071 pg/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antiserum for the Coat-A-Count Aldosterone kit is very specific for aldosterone, with an extremely low crossreactivity to other compounds that might be found in patient samples. Crossreactivities were calculated on a weight-per-weight basis at approximately 50% binding. (See "Specificity" table for representative data.)

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of packed red blood cells in concentrations up to 30 µL/mL has no effect on results, within the precision of the assay.

Lipemia: The 600 pg/mL calibrator was serially diluted with each of two lipemic serum pools, causing the degree of lipemia to increase as the aldosterone concentration (due primarily to the calibrator) decreases. The dilutions were assayed by the Coat-A-Count Aldosterone procedure along with the unspiked lipemic pools. The results show excellent recoveries, even in the presence of severe lipemia.

Alternate Sample Type: To determine whether anticoagulants interfere with the assay, blood was collected from 47 volunteers into plain, heparinized and EDTA vacutainer tubes. All samples were assayed by the Coat-A-Count Aldosterone procedure. The results were subjected to linear regression analysis.

(Heparin) = 0.96 (Serum) – 3.0 pg/mL
 $r = 0.978$

Means:
175 pg/mL (Heparin)
185 pg/mL (Serum)

(EDTA) = 1.04 (Serum) + 22 pg/mL
 $r = 0.985$

Means:
214 pg/mL (EDTA)
185 pg/mL (Serum)

The results show that serum and heparinized plasma yield comparable results, whereas EDTA plasma yields results which are on average about 15% higher.

Protein Effect: To simulate various protein concentrations, experiments were

performed in which 2.0 mL aliquots of a human serum pool were freeze-dried and then reconstituted with various volumes of water (1.0, 1.3, 2.0, 3.0 and 4.0 mL). Each reconstituted aliquot was then assayed by the Coat-A-Count Aldosterone procedure. Note that aliquots reconstituted with half the original volume represent an extremely high protein concentration, on the order of 14 g/dL. (The factor to correct for reconstitution volume is tabulated below.) These results indicate that even wide variations in protein content have relatively little effect on the Coat-A-Count Aldosterone assay. (See "Protein Effect" table.)

Method Comparison 1: Two hundred and thirty-six serum samples, with aldosterone concentrations ranging from 27 to 860 pg/mL, were simultaneously assayed by the Coat-A-Count Aldosterone current procedure (18-hour incubation at room temperature) and by its predecessor (3-hour incubation at 37°C). (see "Method Comparison 1" graph) By linear regression.

(Current) = 1.0 (Previous) – 8.4 pg/mL
r = 0.985

Means:

130 pg/mL (Current)
138 pg/mL (Previous)

Method Comparison 2: In another study, 84 urine samples, with aldosterone concentrations ranging from 26 to 228 pg/mL, were simultaneously assayed by the Coat-A-Count Aldosterone current procedure (18-hour incubation at room temperature) and by its predecessor (3-hour incubation at 37°C). (see "Method Comparison 2" graph) By linear regression.

(Current) = 1.0 (Previous) + 5.9 pg/mL
r = 0.953

Means:

85 pg/mL (Current)
77 pg/mL (Previous)

References

1) Mayes D, et al. Radioimmunoassay for plasma aldosterone. J Clin Endocrinol Metab 1970;30:682-5. 2) Demers LM, et al. Plasma and urinary aldosterone measurement in healthy subjects with a radioimmunoassay kit not requiring chromatography. Clin Biochem 1976;9(5):243-6. 3) Tan SY, et al. Direct non-chromatographic radioimmunoassay of aldosterone: validation of a commercially available kit and observations on age-related changes in concentrations in plasma. Clin Chem 1978;24:1531-3. 4) Walsh P, et al. Non-chromatographic radioimmunoassay procedure for urinary aldosterone. Clin Chem 1979;25:1226-9. 5) Kubasik N, et al. Evaluation of a new commercial radioassay kit for aldosterone using an iodinated tracer. Clin Biochem 1979;12:59-61. 6) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782, 310.645.8200
Fax: 310.645.9999. To place an order:
Tel: 800.372.1782. Fax: 800.234.4372.
Outside the United States, contact your National Distributor.

The Quality System of Diagnostic Products Corporation is registered to ISO 13485:2003.

Tables and Graphs

Example Run

Tube ¹	Duplicate CPM ²	Average CPM ³	Net CPM ⁴	Percent Bound ⁵	Aldosterone pg/mL ⁶
T ⁷	39,951 39,910	39,931			
NSB⁸	498 586	542	0		
A (MB) ⁹	20,576 20,003	20,290	19,748	100%	0
B	17,440 17,750	17,595	17,053	86%	25
C	15,591 16,105	15,848	15,306	78%	50
D	13,317 12,983	13,150	12,608	64%	100
E	9,998 10,177	10,088	9,546	48%	200
F	5,468 5,251	5,360	4,818	24%	600
G	3,613 3,419	3,516	2,974	15%	1,200
Unknowns ¹⁰					
X1	15,808 15,640	15,724	15,182	77%	51
X2	11,370 11,826	11,598	11,056	56%	143
X3	6,697 6,640	6,669	6,127	31%	433

Quality Control Parameters:¹¹

T⁷ = 39,931 cpm

%NSB⁸ = 1.4%

%MB⁹ = 49%

20% Intercept¹² = 810 pg/mL

50% Intercept = 184 pg/mL

80% Intercept = 42 pg/mL

Intraassay Precision (pg/mL)

	Mean ¹	SD ²	CV ³
1	65	3.5	5.4%
2	174	5.5	3.2%
3	243	8.2	3.4%
4	448	11.7	2.6%
5	575	16.8	2.9%
6	813	19.0	2.3%

Interassay Precision (pg/mL)

	Mean ¹	SD ²	CV ³
1	58	9.1	15.7%
2	59	7.7	13.1%
3	162	10.5	6.5%
4	226	15.3	6.8%
5	447	19.4	4.3%
6	548	21.0	3.8%

Linearity (pg/mL)



	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	16 in 16 ⁵	635	—	—
	8 in 16	319	318	100%
	4 in 16	167	159	105%
	2 in 16	92	79	116%
	1 in 16	40	40	100%
2	16 in 16	642	—	—
	8 in 16	331	321	103%
	4 in 16	184	161	114%
	2 in 16	95	80	119%
	1 in 16	44	40	110%
3	16 in 16	750	—	—
	8 in 16	378	375	101%
	4 in 16	210	188	112%
	2 in 16	104	94	111%
	1 in 16	51	47	109%
4	16 in 16	1,008	—	—
	8 in 16	540	504	107%
	4 in 16	292	252	116%
	2 in 16	135	126	107%
	1 in 16	73	63	116%

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O Sistema de Qualidade da Diagnostic Products Corporation está registado sob ISO 13485:2003.
EC REP DPC Biemann GmbH