Cortisol

DPC®

Coat-A-Count® Cortisol

English

Intended Use: Coat-A-Count Cortisol is a radioimmunoassay designed for the quantitative measurement of cortisol (hydrocortisone, Compound F) in serum, urine and heparinized plasma. It is intended strictly for *in vitro* diagnostic use as an aid in the clinical assessment of adrenal status.

Catalog numbers: **TKCO1** (100 tubes), **TKCO2** (200 tubes), **TKCO5** (500 tubes)

The 100-tube kit contains less than 5.5 microcuries (204 kilobecquerels) of radioactive ^{125}I cortisol; the 200-tube kit contains less than 11 microcuries (407 kilobecquerels); and the 500-tube kit contains less than 27.5 microcuries (1018 kilobecquerels).

Summary and Explanation

Cortisol (hydrocortisone, Compound F) is the most abundant circulating steroid and the major glucocorticoid secreted by the adrenal cortex. Physiologically effective in anti-inflammatory activity and blood pressure maintenance, cortisol is also involved in gluconeogenesis, calcium absorption and the secretion of gastric acid and pepsin.

As an indicator of adrenocortical function, measurement of blood cortisol levels is useful in the differential diagnosis of Addison's and Cushing's disease, hypopituitarism, and adrenal hyperplasia and carcinoma.⁴

A variety of stimulation and suppression tests — ACTH Stimulation, ACTH Reserve, Dexamethasone Suppression — can supply supportive information on adrenocortical function.^{4,5,6} Anomalous cortisol concentrations have been shown to exist in patients with acute infections, severe pain, diabetes mellitus or heart failure, and in women either pregnant or on estrogen therapy. In addition, certain virilizing syndromes and iatrogenic conditions raise other naturally occurring steroid levels to high (nonphysiologic) concentrations that may interfere with

assaying cortisol directly, unless a highly specific antiserum is used.

Principle of the Procedure

The Coat-A-Count Cortisol procedure is a solid-phase radioimmunoassay, wherein ^{125}I -labeled cortisol competes for a fixed time with cortisol in the patient sample for antibody sites. Because the antibody is immobilized to the wall of a polypropylene tube, simply decanting the supernatant suffices to terminate the competition and to isolate the antibody-bound fraction of the radiolabeled cortisol. Counting the tube in a gamma counter then yields a number, which converts by way of a calibration curve to a measure of the cortisol present in the patient sample.

Reagents to Pipet: 1

Total Incubation Time: 45 minutes.

Total Counts at Iodination: approximately 90,000 cpm.

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C in a refrigerator designated for incoming radioactive materials. Dispose of in accordance with applicable laws.

Do not use reagents beyond their expiration dates.

Some components supplied in this kit may contain human source material and/or other potentially hazardous ingredients which necessitate certain precautions:

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Water: Use distilled or deionized water.

Radioactivity

A copy of any radioisotope license certificate (Specific or General) issued to a US customer must be on file with Diagnostic Products Corporation before kits or components containing radioactive material can be shipped. These radioactive materials may be acquired by any customer with the appropriate Specific license. Under a General license these radioactive materials may be acquired only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories and hospitals — and strictly for *in vitro* clinical or laboratory tests not involving external or internal administration of the radioactive material or its radiation to human beings or other animals. Its acquisition, receipt, storage, use, transfer and disposal are all subject to the regulations and a (General or Specific) license of the U.S. Nuclear Regulatory Commission or a State with which the NRC has entered into an agreement for the exercise of regulatory control.

Handle radioactive materials according to the requirements of your General or Specific license. To minimize exposure to radiation, the user should adhere to guidelines set forth in the National Bureau of Standards publication on the *Safe Handling of Radioactive Materials* (Handbook No. 92, issued March 9, 1964) and in subsequent publications issued by State and Federal authorities.

Wipe up spills promptly and decontaminate affected surfaces. Avoid generation of aerosols. Dispose of solid radioactive waste according to license requirements. General licensees (holders of NRC Form 483) may dispose of solid radioactive waste as nonradioactive waste, after removing labeling. Specific licensees (NRC Form 313) should refer to Title 10, Code of Federal Regulations, Part 20. Licensees in Agreement States should refer to the appropriate regulations of their own state. General licensees may dispose of liquid radioactive waste of the type contained in this product through a laboratory sink drain. Licensees must remove or deface labels from empty containers of radioactive materials before disposal of solid waste. Specific licensees may dispose of small quantities of liquid radioactive waste of the type used in this

product through a laboratory sink drain. Refer to the appropriate regulations applicable to your laboratory.

Materials Supplied: Initial Preparation

Cortisol Ab-Coated Tubes (TKCO1)

Polypropylene tubes coated with antibodies to cortisol and packaged in zip-lock bags. Protect from moisture, carefully resealing the bags after opening: stable at 2–8°C until the expiration date marked on the bag. Color: light blue.

TKCO1: 100 tubes. **TKCO2:** 200 tubes.

TKCO5: 500 tubes.

¹²⁵I Cortisol (TKCO2)

Iodinated cortisol with blocking agents for steroid-binding proteins, supplied in liquid form, ready to use. Each vial contains 105 mL. Stable at 2–8°C for at least 30 days after opening, or until the expiration date marked on the vial.

TKCO1: 1 vial. **TKCO2:** 2 vials.

TKCO5: 5 vials.

Cortisol Calibrators (TKCO3–8)

Cortisol calibrators, ready to use, in processed human serum. The zero calibrator vial A contains 2.0 mL, while the remaining calibrator vials B through F contain 1.0 mL each. Stable at 2–8°C for at least 30 days after opening, or until the expiration date marked on the vial. The life of the calibrators may be extended by freezing. Aliquot, if necessary, to avoid repeated freezing and thawing.

TKCO1: 1 set. **TKCO2:** 1 set.

TKCO5: 2 sets.

The calibrators contain, respectively, 0, 1, 5, 10, 20 and 50 micrograms of cortisol per deciliter (µg/dL) in processed human serum; equivalently, 0, 27.6, 138, 276, 552 and 1,380 nanomoles per liter (nmol/L). Intermediate calibration points can be obtained by mixing the calibrators in suitable proportions.

Materials Required But Not Provided

Gamma counter – compatible with standard 12×75 mm tubes

Vortex mixer

Urine Extraction

Tubes with caps – suitable for extraction with dichloromethane

Dichloromethane – spectrophotometric grade (also called methylene chloride)

Micropipets: 50 μ L and 500 μ L. The 50 μ L pipet should be of the positive-displacement type.

Evaporator with compressed nitrogen

Radioimmunoassay

Plain 12 $\times$ 75 mm polypropylene tubes — for use as NSB tubes, available from DPC

Micropipets: 25 μ L and 1.0 mL.

Waterbath — capable of maintaining 37°C. Neither an oven nor a heat block is suitable.

Foam decanting rack – available from DPC

Logit-log graph paper – available from DPC (catalog number: ZP797)

A tri-level, human serum-based immunoassay control, containing cortisol as one of over 25 assayed constituents, is available from DPC (catalog number: CON6).

Specimen Collection

Serum and Plasma

The patient need not be fasting, and no special preparations are necessary. Collect blood by venipuncture¹⁰ into plain (without anticoagulant) or heparinized tubes, and separate the serum or plasma from the cells. EDTA plasma is *not* acceptable. (See Alternate Sample Type section). Since cortisol exhibits a circadian rhythm, the time of collection should be noted.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. Coat-A-Count Cortisol has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 25 μ L of serum or plasma per tube.

Storage: 7 days at 2–8°C or 2 months at –20°C.⁸ Before assay, allow samples to come to room temperature (15–28°C) and mix by *gentle* swirling or inversion. Aliquot, if necessary, to avoid repeated thawing and freezing. Do *not* attempt to thaw specimens by heating them in a waterbath. Patient samples expected to contain cortisol concentrations greater than the highest calibrator (50 μ g/dL) should be diluted in the zero calibrator before assay.

Urine

Collect a 24-hour urinary sample without preservative, keeping the specimen refrigerated during collection. Record the total volume in milliliters.

Volume Required (extraction procedure): 0.5 mL of urine.

Storage: 7 days at 2–8°C or one month at –20°C.⁸

Preparation of Urine Samples

Warning – Dichloromethane: Do not pipet by mouth. Conduct the evaporation step under a fume hood. Avoid open flames.

- 1 Label one 12 $\times$ 75 mm polypropylene tube for each urine sample.

Suitable tubes with tightly fitting caps are available from DPC. Do not use antibody-coated tubes at this stage.

- 2 Pipet 500 μ L of each urine sample into the tube prepared.

If the sample is cloudy, or if a precipitate has formed, first centrifuge the urine and work with the supernatant.

- 3 **Extraction:** Add 1.0 mL dichloromethane to every tube. Cap securely.

Do not pipet dichloromethane by mouth.

- 4 Mix by *gentle* inversion for 5 to 10 minutes.

Use a mechanical rotator set at 15 to 20 revolutions per minute.

- 5 Centrifuge for 5 minutes at about 1500 $\times$ g to separate the two phases.

Any sample partially emulsified should be shaken vigorously and centrifuged again.

- 6 Remove most of the top layer, being careful not to disturb the interface.

Use a water-tap aspirator. (Aspirating all the way to the interface may cause the dichloromethane layer to become contaminated with material from the aqueous phase on contact.)

- 7 **Evaporation:** Transfer **50 µL** of the lower (dichloromethane) phase into each of two Cortisol Ab-Coated Tubes.

Pipet all extracts directly to the bottom using a positive-displacement micropipet. Be careful that none of the upper phase enters the assay tube.

- 8 Evaporate to complete dryness under a gentle stream of nitrogen at room temperature (15–28°C).

Conduct the evaporation step under a fume hood. Avoid open flames.

Once dry, let the tubes (suitably labeled) take their place in an assay that includes a set of (serum) calibrator tubes. Continue – step 2 of the radioimmunoassay procedure – by adding 25 µL of the zero calibrator A, and so on. Results in micrograms of (free) cortisol per 100 mL urine (µg/dL) can be read off directly from the calibration curve. No correction factor is necessary.

No more than 30 minutes should elapse between transfer of the extract into the coated tube (step 7 of the extraction procedure) and the addition of tracer (step 3 of the radioimmunoassay procedure). The extraction procedure may be interrupted after step 6 by storing the infranatant layer tightly capped in a –20°C freezer.

Radioimmunoassay Procedure

All components must be at room temperature (15–28°C) before use.

- 1 **Plain Tubes:** Label four plain (uncoated) 12×75 mm polypropylene tubes T (total counts) and NSB (nonspecific binding) in duplicate.

Because nonspecific binding in the Coat-A-Count procedure is characteristically low, the NSB tubes may be safely omitted without

compromising accuracy or quality control.

Coated Tubes: Label twelve Cortisol Ab-Coated Tubes A (maximum binding) and B through F in duplicate. Label additional Cortisol Ab-Coated Tubes, also in duplicate, for controls and patient samples.

Calibrators	µg/dL	nmol/L
A (MB)	0	0
B	1	27.6
C	5	138
D	10	276
E	20	552
F	50	1,380

Urine samples: Once dry, the coated tubes used for evaporating the extracts should be suitably labeled and incorporated into the sequence of patient sample tubes.

- 2 Pipet **25 µL** of the zero calibrator A into the NSB and A tubes, *and also into each of the urine sample evaporation tubes*. Pipet **25 µL** of each remaining calibrator, control and patient sample into the tubes prepared. **Pipet directly to the bottom.**

It is good practice to use a disposable-tip micropipet, changing the tip between samples, in order to avoid carryover contamination.

- 3 Add **1.0 mL** of ¹²⁵I Cortisol to every tube. Vortex.

Laboratories equipped with a reliable pipettor-diluter may handle steps 2 and 3 simultaneously. No more than 10 minutes should elapse during the dispensing of the tracer. Set the T tubes aside for counting at step 6; they require no further processing.

- 4 Incubate for 45 minutes at 37°C.

Use a waterbath; neither an oven nor a heat block is suitable. Longer incubation periods will not significantly affect the assay.

- 5 Decant thoroughly.

Removing all visible moisture will greatly enhance precision. Decant the contents of all tubes (except the T

tubes) using a foam decanting rack, and allow them to drain for 2 or 3 minutes. Then strike the tubes sharply on absorbant paper to shake off all residual droplets.

- 6 Count for **1 minute** in a gamma counter.

Calculation of Results

To obtain results in terms of concentration from a logit-log representation of the calibration curve, first calculate for each pair of tubes the average NSB-corrected counts per minute:

Net Counts = (Average CPM) minus (Average NSB CPM)

Then determine the binding of each pair of tubes as a percent of maximum binding (MB), with the NSB-corrected counts of the A tubes taken as 100%:

Percent Bound = (Net Counts / Net MB Counts) × 100

The calculation can be simplified by omitting the correction for nonspecific binding; samples within range of the calibrators yield virtually the same results when Percent Bound is calculated directly from Average CPM.)

Using logit-log graph paper, plot Percent Bound on the vertical (probability) axis against Concentration on the horizontal (logarithmic) axis for each of the nonzero calibrators, and draw a straight line approximating the path of these points. Results for the unknowns may then be read from the line by interpolation.

It is good practice to inspect results for agreement within replicates, and to construct a graph of the calibration curve (even if the calculations are handled by computer) as a visual check on the appropriateness of the transformation used and as a way to detect deviant calibration points. We also recommend keeping track of these data reduction parameters:

T = Total Counts (as counts per minute)

%NSB = $100 \times (\text{Average NSB Counts} / \text{Total Counts})$

%MB = $100 \times (\text{Net Counts} / \text{Total Counts})$

And the 20, 50 and 80 percent "intercepts," where

20% = Concentration at 20 Percent Bound, etc.

Note that other approaches, e.g. a sound implementation of the 4-parameter logistic, are also acceptable. See Dudley RA, et al. Guidelines for immunoassay data reduction. Clin Chem 1985;31:1264-71.

Urine Samples:

Urinary free cortisol concentrations may be read off directly from the calibration curve. No correction factor is necessary. (Introducing *twice* the usual 25 µL sample into the assay tube at step 7 of the extraction procedure exactly compensates for the 1-in-2 dilution occasioned by extracting 0.5 mL urine into 1.0 mL dichloromethane.)

To determine the urinary free cortisol output in micrograms per day, simply multiply the result in µg/dL by the 24-hour urine volume in deciliters:

$\mu\text{g}/24 \text{ hours} = \mu\text{g}/\text{dL} \times (\text{Total Volume in milliliters}) / 100$

Example Run: For illustration only, not for calculating results from another run. (See "Example Run" table.)

Quality Control

Controls (or serum pools) with at least two cortisol concentration levels — low and high — should routinely be assayed as unknowns.

Report patient results only if the control results for that assay meet your laboratory's established criteria for acceptability.

It is good laboratory practice to record for each assay the lot numbers of the components used, as well as the dates when they were first reconstituted or opened. We also recommend charting control results from day to day — as described, for example, in Westgard JO, et al. A multi-rule chart for quality control. Clin Chem 1981;27:493-501. Note that repeat samples can serve as a valuable additional tool for monitoring interassay precision, and that pairs of control tubes can be spaced throughout the assay to help verify the absence of significant drift.

Solvent Blanks:

For reliable results, the dichloromethane must be of at least spectrophotometric grade and the nitrogen source, if used, must not contaminate the evaporation tubes with oil or moisture. (The use of a trap containing calcium sulfate will help prevent such contamination.) As a check on the purity of these reagents, which are not supplied with the kit, we recommend periodically running a "solvent blank." This involves extracting and assaying distilled or sterile water, in the manner prescribed for urine samples, to verify that it yields an apparent cortisol concentration of essentially zero, that is, below the detection limit of the assay.

Expected Values

Serum

The following expected values for circulating cortisol have been reported.

Diurnal Rhythm Variation	<i>a.m.</i> : 5 – 25 µg/dL (138 – 690 nmol/L) ^{1,8,9} <i>p.m.</i> : Approximately half of <i>a.m.</i> values
ACTH Stimulation	Over twice (usually 3 to 5 times) basal values ^{4,5}
ACTH Reserve (Metyrapone Suppression)	Below basal (control) values ^{5,6}
Dexamethasone Suppression	Below basal (control) values for the screening, low-dose and high-dose tests ^{4,5}

Urine

A study of 24-hour urine samples from 33 healthy laboratory volunteers yielded an absolute range of 21 to 85 µg/24 hours, with a median of 48 µg/24 hours. This is consistent with a reference range of 20 to 90 µg/24 hours for urinary free cortisol, as reported in the literature.⁷

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitation

Serum and urinary cortisol results may be falsely elevated in samples obtained from patients being treated with prednisolone or prednisone (converted to prednisolone *in*

vivo). Caution must therefore be exercised with cortisol determinations for patients undergoing therapy with these and structurally related synthetic corticosteroids.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed as micrograms of cortisol per deciliter (µg/dL). All results were obtained on serum samples.

Conversion Factor:

µg/dL × 27.59 → nmol/L

Calibration Range: 1 – 50 µg/dL
(27.6 – 1,380 nmol/L)

Analytical Sensitivity: 0.2 µg/dL
(5.5 nmol/L)

Intraassay Precision (Within-Run):

Statistics were calculated for each of seven samples from the results of 20 pairs of tubes in a single assay. (See "Intraassay Precision" table.)

Interassay Precision (Run-to-Run):

Statistics were calculated for each of three samples from the results of 20 pairs of tubes. (See "Interassay Precision" table.)

End-of-Run Effect: none up to approximately 700 tubes. (See "End-of-Run Effect" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table.)

Recovery: Samples spiked 1 to 19 with three cortisol solutions (130, 250 and 475 µg/dL) were assayed. (See "Recovery" table.)

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of packed red blood cells in concentrations up to 30 µL/mL has no effect on results, within the precision of the assay.

Lipemia: The 50 µg/dL calibrator was serially diluted with a lipemic serum pool. This causes the degree of lipemia to increase as the cortisol concentration (due primarily to the calibrator) decreases. All dilutions were assayed along with the unspiked lipemic pool. The results show good recoveries, even in the presence of severe lipemia.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 40 volunteers into plain, heparinized, EDTA and Becton Dickinson SST[®] vacutainer tubes. All samples were assayed by the Coat-A-Count Cortisol procedure. By linear regression:

(Heparin) = 0.93 (Serum) + 0.59 µg/dL
r = 0.988

(EDTA) = 1.17 (Serum) – 0.03 µg/dL
r = 0.998

(SST) = 0.94 (Plain Tubes) + 0.70 µg/dL
r = 0.995

Means:

13.8 µg/dL (Serum)
13.5 µg/dL (Heparin)
16.2 µg/dL (EDTA)
13.7 µg/dL (SST)

The results show that EDTA plasma is *not* acceptable.

Specificity: The Coat-A-Count Cortisol antiserum is highly specific for cortisol, with an extremely low crossreactivity to other naturally occurring steroids or therapeutic drugs that may be present in patient samples. (See "Specificity" table.)

Although some steroids exhibit slight crossreactivity, their normal physiological concentrations are low compared to cortisol, hence they will not significantly interfere in the Coat-A-Count Cortisol procedure. There is, however, an observed crossreactivity of approximately 76% with prednisolone, which must be considered for patients undergoing this therapy. Since prednisone is converted to prednisolone *in vivo*, caution must be exercised with cortisol determinations for patients undergoing either therapy.

Method Comparison: The assay was compared to DPC's Double Antibody Cortisol on 89 patient samples. By linear regression:

(CAC) = 1.02 (DA) – 0.25 µg/dL
r = 0.970

Means:

18.8 µg/dL (Coat-A-Count)
18.7 µg/dL (Double Antibody)

Protein Effect: To simulate various protein concentrations, experiments were performed in which 2.0 mL aliquots of a human serum pool were freeze-dried and then reconstituted with various volumes of water (1, 2, 3 and 4 mL). Each

reconstituted aliquot was then assayed by the Coat-A-Count Cortisol procedure. Note that aliquots reconstituted with half the original volume represent an extremely high protein concentration, on the order of 14 g/dL. These results indicate that even wide variations in total protein have no clinically significant effect on the Coat-A-Count Cortisol assay. (See "Protein Effect" table for representative data.)

References

Foster L, Dunn R. Single-antibody technique for radioimmunoassay of cortisol in unextracted serum or plasma. Clin Chem 1974;20:365. 2) Farmer R, Pierce C. Plasma cortisol determination: radioimmunoassay and competitive protein binding compared. Clin Chem 1974;20:411. 3) Ruder H, et al. A radioimmunoassay for cortisol in plasma and urine. J Clin Endo Metab 1972;35:219. 4) Rothfeld B, ed. Plasma cortisol. In: Nuclear medicine *in vitro*. 1974: 120. 5) Murphy B, et al. Clinical studies utilizing a new method for the serial determination of plasma corticoids. Can Med Assoc J 1964;90:775. 6) Spark R. Measurement of serum 11-deoxycortisol and cortisol after metyrapone. Ann Intern Med 1971;75:717. 7) Kowalski A, Paul W. A simple extraction procedure for the determination of free (unconjugated) cortisol in urine by radioimmunoassay. Clin Chem 1979;25:1152. 8) Burtis CA, Ashwood ER, editors. Tietz textbook of clinical chemistry. 2nd ed. Philadelphia: W.B. Saunders, 1994. 9) Wilson JD, Foster DW, editors. Williams textbook of endocrinology. 7th ed. Philadelphia: W.B. Saunders, 1985. 10) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782, 310.645.8200
Fax: 310.645.9999. To place an order:
Tel: 800.372.1782. Fax: 800.234.4372.
Outside the United States, contact your National Distributor.

The Quality System of Diagnostic Products Corporation is registered to ISO 13485:2003.

Tables and Graphs

Example Run

Tube ¹	Duplicate CPM ²	Average CPM ³	Net CPM ⁴	Percent Bound ⁵	Cortisol µg/dL ⁶
T ⁸	84,250 83,899	84,075			
(NSB) ⁹	519 415	467	0		
A	46,585				
(MB) ¹⁰	46,351	46,468	46,001	100%	0
B	38,113 38,034	38,074	37,607	82%	1
C	26,346 26,661	26,504	26,037	57%	5
D	19,778 20,113	19,946	19,479	42%	10
E	14,049 13,851	13,950	13,483	29%	20
F	7,548 7,739	7,644	7,177	16%	50
Unknowns ¹¹					
X1	26,021 26,699	26,360	25,893	56%	4.9
X2	18,879 18,896	18,888	18,421	40%	11
X3	13,402 14,483	13,943	13,476	29%	20

Quality Control Parameters:⁷

T⁸ = 84,075 cpm

%NSB⁹ = 0.56%

%MB¹⁰ = 55%

20% Intercept¹¹ = 37 µg/dL

50% Intercept = 6.7 µg/dL

80% Intercept = 1.2 µg/dL

Intraassay Precision (µg/dL)

	Mean	SD	CV
1	3.1	0.15	4.8%
2	4.5	0.19	4.2%
3	9.9	0.47	4.7%
4	12	0.61	5.1%
5	20	0.60	3.0%
6	25	1.0	4.0%
7	34	1.5	4.4%

Interassay Precision (µg/dL)

	Mean ¹	SD ²	CV ³
1	3.3	0.17	5.2%
2	21	0.83	4.0%
3	36	2.3	6.4%

End-of-run Effect (µg/dL)

	Tubes ¹ 84–89	Tubes 420–425	Tubes 700–705
1	3.1	3.0	3.0
2	12.1	11.9	12.0
3	26.2	25.3	25.5

Linearity (µg/dL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	16	—	—
	4 in 8	8.0	8.0	100%
	2 in 8	4.0	4.0	100%
	1 in 8	1.9	2.0	95%
2	8 in 8	19	—	—
	4 in 8	9.4	9.5	99%
	2 in 8	4.7	4.8	98%
	1 in 8	2.3	2.4	96%
3	8 in 8	20	—	—
	4 in 8	9.3	10	93%
	2 in 8	4.8	5.0	96%
	1 in 8	2.5	2.5	100%
4	8 in 8	27	—	—
	4 in 8	14	14	100%
	2 in 8	6.7	6.8	99%
	1 in 8	3.3	3.4	97%
5	8 in 8	31	—	—
	4 in 8	16	16	100%
	2 in 8	7.9	7.8	101%
	1 in 8	3.8	3.9	97%
6	8 in 8	38	—	—
	4 in 8	19	19	100%
	2 in 8	8.8	9.5	93%
	1 in 8	4.4	4.8	92%

Recovery (µg/dL)

	Spiking Solution ¹	Observed ²	Expected ³	% O/E ⁴
1	—	11	—	—
	A	17	17	100%
	B	23	23	100%
	C	33	34	97%
2	—	18	—	—
	A	22	24	92%
	B	29	30	97%
	C	39	41	95%
3	—	13	—	—
	A	19	19	100%
	B	25	25	100%
	C	36	36	100%
4	—	9.1	—	—
	A	15	15	100%
	B	20	22	91%
	C	30	33	91%

Specificity


Compound ¹	µg/dL Added ²	Apparent Conc. ³ µg/dL	Percent Cross- reactivity ⁴
Aldosterone	2,000	0.6	0.03%
Betamethasone	1,000	1.6	1.6%
Corticosterone	500	4.7	0.94%
Cortisone	500	4.9	0.98%
Danazol	2,000	0.2	0.01%
11-Deoxycorticosterone	2,000	5.1	0.26%
11-Deoxycortisol	50	5.7	11.4%
Dexamethasone	10,000	4.2	0.04%
Estradiol	2,000	0.2	0.01%
Estrone	4,000	0.3	0.007%
Fluticasone	22	ND	ND
Flumethasone	10,000	1.7	0.017%
Methotrexate	10,000	0.4	0.004%
Methylprednisolone	100	12	12%
	200	21	11%
Prednisolone	8	6.1	76%
Prednisone	250	5.8	2.3%
Pregnenolone	2,000	0.4	0.02%
Pregnanediol	2,000	ND	—
Progesterone	10,000	1.8	0.02%
Tetrahydrocortisol	1,000	3.4	0.34%
Triamcinolone	500	0.65	0.13%
	5,000	2.26	0.045%

ND: not detectable.⁵

Protein Effect (µg/dL)

	Recon- stitution ¹	Approx. Conc. ² g/dL	Observed ³	Expected ⁴	% O/E ⁵
1	2.00×	14.0	27.3	29.2	93%
	1.00×	7.0	14.6	—	(100%)
	0.67×	4.7	10.4	9.7	107%
	0.50×	3.5	8.1	7.3	111%
2	2.00×	14.0	32.9	34.8	95%
	1.00×	7.0	17.4	—	(100%)
	0.67×	4.7	10.8	11.6	93%
	0.50×	3.5	8.2	8.7	94%

Procedure for Serum Samples
Refer to Radioimmunoassay Procedure for details.

Add 25 µL Calibrator, Control or Sample. Add 1.0 mL ¹²⁵I Cortisol. Vortex. Incubate 45 minutes at 37°C. Decant and count.

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O Sistema de Qualidade da Diagnostic Products Corporation está registado sob ISO 13485:2003.
DPC Biemann GmbH