DHEA

DPC®

Coat-A-Count® DHEA

English

Intended Use: Coat-A-Count DHEA is a solid-phase ^{125}I radioimmunoassay designed for the quantitative measurement of dehydroepiandrosterone (DHEA) in serum. It is intended strictly for *in vitro* diagnostic use.

Catalog Numbers: **TKDH1** (75 tubes), **TKDH2** (150 tubes).

CDC Analyte Identifier Code: **1309**

CDC Test System Identifier Code: **13030**

CLIA Complexity Category: **High**

The 75-tube kit contains less than 5 microcuries (185 kilobecquerels) of radioactive ^{125}I DHEA, and the 150-tube kit contains less than 10 microcuries (370 kilobecquerels).

Summary and Explanation of the Test

Dehydroepiandrosterone (DHEA) is an adrenal steroid which often circulates at an elevated level in cases of abnormal hair growth (hirsutism) and virilization. Its measurement may also be of value in the assessment of adrenarche and delayed puberty.

The DHEA in circulation originates almost entirely from the adrenals. Although produced at a rate only somewhat lower than dehydroepiandrosterone sulfate (DHEA- SO_4), with which it is peripherally interconvertible, DHEA has a much shorter half-life and so maintains a plasma concentration almost one thousandfold lower. In itself, DHEA is only weakly androgenic, but it can serve as a precursor to more potent androgens.

Plasma levels of DHEA increase steadily from about the seventh year, then gradually decline after the third decade. According to the literature, the upper reference limit for young adults is approximately 10 ng/mL. Like cortisol, DHEA exhibits a diurnal variation, being at its highest in the morning. On the other hand, it shows no significant variation throughout the menstrual cycle. There is a

moderate decrease during pregnancy. Unlike testosterone, DHEA does not circulate bound to sex hormone-binding globulin and therefore is not influenced by alterations in the level of this carrier protein.

Principle of the Procedure

The Coat-A-Count DHEA procedure is a competitive immunoassay based on antibody-coated tubes. Patient samples are first extracted into dichloromethane, evaporated to dryness and resuspended in the zero calibrator. ^{125}I -labeled DHEA then competes with DHEA in the patient extracts for antibody sites. After the tubes are decanted and counted, DHEA concentrations are read directly from a calibration curve.

Reagents to Pipet: 1

Total Incubation Time: 2 Hours

Total Counts at lodination:
approximately 100,000 cpm

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C in a refrigerator designated for incoming radioactive materials. Dispose of in accordance with applicable laws.

Do not use reagents beyond their expiration dates.

Some components supplied in this kit may contain human source material and/or other potentially hazardous ingredients which necessitate certain precautions.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Water: Use distilled or deionized water.

Radioactivity

A copy of any radioisotope license certificate (Specific or General) issued to a US customer must be on file with Diagnostic Products Corporation before kits or components containing radioactive material can be shipped. These radioactive materials may be acquired by any customer with the appropriate Specific license. Under a General license these radioactive materials may be acquired only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories and hospitals — and strictly for *in vitro* clinical or laboratory tests not involving external or internal administration of the radioactive material or its radiation to human beings or other animals. Its acquisition, receipt, storage, use, transfer and disposal are all subject to the regulations and a (General or Specific) license of the U.S. Nuclear Regulatory Commission or a State with which the NRC has entered into an agreement for the exercise of regulatory control.

Handle radioactive materials according to the requirements of your General or Specific license. To minimize exposure to radiation, the user should adhere to guidelines set forth in the National Bureau of Standards publication on the *Safe Handling of Radioactive Materials* (Handbook No. 92, issued March 9, 1964) and in subsequent publications issued by State and Federal authorities.

Wipe up spills promptly and decontaminate affected surfaces. Avoid generation of aerosols. Dispose of solid radioactive waste according to license requirements. General licensees (holders of NRC Form 483) may dispose of solid radioactive waste as nonradioactive waste, after removing labeling. Specific licensees (NRC Form 313) should refer to Title 10, Code of Federal Regulations, Part 20. Licensees in Agreement States should refer to the appropriate regulations of their own state. General licensees may dispose of liquid radioactive waste of the type contained in this product through a laboratory sink drain. Licensees must remove or deface labels from empty containers of radioactive materials before disposal of solid waste. Specific licensees

may dispose of small quantities of liquid radioactive waste of the type used in this product through a laboratory sink drain. Refer to the appropriate regulations applicable to your laboratory.

Materials Supplied: Initial Preparation

DHEA Ab-Coated Tubes (TDH1)

Polypropylene tubes coated with antibodies to DHEA and packaged in zip-lock bags. Store refrigerated and protected from moisture, carefully resealing the bags after opening. Stable at 2–8°C until the expiration date marked on the bag. Color: buff.

TKDH1: 75 tubes. **TKDH2:** 150 tubes.

¹²⁵I DHEA (TDH2)

Lyophilized iodinated DHEA. Reconstitute each vial with **80 mL** distilled water. Stable at 2–8°C for 30 days after reconstitution, or until the expiration date marked on the vial. Do not transfer the tracer to a plastic container.

TKDH1: 1 vial. **TKDH2:** 2 vials.

DHEA Buffer (DHC3)

100 mL of buffer, supplied in liquid form, ready to use. It is employed for resuspending the extracts, and also serves as the zero calibrator (Calibrator A). Store the buffer capped. Stable at 2–8°C for 30 days after opening.

TKDH1: 1 vial. **TKDH2:** 1 vial.

DHEA Calibrators (DHC4–9)

One set of six vials, labeled B through G, containing different levels of DHEA in the same matrix used for the DHEA Buffer. Each vial contains 3.0 mL. Stable at 2–8°C for 30 days after opening.

TKDH1: 1 set. **TKDH2:** 1 set.

The calibrators represent, respectively, 0.5, 1.0, 2.5, 5.0, 10 and 30 nanograms of DHEA per milliliter (ng/mL); equivalently, 1.7, 3.5, 8.7, 17.3, 34.7 and 104 nanomoles per liter (nmol/L). Intermediate calibration points may be obtained by mixing the calibrators in suitable proportions. The life of the calibrators can be extended by freezing. Aliquot, if necessary, in tightly capped silanized glass tubes, to avoid repeated thawing and freezing.

Materials Required But Not Provided

Gamma counter — compatible with standard 12x75 mm tubes

Vortex mixer

Reagent Preparation

Graduated cylinder: 100 mL

Extraction

Glass pipets: 1.0 mL and 6.0 mL

16x125 mm tubes with Teflon-lined caps. The tubes should be made of high-quality borosilicate glass. (Do not use plastic.) To avoid the risk of contamination, use "precleaned" glassware or employ a final rinse with 0.1 N HCl followed by dichloromethane.

Mechanical rotator, capable of 15 to 20 revolutions per minute

Aspirator

Evaporator with fume hood and 37°C heat block or water bath

Compressed nitrogen

Ultrasonic bath

Dichloromethane (methylene chloride) of at least spectrophotometric grade

Radioimmunoassay

Plain 12x75 mm polypropylene tubes — for use as NSB tubes, available from DPC

Micropipets with disposable tips: 200 µL and 1,000 µL.

Waterbath capable of maintaining 37°C — neither an oven nor a heat block is suitable

Foam decanting rack — available from DPC

Logit-log graph paper — available from DPC (catalog number: ZP797)

A tri-level, human serum-based immunoassay control, containing DHEA as one of over 25 assayed constituents, is available from DPC (catalog number: CON6).

Specimen Collection

The patient need not be fasting, and no special preparations are necessary. Collect blood by venipuncture⁶ into plain tubes, and separate the serum from the cells, noting the time of collection.

Before extraction, allow the samples to come to room temperature (15–28°C) and mix by *gentle* swirling or inversion. Aliquot, if necessary, to avoid repeated thawing and freezing. Do *not* attempt to thaw frozen specimens by heating them in a waterbath.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Lipemic samples are unsuitable for use.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. Coat-A-Count DHEA has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 1.0 mL of serum per tube.

Storage: 2 days at 2–8°C or 2 months frozen at –20°C.

Sample Preparation

Caution: Dichloromethane (methylene chloride) is a volatile, flammable organic solvent. Conduct the evaporation step under a fume hood. Avoid open flames and do not pipet by mouth. The dichloromethane must be of at least spectrophotometric grade.

Note that only patient samples and controls require extraction before assay. *Do not extract the calibrators.*

- 1 Label one 16x125 mm glass tube for each patient sample and control.

The tubes should be made of high-quality borosilicate glass. (Do not use plastic.) To avoid the risk of contamination, use "precleaned" glassware or employ a final rinse with 0.1 N HCl followed by dichloromethane.

- 2 Pipet **1.0 mL** of each patient sample and control into the tubes prepared.
- 3 **Extraction:** Add **10 mL** dichloromethane to every tube. Cap securely.
- 4 Mix by inversion using a mechanical rotator for **60 minutes**.

Gentle end-over-end rotation at 15 to 20 revolutions per minute is the preferred method.

- 5** Centrifuge for **5 minutes** at approximately 1500xg to separate the two phases.

Any sample partially emulsified should be shaken vigorously and centrifuged again.

- 6** Remove *most* of the top layer, being careful not to disturb the interface.

Use a water-tap aspirator. (Aspirating all the way to the interface may cause the dichloromethane layer to become contaminated with remains of the serum phase on contact.)

- 7** **Evaporation:** Transfer **6.0 mL** of the lower (dichloromethane) phase to a clean glass tube and evaporate to complete dryness under a gentle stream of nitrogen at 37°C.

- 8** Once dry, resuspend the extracts with **0.6 mL** of DHEA Buffer (Calibrator A). Vortex or sonicate thoroughly. Continue with the radioimmunoassay procedure described below.

The reconstituted extracts may be stored at -20°C. Reconstituted extracts prepared from samples with DHEA concentrations greater than the highest calibrator (30 ng/mL) may be further diluted with the DHEA Buffer before assay.

Radioimmunoassay Procedure

All components must be at room temperature (15–28°C) before use. Calibrators B through G and the DHEA Buffer (Calibrator A) are ready to use. Only patient samples and controls must be subjected to the extraction procedure described above.

- 1** **Plain Tubes:** Label four plain (uncoated) 12×75 mm polypropylene tubes T (total counts) and NSB (nonspecific binding) in duplicate.

Because nonspecific binding in the Coat-A-Count procedure is characteristically low, the NSB tubes may be safely omitted without compromising accuracy or quality control.

Coated Tubes: Label fourteen DHEA Ab-Coated Tubes A (maximum binding) and B

through G in duplicate. Label additional antibody-coated tubes, also in duplicate, for controls and patient samples.

Calibrator	ng/mL	nmol/L
A (MB)	0	0
B	0.5	1.7
C	1.0	3.5
D	2.5	8.7
E	5.0	17.3
F	10	34.7
G	30	104

- 2** Pipet **200 µL** of the zero calibrator A (DHEA Buffer) into the NSB and A tubes, and **200 µL** of each nonzero calibrator into the tubes marked B through G. Pipet **200 µL** of each control and patient sample into the tubes prepared. **Pipet directly to the bottom.**

Reconstituted extracts prepared from samples with DHEA concentrations greater than the highest calibrator (30 ng/mL) may be further diluted with the DHEA Buffer before assay.

It is good practice to use a disposable-tip micropipet, changing the tip between samples, in order to avoid carryover contamination.

- 3** Add **1.0 mL** of ¹²⁵I DHEA to every tube. Vortex.

No more than 10 minutes should elapse during the dispensing of the tracer. Set the T tubes aside for counting at step 6; they require no further processing.

- 4** Incubate for **2 hours at 37°C**.

- 5** Decant thoroughly.

Removing all visible moisture will greatly enhance precision. Using a foam decanting rack, decant the contents of all tubes (except the T tubes) and allow them to drain for 2 or 3 minutes. Then strike the tubes sharply on absorbant paper to shake off all residual droplets.

- 6** Count for **1 minute** in a gamma counter.

Calculation of Results

To obtain results in terms of concentration from a logit-log representation of the calibration curve, first calculate for each pair of tubes the average NSB-corrected counts per minute:

Net Counts = (Average CPM) minus (Average NSB CPM)

Then determine the binding of each pair of tubes as a percent of maximum binding (MB), with the NSB-corrected counts of the A tubes taken as 100%:

Percent Bound = (Net Counts / Net MB Counts) × 100

(The calculation can be simplified by omitting the correction for nonspecific binding; samples within range of the calibrators yield virtually the same results when Percent Bound is calculated directly from Average CPM.)

Using logit-log graph paper, plot Percent Bound on the vertical (probability) axis against Concentration on the horizontal (logarithmic) axis for each of the nonzero calibrators, and draw a straight line approximating the path of these points. Results for the unknowns may then be read from the line by interpolation.

It is good practice to inspect results for agreement within replicates, and to construct a graph of the calibration curve (even if the calculations are handled by computer) as a visual check on the appropriateness of the transformation used and as a way to detect deviant calibration points. We also recommend keeping track of these data reduction parameters:

T = Total Counts (as counts per minute)

%NSB = $100 \times (\text{Average NSB Counts} / \text{Total Counts})$

%MB = $100 \times (\text{Net Counts} / \text{Total Counts})$

And the 20, 50 and 80 percent "intercepts," where

20% = Concentration at 20 Percent Bound, etc.

Note that other approaches, e.g. a sound implementation of the 4-parameter logistic, are also acceptable. See Dudley RA, et al. Guidelines for immunoassay data reduction. Clin Chem 1985;31:1264-71.

Example Run: For illustration only, not for calculating results from another run. (See "Example Run" table.)

Quality Control

Controls or serum pools with at least two DHEA concentration levels — low and high — should routinely be assayed as unknowns.

Report patient results only if the control results for that assay meet your laboratory's established criteria for acceptability.

It is good laboratory practice to record for each assay the lot numbers of the components used, as well as the dates when they were first reconstituted or opened. We also recommend charting control results from day to day — as described, for example, in Westgard JO, et al. A multi-rule chart for quality control. Clin Chem 1981;27:493-501. Note that repeat samples can serve as a valuable additional tool for monitoring interassay precision, and that pairs of control tubes can be spaced throughout the assay to help verify the absence of significant drift.

Solvent Blanks:

For reliable results, the dichloromethane must be of at least spectrophotometric grade and the nitrogen source, if used, must not contaminate the evaporation tubes with oil or moisture. (The use of a trap containing calcium sulfate will help prevent such contamination.) As a check on the purity of these reagents, we recommend periodic use of a "solvent blank." This involves extracting and assaying distilled or sterile water, in the manner prescribed for patient samples, to verify that it yields an apparent DHEA concentration at or below the detection limit.

Expected Values


A preliminary reference range study was performed with the Coat-A-Count DHEA procedure on serum samples collected from 33 healthy adult volunteers, including 19 females and 14 males. The results ranged from 0.2 to 9.8 ng/mL with a median value of 2.5 ng/mL. Ninety-one percent of the results were less than 6.0 ng/mL. Note that these results were

obtained for the recommended sample type, namely serum.

A further reference range study was conducted at a large northeastern hospital on 136 women and 155 men, all of known age. The results, broken down by sex and by age, are summarized in the table and graphs below. In the graphs, solid circles depict medians (50th percentiles), while "T" bars depict 95th percentiles. Shading represents a smoothing of the observed 5th and 95th percentiles according to a simple model of the behavior of circulating DHEA as a function of age.

Females				
Age Group	n	5th %ile	50th %ile	95th %ile
10–19	2	—	3.2	—
20–29	23	1.1	1.6	7.2
30–39	29	0.55	1.3	4.0
40–49	29	0.40	0.85	3.6
50–59	18	0.3	0.7	2.8
60–69	21	0.1	0.3	1.8
70–79	8	0.1	0.25	1.5
80–89	6	—	—	—

Males				
Age Group	n	5th %ile	50th %ile	95th %ile
10–19	6	—	1.1	—
20–29	17	1.0	3.9	7.3
30–39	21	0.75	2.6	5.1
40–49	23	0.7	1.9	4.6
50–59	29	0.55	1.7	4.6
60–69	37	0.4	1.0	2.1
70–79	18	0.1	0.6	1.7
80–89	4	—	0.3	—

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Performance Data

See Tables and Graphs for data *representative* of the Coat-A-Count DHEA kit's performance. DHEA results in the sections are expressed as nanograms per milliliter (ng/mL). Since the molecular mass of DHEA is 288.4 daltons, to convert to S.I. units — nanomoles per liter (nmol/L) — multiply by 3.467:

Conversion Factor:

ng/mL × 3.467 → nmol/L

Calibration Range: 0.5 – 30 ng/mL
(1.7 – 104 nmol/L.)

Analytical Sensitivity: 0.04 ng/mL.

Intraassay Precision (Within-Run)

Statistics were calculated for samples from the results of 20 pairs of tubes in a single assay. (See "Intraassay Precision" table.)

Interextraction: Statistics were calculated for three samples that were individually extracted twenty times, then assayed by the Coat-A-Count DHEA kit. (See "Interextraction" table.)

Specificity: The kit is highly specific for DHEA, with a relatively low crossreactivity to other naturally occurring steroids or therapeutic drugs that may be present in patient samples. Extracted samples dissolved in the zero calibrator were used for this experiment, with the exception of the starred compounds(*), which were assayed unextracted and dissolved in the zero calibrator. (See "Specificity" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table.)

Recovery 1: Samples spiked 1 to 19 with three DHEA solutions (50, 100 and 200 ng/mL) were assayed. (See "Recovery 1" table.)

Recovery 2: In a similar experiment, three more spiking solutions (D, E and F) were prepared at levels of 100, 200 and 400 ng/mL. Expected values were calculated by adding 95% of the unspiked value to 5% of the spiking solution concentration (5, 10 and 20 ng/mL, respectively). (See "Recovery 2" table.)

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Lipemia: It has been established that gross hyperchylomicronemia can interfere with some immunoassays. This interference can be avoided by isolating the clear infranatant layer before assay. To accomplish this manually, draw the lipemic sample into a pipet and seal the tip with Parafilm. Let the pipet stand until the lipids float to the top. Then drain and assay the clear layer of serum. Alternatively, an airfuge may be used. (It has not been established whether other types of lipemia also interfere.)

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 33 volunteers into plain, heparinized and EDTA vacutainer tubes. All samples were assayed by the Coat-A-Count DHEA procedure, with the following results.

(Heparin) = $0.84 (\text{Serum}) + 1.0 \text{ ng/mL}$
 $r = 0.90$

(EDTA) = $0.92 (\text{Serum}) + 0.2 \text{ ng/mL}$
 $r = 0.95$

Means:
3.8 ng/mL (Heparin)

3.4 ng/mL (Serum)
3.4 ng/mL (EDTA)

The serum values ranged from 0.2 to 13.7 ng/mL. Both heparinized and EDTA plasma showed a variability — that is, discrepancies relative to results on serum samples from the same individuals — greater than the variability observed in the interextraction precision study. On this basis, it is recommended that samples be collected as serum.

Method Comparison: The Coat-A-Count DHEA kit was compared to another commercially available DHEA radioimmunoassay on 90 patient serum samples. The statistics displayed below show good correlation between these two methodologies.

(CAC) = $1.17 (\text{Kit A}) + 0.26 \text{ ng/mL}$
 $r = 0.94$

Means:
2.0 ng/mL (Coat-A-Count)
1.5 ng/mL (Kit A)

Effect of Free Fatty Acid: To investigate the effect of nonesterified fatty acid (NEFA), aliquots of three samples were lyophilized, reconstituted and then assayed by the Coat-A-Count DHEA procedure. Aqueous solutions containing 0, 2.5, 5, and 10 mmol/L of oleic acid were used for the reconstitution step. The results demonstrate that oleic acid — at physiologically high concentrations — does have a significant *in vitro* effect on the procedure.

Protein Effect: To simulate various protein concentrations, experiments were performed in which 1.0 mL aliquots of a stripped healthy human serum pool were freeze-dried and then reconstituted with various volumes of water (0.5, 0.67, 1.0, 1.5, 2.0 mL). Each reconstituted aliquot was then assayed by the Coat-A-Count DHEA procedure. Note that aliquots reconstituted with half the original volume represent an extremely high protein concentration, in the order of 14 g/dL. (The factor to correct for reconstitution volume is tabulated below.) The results indicate that total protein concentrations which are extremely low may have a clinically significant effect on the Coat-A-Count DHEA assay. (See "Protein Effect" table.)

References

1) Cohen HN, et al. Clinical value of adrenal androgen measurement in the diagnosis of delayed puberty. *Lancet* 1981;i:689-92. 2) Erb J, et al. Discrimination between schizophrenic and control subjects by means of plasma dehydroepiandrosterone measurements. *J Clin Endocrinol Metab* 1981;52:181-6. 3) Fiet J, et al. Simultaneous radioimmunoassay of androstenedione, dehydroepiandrosterone and 11- β hydroxyandrostenedione in plasma. *Horm Res* 1980;13:133-49. 4) de Peretti E, Forest M. Unconjugated dehydroepiandrosterone plasma levels in normal subjects from birth to adolescence. *J Clin Endocrinol Metab* 1976;43:982-7. 5) Tourniaire J, Pugeat M. Strategic approach of hyperandrogenism in women. *Horm Res* 1983;18:125-34. 6) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture: approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782, 310.645.8200
Fax: 310.645.9999. To place an order:
Tel: 800.372.1782. Fax: 800.234.4372.
Outside the United States, contact your National Distributor.

The Quality System of Diagnostic Products Corporation is registered to ISO 13485:2003.

Tables and Graphs

Example Run

Tube ¹	Duplicate CPM ²	Average CPM ³	Net CPM ⁴	Percent Bound ⁵	DHEA ng/mL ⁶
T ⁷	73,741 72,920	73,331			
(NSB) ⁸	431 411	421			
A (MB) ⁹	34,646 34,588	34,617	34,196	100%	0
B	26,323 26,127	26,225	25,804	75.5%	0.5
C	23,907 23,831	23,869	23,448	68.6%	1
D	19,613 19,552	19,583	19,162	56.0%	2.5
E	16,557 16,556	16,557	16,136	47.2%	5
F	12,823 12,805	12,814	12,393	36.2%	10
G	8,270 8,240	8,255	7,834	22.9%	30
Unknowns: ¹⁰					
X1	23,522 23,019	23,271	22,850	66.8%	1.1
X2	17,590 17,461	17,526	17,105	50.0%	3.8
X3	13,439 13,410	13,425	13,004	38.0%	9.0

Quality Control Parameters:¹¹

T⁷ = 73,331 cpm

%NSB⁸ = 0.6%

%MB⁹ = 47%

20% Intercept¹² = 42.84 ng/mL

50% Intercept = 3.82 ng/mL

80% Intercept = 0.34 ng/mL

Intraassay Precision (ng/mL)

	Mean ¹	SD ²	CV ³
1	1.3	0.14	10.8%
2	2.6	0.17	6.5%
3	9.6	0.50	5.2%

Interextraction (ng/mL)

	Mean ¹	SD ²	CV ³
1	1.2	0.14	11.7%
2	2.8	0.22	7.9%
3	12.6	0.71	5.6%

Specificity

Compound ¹	% Cross reactivity ²
DHEA	(100%)
DHEA-SO ₄	0.066%
Aldosterone	ND
5 α -Androstan-3 α ,17 β -diol*	0.43%
5 β -Androstan-3 α ,17-diol	0.046%
5 α -Androstan-3,17-dione*	0.97%
5 β -Androstan-3,17-dione	0.13%
Δ^5 -Androsten-3 β ,17-diol	6.13%
Androstenedione*	2.18%
Androsterone	0.70%
Corticosterone	ND
Cortisol	ND
Cortisone	ND
Danazol	0.003%
11-Deoxycortisol*	ND
Dexamethasone	0.005%
5 α -Dihydrotestosterone*	0.26%
5(10)-Estron-17 α -ethynyl-17 β -ol-3-one	0.009%
Δ^4 -Estron-17 α -methyl-17 β -ol-3-one	0.17%
Estrone	0.88%
Ethisterone	0.01%
11 β -Hydroxyandrostenedione	0.02%
19-Hydroxyandrostenedione	0.16%
11-Ketotestosterone	0.003%
Norethindrone	0.01%
Norethinodrel	0.009%
Prednisone	ND
Progesterone	0.05%
Spironolactone	0.016%
Triamcinolone	ND
Testosterone	1.07%

ND: not detectable³**Linearity (ng/mL)**

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	4 in 4 ⁵	6.0	—	—
	3 in 4	4.1	4.5	91%
	2 in 4	2.7	3.0	90%
	1 in 4	1.2	1.5	80%
2	4 in 4	11.3	—	—
	3 in 4	8.8	8.5	104%
	2 in 4	5.1	5.7	89%
3	4 in 4	11.8	—	—
	3 in 4	9.5	8.9	107%
	2 in 4	5.6	5.9	95%
4	4 in 4	12.7	—	—
	3 in 4	10.2	9.5	107%
	2 in 4	5.5	6.3	87%
	1 in 4	3.1	3.2	97%
5	4 in 4	19.0	—	—
	3 in 4	13.8	14.3	97%
	2 in 4	9.1	9.5	96%
6	4 in 4	21.4	—	—
	3 in 4	15.8	16.1	98%
	2 in 4	10.7	10.7	100%

Recovery 1

	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	0.0	—	—
	A	2.5	2.5	100%
	B	4.9	5.0	98%
	C	10.5	10.0	105%
2	—	0.0	—	—
	A	2.5	2.5	100%
	B	5.2	5.0	104%
	C	10.5	10.0	105%

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O Sistema de Qualidade da Diagnostic Products Corporation está registado sob ISO 13485:2003.
EC REP DPC Biermann GmbH