COAT-A-COUNT®

DHEA-SO₄

DPC®

Coat-A-Count® DHEA-SO₄

English

Intended Use: Coat-A-Count DHEA-SO₄ is a radioimmunoassay designed for the quantitative measurement of dehydroepiandrosterone sulfate (DHEA-SO₄) in serum and plasma. It is intended strictly for *in vitro* diagnostic use.

Catalog numbers: **TKDS1** (100 tubes), **TKDS2** (200 tubes), **TKDS5** (500 tubes), **TKDSX** (1000 tubes)

The 100-tube kit contains less than 4 microcuries (148 kilobecquerels) of radioactive ¹²⁵I DHEA-SO₄; the 200-tube kit contains less than 8 microcuries (296 kilobecquerels); the 500-tube kit contains less than 20 microcuries (740 kilobecquerels); and the 1000-tube kit contains less than 40 microcuries (1480 kilobecquerels).

Summary and Explanation

Measurement of dehydroepiandrosterone sulfate (DHEA-SO₄, DHEAS), an adrenal steroid, is important to investigations of abnormal hair growth (hirsutism) and balding (alopecia) in women.^{1,4,7} It is also of value in the assessment of adrenarche and delayed puberty.^{3,4}

The DHEA-SO₄ in circulation originates almost entirely from the adrenals, though in men some may also derive from the testes, partly accounting for the sex difference which emerges at about age 15.^{11,17} On the other hand, this hormone is not produced by the ovaries even under pathological conditions.^{2,16} In itself, DHEA-SO₄ is only weakly androgenic, but it can metabolize to more potent androgens like androstenedione and testosterone, and thus be indirectly a cause of hirsutism or virilization.^{1,16}

Plasma levels of DHEA-SO₄ increase steadily from about the seventh year of life, then gradually decline after the third decade.^{11,16,17} According to the literature, the upper limit of normal for young adults is approximately 300 µg/dL for women, and 500 µg/dL for men.⁸ Pregnancy and

oral contraceptives induce a moderate decrease.^{4,16}

DHEA-SO₄ is secreted into the bloodstream at a rate only somewhat greater than DHEA, but because of its much slower turnover — DHEA-SO₄ has a half-life of nearly a full day — it maintains a plasma level almost a thousandfold higher.^{4,13,16} Unlike cortisol, DHEA-SO₄ does not exhibit significant diurnal variation.^{4,15} Unlike testosterone, it does not circulate bound to sex hormone-binding globulin and hence is not influenced by alterations in the level of this carrier protein.¹⁶ Its abundance, together with its within-day and day-to-day stability, makes it an excellent direct indicator of adrenal androgen output — superior, certainly, to the measurement of urinary 17-ketosteroids in this context.^{1,2,4,8}

Accordingly, DHEA-SO₄ is often assayed in conjunction with free testosterone as an initial screen for hyperandrogenism in hirsutism.^{5,9} At least one of these two hormones is likely to be elevated in the great majority of cases — reportedly over 80 percent of the time.⁹ Sometimes DHEA-SO₄ is the only hormone circulating at a level above normal, and is apparently more likely to be elevated during the early stages of hirsutism than most other androgens.⁹

High DHEA-SO₄ levels are often encountered in the polycystic ovary syndrome, showing that adrenal hyperandrogenism is a fairly typical facet of this syndrome.² Elevated plasma levels which — over the course of two weeks or so — are dexamethasone-suppressible may also result from adrenal hyperplasia.^{4,6,9} Extremely high levels (greater than 700 or 800 µg/dL) in women are suggestive of a hormone-secreting adrenal tumor.¹³ (By contrast, DHEA-SO₄ levels are typically normal in the presence of ovarian tumors.¹)

Principle of the Procedure

In the Coat-A-Count DHEA-SO₄ procedure, ¹²⁵I-labeled DHEA-SO₄ competes with DHEA-SO₄ in the patient sample for sites on antibody-coated tubes. After incubation, separation of bound from

free is achieved by simply decanting. The tube is then counted in a gamma counter, the counts being inversely related to the amount of DHEA-SO₄ present in the patient sample. The quantity of DHEA-SO₄ in the sample is determined by comparing the counts to a calibration curve.

Reagents to Pipet: 1

Total Incubation Time: 30 minutes

Total Counts at Iodination:
approximately 50,000 cpm

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C in a refrigerator designated for incoming radioactive materials. Dispose of in accordance with applicable laws.

Do not use reagents beyond their expiration dates.

Some components supplied in this kit may contain human source material and/or other potentially hazardous ingredients which necessitate certain precautions.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Water: Use distilled or deionized water.

Radioactivity

A copy of any radioisotope license certificate (Specific or General) issued to a US customer must be on file with Diagnostic Products Corporation before kits or components containing radioactive material can be shipped. These radioactive materials may be acquired by any customer with the appropriate Specific license. Under a General license these radioactive materials may be acquired only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories and hospitals — and strictly

for *in vitro* clinical or laboratory tests not involving external or internal administration of the radioactive material or its radiation to human beings or other animals. Its acquisition, receipt, storage, use, transfer and disposal are all subject to the regulations and a (General or Specific) license of the U.S. Nuclear Regulatory Commission or a State with which the NRC has entered into an agreement for the exercise of regulatory control.

Handle radioactive materials according to the requirements of your General or Specific license. To minimize exposure to radiation, the user should adhere to guidelines set forth in the National Bureau of Standards publication on the *Safe Handling of Radioactive Materials* (Handbook No. 92, issued March 9, 1964) and in subsequent publications issued by State and Federal authorities.

Wipe up spills promptly and decontaminate affected surfaces. Avoid generation of aerosols. Dispose of solid radioactive waste according to license requirements. General licensees (holders of NRC Form 483) may dispose of solid radioactive waste as nonradioactive waste, after removing labeling. Specific licensees (NRC Form 313) should refer to Title 10, Code of Federal Regulations, Part 20. Licensees in Agreement States should refer to the appropriate regulations of their own state. General licensees may dispose of liquid radioactive waste of the type contained in this product through a laboratory sink drain. Licensees must remove or deface labels from empty containers of radioactive materials before disposal of solid waste. Specific licensees may dispose of small quantities of liquid radioactive waste of the type used in this product through a laboratory sink drain. Refer to the appropriate regulations applicable to your laboratory.

Materials Supplied

DHEA-SO₄ Ab-Coated Tubes (TDS1)

Polypropylene tubes coated with antibodies to DHEA-SO₄ and packaged in zip-lock bags. Store refrigerated and protected from moisture, carefully resealing the bags after opening. Stable at 2–8°C until the expiration date marked on the bag. Color: light pink.

TKDS1: 100 tubes. **TKDS2:** 200 tubes.
TKDS5: 500 tubes. **TKDSX:** 1000 tubes.

125I DHEA-SO₄ (TDS2)

Iodinated DHEA-SO₄ in lyophilized form, with preservative. Reconstitute each vial with **110 mL** of distilled or deionized water. Cap securely and mix by *gentle* inversion. Stable at 2–8°C for 30 days after opening, or until the expiration date marked on the vial. *Do not transfer the tracer to a plastic container.*

TKDS1: 1 vial. **TKDS2:** 2 vials.

TKDS5: 5 vials. **TKDSX:** 10 vials.

DHEA-SO₄ Calibrators (TDS3–9)

Seven vials, labeled A through G, of protein-based calibrators in lyophilized form, with preservative. Reconstitute the zero calibrator A with **3.0 mL** distilled or deionized water and the remaining calibrators B through G with **1.0 mL**. Stable at 2–8°C for 30 days after opening. For longer storage, aliquot into tightly capped glass or polypropylene tubes and freeze at –20°C.

TKDS1: 1 set. **TKDS2:** 1 set.

TKDS5: 2 sets. **TKDSX:** 3 sets.

The calibrators contain, respectively, 0, 5, 20, 50, 200, 500 and 1,000 micrograms of DHEA-SO₄ per deciliter (µg/dL); equivalently: 0, 0.14, 0.54, 1.36, 5.4, 13.6 and 27.1 micromoles per liter (µmol/L). Intermediate calibration points may be obtained by mixing calibrators in suitable proportions.

Materials Required But Not Provided

Gamma counter – compatible with standard 12×75 mm tubes

Vortex mixer

Reagent Preparation

Graduated cylinder: 110 mL

Volumetric pipets: 1.0 mL and 3.0 mL

Distilled or deionized water

Radioimmunoassay

Plain 12×75 mm polypropylene tubes — for use as NSB tubes, available from DPC

Micropipets: 50 µL and 1 000 µL.

Waterbath — capable of maintaining 37°C. Neither an oven nor a heat block is suitable.

Foam decanting rack — available from DPC

Logit-log graph paper — available from DPC (catalog number: ZP797)

A tri-level, human serum-based immunoassay control, containing DHEA-SO₄ as one of over 25 assayed constituents, is available from DPC (catalog number: CON6).

Specimen Collection

The patient need not be fasting, and no special preparations are necessary. Collect blood by venipuncture¹⁹ into plain tubes, and separate the serum from the cells, noting the time of collection. EDTA plasma is also suitable for this assay. (See the section on Alternate Sample Type.)

Studies of nonesterified fatty acid (NEFA) show that oleic acid has a significant *in vitro* effect on the assay at physiologically high concentrations.

An ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. Coat-A-Count DHEA-SO₄ has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 50 µL per tube.

Storage: 2 days at 2–8°C in glass or polypropylene tubes, or 2 months at –20°C.¹⁸ Before assay, allow the samples to come to room temperature (15–28°C) and mix by *gentle* swirling or inversion. Aliquot, if necessary, to avoid repeated thawing and freezing. Do *not* attempt to thaw frozen specimens by heating them in a waterbath.

Radioimmunoassay Procedure

All components must be at room temperature (15–28°C) before use.

- 1 Plain Tubes:** Label four plain (uncoated) 12×75 mm polypropylene

Coat-A-Count DHEA-SO₄ (PITKDS-4, 2005-03-18)

tubes T (total counts) and NSB (nonspecific binding) in duplicate.

Because nonspecific binding in the Coat-A-Count procedure is characteristically low, the NSB tubes may be safely omitted without compromising accuracy or quality control.

Coated Tubes: Label fourteen DHEA-SO₄ Ab-Coated Tubes A (maximum binding) and B through G in duplicate. Label additional DHEA-SO₄ Ab-Coated Tubes, also in duplicate, for controls and patient samples.

Calibrators	µg/dL	µg/mL	µmol/L
A (MB)	0	0	0
B	5	0.05	0.14
C	20	0.2	0.54
D	50	0.5	1.36
E	200	2	5.4
F	500	5	13.6
G	1,000	10	27.1

- 2 Pipet **50 µL** of the zero calibrator A into the NSB and A tubes, and **50 µL** of each remaining calibrator, control and patient sample into the tubes prepared. **Pipet directly to the bottom.**

Patient samples expected to contain DHEA-SO₄ concentrations greater than the highest calibrator (1,000 µg/dL) should be diluted in the zero calibrator before assay.

It is good practice to use a disposable-tip micropipet, changing the tip between samples, in order to avoid carryover contamination.

- 3 Add **1.0 mL** of ¹²⁵I DHEA-SO₄ to every tube. Vortex briefly and gently.

No more than ten minutes should elapse during the dispensing of the tracer. If a repeating pipet with a plastic syringe is used, rinse it twice with the tracer before dispensing. Set the T tubes aside for counting (at step 6); they require no further processing.

- 4 Incubate for 30 minutes at 37°C.

Use a waterbath; neither an oven nor a heat block is suitable.

- 5 Decant thoroughly.

Removing all visible moisture will greatly enhance precision. Decant the contents of all tubes (except the T tubes) using a foam decanting rack, and allow them to drain for 2 or 3 minutes. Then strike the tubes sharply on absorbant paper to shake off all residual droplets.

- 6 Count for **1 minute** in a gamma counter.

Calculation of Results

To obtain results in terms of concentration from a logit-log representation of the calibration curve, first calculate for each pair of tubes the average NSB-corrected counts per minute:

Net Counts = (Average CPM) minus (Average **NSB** CPM)

Then determine the binding of each pair of tubes as a percent of maximum binding (MB), with the NSB-corrected counts of the A tubes taken as 100%:

Percent Bound = (Net Counts / Net **MB** Counts) × 100

The calculation can be simplified by omitting the correction for nonspecific binding; samples within range of the calibrators yield virtually the same results when Percent Bound is calculated directly from Average CPM.)

Using logit-log graph paper, plot Percent Bound on the vertical (probability) axis against Concentration on the horizontal (logarithmic) axis for each of the nonzero calibrators, and draw a straight line approximating the path of these points. Results for the unknowns may then be read from the line by interpolation.

It is good practice to inspect results for agreement within replicates, and to construct a graph of the calibration curve (even if the calculations are handled by computer) as a visual check on the appropriateness of the transformation used and as a way to detect deviant calibration points. We also recommend keeping track of these data reduction parameters:

T = Total Counts (as counts per minute)

%NSB = 100 × (Average NSB Counts / Total Counts)

%MB = 100 × (Net Counts / Total Counts)

And the 20, 50 and 80 percent "intercepts," where

20% = Concentration at 20 Percent Bound, etc.

Note that other approaches, e.g. a sound implementation of the 4-parameter logistic, are also acceptable. See Dudley RA, et al. Guidelines for immunoassay data reduction. Clin Chem 1985;31:1264-71.

Example Run: For illustration only, not for calculating results from another run. (See "Example Run" table.)

Quality Control

Controls (or serum pools) with at least two DHEA-SO₄ concentration levels — low and high — should routinely be assayed as unknowns.

Report patient results only if the control results for that assay meet your laboratory's established criteria for acceptability.

It is good laboratory practice to record for each assay the lot numbers of the components used, as well as the dates when they were first reconstituted or opened. We also recommend charting control results from day to day — as described, for example, in Westgard JO, et al. A multi-rule chart for quality control. Clin Chem 1981;27:493-501. Note that repeat samples can serve as a valuable additional tool for monitoring interassay precision, and that pairs of control tubes can be spaced throughout the assay to help verify the absence of significant drift.

Expected Values

Reference range studies were conducted at two institutions using the Coat-A-Count DHEA-SO₄ kit on a total of 202 adults. Tabulated below are the median and range for the combined data, the range representing the central 95 percent of the observations.

	Median	Range (Central 95%)	<i>n</i>
Females	170 µg/dL	35 – 430 µg/dL	132
Males	280 µg/dL	80 – 560 µg/dL	70

In addition, a total of 74 pregnancy and postpartum samples were analyzed. The results, as expected, showed a progressive decline: median values for the

first and third trimesters and for one day postpartum were 230, 105 and 65 µg/dL, respectively.

A further reference range study was conducted at a large northeastern hospital on 136 women and 155 men, all of known age. The results, broken down by sex and by age, are summarized in the following table and graphs. In the graphs, solid circles depict medians (50th percentiles) while "T" bars depict 95th percentiles. Shading represents a smoothing of the observed 5th and 95th percentiles according to a simple model of the behavior of circulating DHEA-SO₄ as a function of age.


Results of a previous study by the same laboratory on 60 postmenopausal samples are also summarized in the following table and graphs.

Age Group	<i>n</i>	5 th %ile	50 th %ile	95 th %ile
Females				
10–19	2	—	140	—
20–29	23	65	185	380
30–39	29	45	150	270
40–49	29	32	120	240
50–59	18	26	85	200
60–69	21	13	50	130
70–79	8	17	40	90
80–89	6	—	26	—
Post-menopausal	60	10	55	190
Males				
10–19	6	—	215	—
20–29	17	280	420	640
30–39	21	120	300	520
40–49	23	95	250	530
50–59	29	70	160	310
60–69	37	42	130	290
70–79	18	28	80	175
80–89	4	—	36	—

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

See also "DHEA-SO₄ by Age: Males" and "DHEA-SO₄ by Age: Females" graphs. Solid circles depict medians (50th percentiles) while "T" bars depict 95th percentiles. Shading represents a

smoothing of the observed 5th and 95th percentiles according to a simple model of the behavior of circulating DHEA-SO₄ as a function of age.

Performance Data

See Tables and Graphs for data representative of the Coat-A-Count DHEA-SO₄ kit's performance. Results are expressed as micrograms of DHEA-SO₄ per deciliter (µg/dL).

Conversion Factor:

µg/dL × 0.02714 → µmol/L

Calibration Range: 5 – 1,000 µg/dL
(0.14 – 27 µmol/L)

Analytical Sensitivity: 1.1 µg/dL.

Intraassay Precision (Within-Run):

Statistics were calculated for samples from the results of 20 pairs of tubes in a single assay. (See "Intraassay Precision" table.)

Interassay Precision (Run-to-Run):

Statistics were calculated for samples from the results of pairs of tubes in 24 different assays. (See "Interassay Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three DHEA-SO₄ solutions (1,452, 3,046

and 4,820 µg/dL) were assayed. (See "Recovery" table.)

Specificity: The Coat-A-Count DHEA-SO₄ antiserum is highly specific for DHEA-SO₄, with a relatively low crossreactivity to other naturally occurring steroids and drugs that may be present in patient samples. (See "Specificity" table.)

End-of-Run Effect: None up to approximately 250 tubes. (See "End-of-Run Effect" table.)

Protein Effect: To simulate various protein concentrations, experiments were performed in which 2.0 mL aliquots of a human serum pool were freeze-dried and then reconstituted with various volumes of water. Each reconstituted aliquot was then assayed by the Coat-A-Count DHEA-SO₄ assay procedure. Note that aliquots reconstituted with half the original volume represent an extremely high protein concentration, on the order of 14 g/dL. (The factor to correct for reconstitution volume is tabulated.) These results indicate that even wide variations in protein have relatively little effect on the Coat-A-Count DHEA-SO₄ assay. (See "Protein Effect" table.)

Effect of Free Fatty Acid: To investigate the effect of nonesterified fatty acid (NEFA), aliquots of three samples were lyophilized, reconstituted and then assayed by the Coat-A-Count DHEA-SO₄ procedure. Aqueous solutions containing 0, 2.5, 5, and 10 mmol/L of oleic acid were used for the reconstitution step. The results demonstrate that oleic acid — at physiologically high concentrations — does have a significant *in vitro* effect on the procedure.

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of packed red blood cells in concentrations up to 30 µL/mL has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 22 volunteers into plain, heparinized, EDTA and Becton Dickinson SST[®] vacutainer tubes. All samples were assayed by the

Coat-A-Count DHEA-SO₄ procedure. By linear regression:

(Heparin) = 1.03 (Serum) – 3.4 µg/dL
r = 0.991

(EDTA) = 1.02 (Serum) – 7.8 µg/dL
r = 0.990

(SST) = 1.09 (Plain Tubes) – 19 µg/dL
r = 0.987

Means:

208 µg/dL (Serum)

211 µg/dL (Heparin)

204 µg/dL (EDTA)

208 µg/dL (SST)

Method Comparison: The Coat-A-Count DHEA-SO₄ procedure was compared to a commercially available radioimmunoassay (Kit A) on 112 patient serum samples. By linear regression:

(Coat-A-Count) = 1.14 (Kit A) – 6.9 µg/dL
r = 0.975

Means:

183 µg/dL (Coat-A-Count)

168 µg/dL (Kit A)

In another study, 99 patient serum samples were simultaneously assayed by the Coat-A-Count DHEA-SO₄ procedure using the new antibody and previous antibody. By linear regression:

(New) = 0.99 (Previous) – 0.46 µg/dL
r = 0.980

Means:

196 µg/dL (New)

197 µg/dL (Previous)

Curve Displacement: Two samples were each processed by the Coat-A-Count DHEA-SO₄ procedure in assays which included one-to-one dilutions of the sample with each of the calibrators A through G (0, 5, 20, 50, 200, 500 and 1,000 µg/dL). (See "Curve Displacement" table for representative data.)

References

1) Braunstein GD. Coping with hirsutism – management strategies. *Contemp Ob/Gyn* 1981 May;17(5):189-93. 2) Buvat J, et al. Investigative strategy of hyperandrogenism in women. *Horm Res* 1983;18:106-16. 3) Cohen HN, et al. Clinical value of adrenal androgen measurement in the diagnosis of delayed puberty. *Lancet* 1981;i:689-92. 4) Haning RV. Using DHEAS to monitor androgen disorders. *Contemp Ob/Gyn* 1981 Sept;18(9):117-32. 5) Haning RV, et al. Testosterone free index correlates best with dehydroepiandrosterone sulfate. *Fertil Steril* 1981;36:757-65. 6) Hatch R, et al. Hirsutism: implications, etiology, and management. *Am J*

Obstet Gynecol 1981;140:815-28. 7) Kasick JM, et al. Adrenal androgenic female-pattern alopecia: sex hormones and the balding woman. *Clev Clin Quarterly* 1983 Summer;50(2):111-22. 8) Lobo RA, Paul WL, Goebelsmann U. Dehydroepiandrosterone sulfate as an indicator of adrenal androgen function. *Obstet Gynecol* 1981;57:69-73. 9) Lobo RA, Paul WL, Goebelsmann U. Serum levels of DHEAS in gynecologic endocrinopathy and infertility. *Obstet Gynecol* 1981;57:607-12. 10) Maroulis GB, Abraham GE. Diagnostic tests in hirsutism. In: Mahesh V, Greenblatt RB, editors. *Hirsutism and virilism*. Boston: John Wright-PSG, 1983;295-307. 11) De Peretti E, Forest M. Pattern of plasma dehydroepiandrosterone sulfate levels in humans from birth to adulthood: evidence for testicular production. *J Clin Endocrinol Metab* 1978;47:572-77. 12) Rüedi B, et al. Diagnostic strategy in hyperandrogenic syndrome. *Horm Res* 1983;18:117-24. 13) Sciarra F. Diagnosis of virilizing syndromes: endocrinological parameters. In: Molinatti G, et al, editors. *Androgenization in women*. New York: Raven Press, 1983;85-113. 14) Toscano V, et al. Is hirsutism an evolving syndrome? *J Endocrinol* 1983;97:379-87. 15) Tourniaire J, Pugeat M. Strategic approach of hyperandrogenism in women. *Horm Res* 1983;18:125-34. 16) Vermeulen A. Androgen secretion by adrenals and gonads. In: Mahesh V, Greenblatt RB, editors. *Hirsutism and virilism*. Boston: John Wright-PSG, 1983;17-34. 17) Zappulla F, et al. Gonadal and adrenal secretion of dehydroepiandrosterone sulfate in prepubertal and pubertal subjects. *J Endocrinol Invest* 1981;4:197-202. 18) Burtis CA, Ashwood ER, editors. *Tietz textbook of clinical chemistry*. 2nd ed. Philadelphia: W.B. Saunders, 1994. 19) National Committee for Clinical Laboratory Standards, *Procedures for the collection of diagnostic blood specimens by venipuncture*, 4th ed, NCCLS Document H3-A4, 1998.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782, 310.645.8200
Fax: 310.645.9999. To place an order:
Tel: 800.372.1782. Fax: 800.234.4372.
Outside the United States, contact your National Distributor.

The Quality System of Diagnostic Products Corporation is registered to ISO 13485:2003.

Tables and Graphs

Example Run

Tube ¹	Duplicate CPM ²	Average CPM ³	Net CPM ⁴	Percent Bound ⁵	DHEA-SO ₄ µg/dL ⁶
T ⁷	53,225 53,499	53,362			
NSB ⁸	217 242	230	0		
A (MB) ⁹	24,552 24,687	24,620	24,390	100%	0
B	21,704 21,914	21,809	21,579	88%	5
C	16,937 17,449	17,193	16,963	70%	20
D	13,775 12,961	13,368	13,138	54%	50
E	6,189 5,852	6,021	5,791	24%	200
F	3,296 3,257	3,277	3,047	12%	500
G	2,092 2,159	2,126	1,896	8.0%	1,000
Unknowns ¹⁰					
X1	13,053 13,000	13,027	12,797	52%	49
X2	6,227 6,004	6,116	5,886	24%	206
X3	2,990 2,928	2,959	2,729	11%	601

Quality Control Parameters:¹¹

T⁷ = 53,362 cpm

%NSB⁸ = 0.4%

%MB⁹ = 46%

20% Intercept¹² = 273 µg/dL

50% Intercept = 55 µg/dL

80% Intercept = 11 µg/dL

Intraassay Precision (µg/dL)

	Mean ¹	SD ²	CV ³
1	19	0.97	5.1%
2	47	1.8	3.8%
3	190	8.3	4.4%
4	640	34	5.3%

Interassay Precision (µg/dL)

	Mean ¹	SD ²	CV ³
1	18	2.0	11%
2	47	3.6	7.7%
3	206	13	6.3%
4	669	54	8.1%

End-of-Run Effect (µg/dL)

	Tubes ¹ 31–44	Tubes 85–98	Tubes 139–152	Tubes 193–206	Tubes 247–260
1	20	19	18	18	18
2	43	48	44	47	43
3	168	167	169	182	152
4	614	607	618	584	671

Linearity (µg/dL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	16 in 16 ⁵	188	—	—
	8 in 16	91	94	97%
	4 in 16	41	47	87%
	2 in 16	22	24	92%
	1 in 16	13	12	108%
2	16 in 16	202	—	—
	8 in 16	92	101	91%
	4 in 16	48	51	94%
	2 in 16	25	25	100%
	1 in 16	13	13	100%
3	16 in 16	281	—	—
	8 in 16	138	141	98%
	4 in 16	63	70	90%
	2 in 16	31	35	89%
	1 in 16	16	18	89%
4	16 in 16	370	—	—
	8 in 16	170	185	92%
	4 in 16	87	93	94%
	2 in 16	43	46	93%
	1 in 16	24	23	104%
5	16 in 16	528	—	—
	8 in 16	246	264	93%
	4 in 16	138	132	105%
	2 in 16	62	66	94%
	1 in 16	31	33	94%

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EC REP DPC Biermann GmbH
O Sistema de Qualidade da Diagnostic Products Corporation está registado sob ISO 13485:2003.