DUALCOUNT

SOLID PHASE NO BOIL ASSAY
For Vitamin B12 / Folic Acid

DPC®

Solid Phase No Boil Dualcount®

English

Intended Use: Solid Phase No Boil Dualcount is a radioassay designed for the simultaneous measurement of vitamin B12 in serum or plasma, and folic acid in serum, plasma or whole blood. The Solid Phase No Boil Vitamin B12 and Solid Phase No Boil Folic Acid kits are designed for single-analyte determinations of vitamin B12 and folic acid, respectively. All three kits are intended strictly for *in vitro* diagnostic use as aids in clinical diagnosis.

	$\mu\text{Ci} / \text{kBq}$		
	Dualcount kit	Vitamin B12 kit	Folic Acid kit
100 Tubes	KDSP1	KVSP1	KFSP1
^{57}Co Vitamin B12	2/74	2/74	—
^{125}I Folic Acid	6/222	—	6/222
200 Tubes	KDSP2	KVSP2	KFSP2
^{57}Co Vitamin B12	4/148	4/148	—
^{125}I Folic Acid	12/444	—	12/444
500 Tubes	KDSP5	KVSP5	KFSP5
^{57}Co Vitamin B12	10/370	10/370	—
^{125}I Folic Acid	30/1110	—	30/1110

Principle of the Procedure

1. Alkaline Denaturation of Endogenous Proteins

2. Competition for Purified Binder at pH 9.3: B12 analogs do not interfere, since the hog intrinsic factor binder is free of R protein. The folic acid binder has equal affinity for MTHF and PGA.

3. Solid Phase Separation

Total Incubation Time: 2 Hours.

Total counts:

^{125}I -folic acid: approximately 80,000 cpm at lodination.

^{57}Co -vitamin B12: approximately 25,000 cpm.

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C in a refrigerator designated for incoming radioactive materials. Dispose of in accordance with applicable laws.

Do not use reagents beyond their expiration dates.

Some components supplied in this kit may contain human source material and/or other potentially hazardous ingredients which necessitate certain precautions.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chloramphenicol, at concentrations less than 0.1 g/dL has been added as a preservative. Chloramphenicol is known to cause cancer; this disclosure is required by the state of California.

A solution containing sodium hydroxide (NaOH) and potassium cyanide (KCN) is supplied with the kit. Extreme care must be observed to avoid all bodily contact with this reagent, since it contains cyanide and strong alkali.

Water: Use distilled or deionized water.

Radioactivity

A copy of any radioisotope license certificate (Specific or General) issued to a US customer must be on file with Diagnostic Products Corporation before kits or components containing radioactive material can be shipped. These radioactive materials may be acquired by any customer with the appropriate Specific license. Under a General license these radioactive materials may be acquired only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories and hospitals — and strictly for *in vitro* clinical or laboratory tests not

involving external or internal administration of the radioactive material or its radiation to human beings or other animals. Its acquisition, receipt, storage, use, transfer and disposal are all subject to the regulations and a (General or Specific) license of the U.S. Nuclear Regulatory Commission or a State with which the NRC has entered into an agreement for the exercise of regulatory control.

Handle radioactive materials according to the requirements of your General or Specific license. To minimize exposure to radiation, the user should adhere to guidelines set forth in the National Bureau of Standards publication on the *Safe Handling of Radioactive Materials* (Handbook No. 92, issued March 9, 1964) and in subsequent publications issued by State and Federal authorities.

Wipe up spills promptly and decontaminate affected surfaces. Avoid generation of aerosols. Dispose of solid radioactive waste according to license requirements. General licensees (holders of NRC Form 483) may dispose of solid radioactive waste as nonradioactive waste, after removing labeling. Specific licensees (NRC Form 313) should refer to Title 10, Code of Federal Regulations, Part 20. Licensees in Agreement States should refer to the appropriate regulations of their own state. General licensees may dispose of liquid radioactive waste of the type contained in this product through a laboratory sink drain. Licensees must remove or deface labels from empty containers of radioactive materials before disposal of solid waste. Specific licensees may dispose of small quantities of liquid radioactive waste of the type used in this product through a laboratory sink drain. Refer to the appropriate regulations applicable to your laboratory.

Materials Supplied: Initial Preparation

The two components — Vitamin B12/Folic Acid Binder (DSP1) and NaOH/KCN reagent (NHCN) — represent a matched set and are to be used exclusively in combination.

Vitamin B12/Folic Acid Binder and NaOH/KCN reagent supplied in this kit are

not interchangeable with other lots of the same components.

Vitamin B12/Folic Acid Binder (DSP1)

120 mL of a suspension of purified vitamin B12 and folic acid binders immobilized on microcrystalline cellulose particles, with preservative. Stable at 2–8°C for 30 days after opening. Color: orange.

KDSP1, KVSP1, KFSP1: 1 vial

KDSP2, KVSP2, KFSP2: 2 vials

KDSP5, KVSP5, KFSP5: 5 vials

This reagent must be thoroughly suspended before and during use; a magnetic stir bar may be introduced into the vial and stored there. The Binder should be kept in the original wide-mouth vial, which has been treated to prevent adsorption.

NaOH/KCN (NHCN)

15 mL of a solution containing sodium hydroxide (NaOH) and potassium cyanide (KCN), with preservative. Stable at 2–8°C for 30 days after opening. To avoid unnecessary exposure to the atmosphere, recap the bottle immediately after use; otherwise, absorption of CO₂ will gradually result in a change of pH. *Extreme care must be observed to avoid all bodily contact with this reagent.*

KDSP1, KVSP1, KFSP1: 1 vial

KDSP2, KVSP2, KFSP2: 1 vial

KDSP5, KVSP5, KFSP5: 2 vials

Caution! Contains cyanide and strong alkali.

⁵⁷Co Vitamin B12/¹²⁵I Folic Acid Master Tracer (DSP2)

Supplied only in the Dualcount kits.

110 mL of a ready-to-use solution of ⁵⁷Co-vitamin B12 and ¹²⁵I-folic acid. Stable at 2–8°C for 30 days after opening, or until the expiration date marked on the label.

Note: the Master Tracer contains *both* isotopes, whereas the tracers supplied in the single-analyte kits each contain only one isotope.

KDSP1: 1 vial

KDSP2: 2 vials

KDSP5: 5 vials

⁵⁷Co Vitamin B12 (VSP2)

Supplied only in the Vitamin B12 kits.

110 mL of a ready-to-use solution of ⁵⁷Co-vitamin B12. Stable at 2–8°C for 30 days

after opening, or until the expiration date marked on the label.

KVSP1: 1 vial
KVSP2: 2 vials
KVSP5: 5 vials

¹²⁵I Folic Acid (FSP2)

Supplied only in the Folic Acid kits.

110 mL of a ready-to-use solution of ¹²⁵I folic acid, with preservative. Stable at 2–8°C for 30 days after opening, or until the expiration date marked on the label.

KFSP1: 1 vial
KFSP2: 2 vials
KFSP5: 5 vials

Vitamin B12/Folic Acid Calibrators (DSP3–9)

Seven vials, labeled A through G, containing different levels of lyophilized vitamin B12 and folic acid in a protein matrix. Reconstitute the zero calibrator A with **6.0 mL** distilled or deionized water, and the remaining calibrators B through G with **3.0 mL** distilled or deionized water each. Stable at 2–8°C for 30 days after reconstitution, or for 2 months (aliquotted) at –20°C. Avoid excessive exposure to direct light.

KDSP1, KVSP1, KFSP1: 1 set
KDSP2, KVSP2, KFSP2: 1 set
KDSP5, KVSP5, KFSP5: 2 sets

The *lot-specific* values of the calibrators are tabulated in the Radioassay Procedure section. Note that the calibrators supplied with the Solid Phase No Boil Dualcount Kit are *not* interchangeable with those in other DPC kits for vitamin B12 and folic acid. Refer to the vial labels for exact values in appropriate units.

Dithiothreitol Solution (DTT1, 5DTT1)

A solution of dithiothreitol. Stable at 2–8°C for 30 days after opening.

KDSP1, KVSP1, KFSP1: 1 vial x 6 mL (DTT1)
KDSP2, KVSP2, KFSP2: 2 vials x 6 mL (DTT1)
KDSP5, KVSP5, KFSP5: 1 vial x 30 mL (5DTT1)

Anemia Control (ANC)

Lyophilized protein-based control with *lot-specific* vitamin B12 and folic acid concentrations in the deficient range. Reconstitute with **2.0 mL** distilled water.

Stable at 2–8°C for 7 days after reconstitution, or for 2 months (aliquotted) at –20°C. Avoid excessive exposure to direct light.

KDSP1, KVSP1, KFSP1: 1 vial
KDSP2, KVSP2, KFSP2: 1 vial
KDSP5, KVSP5, KFSP5: 1 vial

Materials Required But Not Provided

Gamma Counter: Must be able to discriminate ¹²⁵I from ⁵⁷Co with acceptably low spillover, for simultaneous determinations. A single-channel instrument, set for ⁵⁷Co or ¹²⁵I, is suitable for the single-analyte kits.

Centrifuge: Capable of at least 2000xg. A high-speed, refrigerated model is desirable but not essential. Good results can also be obtained with low-speed, table-top devices, providing the centrifugation time is suitably extended.

Vortex Mixer: The mixer should have the strength to create a vortex in the assay tubes after binder addition. Multiple-tube rack mixers have been found unsuitable for use.

Sterile Distilled Water: For use in the NSB tubes at step 8 of the Radioassay Procedure.

Plastic Containers, with Lids: For preparing the Working Solution.

12x75 mm tubes: Polypropylene tubes (available from DPC) are recommended, but glass and polystyrene tubes have also proved suitable.

Volumetric pipets, 6.0 mL and 3.0 mL: For reconstitution of the calibrators.

Distilled or deionized water: For reconstitution of the calibrators.

Controls: A control with lot-specific values in the deficient range is supplied. Additional materials are available from DPC: CON6 is a tri-level, human serum-based radioassay control containing vitamin B12 and folic acid as two of over 25 assayed constituents.

Micropipet, 200 µL: For sample addition. Disposable-tip pipets with an air interface which prevents the sample from coming into contact with the pipet itself are recommended. The use of positive-

displacement pipets should be avoided because of the risk of carryover.

Micropipet or Dispenser, 50 μ L: For adding the NaOH/KCN. A repeating dispenser has proved reliable and convenient for this reagent addition.

Micropipet or Dispenser, 1,000 μ L: For adding the Working Solution and Vitamin B12/Folic Acid Binder. A repeating dispenser has proved reliable and convenient for these reagent additions.

37°C Waterbath: Do *not* use an oven or a heat block.

Magnetic Stirrer and Stir Bar: For suspending the Vitamin B12/Folic Acid Binder just before and during the addition of this reagent.

Decanting Rack: Foam (or metal) rack capable of holding tubes firmly for simultaneous decanting and draining. Foam racks are available from DPC.

Logit-log graph paper — available from DPC (catalog number: ZP797)

Ascorbic Acid: Required only for whole blood and red cell folic acid determinations. A 1% solution *must be prepared fresh*, on a daily basis. The ascorbic acid crystals should be protected from moisture so that they remain free-flowing, instead of becoming caked from excessive humidity. (Sodium ascorbate is *not* a suitable alternative.)

Foil: Or other means for protecting the hemolysates from light.

Specimen Handling

For folic acid determinations, patient must be in a fasting state.⁴⁰

Sample Type: For vitamin B12 and folic acid determinations, fresh or fresh-frozen serum, heparinized plasma or EDTA plasma may be used. There is no clinically significant interference from the anticoagulants tested.

For whole blood and red cell folic acid determinations, use fresh heparinized or EDTA whole blood. The patient's hematocrit must be known in order to translate whole blood into red cell folic acid results, which are expressed as nanograms per milliliter (ng/mL) of packed red blood cells. Moreover, to calculate the red cell folic acid result exactly, the

patient's serum folic acid level must also be known.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. Solid Phase No Boil Dualcount has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 200 μ L serum, plasma or hemolysate.

Storage and Stability:

Serum or Heparinized Plasma: Storage under suitable conditions is critical for reliable folic acid results. If the sample will *not* be assayed within 8 hours, aliquot and freeze at -20°C : stable for 6–8 weeks.³⁸ Thaw only once (at ambient temperature). Avoid excessive exposure of samples to direct light.

Anti-Folate Drugs: As indicated in the Specificity section, the folate binding protein in the Solid Phase No Boil kit crossreacts to some extent with the "anti-folate" drug methotrexate (MTX), a compound used in the chemotherapy of cancer. Since high-dose infusions of methotrexate may result in circulating levels of 45,000 ng/mL initially and levels in the order of 2,000 ng/mL forty-eight hours later, the kit should *not* be used for folic acid measurements on patients currently receiving the drug. Moreover, great caution must be exercised when interpreting folic acid levels on patients who have recently undergone therapy with methotrexate or any other drug with a similar structure.

Hemolysis: Spuriously elevated serum folic acid results will be obtained to the degree that the sample is hemolyzed, since this vitamin is present in erythrocytes at concentrations many times higher than in plasma. Vitamin B12 results, on the other hand, as measured by the Solid Phase No Boil procedure, are not affected by hemolysis.

Lipemia: The use of an ultracentrifuge is recommended to clear lipemic samples.

High Patient Samples: Patient specimens with vitamin B12 or folic acid

concentrations greater than those of the highest calibrator should be diluted with the zero calibrator and reassayed. Exceptionally high vitamin B12 concentrations — as high as 20,000 pg/mL — may be encountered in a few patients with chronic myelogenous leukemia (CML).^{5,35} Since a carryover of even one microliter from such a sample would cause a spurious increase of 100 pg/mL in the next assay tube, standard precautions should be taken: assaying samples in duplicate, changing pipet tips between samples and avoiding the use of positive-displacement micropipets.

Preparation of Hemolysates: This is required only for folic acid determinations on whole blood. Lyse the cells in freshly obtained heparinized or EDTA whole blood by preparing a 1-in-21 dilution of the sample (fully suspended) in freshly prepared 1% ascorbic acid solution; for example, mix **100 µL** of sample with **2.0 mL** of this solution, then vortex. It is important to use ascorbic acid, rather than sodium ascorbate, and to prepare the 1% solution fresh, on a daily basis.³³

The hemolysate may be assayed immediately; no incubation is required. If not assayed immediately, the hemolysate may be stored at room temperature (15–28°C) in the dark for up to 90 minutes. Hemolysates may be stored at 2–8°C for up to 3 hours, or frozen at –20°C for 2 weeks.^{40,41} Protect from light. Note that the radioassay procedure will not yield meaningful vitamin B12 results for red cells or whole blood, since vitamin B12 is destroyed by ascorbic acid.

Radioassay Procedure

All components must be at room temperature (15–28°C) before use.

- 1 Label eighteen tubes in duplicate: T (total counts), NSB (nonspecific binding), A (maximum binding) and B through G. Label additional tubes, also in duplicate, for controls and patient samples.

Vitamin B12

Calibrator	Approximate pg/mL	Approximate pmol/L
A	0	0
B	50	37
C	100	74
D	300	221
E	600	443
F	1,200	885
G	2,400	1,771

Folic Acid

Calibrator	Approximate ng/mL	Approximate nmol/L
A	0	0
B	0.5	1.1
C	1	2.3
D	3	6.8
E	6	13.6
F	12	27.2
G	24	54.5

- 2 Add **200 µL** of the zero calibrator A into the NSB and A tubes, and **200 µL** of the remaining calibrators B through G into correspondingly labeled tubes. Pipet **200 µL** of each control, serum or plasma sample and hemolysate into the tubes prepared. **Pipet directly to the bottom.**

Patient samples expected to contain concentrations greater than the highest calibrator should be diluted in the zero calibrator before assay.

- 3 **Preparation of the Working Solution:** Prepare the Working Solution *no more than 30 minutes before use*. The amount of each component depends on the number of tubes to be assayed. The volumes required, in *milliliters per assay tube*, are tabulated below. Be sure to multiply these volumes by a number slightly greater than the number of tubes to be assayed. (The components are supplied in volumes sufficient for making up the Working Solution in approximately a 10% excess.)

	Milliliters Per Tube
Dithiothreitol	0.05
Tracer	1.0

- 4** Add **1,000 µL** of freshly prepared Working Solution to all tubes. Vortex.

A repeating dispenser is recommended for the 1,000 µL addition here and at step 8 below. Set the T tubes aside for counting at step 12; they require no further processing.

- 5** Incubate for **30 minutes** at room temperature (15–28°C).

The incubation may be extended up to 45 minutes at room temperature.

- 6** Add **50 µL** of NaOH/KCN to all tubes. Vortex.

Use extreme caution in handling this reagent. A repeating dispenser is recommended for this reagent addition. The tubes should now be blue.

- 7** Incubate for 30 minutes at 37°C.

The incubation may be extended up to *at most* 45 minutes at 37°C. (Longer incubations will result in lower vitamin B12 recoveries.)

- 8** Add **1,000 µL** of Vitamin B12/Folic Acid Binder to all tubes *except the NSB tubes*. To the NSB tubes, add **1,000 µL** of sterile distilled water. **Vortex vigorously.**

The tubes should be vortexed *individually* until a vortex appears in the tube. Multiple-tube rack mixers have been found to be inadequate.

The Binder must be kept *well suspended* while being dispensed: mechanical mixing on a magnetic stirrer is recommended. All tubes should now be purple, except for the T and NSB tubes. Just as much of the Binder as is needed for the assay may be poured off into a glass (not plastic!) beaker, equipped with a stir bar, but the remainder should be stored in the original, silanized vial.

- 9** Incubate for **60 minutes** at room temperature (15–28°C).

- 10** Centrifuge for **at least 15 minutes** at 2000xg or higher.

Higher accelerations are acceptable, but glass or polystyrene tubes may break at higher speeds. Longer centrifugation may be necessary with low-speed, table-top centrifuges.

- 11** Decant the supernatants simultaneously, using a rack, and retain the precipitates for counting.

Thorough, uniform decanting, with removal of all visible moisture, will greatly enhance precision. Using a foam decanting rack, decant the contents of all tubes (except the T tubes) and allow them to drain inverted for 5 minutes. Then tap the tubes lightly to remove all residual droplets.

- 12** Count the precipitate for **1 minute** in a gamma counter.

Note that the precipitate contains the bound fraction, not the free.

Calculation and Quality Control

Serum Samples: To obtain vitamin B12 or folic acid concentrations from a logit-log representation of the calibration curve, first calculate for each pair of tubes the average NSB-corrected counts per minute:

Net Counts = (Average CPM) minus (Average **NSB** CPM)

Then determine the binding of each pair of tubes as a percent of maximum binding (MB), with the NSB-corrected counts of the A tubes taken as 100%:

Percent Bound = (Net Counts / Net **MB** Counts) × 100

The calculation can be simplified by omitting the correction for nonspecific binding; samples within range of the calibrators yield virtually the same results when Percent Bound is calculated directly from Average CPM.

Using logit-log graph paper, plot Percent Bound on the vertical (probability) axis against Concentration on the horizontal (logarithmic) axis for each of the nonzero calibrators, and draw a straight line approximating the path of these points. Results for the unknowns may then be read from the line by interpolation.

Data Reduction: It is good practice to construct a graph of the calibration curve

as a visual check on the appropriateness of the transformation used, even where the calculation of results is handled by computer. See further Davis SE et al. Radioimmunoassay data processing with a small programmable calculator. J Immunoassay 1980;1:15-25; and Dudley RA, et al. Guidelines for immunoassay data reduction. Clin Chem 1985;31:1264-71.

Hemolysates: First calculate the result *R* in nanograms per milliliter of hemolysate, using the folic acid calibration curve. Then multiply by the dilution factor 21 to obtain the whole blood folic acid concentration, in ng/mL. For an approximate measure of the packed red cell concentration, again in ng/mL, multiply the whole blood folic acid concentration by $100/H$, where *H* is the hematocrit *in percent*:

$$\text{Red Cell Folic Acid} \approx 21R \times (100 / H)$$

Strictly speaking, the serum folic acid contribution should be subtracted from the whole blood folic acid concentration before multiplying by $100/H$. Using the patient's serum folic acid level *S*, the exact equation is:

$$\text{Red Cell Folic Acid} = (21R - [S \times ((100 - H) / 100)]) \times (100 / H)$$

The term in the square brackets is, however, in most cases so small compared to the term $21R$ that it may be justifiably neglected.

Controls: Use the control supplied with the kit, and chart results from day to day as described in Westgard JO, et al. A multi-rule chart for quality control. Clin Chem 1981;27:493-501. Repeat samples are a valuable additional tool for monitoring interassay precision.

Q. C. Parameters: We recommend keeping track of these performance measures:

T = Total Counts (as counts per minute)

$\%NSB = 100 \times (\text{Average NSB Counts} / \text{Total Counts})$

$\%MB = 100 \times (\text{Net Counts} / \text{Total Counts})$

And the 20, 50 and 80 percent "intercepts," where

20% = Concentration at 20 Percent Bound, etc.

Centrifugation: Use the formulas below to calculate the acceleration of your centrifuge at a given speed, or the speed

(in revolutions per minute) required to achieve the desired *g* force.

$$g \text{ force} = 28.38 \times (\text{rpm} / 1,000)^2 \times \text{radius [inches]}$$

$$\text{rpm} = 187.7 \times \sqrt{(g \text{ force} / \text{radius [inches]})}$$

or

$$g \text{ force} = 11.18 \times (\text{rpm} / 1,000)^2 \times \text{radius [cm]}$$

$$\text{rpm} = 299.07 \times \sqrt{(g \text{ force} / \text{radius [cm]})}$$

Acid: For illustration only, not for calculating results from another run. Because the calibrator values are *lot-specific*, concentrations listed in the right-most column may not match the values of the calibrators supplied in your shipment. (See "Example Run, Vitamin B12" and "Example Run, Folic Acid" tables.)

Expected Values

Reference Range Study

Fresh serum samples from a total of one hundred healthy individuals were assayed for vitamin B12 and folic acid by the DPC Solid Phase No Boil Dualcount kit. The kit was also used for folic acid determinations on whole blood samples from eighty-eight healthy individuals. The reference groups consisted primarily of laboratory personnel, both men and women, in apparent good health. In each case, the reference range was determined nonparametrically as the central region encompassing 95% of the observations. The reference ranges for the Solid Phase No Boil Dualcount kit suggested by this study are tabulated below.

Serum Vitamin B12:	200 – 950 pg/mL	150 – 700 pmol/L
Serum Folic Acid:	3.0 – 17 ng/mL	7 – 39 nmol/L
Whole Blood Folic Acid:	80 – 320 ng/mL	180 – 725 nmol/L
Red Cell Folic Acid:	175 – 700 ng/mL	395 – 1,585 nmol/L

The folic acid reference ranges were confirmed by further studies at two locations on a total of 162 individuals. Note that results near the lower limit of normal — e.g. red cell folic acid results in the 100 to 200 ng/mL range — should be interpreted with caution, in the light of other clinical and laboratory data. Also,

serum folic acid results should be interpreted in the light of *red cell* determinations.

The reference ranges tabulated above were obtained in the United States. Studies performed in Sweden and Norway suggest reference ranges of 100–350 ng/mL and 130–350 ng/mL, respectively, for red cell folic acid. Consider these limits as guidelines only. Each laboratory should establish its own reference ranges.

Vitamin B12: There is no overlap between the 200–950 pg/mL vitamin B12 range for healthy individuals and the absolute range encountered for the confirmed vitamin B12-deficient samples. (See the "Expected Values: Comparative Discrimination Plot" in Tables and Graphs section.) 200 pg/mL may thus be taken as the lower limit of normal for vitamin B12 values measured by the Solid Phase No Boil kit. Note that values less than 200 pg/mL and greater than 124 pg/mL — the highest Solid Phase No Boil value observed for the vitamin B12-deficient samples — should be treated as statistically subnormal; further investigation may be warranted in such cases.

Common causes of high vitamin B12 levels include liver disease, myeloproliferative disease (with chronic myelogenous leukemia as a special case) and the use of multivitamin supplements.

Confirmed Vitamin B12-Deficient Samples

DPC has studied over one hundred vitamin B12-deficient samples confirmed as vitamin B12-deficient by independent clinical and hematological criteria.¹⁴ Two reference methods were used: a microbiological method based on *Lactobacillus leichmannii* and DPC's Charcoal Boil Dualcount kit, a radioassay based on heat denaturation, charcoal separation and purified intrinsic factor. On most, but not all of the samples, vitamin B12 values by all three assay methods were available. These clinical results are displayed in a comparative discrimination plot showing the vitamin B12-deficient range for all three methods side by side, and the extent of overlap between normals and deficient for the two

radioassays. (See the "Expected Values: Comparative Discrimination Plot" in the Tables and Graphs section.)

Performance Data

See Tables and Graphs for data representative of the Solid Phase No Boil Dualcount kit's performance. Vitamin B12 results are expressed in pg/mL. Folic Acid results are expressed in ng/mL. Except as noted, all results were generated on serum or on heparinized whole blood.

Vitamin B12 Conversion Factor:

pg/mL $\times$ 0.7378 $\rightarrow$ pmol/L

Folic Acid Conversion Factor:

ng/mL $\times$ 2.266 $\rightarrow$ nmol/L

Calibration Range:

Vitamin B12:
approximately 50 – 2,400 pg/mL
The calibrator values are *lot-specific*.

Folic Acid:
approximately 0.5 – 24 ng/mL.
The calibrator values are *lot-specific*.

Analytical Sensitivity:

Vitamin B12: 34 pg/mL.

Folic Acid: 0.3 ng/mL.

Intraassay Precision (Within-Run):

Statistics were calculated for samples from the results of twenty pairs of tubes in a single assay. (See "Intraassay Precision" tables.)

Interassay Precision (Run-to-Run):

Statistics were calculated from the results of pairs of tubes in eighteen different assays for 14 samples. (See "Interassay Precision" tables.)

Cumulative Interassay Precision:

Cumulative interassay statistics were generated for nine controls over the course of four months. For each control, one pair of tubes was processed in each of approximately 100 assays. The statistics reflect a number of lot changes and the benchwork of several technicians. (See "Cumulative Interassay Precision" tables.)

End-of-Run Effect: Pairs of tubes were spaced throughout an assay involving over 350 tubes. Using a repeating pipet, addition of the Vitamin B12/Folic Acid Binder was completed in eight minutes. The results indicate no position effect due

to delays in the addition of reagents, even in assays of this magnitude.

Specificity

Binder Purity: Two experiments were performed to demonstrate the purity of the intrinsic factor (IF) binder employed in the Solid Phase No Boil Vitamin B12 radioassay.

- Overnight incubation of the Vitamin B12/Folic Acid Binder with serum containing antibodies specific for intrinsic factor resulted in complete destruction of this reagent's ability to bind vitamin B12. Had R protein been present, it would have retained at least some of its vitamin B12-binding capacity.
- Substitution of the Vitamin B12/Folic Acid Binder in a radioassay system based on heat denaturation and coated charcoal, and operating at an assay pH of approximately 2, again produced complete inactivation of the Binder. Had it contained R protein, it would have retained at least some of its activity under these acidic conditions.

Crossreactivity: The specificity of the purified intrinsic factor binder was assessed by assaying a cobinamide preparation. The results show that, even in concentrations as high as 50,000,000 pg/mL, cobinamide causes virtually no reduction in binding of the radiolabeled vitamin B12, i.e. it yields an apparent vitamin B12 concentration at the detection limit; crossreactivity with cobinamide is therefore less than 0.0001%. (See "Specificity 1" table.)

Also studied was the crossreaction between methotrexate (MTX), an antifolate drug, and the folate-binding protein employed in the Solid Phase No Boil kit. The percent crossreactivity, even though numerically small, is clinically significant, since methotrexate is often administered in extremely high doses. Accordingly, folic acid measurements are contraindicated in patients who are undergoing, or who have recently undergone, therapy with this or any structurally related antifolate drug. (See "Specificity 2" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" tables for representative data.)

Recovery: Six spiking solutions were made using the zero calibrator as diluent. Three solutions (A, B and C) were made to represent 10,000, 20,000 and 30,000 pg/mL of vitamin B12. The other three (D, E and F) were made to represent 60, 120 and 180 ng/mL of folic acid. A 50 µL aliquot of each solution was spiked into 950 µL aliquots of six patient samples, for a spiking ratio of 1-in-20, leaving the matrix of the sample relatively intact. Values tabulated represent recovery of the spiking solution only: the contribution of the serum sample has been subtracted out. Numbers in parentheses are the concentrations of the original unspiked serum samples. (See "Recovery" tables.)

Protein Effect: To determine the effects of variations in protein content on the Solid Phase No Boil procedure, two patient serum samples were spiked with varying amounts of human serum albumin (HSA) and normal human IgG. All samples, spiked and unspiked, were then processed in a single assay. The results show that even large variations in protein content produce only small variations in the observed result. We may conclude that there is no clinically significant protein effect in the Solid Phase No Boil Dualcount system.

Alternate Sample Type: Samples from healthy individuals were collected simultaneously into plain, heparinized and EDTA tubes, then assayed by the Solid Phase No Boil Dualcount procedure. The results indicate that there is no clinically significant interference from the anticoagulants tested. (See "Alternate Sample Type" table.)

Patient Blanks: To verify that DPC's Solid Phase No Boil Dualcount does not suffer from nonspecific binding problems of the kind described in the literature,^{6,24,27,35} an assay was performed which included a calibration curve, calibrator blanks (Binder omitted) and two sets of patient blank tubes. From one set of the patient blank tubes, Solid Phase Binder was simply omitted during the course of the assay. The other set was assayed substituting for the Binder a "Cellulose Blank" reagent: a suspension of cellulose particles identical

to the Solid Phase Binder but for the absence of purified intrinsic factor and folate binding protein. (For purposes of kit evaluation, cellulose blank reagent may be obtained on request from DPC.)

The results indicate that (a) omission of the binder is a suitable alternative to the use of cellulose blank reagent, and (b) nonspecific binding is *uniformly low* for both patient samples and calibrators in the Solid Phase No Boil Dualcount system. Hence a single set of NSB tubes — zero calibrator tubes with binder omitted — suffices for the routine determination of nonspecific binding. Indeed, the NSB tubes may be omitted altogether: calculations show that correction for nonspecific binding can be neglected without any significant effect on results for samples within range of the calibrators. (See "Patient Blanks" tables.)

The Denaturation Step

To test the effectiveness of the 37°C alkaline denaturation step in the Solid Phase No Boil Dualcount system, independent studies were performed with the kit on samples from patients with chronic myelogenous leukemia (CML) and other myeloproliferative disorders, and on pernicious anemia (PA) samples which had been tested for the presence of anti-intrinsic factor blocking antibodies.

CML Samples: In an experiment with the Solid Phase No Boil kit, charcoal precipitation was used to assess nonspecific binding in patient blank tubes at the end of the second incubation. The results showed virtually complete precipitation of the tracer in all tubes; the behavior of CML patient blank tubes in this experiment was similar to that of blanks generated with healthy patient samples. (The two CML samples had vitamin B12 levels on the order of 3,000 and 5,000 pg/mL.) This experiment shows that the 37°C alkaline denaturation step employed in the Solid Phase No Boil system effectively inactivates even the high levels of vitamin B12-binding protein sometimes found in CML patients; for if endogenous carrier protein had escaped denaturation, it would have bound to the tracer, causing a significant amount of radioactivity to remain in the supernatant after charcoal precipitation.

PA Samples: Blocking antibodies specific for intrinsic factor are present in more than half of all patients with pernicious anemia. Fourteen PA samples, seven of them positive for the presence of such antibodies, were assayed both by the Solid Phase No Boil kit and the Charcoal Boil Dualcount kit.⁸ The results show that the 37°C alkaline denaturation step in the Solid Phase No Boil kit is as effective as a boiling step in denaturing anti-intrinsic factor antibodies which would otherwise cause spuriously elevated (false normal) results in many cases of PA.

As further confirmation, four antibody-positive and four antibody-negative PA samples were each subjected to the 37°C alkaline denaturation procedure of the Solid Phase No Boil kit and to a control procedure involving substitution of water for the sodium hydroxide.⁸ Antibody-positive samples subjected to the control procedure retained their ability to block the activity of intrinsic factor, whereas antibody-positive samples subjected to the kit's alkaline denaturation procedure did not: the latter showed behavior indistinguishable from that of the antibody-negative samples, indicating that the kit's "no boil" denaturation step had effectively inactivated the anti-intrinsic factor antibodies present in the positive samples. Note that while the experiments reported here involved many samples with extremely high antibody titers, they cannot rule out the possibility that antibodies in a rare sample with a still higher titer might not be fully denatured by the "no boil" approach. Where such a doubt exists, the "no boil" result can be confirmed by reassaying the sample by the alternate "boil" procedure.

Effect of Incubation Time on

Hemolysates: Blood samples were collected from four individuals into heparinized and EDTA tubes. Hemolysates were prepared by incubating aliquots of each sample with the 1% ascorbic acid solution for increasing lengths of time (0 to 90 minutes in the dark at room temperature (15–28°C) prior to assay. The table shows no significant difference in red cell folic acid results for lysis times ranging from zero to 90 minutes. Accordingly, hemolysates may be assayed immediately following the ascorbic acid dilution step. The table also

shows that results are not influenced by the choice of anticoagulant: either heparinized or EDTA whole blood may be used. (See "Effect of Incubation Time on Hemolysates" table.)

Method Comparison for B12: The assay was compared to DPC's IMMULITE Vitamin B12 assay on 98 patient serum samples. (Concentration range: approximately 128 to 1,037 pg/mL. See Method Comparison graph.) By linear regression:

(SPNB) = 0.98 (IML) – 5.3 pg/mL
r = 0.982

Means:

371 pg/mL (Solid Phase No Boil)
385 pg/mL (IMMULITE – LKVB)

Method Comparison for Folic Acid: The assay was compared to DPC's IMMULITE Folic Acid assay on 71 patient serum samples. (Concentration range: approximately 1.4 to 20 ng/mL. See Method Comparison graph.) By linear regression:

(SPNB) = 0.92 (IML) + 0.71 ng/mL
r = 0.985

Means:

6.2 ng/mL (Solid Phase No Boil)
5.9 ng/mL (IMMULITE – LKFO)

References

1) Allen RH. Clinical role and current status of serum cobalamin (Vitamin B12) assays. *Ligand Quarterly* 1981 Fall;4(3):37-44, 67. 2) Allen RH. Quantitative testing for vitamin B12. *United States Patent* 4,188,189. 1982 Feb 12. 3) Allen RH. Personal communication. 4) Beck WS. More on the vitamin B12 problem. *Ligand Quarterly* 1979 Spring;2(12):10. 5) Chanarin I. *The megaloblastic anaemias*. 2nd ed. Oxford: Blackwell Scientific, 1979. 6) Chen I-W, et al. Clinical significance of serum vitamin B12 measured by radioassay using pure intrinsic factor. *J Nucl Med* 1981;22:447-51. 7) Chen I-W, et al. Semi-automated system for simultaneous assays of serum vitamin B12 and folic acid in serum evaluated. *Clin Chem* 1982;28:2161-5. See also *J Clin Chem Clin Biochem* 1981;19:635. 8) Chen I-W. Personal communication. 9) Cohen KL, Donaldson RM. Unreliability of radiodilution assays as screening tests for cobalamin (vitamin B12) deficiency. *J Am Med Assoc* 1980;244:1942-5. 10) Colman N. Technological factors affecting the use of radioassays in the assessment of nutritional anemias. *Ligand Quarterly* 1981 Fall;4(3):31-6, 67. 11) Conrad ME, et al. Iron, folic acid and vitamin B12. In: Fairbanks VF, editor. *Current hematology*. Vol 1. New York: John Wiley, 1981:

123-90. 12) Cooper BA, Whitehead VM. Evidence that some patients with pernicious anemia are not recognized by radiodilution assay for cobalamin in serum. *New Engl J Med* 1978;299:816-8. See further: editorial by Donaldson RM, *ibid*, 827-8. 13) Cooper BA, et al. Vitamin B12 assay: an evaluation of radiodilution assay using cobinamide to increase specificity. In: Zagalak B, Friedrich W, editors. *Vitamin B12*. Berlin: Walter de Gruyter, 1979: 649-54. See also *Clin Chem* 1979;25(6):abstract no. 369. 14) DPC, Data on file. 15) Food and Drug Administration (FDA), Letters from the Bureau of Medical Devices (BMD) to manufacturers of vitamin B12 radioassay kits (Apr 1980, Jan 1981). 16) Gilbert RB, Mailhot EA. Cobalamin radioassays. *J Am Med Assoc* 1981;246:734-5. 17) Herbert V. The nutritional anemias. *Hosp Pract* 1980 Mar;15(3):65-9, 73-6, 81-3, 87-9. 18) Introduction to folate by radioassay. *Mini-Seminar* 1 no 4. Los Angeles: DPC, 1980. 19) Kolhouse JF, et al. Cobalamin analogues are present in human plasma and can mask cobalamin deficiency because current radioisotope dilution assays are not specific for true cobalamin. *New Engl J Med* 1978;299:785-92. 20) Kondo H, et al. Presence and formation of cobalamin analogues in multivitamin-mineral pills. *J Clin Invest* 1982;70:889-98. 21) Kubasik NP, Ricotta M, Sine HE. Commercially supplied binders for plasma cobalamin (vitamin B12) analysis compared to microbiological assay. *Clin Chem* 1980;26:598-600. 22) Laso F, et al. Intr clinique du dosage de l'acide folique par radio-essai. *Schweiz Med Wochenschr* 1978;108:1393-8. 23) Lee-Own V, et al. Formation of a vitamin B12-serum complex on heating at alkaline pH. *Clin Chim Acta* 1979;93:239-44. 24) LeFebvre RJ, Virji AS, Mertens BF. Erroneously low results due to high nonspecific binding encountered with a radioassay kit that measures "true" serum vitamin B12. *Am J Clin Pathol* 1980;74:209-13. 25) Linnell JC. The value of radioisotopic assays for "serum B12" in the diagnosis of cobalamin deficiency disorders. *Clin Lab Haematol* 1981;3:99-106. 26) Magnus EM, et al. Folate and vitamin B12 (cobalamin) blood levels in elderly persons in geriatric homes. *Scand J Haematol* 1982;28:360-6. 27) Mertens BF. CAP/ASCP presentation. March 1980. Reported in: Byrne M. Caution advised on B12 kits: lower values may not be "truer." *Clin Chem News* 1980 Apr;6(4):1, 5. 28) Mitchell C, et al. Chloral hydrate interferes with radioassay of vitamin B12. *Clin Chem* 1981;27:1480-1. 29) Mollin DL, et al. The serum vitamin B12 level: its assay and significance. *Clinics in Haematol* 1976;5(3):521-46. 30) Mollin DL, Ross GIM. The pathophysiology of vitamin B12 deficiency in the megaloblastic anaemias. In: Heinrich HC, editor. *Vitamin B12 and intrinsic factor*. Stuttgart: Enke, 1961: 413. 31) Nash BJ, Lim-Yin YW. A discrepancy between plasma and serum values for vitamin B12 and folate. *Med Lab Sci* 1979;36:297-8. 32) National Committee for Clinical Laboratory Standards (NCCLS).

Guidelines for evaluating a B12 (cobalamin) assay. (Villanova, PA: NCCLS proposed standard PLSA-12, March 1979). Revised as "NCCLS recommendations with respect to vitamin B12 assays" – presented to the Food and Drug Administration (FDA) panel, 9 July 1979. 33) Netteland B, Bakke OM. Inadequate sample preparation as a source of error in determination of erythrocyte folate by competitive binding radioassay. Clin Chem 1977;23:1505-6. See further: letters by Gutcho S. and by Bakke OM, *ibid*, 1978;24:388-9. 34) R-proteins. Mini-Seminar 2 no. 1. Los Angeles: DPC, 1980. 35) Zucker RM, Podell E, Allen RH. Multiple problems with current no-boil assays for serum cobalamin. Ligand Quarterly 1981 Fall;4(3):52-8. See further: letters by Ziering S, et al, and by Zucker RM, et al, *ibid*, 60-3, and by Allen RH, *ibid*, 1982 Spring;5(1):48-9. 36) El Shami AS, et al. More on vitamin B12 results as measured with boil and no boil kits. Clin Chem 1983;29:2115-6. 37) Carmel R, Sinow RM, Karnaze DS. Atypical cobalamin deficiency: subtle biochemical evidence of deficiency is commonly demonstrable in patients without megaloblastic anemia and is often associated with protein-bound cobalamin malabsorption. J Lab Clin Med 1987;109:454-63. 38) Karnaze DS, Carmel R. Low serum cobalamin levels in primary degenerative dementia: do some patients harbor atypical cobalamin states? Arch Int Med 1987;147:429-31. 39) Tietz NW, editor. Clinical guide to laboratory tests. 3rd ed. Philadelphia: WB Saunders, 1995:246, 636. 40) Burtis CA, Ashwood ER, editors. Tietz textbook of clinical chemistry. 2nd ed. Philadelphia: WB Saunders, 1994:2056. 41) Kubasik NP, Graham M, Sine HE. Storage and stability of folate and vitamin B-12 in plasma and blood samples. Clinica Chimica Acta 1979;95:147-149.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782, 310.645.8200
Fax: 310.645.9999. To place an order:
Tel: 800.372.1782. Fax: 800.234.4372.
Outside the United States, contact your National Distributor.

The Quality System of Diagnostic Products Corporation is registered to ISO 13485:2003.

Tables and Graphs

Example Run, Vitamin B12 (pg/mL)

Tube ¹	Duplicate CPM ²	Average CPM ³	Net CPM ⁴	Percent Bound ⁵	Approx. B12 pg/mL ⁶
T ⁷	22,379 22,214	22,297			
NSB ⁸	112 151	132	0		
A (MB) ⁹	9,312 9,698	9,505	9,373	100%	0
B	8,377 8,341	8,359	8,227	88%	50
C	7,735 7,446	7,591	7,459	80%	100
D	5,748 6,449	6,099	5,967	64%	300
E	4,223 4,152	4,188	4,056	43%	600
F	2,832 2,743	2,788	2,656	28%	1,200
G	1,736 1,757	1,747	1,615	17%	2,400
Unknowns ¹⁰					
X1	7,007 6,985	6,996	6,864	73%	161
X2	4,622 4,629	4,626	4,494	48%	510
X3	3,239 3,239	3,239	3,107	33%	983

Quality Control Parameters:¹¹

T⁷ = 22,297 cpm
%NSB⁸ = 0.59%
%MB⁹ = 42%
20% Intercept¹² = 2,029 pg/mL
50% Intercept = 468 pg/mL
80% Intercept = 108 pg/mL

Example Run, Folic Acid (ng/mL)

Tube ¹	Duplicate CPM ²	Average CPM ³	Net CPM ⁴	Percent Bound ⁵	Approx Folic Acid ng/mL ⁶
T ⁷	68,548 68,625	68,587			
NSB ⁸	279 442	361	0		
A (MB) ⁹	30,232 30,553	30,393	30,032	100%	0
B	25,506 25,907	25,707	25,346	84%	0.5
C	22,556 22,109	22,333	21,972	73%	1.0
D	14,711 14,624	14,668	14,307	48%	3.0
E	9,929 9,551	9,740	9,379	31%	6.0
F	6,606 6,196	6,401	6,040	20%	12
G	3,818 3,894	3,856	3,495	12%	24
Unknowns ¹⁰					
X1	19,129 18,897	19,013	18,652	62%	1.7
X2	11,273 11,309	11,291	10,930	36%	4.9
X3	5,817 5,957	5,887	5,526	18%	13

Quality Control Parameters:¹¹

T⁷ = 68,587 cpm

%NSB⁸ = 0.53%

%MB⁹ = 44%

20% Intercept¹² = 12 ng/mL

50% Intercept = 2.8 ng/mL

80% Intercept = 0.66 ng/mL

Intraassay Precision

Vitamin B12	Mean ¹ pg/mL	SD ²	CV ³
1	211	14	6.6%
2	240	12	5.0%
3	300	13	4.3%
4	522	30	5.7%
5	823	22	2.7%
6	927	18	1.9%

Folic Acid	Mean ¹ ng/mL	SD ²	CV ³
1	2.3	0.12	5.2%
2	2.5	0.09	3.6%
3	5.4	0.23	4.3%
4	5.4	0.33	6.1%
5	11.9	0.29	2.4%
6	13.0	0.33	2.5%

Interassay Precision

Vitamin B12	Mean ¹ pg/mL	SD ²	CV ³
Patients ⁴			
1	222	20	8.9%
2	258	23	5.4%
3	323	25	7.8%
4	439	29	6.5%
5	541	29	5.4%
6	762	29	3.9%
7	805	37	4.6%
8	971	40	4.1%
Controls ⁵			
1	201	19	9.2%
2	242	19	8.0%
3	309	23	7.4%
4	527	33	6.3%
5	835	36	4.3%
6	904	27	3.0%

Folic Acid	Mean ¹ ng/mL	SD ²	CV ³
Patients⁴			
1	2.4	0.18	7.5%
2	3.3	0.26	7.9%
3	3.6	0.24	6.6%
4	4.5	0.24	5.3%
5	4.8	0.37	7.7%
6	4.9	0.22	4.5%
7	6.2	0.27	4.4%
8	9.4	0.44	4.7%
Controls⁵			
1	2.4	0.16	6.7%
2	2.5	0.18	7.2%
3	5.5	0.25	4.6%
4	5.5	0.29	5.3%
5	12.1	0.47	3.9%
6	12.3	0.72	5.8%

Cumulative Interassay Precision

Vitamin B12	Mean ¹ pg/mL	SD ²	CV ³
1	190	21	11.1%
2	221	20	9.0%
3	233	28	12.0%
4	319	29	9.1%
5	365	25	6.8%
6	506	35	6.9%
7	823	43	5.2%
8	858	40	4.7%
9	924	43	5.2%

Folic Acid	Mean ¹ ng/mL	SD ²	CV ³
1	1.9	0.28	15.0%
2	2.3	0.22	9.6%
3	2.4	0.22	9.3%
4	5.1	0.45	8.9%
5	5.3	0.38	7.2%
6	5.6	0.39	7.0%
7	11.5	0.59	5.1%
8	12.0	0.66	5.5%
9	12.2	0.71	5.8%

Specificity 1

Cobinamide ¹	%B/B ₀ ²	Apparent Vitamin B12 Conc. ³	% Cross- reactivity ⁴
50,000 ng/mL	93%	34 pg/mL	0.0001%

Specificity 2


Methotrexate ¹	%B/B ₀ ²	Apparent Folic Acid Conc. ³	% Cross- reactivity ⁴
100 ng/mL	55.6%	3.7 ng/mL	3.7%

Linearity

Vitamin B12 (pg/mL)


	Dilution ¹	Observed ²	Corrected ³
1	4 in 4 ⁴	1,359	1,359
	3 in 4	992	1,323
	2 in 4	623	1,246
	1 in 4	279	1,116
	Mean (CV) ⁵ :		1,261 (8.6%)
2	4 in 4	1,935	1,935
	3 in 4	1,424	1,899
	2 in 4	952	1,904
	1 in 4	472	1,888
	Mean (CV):		1,907 (1.0%)
3	4 in 4	3,399	3,399
	3 in 4	2,640	3,520
	2 in 4	1,843	3,686
	1 in 4	890	3,560
	Mean (CV):		3,541 (3.3%)

Method Comparison (Vitamin B12)

(SPNB) = 0.98 (IML) – 5.3 pg/mL
 $r = 0.982$

Method Comparison (Folic Acid)

(SPNB) = 0.92 (IML) + 0.71 ng/mL
 $r = 0.985$

Linearity

Folic Acid (ng/mL)

	Dilution ¹	Observed ²	Corrected ³
1	4 in 4 ⁴	7.0	7.0
	3 in 4	5.4	7.2
	2 in 4	3.7	7.4
	1 in 4	2.1	8.4
Mean (CV) ⁵ :		7.5 (8.3%)	
2	4 in 4	8.8	8.8
	3 in 4	6.3	8.4
	2 in 4	4.4	8.8
	1 in 4	2.5	10.0
Mean (CV):		9.0 (7.7%)	
3	4 in 4	10.1	10.1
	3 in 4	7.9	10.5
	2 in 4	5.4	10.8
	1 in 4	2.8	11.2
Mean (CV):		10.7 (4.4%)	

Recovery

Vitamin B12 (pg/mL)

Amount added: ¹	500	1,000	1,500
Sample ²	A	B	C
1 (400)	494	926	1,419
2 (444)	457	966	1,408
3 (521)	448	975	1,464
4 (561)	481	914	1,416
5 (569)	481	947	1,413
6 (598)	520	921	1,389
Mean ± SD: ³	480 ± 26	942 ± 25	1,418 ± 25
% Recovery: ⁴	96 ± 5%	94 ± 3%	95 ± 2%

Folic Acid (ng/mL)

Amount added: ¹	3	6	9
Sample ²	D	E	F
1 (2.2)	3.0	6.3	9.1
2 (2.9)	3.5	6.3	9.6
3 (3.0)	3.3	6.2	9.8
4 (4.1)	3.3	6.2	9.6
5 (4.2)	3.1	5.8	8.9
6 (5.8)	2.9	6.0	9.5
Mean ± SD: ³	3.2 ± 0.22	6.1 ± 0.20	9.4 ± 0.34
% Recovery: ⁴	107 ± 7%	102 ± 3%	104 ± 4%

Alternate Sample Type

Analyte ¹	x	y	Slope ²	Inter-cept ³	r	n
Vitamin B12	Serum ⁴	Heparin ⁵	1.09	−34	.991	29
Vitamin B12	Serum	EDTA	0.97	−7	.985	29
Folic Acid	Serum	Heparin	0.98	+ 0.1	.994	30
Folic Acid	Serum	EDTA	0.97	+ 0.2	.993	30
RBC Folate	Heparin	EDTA	0.97	+ 8	.994	30

Patient Blanks

	Patient Blanks ¹	
	Vitamin B12	Folic Acid
Binder Omitted ²	1.13 ± 0.34%	1.64 ± 0.19%
Cellulose Blanks ³	1.14 ± 0.04%	1.19 ± 0.11%
Mean ± SD; Percent of Total Counts ⁴		
	Binder Omitted ²	
	Vitamin B12	Folic Acid
Calibrators ⁵	0.97 ± 0.05%	1.01 ± 0.01%
Patients ⁶	1.13 ± 0.34%	1.64 ± 0.19%
Mean ± SD; Percent of Total Counts ⁴		

Effect of Incubation Time on Hemolysates

	0 min	15 min	30 min	60 min	90 min
Heparin ¹					
1	314	316	308	318	335
2	175	172	184	173	175
3	178	180	184	172	185
4	133	135	130	131	131
Mean: ²	200	201	202	199	207
EDTA					
1	294	306	312	329	312
2	180	165	176	171	174
3	177	176	180	192	178
4	132	132	127	135	129
Mean:	196	195	199	207	198

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
Expected Values: Comparative Discrimination Plot

Deutsch. Solid Phase No Boil: Solid Phase No Boil. Charcoal Boil: Charcoal Boil. Reference Group (Adult Normals): Referenz Gruppe (Normale Erwachsene). Not shown: Nicht gezeigt. Normals: Normale. Confirmed Vitamin B12 Deficients: bestätigter Vitamin B12 Mangel.

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O Sistema de Qualidade da Diagnostic Products Corporation está registado sob ISO 13485:2003.
EC REP DPC Biemann GmbH