IF bAb

SOLID PHASE

DPC[®]

Solid Phase IF bAb

English

Intended Use: IF bAb is a solid phase, ^{57}Co radioassay designed for the detection of intrinsic factor blocking antibody in serum. It is intended strictly for *in vitro* diagnostic use as an aid in the diagnosis of pernicious anemia.

Catalog Number: **KISPZ** (50 tubes), **KISP1** (100 tubes).

The 50-tube kit contains less than 1 microcurie (37 kilobecquerels) of radioactive ^{57}Co -Vitamin B12, and the 100-tube kit contains less than 2 microcuries (74 kilobecquerels).

Summary and Explanation of the Test

Pernicious anemia is one of the major causes of vitamin B12 ("B12") deficiency and megaloblastic anemia.⁶⁻⁸ It is a disease of the stomach in which secretion of intrinsic factor is severely reduced or absent, resulting in malabsorption of B12.

In view of its association with a variety of antibodies, including parietal cell antibody and at least two types of anti-intrinsic factor antibody, pernicious anemia appears to be an autoimmune process.^{8,15,20} Intrinsic factor *blocking* antibody, which prevents the binding of vitamin B12 to the intrinsic factor molecule, is present in more than half of all patients with pernicious anemia and only rarely encountered in other conditions.^{11,17,19} Other types of anti-intrinsic factor antibody are less common. Parietal cell antibody, on the other hand, is present in a larger percentage of pernicious anemia patients; but it is also associated with a number of other disorders, making it a less specific marker for pernicious anemia.

Qualitative assay for intrinsic factor blocking antibody can therefore make an important contribution to the diagnosis and differential diagnosis of B12 deficiency.^{3,4,11,13,18} In the context of a low or borderline B12 result, where other clinical and hematological findings are compatible with a diagnosis of B12

deficiency, the presence of intrinsic factor blocking antibody can be taken as confirmation of this diagnosis and at the same time as an indication of its cause. In fact, the combination of megaloblastic anemia, low serum B12 and serum antibodies to intrinsic factor is essentially diagnostic of pernicious anemia and obviates further investigation, including Schilling tests.¹⁶

A negative result, on the other hand, cannot rule out the possibility of pernicious anemia since blocking antibody is *not* demonstrable in nearly half of all patients with this disorder.^{4,11}

Principle of the Procedure

DPC's IF bAb procedure detects anti-intrinsic factor blocking antibody in serum by its effect on the vitamin B12-binding capacity of a solid-phase intrinsic factor binder. Two steps are involved:

- 1 Serum and binder are incubated together for a fixed time under conditions which leave the equilibrium between free and bound B12 in the patient sample essentially undisturbed. Intrinsic factor blocking antibody, if present, will interact with the solid-phase binder, reducing the number of available B12-binding sites. (The influence of endogenous B12 is minimal, since at least ninety-nine percent circulates in bound form.)^{1,5,14}
- 2 The solid-phase binder is isolated by centrifuging and decanting, and then resuspended in the presence of radiolabeled B12, which binds to sites that are still available. Finally, the bound fraction is isolated and counted; and the counts are compared to those of antibody-positive and antibody-negative reference materials.

Total Incubation Time: 1 Hour.

Total counts: ^{57}Co -Vitamin B12 approximately 25,000 cpm

Reagents to Pipet: 3

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C in a refrigerator designated for incoming radioactive materials. Dispose of in accordance with applicable laws.

Do not use reagents beyond their expiration dates.

Some components supplied in this kit may contain human source material and/or other potentially hazardous ingredients which necessitate certain precautions.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Water: Use distilled or deionized water.

Radioactivity

A copy of any radioisotope license certificate (Specific or General) issued to a US customer must be on file with Diagnostic Products Corporation before kits or components containing radioactive material can be shipped. These radioactive materials may be acquired by any customer with the appropriate Specific license. Under a General license these radioactive materials may be acquired only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories and hospitals — and strictly for *in vitro* clinical or laboratory tests not involving external or internal administration of the radioactive material or its radiation to human beings or other animals. Its acquisition, receipt, storage, use, transfer and disposal are all subject to the regulations and a (General or Specific) license of the U.S. Nuclear Regulatory Commission or a State with which the NRC has entered into an agreement for the exercise of regulatory control.

Handle radioactive materials according to the requirements of your General or Specific license. To minimize exposure to radiation, the user should adhere to

guidelines set forth in the National Bureau of Standards publication on the *Safe Handling of Radioactive Materials* (Handbook No. 92, issued March 9, 1964) and in subsequent publications issued by State and Federal authorities.

Wipe up spills promptly and decontaminate affected surfaces. Avoid generation of aerosols. Dispose of solid radioactive waste according to license requirements. General licensees (holders of NRC Form 483) may dispose of solid radioactive waste as nonradioactive waste, after removing labeling. Specific licensees (NRC Form 313) should refer to Title 10, Code of Federal Regulations, Part 20. Licensees in Agreement States should refer to the appropriate regulations of their own state. General licensees may dispose of liquid radioactive waste of the type contained in this product through a laboratory sink drain. Licensees must remove or deface labels from empty containers of radioactive materials before disposal of solid waste. Specific licensees may dispose of small quantities of liquid radioactive waste of the type used in this product through a laboratory sink drain. Refer to the appropriate regulations applicable to your laboratory.

Materials Supplied: Initial Preparation

The components are *not* interchangeable with those in any of DPC's vitamin B12 and folic acid radioassay kits.

IF bAb Binder (ISP1)

35 mL ready-to-use suspension of purified intrinsic factor immobilized on microcrystalline cellulose particles, with preservative. Stable at 2–8°C for 30 days after opening. *Do not transfer the Binder to a plastic container.*

KISPZ: 1 vial. **KISP1:** 2 vials.

Both during storage and in use, the Binder should be kept in its original wide-mouth container, which has been treated to prevent adsorption: the dimensions of the bottle allow for inserting a dispenser and for storing a magnetic stir bar.

⁵⁷Co Vitamin B12 (ISP2)

Lyophilized cobalt-57-labeled vitamin B12. Reconstitute by adding **3.0 mL** of *sterile distilled* water using a volumetric pipet. Let

stand for 10 minutes, then mix by *gentle* inversion. Stable at 2–8°C for 30 days after reconstitution, or until the expiration date marked on the vial. The reconstituted ⁵⁷Co-Vitamin B12 must be diluted 1-in-10 with IF bAb Buffer just before use (step 6).
KISPZ: 1 vial. **KISP1:** 2 vials.

IF bAb Buffer (ISPB)

35 mL of a ready-to-use buffer. Stable at 2–8°C for 30 days after opening.

KISPZ: 1 vial. **KISP1:** 2 vials.

Positive and Negative IF bAb References and IF bAb Control (ISP3–5)

Three vials containing lyophilized human serum, with preservatives. Reconstitute each vial by adding **1.0 mL sterile distilled** water. Use volumetric pipets, and mix by *gentle* swirling or inversion. *Avoid bacterial contamination and discard if there is any evidence of growth.* Stable at 2–8°C for 30 days after reconstitution, or for 2 months (aliquotted) at –20°C.

KISPZ: 1 set. **KISP1:** 1 set.

The Negative IF bAb Reference consists of pooled human serum negative for intrinsic factor blocking antibody, with a vitamin B12 value in the mid-normal range. The Positive IF bAb Reference is an antibody-positive patient sample diluted in the Negative IF bAb Reference to a low titer; it is used to establish the positive cutoff, hence its titer determines the sensitivity of the IF bAb procedure. The IF bAb Control is an antibody-positive patient sample diluted in the Negative IF bAb Reference and of a higher titer than the Positive Reference.

For the current control ratio ranges, please refer to the Control insert.

Materials Required But Not Provided

Gamma counter: set for cobalt-57.

Centrifuge: capable of at least 1500xg. A high-speed, refrigerated model is desirable but not essential. Good results can also be obtained with low-speed, table-top devices, providing the centrifugation time is suitably extended.

Waterbath: 37°C (optional). Do not use an oven or a heat block. (The assay can also be run at room temperature (15–28°C).)

Magnetic stirring apparatus. The stir bar should fit inside the Solid Phase Binder's original container.

Vortex mixer.

Reagent Preparation

Sterile distilled water. This is essential for good results, since bacterial contamination may lead to an increase in free B12 levels. The use of *commercially bottled* sterile distilled water is strongly recommended.

Volumetric pipets: 1.0 mL and 3.0 mL

Radioimmunoassay

Tubes: 12x75 mm. Polypropylene tubes are available from DPC.

Micropipets: 50 µL and 100 µL. Use disposable-tip pipets with an air interface, rather than positive-displacement pipets, to avoid the risk of carryover.

Micropipet or dispenser: 500 µL.

Dispenser: 1 mL.

Containers and volumetric pipets for preparing the working solution.

Foam decanting rack — available from DPC.

Specimen Collection

The patient need not be fasting. Collect blood by venipuncture²² into plain — preferably sterile — tubes, noting the time of collection. Let the specimen clot, then separate the serum from cells by centrifugation.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. Solid Phase IF bAb has not been tested with all possible variations of tube types.

Volume Required:

Basic procedure: 50 µL serum per tube

Alternate "confirmation" procedure: 100 µL serum per tube.

Storage: 2 days at 2–8°C, or for 2 months frozen at –20°C.

Freezing samples without delay is recommended to prevent bacterial growth. Before assay, allow the samples to come to room temperature (15–28°C), then mix by *gentle* swirling or inversion.

Precautions: Handle all samples in a way that will minimize bacterial contamination. If possible, determine the patient's serum B12 level. In the case of patients undergoing B12 therapy, wait at least 48 hours — preferably one week — after any B12 injection before collecting a sample.^{4,11,21}

Radioassay Procedure

All components must be at room temperature (15–28°C) before use. Assays should not exceed 70 tubes.

- 1 Label six tubes in duplicate **N** (Negative IF bAb Reference), **PR** (Positive IF bAb Reference) and **C** (IF bAb Control). Label additional tubes, also in duplicate, for each patient sample and control.

- 2 **Sample:** Pipet **50 µL** of each reference, control and patient sample into the tubes prepared.

Alternatively, pipet 100 µL of each patient sample and 50 µL of the References and Control into the tubes prepared.

Samples yielding indeterminate results with the basic 50 µL procedure should be reassayed by the alternate 100 µL procedure. To obviate the need for a separate confirmation assay, patient samples may be assayed simultaneously. (**Note:** always assay 50 µL of the references and control.)

- 3 **Binder:** Add **500 µL** of IF bAb Binder to each tube. Vortex vigorously for 5 seconds.

The binder must be kept well suspended with the help of a magnetic stirring device: the stir bar should be stored in the reagent's wide-mouth container. Do not transfer the Binder to a plastic container.

- 4 Incubate for **30 minutes at 37°C** in a waterbath. Then add **1 mL** water to each tube.

Alternatively, incubate for 60 minutes at room temperature (15–28°C). Then add 1 mL water to each tube.

Do not use an oven or a heat block. Use sterile distilled or deionized water.

- 5 Centrifuge for **15 minutes at 3000xg** and decant, retaining the pellet.

Inadequate centrifugation will have a serious effect on the assay. Centrifuge all tubes for 15 minutes at 3000xg or for a suitably longer period of time at a lower g force, e.g. 30 minutes at 1500xg. A high-speed, refrigerated model is desirable but not essential.

Using a foam rack, decant all tubes simultaneously and allow them to drain for 2 or 3 minutes on absorbant paper. Then tap the tubes gently and blot the rims to remove all visible moisture.

- 6 **Tracer/Buffer:** Prepare a working solution by diluting ⁵⁷Co Vitamin B12 **1-in-10** with IF bAb Buffer, in a volume sufficient for adding 500 µL to each tube.

At this point some laboratories may wish to set up a pair of tubes labeled T for total counts.

The working solution should be prepared no more than 30 minutes before use. The procedure calls for 450 µL of buffer and 50 µL of ⁵⁷Co Vitamin B12 per assay tube. Prepare a volume in slight excess to allow for pipetting losses.

- 7 Add **500 µL** of freshly prepared working solution to each tube. Vortex. Set the (optional) total T tubes aside for counting at step 10; they require no further processing.

- 8 Incubate for **30 minutes at room temperature** (15–28°C). Then add **1 mL** water to each tube.

- 9 Centrifuge for **15 minutes at 3000xg** and decant, as before, retaining the pellet.

Inadequate centrifugation will have a serious effect on the assay. Centrifuge all tubes for 15 minutes at 3000xg or for a suitably longer period of time at a lower g force, e.g. 30 minutes at 1500xg. A high-speed, refrigerated model is desirable but not essential.

- 10 Count each tube for **1 minute** in a gamma counter set for ^{57}Co .

Calculation and Quality Control

Divide the average counts of each patient sample, control and reference *into* the average counts of the Negative IF bAb Reference, carrying out the computation to two decimal places. An example worksheet is provided at the back of the package insert.

Ratio = (Negative Ref Average CPM / Average CPM)

The negative cutoff is equal to 1.00, by definition, while the positive cutoff is equal to the observed ratio of the Positive IF bAb Reference.

Record the lot numbers and reconstitution dates of the references and control. The observed ratios of the Positive IF bAb Reference and IF bAb Control should fall within the lot-specific ranges.

Record Keeping: It is good laboratory practice to record for each assay the lot numbers of the components used, as well as the dates when they were first reconstituted or opened. Also record the *lot-specific* ranges for the Positive Reference and Control provided in the box insert. Reject an assay as unacceptable if the ratio observed for either standard falls outside the limits indicated.

Sample Handling: The instructions for handling and storing patient samples and components should be carefully observed. All samples, including the calibrators and controls, should be assayed in duplicate. It is important to use a *disposable-tip* micropipet, changing the tip between samples, to avoid carryover contamination. Inspect the results for agreement within tube pairs. *Assays should be limited to a maximum of 70 tubes*, and it is recommended that pairs of PR tubes be spaced throughout the assay to verify the absence of drift.

Controls: One antibody-positive, human serum-based control is supplied with the kit. Individual patient samples or pooled human serum with negative and intermediate ratios should be employed as additional controls on a routine basis. Their ratios and also those of the Positive Reference and Control should be charted from run to run as described in Westgard

JO, et al. A multi-rule chart for quality control. Clin Chem 1981;27:493-501. Repeat samples are a valuable additional tool for monitoring interassay precision.

Centrifugation: Inadequate centrifugation will have a serious adverse effect on accuracy and reproducibility. The procedure calls for centrifuging at 3000xg for 15 minutes. Lower accelerations have proven satisfactory when the centrifugation time is increased appropriately, for example, to 30 minutes at 1500xg. A high-speed, refrigerated centrifuge is desirable but not essential. Use the formulas below to calculate the acceleration of your centrifuge at a given speed, or the speed (in revolutions per minute) required to achieve a desired g force.

$$g \text{ force} = 28.38 \times (\text{rpm} / 1,000)^2 \times \text{radius [inches]}$$

$$\text{rpm} = 187.7 \times \sqrt{(g \text{ force} / \text{radius [inches]})}$$

Matrix Effect: It is important to know what range of variation to expect from antibody-negative patient samples. The median and absolute range — or 95% confidence limits — for the set of all samples in an assay known to be antibody-negative constitute valuable points of reference in judging the status of unknowns. Assaying a number of antibody-negative patient samples or serum pools in duplicate is more informative than assaying several replicates of the Negative Reference.

Q.C. Parameters: We recommend keeping track of these performance measures:

T = Total Counts (as counts per minute)

$$\%MB = 100 \times (\text{Negative Ref. Average Counts} / \text{Total Counts})$$

Example: For illustration only, not for the assessment of results from another assay. In the table, the indeterminate range is equal to 1.11 – 1.27. (See "Example" table.)

Interpretation of Results

Positive: Samples with ratios equal to or greater than the positive cutoff should be reported as *positive* for intrinsic factor blocking antibody. The positive cutoff is defined as the value actually observed for the Positive Reference minus 0.10.

Negative: Samples with ratios less than the negative cutoff should be reported as *negative*. The negative cutoff is defined as the value observed for the Negative Reference plus 0.10.

Indeterminate: Samples which are neither positive nor negative by the criteria stated above are said to have *indeterminate* ratios. Thus a ratio is indeterminate if it is greater than 1.10 and less than the positive ratio cutoff minus 0.10.

Any 50 µL sample yielding an indeterminate ratio should be reassayed by the alternate 100 µL confirmation procedure. A sample analyzed (simultaneously or consecutively) at both sample volumes yielding an indeterminate ratio with 50 µL but a positive result with 100 µL should be reported as positive for intrinsic factor blocking antibody. Note that the status of many specimens can be read directly from the raw counts without calculation. (See the "Example" table.)

Limitations

A negative or indeterminate result is compatible with the presence of intrinsic factor blocking antibody in an undetectably low titer. The IF bAb procedure can detect the presence of antibody only when the titer is at least half that of the Positive IF bAb Reference. Furthermore, only antibody which has retained its ability to block the binding of intrinsic factor to B12 will be measured in the procedure.

False positive results will occur if the concentration of free B12 in the patient sample is sufficiently high.²¹ This is ordinarily not a cause for concern, as B12 circulates primarily in protein-bound form.^{1,5,14} Even so, care should be exercised in handling specimens and kit components to avoid bacterial contamination which might lead to levels of free B12 high enough to interfere with the assay. Likewise, no specimen should be collected from a patient currently undergoing B12 therapy less than one week after the last B12 injection.^{4,11,21}

The combination of a high B12 result and a strongly positive result by the IF bAb procedure in the same patient should cast doubt on the validity of both measurements. The IF bAb result may

represent a false positive due to high free B12 levels in the sample.²¹ On the other hand, high levels of intrinsic factor blocking antibody may interfere with certain B12 radioassays, causing spuriously high results.² The B12 result should therefore be verified by reassaying several different dilutions of the sample.

A positive result should not be used on its own to support a diagnosis of pernicious anemia, since intrinsic factor blocking antibody has been known to occur in other conditions, though infrequently.^{8,15,20} For example, it is sometimes encountered in certain disorders of the nervous system, such as the Eaton-Lambert syndrome.¹¹ IF bAb results should therefore be interpreted only against the background of other clinical and hematological findings.¹⁶

Performance Data

See Tables and Graphs for data representative of the Solid Phase IF bAb procedure. Results were generated using the basic (50 µL) procedure, and are expressed as ratios of Negative IF bAb Reference CPM to sample CPM, unless otherwise specified.

Intraassay Precision (Within-Run):

Statistics were calculated from the results of twenty pairs of tubes in a single run for the Negative and Positive References and two patient samples, one weakly positive and the other negative for intrinsic factor blocking antibody. The sample designated PR/2 is the Positive Reference diluted 1-in-2 with the Negative Reference. Note that the second patient sample, which was known to be antibody-positive from assay results using 100 µL, has a 50 µL ratio intermediate between those of the Negative and Positive References. Its range, however, shows no overlap with the ranges of either the Negative Reference or the antibody-negative patient sample. (See the "Intraassay Precision" table.)

Interassay Precision (Run-to-Run):

Statistics were calculated from the results of pairs of tubes in at least 16 different assays for the three references and two patient samples, one positive and the other negative for intrinsic factor blocking antibody. The statistics show that even on an assay-to-assay basis, the IF bAb procedure yields good discrimination

between antibody-negative samples and samples with an intrinsic factor blocking antibody titer comparable to that of the Positive Reference. Even better discrimination can be expected within the assay: although the absolute range of values observed for the Positive Reference overlaps somewhat with the range observed for the antibody-positive patient sample, the latter had a higher ratio in every run. (See "Interassay Precision" table.)

End-of-Run Effect: None up to approximately 70 tubes. (See "End-of-Run Effect" table.)

Protein and Matrix Effects: To simulate various protein concentrations, aliquots of two antibody-negative patient sample pools were freeze-dried, reconstituted with various volumes of water and then assayed by the basic 50 μ L IF bAb procedure. The results indicate that even wide variations in protein concentration, ranging from approximately 3.5 to 11 g/dL, have no clinically significant effect on the ratios observed. (See "Protein and Matrix Effects 1" table.)

Since doubling the sample volume effectively doubles the protein concentration in each assay tube, the IF bAb procedure was applied to several patient samples using the alternate 100 μ L sample size in order to confirm the absence of a significant matrix effect. Twenty-four antibody-negative and two antibody-positive samples were assayed in duplicate, using 50 μ L and 100 μ L volumes of each sample. All samples were processed in a single assay. The positive cutoff was 1.27. Mean values and absolute ranges are displayed for the group of antibody-negative samples. The results show that antibody-negative patient samples assayed according to the basic procedure (that is, using 50 μ L) yield essentially the same range of values as the Negative Reference: they have an average ratio of 1.00 and an observed variation comparable to the imprecision of the Negative Reference assayed as an unknown. (See "Protein and Matrix Effects 2" table.)

There is a small but limited volume effect when samples are assayed using twice the volume specified in the basic procedure. Ratios of antibody-negative

samples show a typical increase of about 0.04 when the sample size is increased from 50 μ L to 100 μ L. Virtually all of the antibody-negative samples tested here (23 out of 24) showed increases of at most 0.08, while the majority (20 out of 24) showed increases of 0.05 or less. None of them, when assayed using 100 μ L, yielded a ratio greater than that of the Positive Reference assay using 50 μ L; in fact, their ratios fell entirely within the negative region (1.00 ± 0.10).

Observe that the increments associated with doubling the volume are significantly greater for the two positive samples than for any of the antibody-negative samples, i.e. 0.21 and 0.28 as opposed to 0.08 or less. These represent increases beyond what can be attributed to volume effect, as should be expected since doubling the sample volume effectively doubles the level of intrinsic factor blocking antibody present in the sample.

B12 Effect: Twenty-three antibody-negative patient samples with a wide range of vitamin B12 values were assayed by the basic 50 μ L IF bAb procedure. The results, tabulated in order of increasing B12 concentration, show no apparent trend towards higher ratios as a function of B12 concentration. This was confirmed by linear regression analysis, which yielded a slope of zero to three decimal places. (See "B12 Effect 1" table.)

The assay is affected by variations in *free* B12. To demonstrate the extent of this effect, B12 was spiked at various concentrations into charcoal-absorbed human serum albumin and then assayed by the basic 50 μ L IF bAb procedure. Because B12 carrier proteins (the transcobalamins) are absent, essentially all of the B12 in the spiked samples is present as free B12. Results obtained by two different technicians are displayed. The results show that samples with *free* B12 levels above 300 pg/mL will yield antibody-positive results in the basic 50 μ L IF bAb procedure. On the assumption that protein-bound B12 constitutes roughly 99% of the B12 in circulation, this corresponds to a free-plus-protein-bound B12 level of 30,000 pg/mL. (On the same assumption, the B12 values displayed in the previous table correspond to free B12 values ranging from 3.5 to 19 pg/mL.) (See "B12 Effect 2" table.)

Linearity: Three antibody-positive patient samples were serially diluted with the Negative IF bAb Reference and assayed by the basic 50 µL IF bAb procedure. The positive cutoff was 1.24 in the assay that included samples A and B, and 1.37 in the assay which included sample C. Results for the first two samples show that the ratios of highly positive samples decrease steadily with serial dilution, passing from positive to negative values as expected. (See "Linearity 1" table.)

In another experiment, three antibody-positive 50 µL samples were serially diluted and assayed. The higher dilutions were also assayed using 100 µL. All results were obtained in a single assay which had a positive cutoff of 1.29. As expected, the ratio determined for a given dilution assay using 50 µL is roughly equal to the ratio obtained for the next higher dilution assay using 100 µL. Collectively, the results presented demonstrate that even a strongly positive sample will cease to yield a positive result in the IF bAb procedure when diluted far enough. Serial dilution eventually transforms a positive sample to one with an antibody level that may be *undetectably* low. Hence, a negative result by the IF bAb procedure could mean *either* that intrinsic factor blocking antibody is not present *or* that its level is below the sensitivity of the kit. (See "Linearity 2" tables.)

Hemolysis: Presence of packed red blood cells in concentrations up to 30 µL/mL has no effect on results, within the precision of the assay.

References

1) Allen RH. Human vitamin B12 transport proteins. *Prog Hematol* 1975;9:57-84. 2) Allen RH. More on no-boil assay. *Lig Quarterly* 1982(Spring);5(1): 48-9. 3) Beck WS. Vitamin B12 deficiency. In: Williams WJ et al, editors. *Hematology*, 2nd edition. New York: McGraw-Hill, 1977: 307-34. 4) Carmel R. The laboratory diagnosis of megaloblastic anemias. *Wes J Med* 1978;128:294-304. 5) Carmel R. Cobalamin-binding proteins of man. *Contemp Hematol/Oncol* 1981;2:79-129. 6) Carmel R. Megaloblastic anemia: vitamin B12 and folate. *Curr Hematol* 1983;2:243-80. 7) Castle WB. The conquest of pernicious anemia. In: Wintrobe MM, editor. *Blood: pure and eloquent*. New York: McGraw-Hill, 1980: 283-317. 8) Chanarin I. *The Megaloblastic anaemias*, 2nd ed. Oxford: Blackwell Scientific, 1979. 9) Conrad ME, et al. Iron, folic acid and vitamin B12. *Curr Hematol*

1981;1:123-90. 10) Diagnostic Products Corporation. *Developments in vitamin B12 testing*, 1983. Available on request. 11) Fairbanks VF, et al. Tests for pernicious anemia: serum intrinsic factor blocking antibody. *Proc Mayo Clin* 1983;58:203-4. 12) Gottlieb C, Lau K.-S, Wasserman LR, Herbert V. Rapid charcoal assay for intrinsic factor (IF), gastric juice unsaturated B12 binding capacity, antibody to IF, and serum unsaturated B12 binding capacity. *Blood* 1965;25:875-85. 13) Herbert V. Megaloblastic anemias. In: Beeson PB, McDermott W, Wyngaarden JB, editors. *Cecil's textbook of medicine*, 15th ed. Philadelphia: W.B. Saunders, 1979: 1719-29. 14) Jacob E, Baker SJ, Herbert V. Vitamin B12-binding proteins. *Physiol Rev* 1980;60:918-59. 15) Jeffries GH. Pernicious anemia and atrophic gastritis. In: Samter M, editor. *Immunological diseases*, vol 2, 3rd ed. Boston: Little, Brown, 1978: 1296-1307. 16) Lindenbaum J. Status of laboratory testing in the diagnosis of megaloblastic anemia. *Blood* 1983;61:624-7. 17) Nicolas JP. Heterogeneity of antibodies against B12 binders in pernicious anaemia. In: Zagalak B, Friedrich W, editors. *Vitamin B12*. Berlin: Walter de Gruyter, 1979: 791-802. 18) Rothenberg SP, et al. Autoantibodies to intrinsic factor: their determination and clinical usefulness. *J Lab Clin Med* 1971;77:476-84. 19) Solanki DL, et al. Pernicious anemia in blacks. *Am J Clin Pathol* 1981;75:96-9. 20) Taylor KB. Immune aspects of pernicious anaemia and atrophic gastritis. *Clin Hematol* 1976;5(3):497-519. 21) Ungar B. Detection and titration of antibody to intrinsic factor by a modified charcoal method. *Aus J Exper Biol Med Sci* 1967;45:317-21. 22) National Committee for Clinical Laboratory Standards. *Procedures for the collection of diagnostic blood specimens by venipuncture*; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782, 310.645.8200
Fax: 310.645.9999. To place an order:
Tel: 800.372.1782. Fax: 800.234.4372.
Outside the United States, contact your National Distributor.

The Quality System of Diagnostic Products Corporation is registered to ISO 13485:2003.

Tables and Graphs

Example Run

Tube ¹	Duplicate CPM ²	Average CPM ³	IF bAb Ratio ⁴	Status ⁵
Total Counts ⁶	22,819 22,533			
Neg IF bAb Ref. ⁷	14,800 14,907	14,854	1.00	—
Pos IF bAb Ref. ⁸	10,725 10,848	10,787	1.38	—
X1	14,624 14,549	14,587	1.02	Negative
X2	13,149 13,204	13,177	1.13	Indeterminate
X3	11,039 10,968	11,004	1.35	Positive
X4	5,167 5,331	5,249	2.83	Positive

Quality Control Parameters:⁹

T⁶ = 22,676 cpm

Positive Cutoff¹⁰ = 1.38

%MB¹¹ = 66%

Intraassay Precision (ratio)

	Mean ¹	SD ²	CV ³	Absolute Range ⁴
Negative Ref.	1.01	0.01	1.0%	0.98 – 1.02
Patient Sample 1	1.04	0.01	1.0%	1.02 – 1.05
PR/2	1.16	0.01	0.9%	1.14 – 1.18
Patient Sample 2	1.21	0.01	0.8%	1.19 – 1.23
Positive Ref.	1.31	0.02	1.5%	1.27 – 1.33

Interassay Precision (ratio)

	n	Mean ¹	SD ²	CV ³	Absolute Range ⁴
Neg. Ref.	20	1.00	0.01	1.3%	0.97 – 1.10
Patient Sample 1	20	1.04	0.02	1.6%	1.01 – 1.07
Pos. Ref.	19	1.29	0.03	1.9%	1.24 – 1.34
Patient Sample 2	20	1.34	0.03	2.3%	1.28 – 1.39
Control	16	2.38	0.08	3.3%	2.21 – 2.50

End-of-Run Effect (ratio)

	Tubes ¹ 9–26	Tubes 23–28	Tubes 49–54	Tubes 75–80
1	1.30	1.29	1.32	1.26
2	1.00	1.00	1.07	1.07
3	2.41	2.39	2.38	2.65

Protein and Matrix Effects 1

	Total Protein ¹	Observed Ratio ²
1	14.0 10.8 7.0 4.7 3.5	1.03 1.04 1.06 1.04 1.07
2	14.0 10.8 7.0 4.7 3.5	1.14 1.08 1.10 1.05 1.03

Protein and Matrix Effects 2

	50 µL Assay ¹	100 µL Assay ²
24 Negative Samples ³	1.00 [0.98–1.02]	1.04 [1.00–1.09]
Two Positive Samples ⁴		
1	1.19	1.40
2	1.22	1.50

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O Sistema de Qualidade da Diagnostic Products Corporation está registado sob ISO 13485:2003.
DPC Biemann GmbH