Glucagon

DOUBLE ANTIBODY

DPC[®]

Double Antibody Glucagon

English

Intended Use: Double Antibody Glucagon is an ^{125}I radioimmunoassay designed for the quantitative measurement of glucagon in EDTA plasma. It is intended for *in vitro* diagnostic use to aid in the diagnosis and treatment of patients with disorders of carbohydrate metabolism, including diabetes mellitus, hypoglycemia and hyperglycemia.

Catalog number: **KGND1** (100 tubes)

CDC Analyte Identifier Code: **2206**

CDC Test System Identifier Code: **13031**

CLIA Complexity Category: **High**

The 100-tube kit contains less than 2 microcuries (74 kilobecquerels) of radioactive ^{125}I glucagon.

Summary and Explanation

Glucagon is a single-chain polypeptide, 29 amino acids in length, produced in the pancreatic alpha cells. In its action, it is generally antagonistic to insulin, which is secreted by the pancreatic beta cells. Like insulin, glucagon is involved in the regulation of circulating glucose levels; but whereas insulin typically acts to oppose hyperglycemia, glucagon is part of a counterregulatory system – which includes epinephrine and several other factors, usually of lesser significance – designed to defend against hypoglycemia.^{6,25,26} In response to increased glucose utilization (as during periods of exercise) or depletion induced by insulin, glucagon helps to restore circulating levels by stimulating glucose production in the liver.⁶

Glucagon determinations may thus be of relevance to patients with insulin-dependent diabetes who are undergoing, or are being considered for, intensive insulin therapy. In such a context, defective counterregulation – characterized by a *combined* deficiency in both glucagon and epinephrine responses to decrements in circulating glucose – would entail increased risk for severe, insulin-induced hypoglycemic episodes.^{1,4,5,6,24}

Glucagon is structurally related to other neurohormonal peptides, including secretin, vasoactive intestinal polypeptide (VIP) and gastric inhibitory peptide (GIP).²² The hormone shows considerable species-to-species similarity, especially among Old World mammals.^{7,22} The amino acid sequence of human glucagon is identical to that of the glucagon found in cows, pigs, rabbits, rats and hamsters; guinea pig glucagon, on the other hand, differs at five positions. Turkey and chicken glucagons differ from human glucagon at just one position; duck and alligator glucagons at two.

Glucagon derives from preproglucagon, which is present in the alpha cells of the pancreas and also in the gut. This precursor includes an N-terminal signal peptide (residues –24 to –1), followed by glicentin (residues 1 to 69) and what has been called the major proglucagon fragment (residues 72 to 160). The latter includes two glucagon-like peptides, designated GLP-1 and GLP-2. Glicentin, which is also known as gut glucagon and glucagon-like immunoreactivity 1 (GLI-1), includes glicentin-related pancreatic peptide (GRPP, glicentin-1) followed by glucagon itself (residues 33 to 61) and a hexapeptide. The C-terminal fragment of glicentin, which contains glucagon and the hexapeptide, has been variously referred to as oxyntomodulin, enteroglucagon, GLI-2 and glucagon-37.²²

Whereas the gut secretes mainly the intact N-terminal half of preproglucagon (i.e., glicentin) along with oxyntomodulin and the C-terminal fragments GLP-1 and GLP-2, the pancreas secretes mainly the intact C-terminal half (i.e., the so-called major proglucagon fragment) along with the N-terminal fragments GRPP, glucagon and oxyntomodulin.^{2,22} Since glicentin, which has sometimes been called proglucagon, is not a typical secretory product of the pancreas, it may not be a normal intermediary in the biosynthesis of glucagon.² Even so, the fact that circulating glucagon levels tend to remain normal in human subjects after total pancreatectomy hints at the existence of some mechanism for transforming intestinal products to glucagon.^{2,26}

Principle of the Procedure

DPC's Double Antibody Glucagon procedure is a sequential radioimmunoassay. After preincubation of the patient sample with anti-glucagon antibody, ¹²⁵I-labeled glucagon competes with glucagon in the sample for antibody sites. After incubation for a fixed time, separation of bound from free is achieved by the PEG-accelerated double-antibody method, followed by centrifugation. The precipitate containing the antibody-bound fraction is then counted, and patient concentrations are read from a calibration curve.

Reagents to Pipet: 3

Total Incubation Time: Two 24-hour incubations.

Total Counts at Iodination: approximately 25,000 cpm.

Separation: The ready-to-use Precipitating Solution combines second antibody and dilute PEG.

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C in a refrigerator designated for incoming radioactive materials. Dispose of in accordance with applicable laws.

Do not use reagents beyond their expiration dates.

Some components supplied in this kit may contain human source material and/or other potentially hazardous ingredients which necessitate certain precautions.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Water: Use distilled or deionized water.

Radioactivity: A copy of any radioisotope license certificate (Specific or General) issued to a US customer must be on file with Diagnostic Products Corporation before kits or components containing radioactive material can be shipped. These radioactive materials may be acquired by any customer with the appropriate Specific license. Under a General license these radioactive materials may be acquired only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories and hospitals — and strictly for *in vitro* clinical or laboratory tests not involving external or internal administration of the radioactive material or its radiation to human beings or other animals. Its acquisition, receipt, storage, use, transfer and disposal are all subject to the regulations and a (General or Specific) license of the U.S. Nuclear Regulatory Commission or a State with which the NRC has entered into an agreement for the exercise of regulatory control.

Handle radioactive materials according to the requirements of your General or Specific license. To minimize exposure to radiation, the user should adhere to guidelines set forth in the National Bureau of Standards publication on the *Safe Handling of Radioactive Materials* (Handbook No. 92, issued March 9, 1964) and in subsequent publications issued by State and Federal authorities.

Wipe up spills promptly and decontaminate affected surfaces. Avoid generation of aerosols. Dispose of solid radioactive waste according to license requirements. General licensees (holders of NRC Form 483) may dispose of solid radioactive waste as nonradioactive waste, after removing labeling. Specific licensees (NRC Form 313) should refer to Title 10, Code of Federal Regulations, Part 20. Licensees in Agreement States should refer to the appropriate regulations of their own state. General licensees may dispose of liquid radioactive waste of the type contained in this product through a laboratory sink drain. Licensees must remove or deface labels from empty containers of radioactive materials before disposal of solid waste. Specific licensees may dispose of small quantities of liquid radioactive waste of the type used in this product through a laboratory sink drain.

Refer to the appropriate regulations applicable to your laboratory.

Materials Supplied

Glucagon Antiserum (GND1)

Lyophilized glucagon antiserum, raised in rabbits, with preservative. Reconstitute by adding **10 mL** distilled water. Mix by *gentle* inversion. Stable at 2–8°C for 30 days after reconstitution. Color: blue.

KGND1: 1 vial

¹²⁵I Glucagon (GND2)

Lyophilized, iodinated synthetic human glucagon, with preservative. Just before assay, reconstitute by adding **10 mL** distilled water. Mix by *gentle* inversion. Stable at 2–8°C for 30 days after reconstitution.

KGND1: 1 vial

Glucagon Zero Calibrator (GND3)

Lyophilized zero calibrator in a protein-based matrix, with preservative. *At least 30 minutes before use*, reconstitute each vial with **10 mL** distilled water. Mix by *gentle* inversion. Stable at 2–8°C for 30 days after reconstitution, or for 2 months (aliquotted) at –20°C.

KGND1: 2 vials

Glucagon Calibrator F (GND8)

Lyophilized calibrator containing synthetic human glucagon in a protein-based matrix, with preservative. *At least 30 minutes before use*, reconstitute each vial with **2.0 mL** distilled or deionized water. Use a volumetric pipet and mix by *gentle* inversion. *Discard after use*.

KGND1: 3 vials

The calibrator F has a *lot-specific* value of approximately 500 picograms of synthetic human glucagon per milliliter (pg/mL); equivalently 144 picomoles per liter (pmol/L). Refer to the vial labels for exact values in pg/mL. The assay is standardized in terms of the World Health Organization's International Standard for Glucagon, number 69/194.

Precipitating Solution (GNRG)

110 mL of Precipitating Solution consisting of goat anti-rabbit gamma globulin (GARGG) and dilute polyethylene glycol (PEG) in saline, with gentamicin as a preservative. The Precipitating Solution is

supplied in liquid form, ready to use. Remove the aluminum foil seal completely and store refrigerated: stable at 2–8°C for 30 days after opening. Since a fine (but visible) precipitate may form after refrigeration, the Precipitating Solution should be thoroughly mixed before use, without foaming. Color: red.

KGND1: 1 vial

Glucagon Controls (GNCO1, GNCO2)

Two vials, labeled Glucagon Control 1 and 2, containing synthetic human glucagon in a protein-based matrix, with preservative. The controls are supplied lyophilized. *At least 30 minutes before use*, reconstitute each vial with **1.0 mL** distilled water. Use a volumetric pipet and mix by *gentle* inversion. *Discard after use*. Refer to the control insert for glucagon values in pg/mL.

KGND1: 3 sets.

Required But Not Provided

Gamma counter

Centrifuge – refrigerated and capable of at least 1500xg

Vortex mixer

Reagent Preparation

Distilled or deionized water

Pipets for measuring 5.0 mL and 10 mL

Volumetric pipet: 1.0 mL

Radioimmunoassay

Plain 12x75 mm *glass* tubes – available from DPC. Tubes should be made of high-quality borosilicate glass. Because of glucagon's tendency to adsorb to plastic surfaces, glass tubes (rather than plastic) should be employed for the assay.

Micropipets – 50 μ L, 100 μ L, 200 μ L, 500 μ L, and 1.0 mL. For the 100 μ L additions, a reliable repeating dispenser is recommended. A syringe-style dispenser accurate to within ± 0.05 mL is recommended for the 1.0 mL addition of Precipitating Solution.

Parafilm

Foam decanting rack – available from DPC

Logit-log graph paper – available from DPC (catalog number: ZP797)

Specimen Collection

The patient need not be fasting, and no special preparations are necessary. Collect blood by venipuncture²⁷ into *iced* EDTA tubes, noting the time of collection. The tubes should be immersed in an ice bath at 2–8°C throughout the entire collection and handling process.

Immediately following collection, add 5,000 Kallikrein Inactivating Units (KIU) per 10 mL of whole blood, in the form of Trasylol (500 µL) or Aprotinin (1 mg). Mix thoroughly and keep chilled until separation.

Separate the plasma from the cells by centrifugation in a *refrigerated* centrifuge. Aspirate the plasma with a glass Pasteur pipet, then transfer it into a *glass* tube and freeze immediately. Aliquot, if necessary.

Lipemia may interfere with the assay. The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. Double Antibody Glucagon has not been tested with all possible variations of tube types.

Volume Required: 200 µL of EDTA plasma per tube.

Storage: *Caution:* Do not use plastic tubes or pipets in any part of the preparation. Keep samples chilled at all times and away from direct light. Thawing and refreezing, or prolonged exposure to room temperature (15–28°C), may result in significantly lower glucagon values.

Radioimmunoassay Procedure

All components except the Precipitating Solution must be at room temperature (15–28°C) before use.

Caution: Because of the tendency of glucagon to adsorb to plastic surfaces, *glass* tubes rather than plastic should be employed for the assay. (The use of plastic tubes could result in the loss of all detectable glucagon.)

- 1 Just before assay, dilute the 500 pg/mL Glucagon Calibrator F in the Glucagon Zero Calibrator, using *glass* tubes or *glass* vials, as described in the table below. Vortex each dilution thoroughly.

Add this...		To this...		Yielding this Calibrator		
µL	Calib.	µL	Calib.	pg/mL	pmol/L	
50	Calib. F	950	Zero	B	25	7
100	Calib. F	900	Zero	C	50	14
200	Calib. F	800	Zero	D	100	29
500	Calib. F	500	Zero	E	250	72

- 2 Label sixteen glass tubes in duplicate: T (total counts), NSB (nonspecific binding), A (maximum binding) and B through F. Label additional tubes, also in duplicate, for plasma samples and controls.

Calibrator	Approximate pg/mL	Approximate pmol/L
A (MB)	0	0
B	25	7
C	50	14
D	100	29
E	250	72
F	500	144

Note: The value of the calibrators are *lot-specific*. Refer to the vial labels for exact values in pg/mL.

- 3 Pipet **200 µL** of the zero calibrator A into the NSB and A tubes, and **200 µL** of each of the remaining calibrators B through F into correspondingly labeled tubes. Pipet **200 µL** of each patient plasma sample and each control into the tubes prepared.

It is good practice to use a *disposable-tip* micropipet, changing the tip between samples, in order to avoid carryover contamination.

Samples expected to contain glucagon concentrations greater than the highest calibrator (calibrator F, 500 pg/mL) should be diluted in the zero calibrator before assay.

- 4 Add **100 µL** of Glucagon Antiserum (BLUE) to all tubes except the NSB and T tubes. Vortex.

A repeating dispenser is recommended for this step and for the addition of tracer at step 6.

5 Cover the rack and incubate for 24 hours at 2–8°C.

Cover the entire rack with Parafilm or place the racks in plastic bags to avoid evaporation.

6 Add 100 µL of ¹²⁵I Glucagon (CLEAR) to all tubes. Vortex.

Remove the T tubes for counting at step 11; they require no further processing.

7 Cover the rack and incubate for 24 hours at 2–8°C.

Cover the entire rack with Parafilm or place the racks in plastic bags to avoid evaporation.

8 Add 1.0 mL of cold Precipitating Solution (RED) to all tubes. Vortex.

For the 1.0 mL reagent addition, a repeating dispenser may be employed.

9 Centrifuge for 15 minutes at 1500xg.

10 Using a foam decanting rack, decant (or aspirate) the supernatant, retaining the precipitate for counting.

Let the tubes stand inverted on absorbant paper for 10 minutes. Then tap the tubes gently and blot the rims, to remove residual droplets.

11 Count each tube for 1 minute.

Calculation of Results

To obtain results in terms of concentration from a logit-log representation of the calibration curve, first calculate for each pair of tubes the average NSB-corrected counts per minute:

Net Counts = (Average CPM) minus (Average NSB CPM)

Then determine the binding of each pair of tubes as a percent of maximum binding (MB), with the NSB-corrected counts of the A tubes taken as 100%:

Percent Bound = (Net Counts / Net MB Counts) × 100

(The calculation can be simplified by omitting the correction for nonspecific binding; samples within range of the calibrators yield virtually the same results

when Percent Bound is calculated directly from Average CPM.)

Using logit-log graph paper, plot Percent Bound on the vertical (probability) axis against Concentration on the horizontal (logarithmic) axis for each of the nonzero calibrators, and draw a straight line approximating the path of these points. Results for the unknowns may then be read from the line by interpolation.

It is good practice to inspect results for agreement within replicates, and to construct a graph of the calibration curve (even if the calculations are handled by computer) as a visual check on the appropriateness of the transformation used and as a way to detect deviant calibration points. We also recommend keeping track of these data reduction parameters:

T = Total Counts (as counts per minute)

%NSB = 100 × (Average NSB Counts / Total Counts)

%MB = 100 × (Net Counts / Total Counts)

And the 20, 50 and 80 percent "intercepts," where

20% = Concentration at 20 Percent Bound, etc.

Note that other approaches, e.g. a sound implementation of the 4-parameter logistic, are also acceptable. See Dudley RA, et al. Guidelines for immunoassay data reduction. Clin Chem 1985;31:1264-71.

Example Run: For illustration only, not for calculating results from another run. Because the calibrator values are *lot-specific*, concentrations listed in the right-most column may not match the values of the calibrators supplied in your shipment. (See "Example Run" table.)

Quality Control

The Glucagon control(s) supplied with the kit should be used as quality control material to monitor the performance of the assay.

Report patient results only if the control results for that assay meet your laboratory's established criteria for acceptability.

It is good laboratory practice to record for each assay the lot numbers of the components used, as well as the dates when they were first reconstituted or

opened. We also recommend charting control results from day to day — as described, for example, in Westgard JO, et al. A multi-rule chart for quality control. Clin Chem 1981;27:493-501. Note that repeat samples can serve as a valuable additional tool for monitoring interassay precision, and that pairs of control tubes can be spaced throughout the assay to help verify the absence of significant drift.

Centrifugation: The procedure calls for centrifuging at 1500xg for 15 minutes. A high-speed, refrigerated centrifuge is desirable but not essential. Use the formulas below to calculate the acceleration of your centrifuge at a given speed, or the speed (in revolutions per minute) required to achieve a desired g force.

$$g \text{ force} = 28.38 \times (\text{rpm} / 1,000)^2 \times \text{radius [inches]}$$

$$\text{rpm} = 187.7 \times \sqrt{(g \text{ force} / \text{radius [inches]})}$$

or

$$g \text{ force} = 11.18 \times (\text{rpm} / 1,000)^2 \times \text{radius [cm]}$$

$$\text{rpm} = 299.1 \times \sqrt{(g \text{ force} / \text{radius [cm]})}$$

Expected Values

A reference range study was performed in which plasma samples from 51 laboratory volunteers were assayed by the Double Antibody Glucagon procedure. This study yielded a median of 112 pg/mL (32 pmol/L) and a range — corresponding to the central 95% of the observations — of 59 – 177 pg/mL (17 – 51 pmol/L)

Laboratories should consider these results as *guidelines only*. Each laboratory should establish its own reference ranges.

Limitations

When monitoring a patient for glucagon levels, a baseline must be established for that patient using the same assay system, as differences have been observed in glucagon values measured by different assays.

Performance Data

See Tables and Graphs for data representative of the Double Antibody Glucagon kit's performance. Results are expressed as pg/mL.

Conversion Factor:

$$\text{pg/mL} \times 0.2872 \rightarrow \text{pmol/L}$$

Calibration Range: 25 – 500 pg/mL (WHO IS 69/194)

The calibrator values are *lot-specific*.

Analytical Sensitivity: 13 pg/mL (3.7 pmol/L).

Intraassay Precision (Within-Run):

Statistics were calculated for samples from the results of 20 pairs of tubes in a single assay. (See "Intraassay Precision" table.)

Interassay Precision (Run-to-Run):

Statistics were calculated for samples from the results of pairs of tubes in 20 different assays. (See "Intraassay Precision" table.)

End-of-Run Effect: None up to approximately 180 tubes. (See "End-of-Run" Effect table.)

Specificity: Percent crossreactivity was calculated for intact glucagon, glucagon fragments, and related compounds on a weight-per-weight basis; the table shows that Double Antibody Glucagon is highly specific for glucagon. (See "Specificity" table.)

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of packed red blood cells in concentrations up to 30 $\mu\text{L/mL}$ has no effect on results, within the precision of the assay.

Linearity: Samples were assayed under various dilutions. (See "Linearity" table.)

Effect of Storage in Plastic and Glass Tubes:

To determine the effect of storing patient samples in plastic and in glass tubes, blood was collected from 10 volunteers into EDTA tubes and spiked with aprotinin to a concentration of one gram aprotinin per liter (1.0 g/L). The plasma was aliquoted into glass and plastic tubes, which were stored frozen for 5 days. All samples were then assayed by the Double Antibody Glucagon procedure, with the following results.

$$(\text{Plastic}) = 0.87 (\text{Glass}) - 6.0 \text{ pg/mL}$$

$$r = 0.98$$

Means:

57 pg/mL (Plastic)

72 pg/mL (Glass)

The values obtained for aliquots stored in plastic tubes were consistently lower than

matched values for aliquots stored in glass tubes. Based on this study, it is recommended that patient samples be stored in glass tubes.

Effect of Aprotinin: To determine the effect of aprotinin on patient sample stability, two sets of twelve EDTA plasma samples were collected from 12 volunteers. One set was spiked with aprotinin (1.0 g/L). All samples were then assayed by the Double Antibody Glucagon procedure. By linear regression:

(Without Aprotinin) =
1.07 (With Aprotinin) + 52 pg/mL
 $r = 0.77$

Means:
111 pg/mL (Without Aprotinin)
55 pg/mL (With Aprotinin)

The analysis shows that results are higher in the absence of aprotinin, presumably due to crossreacting glucagon fragments arising from proteolytic degradation. This confirms the need to stabilize samples by adding aprotinin, as indicated in the Specimen Collection section.

Method Comparison: DPC's Double Antibody Glucagon kit was compared to another commercially available glucagon assay (Kit A) on 20 samples. DPC results were generally lower than results for Kit A. By linear regression:

(DPC) = 0.015 (Kit A) + 142 pg/mL
 $r = 0.051$

Means:
145 pg/mL (DPC)
251 pg/mL (Kit A)

The two kits likewise differ in their expected values; the upper limit of the reference range for the DPC kit is 130 pg/mL, but 250 pg/mL for Kit A.

Differences in reference ranges are actually a common feature of immunoassays employing different antibodies and heterogeneous analytes. This may help to explain the lower values generated by the DPC kit; they suggest that Kit A recognizes glucagon precursors and/or fragments to a greater extent. This hypothesis is corroborated by the DPC assay's demonstrated lack of crossreactivity to preproglucagon 72–108, preproglucagon 72–107, preproglucagon 126–159, gastric inhibitory peptide, oxyntomodulin, secretin and vasoactive intestinal peptide (VIP). (See "Method Comparison" table.)

Standardization: To verify the standardization of the Double Antibody Glucagon procedure, the World Health Organization's International Standard for Glucagon, number 69/194, was diluted to concentrations corresponding to the nonzero calibrator levels. The DPC zero calibrator was used as diluent. These dilutions of the WHO preparation, along with a set of the DPC calibrators, were processed by DPC's Double Antibody Glucagon procedure.

The percent bound (%B/B₀) values displayed in the table confirm parallelism with the WHO preparation. The DPC calibrators behave like dilutions of the reference material. (See "Standardization" table.)

References

- 1) Amiel S, et al. Defective glucose counterregulation after strict glycemic control of insulin-dependent diabetes mellitus. *New Eng J Med* 1987;316:1376-83.
- 2) Bataille D, et al. Glucagon and related peptides; molecular structure and biological specificity. *Ann N Y Acad Sci* 1988;527:168-85.
- 3) Bodanszky M, Bodanszky A. Conformation of peptides of secretin-VIP-glucagon family in solution. *Peptides* 1986;7 Sup 1:43-8.
- 4) Bolli GB, et al. A reliable and reproducible test for adequate glucose counterregulation in type I diabetes mellitus. *Diabetes* 1984;33:732-7.
- 5) Cryer PE, Gerich JE. Glucose counterregulation, hypoglycemia, and intensive insulin therapy in diabetes mellitus. *New Eng J Med* 1985;313:232-41.
- 6) Cryer PE, White NH, Santiago JV. The relevance of glucose counterregulatory systems to patients with insulin-dependent diabetes mellitus. *Endocr Rev* 1986;7:131-9.
- 7) Epple A, Brinn JE. The Comparative Physiology of the Pancreatic Islets. Berlin: Springer-Verlag, 1987.
- 8) Gulan M, et al. Counterregulatory hormone responses preserved after long-term intravenous insulin infusion compared to continuous subcutaneous insulin infusion. *Diabetes* 1988;37:526-31.
- 9) Hagopian WA, Tager HS. Hepatic glucagon metabolism. *J Clin Invest* 1987;79:409-17.
- 10) Heller SR, Macdonald IA, Tattersall RB. Counterregulation in type II (non-insulin-dependent) diabetes mellitus; normal endocrine and glycaemic responses, up to ten years after diagnosis. *Diabetologia* 1987;30:924-9.
- 11) Hermansen K, et al. Glucagon responses to intravenous arginine and oral glucose in insulin-dependent diabetic patients during six months conventional or continuous subcutaneous insulin infusion. *Metabolism* 1988;37:640-4.
- 12) Iyengar R, Herberg JT, Rich KA. The glucagon receptor. *Pharmacol Ther* 1988;37:151-65.
- 13) Kumai M, et al. Glucagon secretion in anorexia nervosa.

Am J Clin Nutr 1984;47:239-42. 14) Lefévre PJ. Accumulated experience in the physiology of glucagon. In: Picazo J, editor. Glucagon in 1987. Lancaster: MTP Press, 1987: 3-12. 15) Matschinsky FM, Bedoya FJ. Metabolism of pancreatic islets and regulation of insulin and glucagon secretion. In: DeGroot LJ, editor. Endocrinology. Philadelphia: Saunders, 1989;2:1290-1303. 16) Nishino T, et al. Glucagon radioimmunoassay with use of antiserum to glucagon C-terminal fragment. Clin Chem 1981;27:1690-7. 17) Reaven GM, et al. Documentation of hyperglucagonemia throughout the day in nonobese and obese patients with noninsulin-dependent diabetes mellitus. J Clin Endocrinol Metab 1987;64:106-10. 18) von Schenck H, Grubb AO. Interference of immunoglobulins in two glucagon radioimmunoassays. Clin Chem 1982;28:1103-7. 19) Service FJ. Hypoglycemia. In: DeGroot LJ, editor. Endocrinology. Philadelphia: Saunders, 1989;2:1524-40. 20) Shima K, et al. The normal range of the plasma immunoreactive glucagon level during a 75 g oral glucose tolerance test. Diabetes Res Clin Pract 1986;1:265-71. 21) Smith RJ. Biological actions and interactions of insulin and glucagon. In: DeGroot LJ, editor. Endocrinology. Philadelphia: Saunders, 1989;2:1333-45. 22) Steiner DF, Bell GI, Tager HS. Chemistry and biosynthesis of pancreatic protein hormones. In: DeGroot LJ, editor. Endocrinology. Philadelphia: Saunders, 1989;2:1263-89. 23) Tahara Y, et al. Discrepancy of fasting plasma glucagon levels measured by different glucagon antisera. Diabetes Res 1988;7:179-83. 24) Unger RH. Insulin-glucagon relationships in the defense against hypoglycemia. Diabetes. 1983;32:575-83. 25) Unger RH. Glucagon physiology and pathophysiology in the light of new advances. Diabetologia. 1985;28:574-8. 26) Unger RH, Orci L. Glucagon secretion and metabolism in man. In: DeGroot LJ, editor. Endocrinology. Philadelphia: Saunders, 1989;2:1318-32. 27) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782, 310.645.8200
Fax: 310.645.9999. To place an order:
Tel: 800.372.1782. Fax: 800.234.4372.
Outside the United States, contact your National Distributor.

The Quality System of Diagnostic Products Corporation is registered to ISO 13485:2003.

Tables and Graphs

Example Run

Tube ¹	Duplicate CPM ²	Average CPM ³	Net CPM ⁴	% Bound ⁵	Approx. Glucagon pg/mL ⁶
T ⁷	22,823 22,235	22,529			
NSB ⁸	555 600	578	0		
A (MB) ⁹	9,245 9,367	9,306	8,729	100%	0
B	8,328 8,272	8,300	7,723	88%	25
C	7,399 7,308	7,354	6,776	78%	50
D	5,397 5,824	5,611	5,033	58%	100
E	3,358 3,384	3,371	2,794	32%	250
F	2,272 2,259	2,266	1,688	19%	500
Unknowns: ¹⁰					
X1	6,494 6,982	6,738	6,161	71%	65
X2	2,768 2,637	2,703	2,125	24%	357

Quality Control Parameters:¹¹

T⁷ = 22,529 cpm

%NSB⁸ = 2.6%

%MB⁹ = 39%

20% Intercept¹² = 442 pg/mL

50% Intercept = 136 pg/mL

80% Intercept = 42 pg/mL

Intraassay Precision (pg/mL)

	Mean ¹	SD ²	CV ³
1	52	3.4	6.5%
2	209	6.6	3.2%
3	343	16.6	4.8%

Interassay Precision (pg/mL)

	Mean ¹	SD ²	CV ³
1	59	7.0	11.9%
2	205	16.5	8.0%
3	351	21.1	6.0%

End-of-Run Effect

	Tubes 28 – 33	Tubes 89 – 94	Tubes 135 – 140	Tubes 181 – 186
1	51	52	49	48
2	222	195	208	197
3	330	344	331	359

Linearity (pg/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	4 in 4 ⁵	196	—	—
	2 in 4	97.3	98.0	99%
	1 in 4	49.7	49.0	101%
2	4 in 4	212	—	—
	2 in 4	108	106	102%
	1 in 4	51.5	53.0	97%
3	4 in 4	230	—	—
	2 in 4	123	115	107%
	1 in 4	59.8	57.5	104%
4	8 in 8	432	—	—
	4 in 8	239	216	111%
	2 in 8	105	108	97%
	1 in 8	50.1	54.0	93%

Standardization

Concentration pg/mL	DPC Glucagon Calibrators (%B/B ₀)	WHO Preparation in DPC Zero (%B/B ₀)
25	91%	95%
50	87%	86%
100	72%	74%
250	43%	43%
500	24%	25%

Specificity

Compound ¹	pg/mL Added ²	Apparent pg/mL ³	% Cross- reactivity ⁴
Glucagon (Human) M.W. 3,483			
	1,000	954	95%
	500	609	122%
	250	284	114%
Glucagon (22–29, Human) M.W. 1,039			
	1,000	52	5.2%
	500	41	8.2%
	250	24	9.6%
Glucagon (19–29, Human) M.W. 1,352			
	1,000	22.8	2.3%
	500	21.9	4.4%
	250	11.8	4.7%
Glucagon-like Peptide 1 (Preproglucagon 72–108) (GLP-1, Human) M.W. 4,170			
	1,000,000	ND	ND
Glucagon-like Peptide 1 (Preproglucagon 72–107 amide) (GLP-1, Human) M.W. 4,112			
	1,000,000	ND	ND
Glucagon-like Peptide 2 (Preproglucagon 126–159) (GLP-2, Human) M.W. 3,923			
	1,000,000	ND	ND
Gastric Inhibitory Peptide (GIP, Human) M.W. 4,984			
	20,000	ND	ND
Oxyntomodulin (Glucagon 37) M.W. 4,422			
	40,000	42	0.1%
Secretin (Human) M.W. 2,883			
	100,000	ND	ND
Vasoactive Intestinal Peptide (VIP; Human, Porcine, Rat) M.W. 3,326			
	100,000	ND	ND

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O Sistema de Qualidade da Diagnostic Products Corporation está registado sob ISO 13485:2003.
**EC REP** DPC Biermann GmbH