Insulin

DPC®

Coat-A-Count® Insulin

English

Intended Use: Coat-A-Count Insulin is a solid-phase ^{125}I radioimmunoassay designed for the quantitative measurement of insulin in serum. It is intended strictly for *in vitro* diagnostic use as an aid in clinical diagnosis.

Catalog Numbers: **TKIN1** (100 tubes), **TKIN2** (200 tubes), **TKIN5** (500 tubes), **TKINX** (1000 tubes).

The 100-tube kit contains less than 3 microcuries (111 kilobecquerels) of radioactive ^{125}I insulin; the 200-tube kit contains less than 6 microcuries (222 kilobecquerels); the 500-tube kit contains less than 15 microcuries (555 kilobecquerels); and the 1000-tube kit contains less than 30 microcuries (1110 kilobecquerels).

CDC Analyte Identifier Code: **2812**
CDC Test System Identifier Code: **13030**
CLIA Complexity Category: **High**

Summary and Explanation of the Test

Human insulin is a polypeptide hormone originating in the beta cells of the pancreas and serving as a principal regulator for the storage and production of carbohydrates. Its secretion is normally stimulated by increases in the amount of glucose in circulation. This leads to higher insulin levels and more rapid tissue-assimilation of glucose — followed by a decline in the insulin level as the glucose level subsides.

In a number of conditions, notably insulinoma and diabetes, this relationship is impaired. Insulin tends to circulate at inappropriately high levels in patients with insulin-secreting pancreatic tumors; such tumors can thus be a cause of hypoglycemia. Accordingly, insulin immunoassays — used sometimes in connection with provocative doses of tolbutamide or calcium — play an essential role in the identification (and localization) of insulinomas. The finding of fasting hypoglycemia in association with

an *inappropriately high* serum insulin concentration is considered diagnostic.

Insulin levels do not figure in the subclassification of diabetes worked out by the National Diabetes Data Group. Nevertheless, when obtained in the course of a glucose tolerance test, they appear to be of some prognostic value in predicting the benefits of insulin therapy and the likelihood of progression to insulin-dependence and the complications (such as retinopathy) characteristic of diabetes.

The application of insulin radioimmunoassay to patients already undergoing insulin therapy is complicated by the fact that such therapy typically leads to the formation of anti-insulin antibodies capable of interfering with the assay. Some investigators have sought therefore to measure insulin in urine, or in serum samples subjected to column chromatography or PEG precipitation. But the measurement of "free" insulin remains of limited interest as a technique for monitoring insulin therapy in the absence of statistics establishing therapeutic or toxic ranges. So far it appears that glucose control in diabetics cannot in general be achieved by normalizing the insulin profile. Nor is it known at what point abnormally high insulin levels become dangerous.

Principle of the Procedure

The Coat-A-Count Insulin procedure is a solid-phase radioimmunoassay, wherein ^{125}I -labeled insulin competes for a fixed time with insulin in the patient sample for sites on insulin-specific antibody. Because the antibody is immobilized to the wall of a polypropylene tube, simply decanting the supernatant suffices to terminate the competition and to isolate the antibody-bound fraction of the radiolabeled insulin. Counting the tube in a gamma counter then yields a number, which converts by way of a calibration curve to a measure of the insulin present in the patient sample.

Reagents to Pipet: 1

Total Incubation Time: Overnight (18–24 hours) or 3 hours (with Calibrator B omitted).

Total Counts at Iodination:
approximately 40,000 cpm

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C in a refrigerator designated for incoming radioactive materials. Dispose of in accordance with applicable laws.

Do not use reagents beyond their expiration dates.

Some components supplied in this kit may contain human source material and/or other potentially hazardous ingredients which necessitate certain precautions.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Water: Use distilled or deionized water.

Radioactivity

A copy of any radioisotope license certificate (Specific or General) issued to a US customer must be on file with Diagnostic Products Corporation before kits or components containing radioactive material can be shipped. These radioactive materials may be acquired by any customer with the appropriate Specific license. Under a General license these radioactive materials may be acquired only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories and hospitals — and strictly for *in vitro* clinical or laboratory tests not involving external or internal administration of the radioactive material or its radiation to human beings or other animals. Its acquisition, receipt, storage, use, transfer and disposal are all subject to the regulations and a (General or Specific) license of the U.S. Nuclear Regulatory Commission or a State with which the NRC has entered into an

agreement for the exercise of regulatory control.

Handle radioactive materials according to the requirements of your General or Specific license. To minimize exposure to radiation, the user should adhere to guidelines set forth in the National Bureau of Standards publication on the *Safe Handling of Radioactive Materials* (Handbook No. 92, issued March 9, 1964) and in subsequent publications issued by State and Federal authorities.

Wipe up spills promptly and decontaminate affected surfaces. Avoid generation of aerosols. Dispose of solid radioactive waste according to license requirements. General licensees (holders of NRC Form 483) may dispose of solid radioactive waste as nonradioactive waste, after removing labeling. Specific licensees (NRC Form 313) should refer to Title 10, Code of Federal Regulations, Part 20. Licensees in Agreement States should refer to the appropriate regulations of their own state. General licensees may dispose of liquid radioactive waste of the type contained in this product through a laboratory sink drain. Licensees must remove or deface labels from empty containers of radioactive materials before disposal of solid waste. Specific licensees may dispose of small quantities of liquid radioactive waste of the type used in this product through a laboratory sink drain. Refer to the appropriate regulations applicable to your laboratory.

Materials Supplied: Initial Preparation

Insulin Ab-Coated Tubes (TIN1)

Polypropylene tubes coated with antibodies to insulin and packaged in zip-lock bags. Store refrigerated and protected from moisture, carefully resealing the bags after opening. Stable at 2–8°C until the expiration date marked on the bag. Color: lime green.

TKIN1: 100 tubes. **TKIN2:** 200 tubes.

TKIN5: 500 tubes. **TKINX:** 1000 tubes.

¹²⁵I Insulin (TIN2)

A *concentrate*, consisting of iodinated insulin. To each vial add a measured **100 mL** of distilled or deionized water, and mix by *gentle* inversion. Stable at 2–8°C for 30 days after preparation, or until the

expiration date marked on the label.

TKIN1: 1 vial. **TKIN2:** 2 vials.

TKIN5: 5 vials. **TKINX:** 10 vials.

Insulin Calibrators (INC3–9)

Seven vials, labeled A through G, of lyophilized processed in nonhuman serum. *At least 30 minutes before use*, reconstitute the zero calibrator A with **6.0 mL** of distilled or deionized water, and each of the remaining calibrators B through G with **3.0 mL** of distilled or deionized water. Use volumetric pipets and mix by *gentle* swirling. Stable at -20°C for 30 days after reconstitution. Aliquot, if necessary, to avoid repeated freezing and thawing.

TKIN1: 1 set. **TKIN2:** 1 set.

TKIN5: 2 sets. **TKINX:** 3 sets.

The reconstituted calibrators have *lot-specific* insulin values representing *approximately* 0, 5, 15, 50, 100, 200 and 350 micro-International Units of insulin per milliliter ($\mu\text{IU/mL}$, WHO 1st IRP 66/304) in processed in nonhuman serum. Intermediate calibration points can be obtained by mixing the calibrators in suitable proportions.

Materials Required But Not Provided

Gamma counter — compatible with standard 12x75 mm tubes

Vortex mixer

Reagent Preparation

Distilled or deionized water

Graduated cylinder: 100 mL

Volumetric pipets: 3.0 mL and 6.0 mL

Radioimmunoassay

Plain 12x75 mm polypropylene tubes — for use as NSB tubes, available from DPC

Micropipets: 200 μL and 1,000 μL .

Foam decanting rack — available from DPC

Logit-log graph paper — available from DPC (catalog number: ZP797)

A tri-level, human serum-based immunoassay control, containing insulin as one of over 25 assayed constituents, is available from DPC (catalog number: CON6).

Specimen Collection

No special preparations are necessary. Collect blood by venipuncture²⁰ into plain tubes (without anticoagulant), noting the time of collection, and separate the serum from the cells.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed, icteric or grossly contaminated samples may give erroneous results.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. Coat-A-Count Insulin has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 200 μL of serum per tube.

In collecting samples for insulin determinations, it is important to avoid hemolysis, which can lead to spuriously low results.^{22,23}

Storage: $2-8^{\circ}\text{C}$ for 7 days, or for up to 3 months frozen at -20°C .²¹

Before assay, allow the samples to come to room temperature ($15-28^{\circ}\text{C}$) and mix by *gentle* swirling or inversion. Aliquot, if necessary, to avoid repeated thawing and freezing. Do *not* attempt to thaw frozen specimens by heating them in a waterbath. Samples contaminated from the administration of radioactivity to the patient are unsuitable for use.

Radioimmunoassay Procedure

All components must be at room temperature ($15-28^{\circ}\text{C}$) before use. An automatic pipettor-diluter is *not* recommended for this assay.

- 1 Plain Tubes:** Label four plain (uncoated) 12x75 mm polypropylene tubes T (total counts) and NSB (nonspecific binding) in duplicate.

Because nonspecific binding in the Coat-A-Count procedure is characteristically low, the NSB tubes may be safely omitted without compromising accuracy or quality control.

Coated Tubes: Label fourteen Insulin Ab-Coated Tubes A (maximum binding) and B through G in duplicate. Label additional antibody-coated tubes, also in duplicate, for controls and patient samples.

Calibrator	Approximate $\mu\text{IU/mL}$ 1st IRP [66/304]
A (MB)	0
B*	5
C	15
D	50
E	100
F	200
G	350

* Omit calibrator B if the alternate, 3-hour incubation at room temperature (15–28°C) will be used at step 4.

- Pipet **200 μL** of the zero calibrator A into the NSB and A tubes, and **200 μL** of each remaining calibrator, control and patient sample into the tubes prepared. **Pipet directly to the bottom.**

Dilute patient serum samples expected to contain insulin concentrations greater than the highest calibrator with the zero calibrator before assay.

It is good laboratory practice to use a disposable-tip micropipet, changing the tip between samples, in order to avoid carryover contamination.

- Add **1.0 mL** of ^{125}I Insulin to every tube. Vortex.
Samples should not be left in the tubes for extended periods of time. Following step 2 (sample addition), step 3 (tracer addition) should be completed with minimal delay, with no more than 40 minutes elapsing between the addition of the first sample and the completion of tracer addition. Set the T tubes aside for counting (at step 6); they require no further processing.
- Incubate for **18–24 hours** at room temperature (15–28°C).
Alternatively, incubate for 3 hours at room temperature (15–28°C), omitting the 5 $\mu\text{IU/mL}$ calibrator (Calibrator B). This approach allows for same-day

comparable results (see the Performance Data sections).

- Decant thoroughly.

Removing all visible moisture will greatly enhance precision. Using a foam decanting rack, decant the contents of all tubes (except the T tubes) and allow them to drain for 2 or 3 minutes. Then strike the tubes sharply on absorbant paper to shake off all residual droplets.

- Count for **1 minute** in a gamma counter.

Calculation of Results

To calculate insulin concentrations from a logit-log representation of the calibration curve, first calculate for each pair of tubes the average NSB-corrected counts per minute:

Net Counts = (Average CPM) minus (Average NSB CPM)

Then determine the binding of each pair of tubes as a percent of maximum binding (MB), with the NSB-corrected counts of the A tubes taken as 100%:

Percent Bound = (Net Counts / Net MB Counts) $\times$ 100

(The calculation can be simplified by omitting the correction for nonspecific binding; samples within range of the calibrators yield virtually the same results when Percent Bound is calculated directly from Average CPM.)

Using logit-log graph paper, plot Percent Bound on the vertical (probability) axis against Concentration on the horizontal (logarithmic) axis for each of the nonzero calibrators, and draw a straight line approximating the path of these points. Results for the unknowns may then be read from the line by interpolation. (The calibration curve for the alternate 3-hour, room temperature (15–28°C) procedure should be based on calibrators C through G only.)

It is good practice to inspect results for agreement within replicates, and to construct a graph of the calibration curve (even if the calculations are handled by computer) as a visual check on the appropriateness of the transformation used and as a way to detect deviant calibration points. We also recommend

keeping track of these data reduction parameters:

T = Total Counts (as counts per minute)

$\%NSB = 100 \times (\text{Average NSB Counts} / \text{Total Counts})$

$\%MB = 100 \times (\text{Net Counts} / \text{Total Counts})$

And the 20, 50 and 80 percent "intercepts," where

20% = Concentration at 20 Percent Bound, etc.

Note that other approaches, e.g. a sound implementation of the 4-parameter logistic, are also acceptable. See Dudley RA, et al. Guidelines for immunoassay data reduction. Clin Chem 1985;31:1264-71.

Example Run (Basic Procedure): For illustration only, not for calculating results from another run. The values were generated with the overnight procedure. Because the calibrator values are *lot-specific*, concentrations listed in the right-most column may not match the values of the calibrators supplied in your shipment. (See "Example Run – Basic Procedure" table.)

Example Run (Alternate 3-hour Procedure): For illustration only, not for calculating results from another run. The values were generated with the 3-hour room-temperature procedure. Because the calibrator values are *lot-specific*, concentrations listed in the right-most column may not match the values of the calibrators supplied in your shipment. (See "Example Run – Alternate Procedure" table.)

Quality Control

Controls or serum pools with at least two insulin concentration levels (low and high) should routinely be assayed as unknowns.

Report patient results only if the control results for that assay meet your laboratory's established criteria for acceptability.

It is good laboratory practice to record for each assay the lot numbers of the components used, as well as the dates when they were first reconstituted or opened. We also recommend charting control results from day to day — as described, for example, in Westgard JO, et al. A multi-rule chart for quality control. Clin Chem 1981;27:493-501. Note that

repeat samples can serve as a valuable additional tool for monitoring interassay precision, and that pairs of control tubes can be spaced throughout the assay to help verify the absence of significant drift.

Expected Values

In a study of 103 serum samples from fasting adult volunteers with fasting glucose levels less than or equal to 110 mg/dL, the Coat-A-Count Insulin assay yielded a median of 8.3 $\mu\text{IU/mL}$, with 95 percent of the results being 22 $\mu\text{IU/mL}$ or less.

Laboratories should consider these results *as guidelines only*. Each laboratory should establish its own reference ranges.

Performance Data

See Tables and Graphs for data *representative* of the Coat-A-Count Insulin. Insulin results in the sections below are expressed as $\mu\text{IU/mL}$. Except as noted, all results were obtained using an overnight, room-temperature incubation.

Calibration Range: 5 – 350 $\mu\text{IU/mL}$ (WHO 1st IRP 66/304). The calibrators are *lot-specific*.

Analytical Sensitivity:

Basic Overnight Procedure: 1.2 $\mu\text{IU/mL}$

Alternate 3-hour incubation: 4.6 $\mu\text{IU/mL}$

Precision Profile: A precision profile is displayed, based on approximately 20 degrees of freedom and representing the intraassay CVs to be expected for samples assayed in duplicate using the Basic Procedure. (See "Precision Profile" graph.)

Intraassay Precision (Within-Run):

Statistics were calculated for each of five samples from the results of 20 pairs of tubes in a single assay. The study was done using the Alternate Procedure. (See "Intraassay Precision" table.)

Interassay Precision (Run-to-Run):

Statistics were calculated for each of four samples from the results of pairs of tubes in 20 different assays. The study was done using the Basic Procedure. (See "Interassay Precision" table.)

Specificity: The Coat-A-Count Insulin antiserum is highly specific for insulin, with a particularly low crossreactivity to compounds (other than insulin and pro-insulin) that might be present in patient samples. (See "Specificity" table.)

End-of-Run Effect: None up to approximately 250 tubes. (See "End-of-Run Effect" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three insulin solutions (329, 1,420 and 3,535 $\mu\text{IU/mL}$) were assayed. (See "Recovery" table for representative data.)

Bilirubin: Severe icterus (bilirubin up to 200 mg/L) may cause an increase of values. (See "Bilirubin" table.)

Hemolysis: Presence of packed red blood cells in concentrations up to 30 $\mu\text{L/mL}$ can lead to spuriously low results.^{22,23} (See "Hemolysis" table.)

Lipemia: To simulate various degrees of lipemia, patient samples were spiked with a fat emulsion consisting of a mixture of lipids. The stock solution had a triglycerides concentration of 20,000 mg/dL . The results showed that lipemia has no clinically significant effect on the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 20 volunteers into plain, heparinized, EDTA and Becton Dickinson SST[®] vacutainer tubes. All samples were assayed by the Coat-A-Count Insulin procedure, with the following results.

(Heparin) = $0.863 (\text{Serum}) + 1.05 \mu\text{IU/mL}$
 $r = 0.972$

(EDTA) = $0.774 (\text{Serum}) + 0.54 \mu\text{IU/mL}$
 $r = 0.952$

(SST) = $0.97 (\text{Plain Tubes}) + 1.51 \mu\text{IU/mL}$
 $r = 0.986$

Means:

21.3 $\mu\text{IU/mL}$ (Serum)
19.4 $\mu\text{IU/mL}$ (Heparin)
17.0 $\mu\text{IU/mL}$ (EDTA)
22.2 $\mu\text{IU/mL}$ (SST)

Protein Effect: To simulate various protein concentrations, experiments were performed in which 2.0 mL aliquots of a human serum pool were freeze-dried and

then reconstituted with various volumes of water. Each reconstituted aliquot was then assayed by the Coat-A-Count Insulin procedure. Note that aliquots reconstituted with half the original volume represent an extremely high protein concentration, on the order of 14 g/dL . (The factor to correct for reconstitution volume is tabulated below.) These results indicate that even wide variations in protein have virtually no effect on the Coat-A-Count Insulin assay. (See "Protein Effect" table.)

Method Comparison: The Coat-A-Count assay was compared to DPC's IMMULITE Insulin assay on 172 patient samples, with insulin concentrations ranging from approximately 5.3 to 200 $\mu\text{IU/mL}$. (See graph.) By linear regression:

(CAC) = $0.90 (\text{IML}) + 0.99 \mu\text{IU/mL}$
 $r = 0.973$

Means:

36 $\mu\text{IU/mL}$ (Coat-A-Count)
39 $\mu\text{IU/mL}$ (IMMULITE)

References

- 1) Bates HM. Insulinoma and pheochromocytoma. Lab Management 1983 (May); 21:11–12, 15. 2) Bennett PH. The diagnosis of diabetes; new international classification and diagnostic criteria. Ann Rev Med 1983;34:295–309. 3) Berntorp K, et al. Relation between plasma insulin and blood glucose in a cross-sectional population study of the oral glucose tolerance test. Acta Endocrinol 1983;102:549–56. 4) Boehm TM, Lebovitz HE. Statistical analysis of glucose and insulin responses to intravenous tolbutamide: evaluation of hypoglycemic and hyperinsulinemic states. Diabetes Care 1979;2:479–90. 5) Fajans SS, Floyd JC. Diagnosis and medical management of insulinomas. Ann Rev Med 1979;30:313–29. 6) Hirata Y, et al. Autoimmunity in diabetes mellitus. Gunma Symposia Endocrinol 1980;17:81–91. 7) Kahn CR, Rosenthal AS. Immunologic reactions to insulin: insulin allergy, insulin resistance and the autoimmune insulin syndrome. Diabetes Care 1979;2:283–95. 8) Kaplan EL, Rubinstein AH, et al. Calcium infusion: a new provocative test for insulinomas. Annals Surgery 1979;190:501–7. 9) Malone JI, Root AW. Plasma free insulin concentrations: keystone to effective management of diabetes mellitus in children. J Pediatr 1981;99:862–7. See also *ibid* 1983;102:800–2. 10) Marschner I, et al. Group experiments on the radioimmunological insulin determination. Horm Metabol Res 1974;6:293–6. 11) National Diabetes Data Group. Classification and diagnosis of diabetes mellitus and other categories of glucose intolerance. Diabetes

1979;28:1039–57. 12) Pfeifer MA, et al. Insulin secretion in diabetes mellitus. *Am J Med* 1981;70:579–88. 13) Reeves WG. Insulin antibody determination: theoretical and practical considerations. *Diabetologia* 1983;24:399–403. 14) Sachse G, Willms B. Klinische Bedeutung der Seruminsulinbestimmung: Diagnose und Prognose bei drohendem Sekundrversagen der Sulfonylharnstofftherapie. *Medizin Klinik* 1979;74:1635–40. 15) Schade DS, Eaton RP, Spencer W. Normalization of plasma insulin profiles in diabetic subjects with programmed insulin delivery. *Diabetes Care* 1980;3:9–14. 16) Service FJ, et al. Insulinoma: clinical and diagnostic features of 60 consecutive cases. *Mayo Clinic Proceed* 1976;51:417–29. 17) Turkington RW, Estkowski A, Link M. Secretion of insulin or connecting peptide; a predictor of insulin dependence of obese diabetics. *Arch Intern Med* 1982;142:1102–5. 18) Turkington RW, Weindling HK. Insulin secretion in the diagnosis of adult-onset diabetes mellitus. *JAMA* 1978;240:833–36. See also *ibid* 1979;241:1462–3. 19) Yalow R, Bauman WA. Plasma insulin in health and disease. In: Ellenberg M, Rifkin H, editors. *Diabetes mellitus; theory and practice*. New York: Excerpta Medica, 1983: 119–50. 20) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 21) Burtis CA, Ashwood ER, editors. *Tietz textbook of clinical chemistry*. 2nd ed. Philadelphia: WB Saunders, 1994:943–4. 22) Sapin R, Ongagna JC, Gasser F, Grucker D. Insulin measurements in haemolysed serum; influence of insulinase inhibitors. *Clin Chim Acta* 1998;274:111–7. 23) Chevenne D, Letailleur A, Trivin F, Porquet D. Effect of hemolysis on the concentration of insulin in serum determined by RIA and IRMA. *Clin Chem* 1998;44:354–6.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782, 310.645.8200
Fax: 310.645.9999. To place an order:
Tel: 800.372.1782. Fax: 800.234.4372.
Outside the United States, contact your National Distributor.

The Quality System of Diagnostic Products Corporation is registered to ISO 13485:2003.

Tables and Graphs

Example Run (Basic Procedure)

Tube ¹	Duplicate CPM ²	Average CPM ³	Net CPM ⁴	Percent Bound ⁵	Approx. Insulin μ U/mL ⁶
T ⁷	36,995 37,641	37,318			
(NSB) ⁸	365 344	355			
A (MB) ⁹	13,266 13,326	13,296	12,942	100%	0
B	10,516 10,862	10,689	10,335	80%	8.5
C	8,722 8,480	8,601	8,247	64%	24
D	5,654 5,612	5,633	5,279	41%	64
E	4,082 3,947	4,015	3,660	28%	122
F	2,878 2,848	2,863	2,509	19%	225
G	1,915 1,910	1,913	1,558	12%	420
Unknowns: ¹⁰					
X1	10,554 10,771	10,663	10,308	80%	9.1
X2	7,137 7,204	7,171	6,816	53%	38
X3	4,243 4,085	4,164	3,810	29%	117

Quality Control Parameters:¹¹

T⁷ = 37,318 cpm

%NSB⁸ = 0.95%

%MB⁹ = 35%

20% Intercept¹² = 209 μ U/mL

50% Intercept = 43 μ U/mL

80% Intercept = 8.9 μ U/mL

Example Run (Alternate Procedure)

Tube ¹	Duplicate CPM ²	Average CPM ³	Net CPM ⁴	Percent Bound ⁵	Approx. Insulin μ U/mL ⁶
T ⁷	46,119 46,295	46,207			
(NSB) ⁸	476 449	463			
A (MB) ⁹	9,400 9,569	9,485	9,022	100%	0
C	7,492 7,557	7,525	7,062	78%	24
D	5,695 5,741	5,718	5,256	58%	64
E	4,221 4,164	4,193	3,730	41%	122
F	2,926 2,988	2,957	2,495	28%	225
G	2,183 2,123	2,153	1,691	19%	420

Unknowns:¹⁰

X1	7,069 6,975	7,022	6,560	73%	32
X2	4,505 4,753	4,629	4,167	46%	103

Quality Control Parameters:¹¹

T⁷ = 46,207 cpm

%NSB⁸ = 1.0%


%MB⁹ = 20%

20% Intercept¹² = 359 μ U/mL

50% Intercept = 88 μ U/mL

80% Intercept = 21 μ U/mL

Precision Profile

Intraassay Precision (μ U/mL) (Alternate procedure)

	Mean ¹	SD ²	CV ³
1	17	1.6	9.3%
2	39	2.0	5.1%
3	80	2.8	3.5%
4	117	3.6	3.1%
5	278	14	5%

Interassay Precision (μ U/mL) (Basic Procedure)

	Mean ¹	SD ²	CV ³
1	16	1.6	10.0%
2	35	2.5	7.1%
3	76	5.5	7.2%
4	95	4.7	4.9%

End-of-Run Effect (μ U/mL)

	Tubes ¹ 31-44	Tubes 85-98	Tubes 139-152	Tubes 193-206	Tubes 247-260
1	8.6	7.6	8.3	7.7	9.3
2	12	10	11	9	11
3	13	14	14	15	16
4	30	32	33	29	30
5	40	41	37	39	41
6	110	116	122	116	112
7	116	117	116	123	124

Linearity (μIU/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	32 in 32 ⁵	126	—	—
	16 in 32	65	63	103%
	8 in 32	31	31	100%
	4 in 32	14	16	88%
	2 in 32	7.4	7.8	95%
	1 in 32	3.9	3.9	100%
2	32 in 32	148	—	—
	8 in 16	79	74	107%
	4 in 16	43	37	116%
	2 in 16	18	19	95%
	1 in 16	8.6	9.2	93%
3	32 in 32	169	—	—
	8 in 16	86	85	101%
	4 in 16	45	42	107%
	2 in 16	22	21	105%
	1 in 16	9.4	10.5	90%

Specificity


Compound ¹	Amount Added ² ng/mL	% Cross-reactivity ³
Pro-insulin	5	32%
C-Peptide	1,000	ND
Glucagon	1,000	ND

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O Sistema de Qualidade da Diagnostic Products Corporation está registado sob ISO 13485:2003.
EC REP DPC Biermann GmbH