Buprenorphine

DOUBLE ANTIBODY

DPC[®]

Double Antibody Buprenorphine

English

Intended Use: Buprenorphine is a double-antibody ^{125}I radioimmunoassay designed for the quantitative and qualitative measurement of buprenorphine in urine. It is intended strictly for *in vitro* diagnostic use in the context of a program involving an established confirmatory test for buprenorphine and its metabolites.

Catalog Numbers: **KBPD1** (100 tubes), **KBPD2** (200 tubes).

The 100-tube kit contains less than 4 microcuries (148 kilobecquerels) of radioactive ^{125}I buprenorphine, and the 200-tube kit contains less than 8 microcuries (296 kilobecquerels).

The Double Antibody Buprenorphine kit provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method.¹ Clinical consideration and professional judgement should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

Summary and Explanation of the Test

Buprenorphine is an oripavine derivative and is a potent mixed opiate agonist/antagonist used for pain relief. It is administered clinically by intravenous, intramuscular or sublingual routes.²

Buprenorphine is metabolized by *N*-dealkylation and 3-*O*-glucuronidation.³⁻⁷ Plasma levels of metabolites are undetectable after single doses, but in chronic dosing they are of the same order as unchanged buprenorphine.⁸ The major part (about 70%) of a dose of buprenorphine appears in the feces,⁶ and only small amounts are found in the urine, mainly as buprenorphine glucuronide,

N-dealkylbuprenorphine glucuronide and *N*-dealkylbuprenorphine.^{6,7} There is little unchanged buprenorphine in urine.

Effective doses of buprenorphine are of the order of a few hundred micrograms. Combined with extensive metabolism and low urine excretion, levels of buprenorphine equivalents in urine will be low. In a study of patients taking 0.4 to 2.0 mg buprenorphine daily for the relief of chronic pain, urine concentrations measured by the DPC Double Antibody Buprenorphine kit ranged from 5.6 to more than 100 ng/mL. The effect of enzymatic hydrolysis with β -glucuronidase and sulfatase gave a mean increase in observed concentration of 48%.

Because analysis of buprenorphine by conventional techniques is difficult, radioimmunoassay has been widely used.^{8,9} HPLC methods have been reported, but have limited sensitivity for use as confirmatory tests.^{10,11} GC/MS, by comparison, has better sensitivity.^{6,7}

Principle of the Procedure

DPC's Double Antibody Buprenorphine procedure is a competitive radioimmunoassay in which ^{125}I -labeled buprenorphine competes with buprenorphine in the sample for antibody binding sites. After incubation for a fixed time, separation of bound from free is achieved by the PEG-accelerated double-antibody method. Finally, the antibody-bound fraction is precipitated and counted. Patient sample concentrations are read from a calibration curve.

Reagents to pipet: 3

Total Incubation Time: 70 minutes

Separation: The ready-to-use Precipitating Solution combines second antibody and dilute PEG.

Total counts at iodination: approximately 100,000 cpm

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C in a refrigerator designated for incoming radioactive

materials. Dispose of in accordance with applicable laws.

Do not use reagents beyond their expiration dates.

Some components supplied in this kit may contain human source material and/or other potentially hazardous ingredients which necessitate certain precautions.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Water: Use distilled or deionized water.

Radioactivity

A copy of any radioisotope license certificate (Specific or General) issued to a US customer must be on file with Diagnostic Products Corporation before kits or components containing radioactive material can be shipped. These radioactive materials may be acquired by any customer with the appropriate Specific license. Under a General license these radioactive materials may be acquired only by physicians, veterinarians in the practice of veterinary medicine, clinical laboratories and hospitals — and strictly for *in vitro* clinical or laboratory tests not involving external or internal administration of the radioactive material or its radiation to human beings or other animals. Its acquisition, receipt, storage, use, transfer and disposal are all subject to the regulations and a (General or Specific) license of the U.S. Nuclear Regulatory Commission or a State with which the NRC has entered into an agreement for the exercise of regulatory control.

Handle radioactive materials according to the requirements of your General or Specific license. To minimize exposure to radiation, the user should adhere to guidelines set forth in the National Bureau of Standards publication on the *Safe*

Handling of Radioactive Materials (Handbook No. 92, issued March 9, 1964) and in subsequent publications issued by State and Federal authorities.

Wipe up spills promptly and decontaminate affected surfaces. Avoid generation of aerosols. Dispose of solid radioactive waste according to license requirements. General licensees (holders of NRC Form 483) may dispose of solid radioactive waste as nonradioactive waste, after removing labeling. Specific licensees (NRC Form 313) should refer to Title 10, Code of Federal Regulations, Part 20. Licensees in Agreement States should refer to the appropriate regulations of their own state. General licensees may dispose of liquid radioactive waste of the type contained in this product through a laboratory sink drain. Licensees must remove or deface labels from empty containers of radioactive materials before disposal of solid waste. Specific licensees may dispose of small quantities of liquid radioactive waste of the type used in this product through a laboratory sink drain. Refer to the appropriate regulations applicable to your laboratory.

Materials Supplied: Initial Preparation

Prepare all components at least 10 minutes before use.

Buprenorphine Antiserum (BPD1)

Lyophilized buprenorphine antiserum, with preservative. Reconstitute each vial by adding **11 mL** distilled or deionized water. Mix by *gentle* inversion. Stable at 2–8°C for 30 days after reconstitution.

KBPD1: 1 vial **KBPD2:** 2 vials

¹²⁵I Buprenorphine (BPD2)

Lyophilized iodinated buprenorphine, with preservative. Reconstitute each vial by adding **11 mL** distilled or deionized water. Mix by *gentle* inversion. Stable at 2–8°C for 30 days after reconstitution, or until the expiration date marked on the label.

KBPD1: 1 vial **KBPD2:** 2 vials

Buprenorphine Calibrators (BPD3–8)

Six vials, labeled A (3 mL) and B through F (1 mL each), of buprenorphine calibrators supplied in liquid form, ready to use with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months

(aliquotted) at –20°C.

KBPD1: 1 set **KBPD2:** 1 set

The calibrators have *lot-specific* values of *approximately* 0, 0.5, 1, 2, 5 and 10 ng/mL in processed human urine. Refer to the vial labels for exact concentrations in ng/mL. Intermediate calibration points can be obtained by mixing calibrators in suitable proportions.

Precipitating Solution (BPGG)

Precipitating Solution consisting of goat anti-rabbit gamma globulin (GARGG) and dilute polyethylene glycol (PEG) in saline, with preservative. The Precipitating Solution is supplied in liquid form, ready to use. Stable at 2–8°C for 30 days after opening. Since a fine precipitate may form after refrigeration, the Precipitating Solution should be thoroughly mixed before use, without foaming. Color: blue.

KBPD1: 1 vial **KBPD2:** 2 vials

Materials Required But Not Provided

Gamma counter – compatible with standard 12x75 mm tubes

Centrifuge – preferably refrigerated and capable of at least 1500xg

Vortex mixer

Reagent Preparation

Distilled or deionized water

Pipets – for delivering 11 mL

Radioimmunoassay

Plain 12x75 mm polypropylene tubes

Micropipets – 25 µL, 100 µL, and 1.0 mL

Foam decanting rack – available from DPC

Logit-log graph paper – available from DPC (catalog number: ZP397)

Specimen Collection

Collect urine without preservative. The specimen can be refrigerated at 2–8°C or frozen at –20°C. If cloudy, it should be cleared by filtration or centrifugation before use, and mixed by gentle swirling.

Volume Required: 25 µL of urine per tube.

Samples with buprenorphine concentrations greater than the highest calibrator (10 ng/mL) should be diluted

either with the zero calibrator or with buprenorphine-free urine.

If adulteration of the specimen is suspected, do *not* accept for analysis.

Radioimmunoassay Procedure

All components except the Precipitating Solution must be at room temperature (15–28°C) before use.

- 1 Label sixteen tubes in duplicate: T (total counts), NSB (nonspecific binding), A (maximum binding) and B through F. Label additional tubes, also in duplicate, for urine samples and controls.

Calibrator	Approximate ng/mL
A (MB)	0
B	0.5
C	1
D	2
E	5
F	10

Note: The calibrators are *lot-specific*. Refer to the vial labels for exact concentrations in ng/mL.

- 2 Pipet **25 µL** of the zero calibrator A into the NSB and A tubes, and **25 µL** of each of the remaining calibrators B through F into correspondingly labeled tubes. Pipet **25 µL** of each urine sample and each control into the tubes prepared.

Patient samples expected to contain buprenorphine concentrations greater than the highest calibrator (10 ng/mL) should be diluted in the zero calibrator before assay.

- 3 Add **100 µL** of ¹²⁵I Buprenorphine to all tubes.

A repeating dispenser is recommended for this step and for the addition of antiserum at step 4. Remove the T tubes for counting at step 10; they require no further processing.

- 4 Add **100 µL** of Buprenorphine Antiserum to all tubes except the NSB and T tubes. Vortex.

- 5 Incubate for **60 minutes** at room temperature (15–28°C).

- 6 Add **1.0 mL** of *cold* Precipitating Solution to all tubes. Vortex.

Before dispensing, the Precipitating Solution should be thoroughly mixed (without foaming) until all visible precipitate has disappeared.
- 7 Incubate for **10 minutes** at room temperature (15–28°C).
- 8 Centrifuge for **15 minutes** at 3000xg – or for a suitably longer period of time at a lower acceleration, e.g. 30 minutes at 1500xg.
- 9 Using a foam decanting rack, decant (or aspirate) the supernatant, retaining the precipitate for counting.

Let the tubes stand inverted on absorbant paper for 10 minutes. Then tap the tubes gently and blot the rims, to remove residual droplets.
- 10 Count each tube for **1 minute**.

Calculation of Results

To calculate buprenorphine concentrations from a logit-log representation of the calibration curve, first calculate for each pair of tubes the average NSB-corrected counts per minute:

Net Counts = (Average CPM) minus (Average NSB CPM)

Then determine the binding of each pair of tubes as a percent of maximum binding (MB), with the NSB-corrected counts of the A tubes taken as 100%:

Percent Bound = (Net Counts / Net **MB** Counts) × 100

Using logit-log graph paper, plot Percent Bound on the vertical (probability) axis against Concentration on the horizontal (logarithmic) axis for each of the nonzero calibrators, and draw a straight line approximating the path of these points. Results for the unknowns may then be read from the line by interpolation.

It is good practice to inspect results for agreement within replicates, and to construct a graph of the calibration curve (even if the calculations are handled by computer) as a visual check on the appropriateness of the transformation used and as a way to detect deviant

calibration points. We also recommend keeping track of these data reduction parameters:

T = Total Counts (as counts per minute)

%NSB = $100 \times (\text{Average NSB Counts} / \text{Total Counts})$

%MB = $100 \times (\text{Net Counts} / \text{Total Counts})$

And the 20, 50 and 80 percent "intercepts," where

20% = Concentration at 20 Percent Bound, etc.

Note that other approaches, e.g. a sound implementation of the 4-parameter logistic, are also acceptable. See Dudley RA, et al. Guidelines for immunoassay data reduction. Clin Chem 1985;31:1264-71.

Example Run: For illustration only, not for calculating results from another run. Because the calibrator values are *lot-specific*, concentrations listed in the right-most column may not match the values of the calibrators supplied in your shipment. (See "Example Run" table.)

Quality Control

Controls or urine pools with at least two buprenorphine concentration level, low and high, should routinely be assayed as unknowns.

Report patient results only if the control results for that assay meet your laboratory's established criteria for acceptability.

It is good laboratory practice to record for each assay the lot numbers of the components used, as well as the dates when they were first reconstituted or opened. We also recommend charting control results from day to day — as described, for example, in Westgard JO, et al. A multi-rule chart for quality control. Clin Chem 1981;27:493-501. Note that repeat samples can serve as a valuable additional tool for monitoring interassay precision, and that pairs of control tubes can be spaced throughout the assay to help verify the absence of significant drift.

Centrifugation: The procedure calls for centrifuging at 3000xg for 15 minutes. Lower accelerations are satisfactory only if the centrifugation time is increased appropriately — for example, to 30 minutes at 1500xg. A high-speed, refrigerated centrifuge is desirable but not essential.

Use the formulas below to calculate the acceleration of your centrifuge at a given speed, or the speed (in revolutions per minute) required to achieve a desired g force.

$$g \text{ force} = 28.38 \times (\text{rpm} / 1,000)^2 \times \text{radius [inches]}$$

$$\text{rpm} = 187.7 \times \sqrt{(g \text{ force} / \text{radius [inches]})}$$

or

$$g \text{ force} = 11.18 \times (\text{rpm} / 1,000)^2 \times \text{radius [cm]}$$

$$\text{rpm} = 299.07 \times \sqrt{(g \text{ force} / \text{radius [cm]})}$$

Limitations

Based on a review of the literature and on studies summarized in the Performance Data section, there do not appear to be any substances that could cause a false positive reaction. Nevertheless, there is the possibility that substances and/or factors not listed in the Specificity section — e.g. technical or procedural errors — may interfere with the test and cause false positive results.

The kit measures buprenorphine and crossreacts extensively with buprenorphine metabolites. Accordingly, the kit should not be used for pharmacokinetic analysis.

Performance Data

See Tables and Graphs for data *representative* of the Double Antibody Buprenorphine kit's performance. Results are expressed in ng/mL.

Calibration Range:

Approximately 0.5 – 10 ng/mL.

The calibrator values are lot-specific.

Analytical Sensitivity: 0.1 ng/mL.

Intraassay Precision (Within-Run):

Statistics were calculated for samples from the results of 20 pairs of tubes in a single assay. (See "Intraassay Precision" table.)

Interassay Precision (Run-to-Run):

Statistics were calculated for samples from the results of pairs of tubes in 20 different assays. (See "Interassay Precision" table.)

End-of-Run Effect: none up to 100 tubes. (See "End-of-run Effect" table.)

Specificity: The antiserum is highly specific for buprenorphine, with an extremely low crossreactivity to other compounds that may be present in equine

urine samples. All compounds tested were weighed out and aliquoted in urine, and then frozen at –20°C until assay. (See "Specificity" table.)

The following compounds were found to be nondetectable by the Double Antibody Buprenorphine procedure at a level of 100,000 ng/mL.

Acetaminophen
Acetylsalicylic acid
Allylcyclopentylbarbituric acid
Amobarbital
d-Amphetamine
d,l-Amphetamine
Aprobarbital
Barbituric acid
Benzoylcegonine (cocaine metab.)
Butabarbital
Butethal
Caffeine
Chlorazepate
Chlorpromazine
Clonazepam
Cocaine
Codeine
Cotinine (nicotine metabolite)
5,5-Diallylbarbituric acid
Diazepam
Ecgonine (cocaine metabolite)
Ethylmorphine
Fenfluramine
Fentanyl
Lidocaine
Methadone
Methamphetamine
Methaqualone
Morphine
Morphine-3-glucuronide
Nalorphine
Normorphine
Phencyclidine
Phenobarbital
d-Propoxyphene
Secobarbital

Linearity: Three urine samples from subjects who had buprenorphine administered for pain relief were assayed, both undiluted and diluted. Some loss in linearity should be expected with actual samples due to the presence of buprenorphine metabolites.¹² (See "Linearity" Table.)

Recovery: Samples spiked 1 to 19 with three Buprenorphine solutions (40, 60 and 120 ng/mL) were assayed. (See "Recovery" table for representative data.)

References

1) National Institute on Drug Abuse (NIDA).
Urine testing for drugs of abuse. Research

Monograph 1986; 73. 2) McQuay HJ, Moore RA, Bullingham RES. Buprenorphine kinetics. In: Foley K, Inturrisi C, editors. Advances in pain research and therapy. New York: Raven Press, 1986: 271-8. 3) Brewster D, Humphrey MJ, McLeavy MA. Biliary excretion, metabolism and enterohepatic circulation of buprenorphine. *Xenobiotica* 1981;11:189-96. 4) Brewster D, Humphrey MJ, McLeavy MA. The systemic bioavailability of buprenorphine by various routes of administration. *J Pharm Pharmacol* 1981;33:500-6. 5) Rance MJ, Shillingford JS. The metabolism of phenolic opiates by rat intestine. *Xenobiotica* 1977;7:529-36. 6) Cone EJ, et al. The metabolism and excretion of buprenorphine in humans. *Drug Metab Dispos* 1984;12:577-81. 7) Blom Y, Bondesson U, Anggard E. Analysis of buprenorphine and its *N*-dealkylated metabolite in plasma and urine by selected-ion monitoring. *J Chromatogr* 1985;338:89-98. 8) Hand CW, et al. Radioimmunoassay of buprenorphine with iodine label: analysis of buprenorphine and metabolites in human plasma. *Ann Clin Biochem* 1986;23:47-53. 9) Bartlett AJ, et al. The radioimmunoassay of buprenorphine. *Eur J Clin Pharmacol* 1980;18:339-45. 10) Tebbett IR. Analysis of buprenorphine by high-performance liquid chromatography. *J Chromatogr* 1985;347:411-3. 11) Noda J, Kojima T. Quantitative determination of buprenorphine in human plasma by high-performance liquid chromatography. *Hiroshima J Med Sci* 1985;43:261-5. 12) Lee JW, Shaver JM. Sensitive and specific radioimmunoassays for opiates using commercially available materials. III. Buprenorphine assay. *J Clin Immunoassay* 1991;14:49. 13) Hand CW, et al. Radioimmunoassay of buprenorphine in urine: studies in patients and in a drug clinic. *J Anal Toxicol* 1989;13:100-4.

Technical Assistance

Available outside the United States only.
For technical assistance, contact your National Distributor.

The Quality System of Diagnostic Products Corporation is registered to ISO 13485:2003.

Tables and Graphs

Example Run

Tube	Duplicate CPM	Average CPM	Net CPM	Percent Bound	Approx. Bupren. ng/mL
T	64,577 64,114	64,346			
(NSB)	905 882	894			

Tube	Duplicate CPM	Average CPM	Net CPM	Percent Bound	Approx. Bupren. ng/mL
A(MB)	38,494 37,940	38,217	37,323	100%	0
B	31,966 31,684	31,825	30,931	82.9%	0.5
C	26,396 26,262	26,329	25,435	68.2%	1
D	19,187 18,277	18,732	17,838	47.8%	2
E	10,254 10,056	10,155	9,261	24.8%	5
F	5,985 5,587	5,786	4,892	13.1%	10

Unknowns:

X1	16,858 16,463	16,660	15,766	42.2%	2.5
X2	8,948 8,827	8,888	7,994	21.4%	5.9
X3	7,448 7,103	7,276	6,383	17.1%	7.5

Quality Control Parameters:

T = 64,346 cpm
%NSB = 1.4%
%MB = 58%
20% Intercept = 6.4 ng/mL
50% Intercept = 1.9 ng/mL
80% Intercept = 0.58 ng/mL

Intraassay Precision (ng/mL)

	Mean ¹	SD ²	CV ³
1	1.4	0.09	6.4%
2	2.4	0.17	7.1%
3	4.9	0.29	5.9%

Interassay Precision (ng/mL)

	Mean ¹	SD ²	CV ³
1	1.5	0.11	7.3%
2	2.8	0.15	5.4%
3	6.2	0.41	6.6%

End-of-Run Effect (ng/mL)

	Tubes ¹ 17–22	Tubes 73–78	Tubes 145–150
1	1.9	1.8	1.9
2	4.0	4.1	4.2
3	7.6	7.5	8.4

Specificity

Compound ¹	ng/mL Added ²	Apparent ng/mL ³	% Cross- reactivity ⁴
N-dealkyl- buprenorphine	1	2.3	234%
	10	9.8	98.1%
	100	22.1	22.1%
Etorphine	100	1.8	1.8%
	1,000	3.8	0.38%
	10,000	7.1	0.07%
Diprenorphine	100,000	6.56	0.0066%
	50,000	4.20	0.0084%
	10,000	1.69	0.0169%
	5,000	1.22	0.0244%
	1,000	0.64	0.064%
Nalbuphine	500	0.48	0.096%
	100,000	1.52	0.0015%
	50,000	0.56	0.0011%

Linearity (ng/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8	8.5	—	—
	4 in 8	5.1	4.3	119%
	2 in 8	2.8	2.1	133%
	1 in 8	1.4	1.1	127%
2	8 in 8	8.7	—	—
	4 in 8	5.2	4.3	121%
	2 in 8	2.7	2.2	123%
	1 in 8	1.4	1.1	127%
3	8 in 8	10.4	—	—
	4 in 8	6.4	5.2	123%
	2 in 8	3.3	2.6	127%
	1 in 8	1.5	1.3	115%

Recovery (ng/mL)

	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	0	—	—
	A	1.6	2.0	80%
	B	2.6	3.0	87%
	C	5.5	6.0	92%
2	—	0.5	—	—
	A	2.6	2.5	104%
	B	3.7	3.5	106%
	C	6.0	6.5	92%
3	—	2.9	—	—
	A	5.2	4.8	108%
	B	5.5	5.8	95%
	C	7.8	8.8	89%
4	—	4.4	—	—
	A	6.5	6.2	105%
	B	7.0	7.2	97%
	C	11.1	10.2	109%

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Il Sistema Qualità della Diagnostic Products Corporation è certificato secondo le norme ISO 13485:2003.
EC REP DPC Biemann GmbH