# **AFP Trisomy 21 Application Sheet**

**DPC®**

# AFP Trisomy 21 Application Sheet

## English

### Intended Use:

IMMULITE/IMMULITE 1000 AFP (LKAP), IMMULITE 2000 AFP (L2KAP) and IMMULITE 2500 AFP (L5KAP) may be used for the risk assessment for Down's syndrome (trisomy 21), Edwards syndrome (trisomy 18), as well as for neural tube defects (NTD), in the second trimester of pregnancy. A pregnant woman's risk, as calculated by the PRISCA 4.0 software, is not confirmatory for these chromosomal abnormalities and for NTD; rather, it is intended for *in vitro* diagnostic use as additional support for her decision whether or not to undergo a diagnostic procedure.

To calculate the risk assessment during the second trimester of pregnancy, AFP test (LKAP, L2KAP, L5KAP) should be used in conjunction with HCG test (LKCG, L2KCG, L5KCG) – **Double Test**; or in conjunction with HCG test (LKCG, L2KCG, L5KCG) and with Unconjugated Estriol test (LKEF, L2KEF, L5KEF) – **Triple Test**.

The concentration of AFP in a given specimen determined with assays from different manufacturers can vary due to differences in assay methods and reagent specificity. The results reported by the laboratory to the physician must include the identity of the assay used. Values obtained with different AFP assays cannot be used interchangeably. Before changing assays, the laboratory must establish a range of median values from non-affected pregnancies for the new assay based on normal sera from pregnant women with confirmed gestational age.

**WARNING:** This assay has not been FDA approved for the Down's syndrome (trisomy 21) or for the Edwards syndrome (trisomy 18) risk assessment in the US.

## Summary and Explanation

### Down's Syndrome (Trisomy 21)

Down's syndrome (trisomy 21) is the most common chromosomal abnormality in newborns, with an incidence of approximately 1:800. The only prenatal test diagnostic for chromosomal abnormalities is cytogenetic examination (karyotyping) of fetal cells using invasive procedures such as amniocentesis or chorion villus sampling (CVS).

Bases for the selection of pregnant women to undergo diagnostic procedures of this kind include: advanced maternal age, the subject's medical history, or a statistical risk assessment derived from maternal age plus some combination of biochemical screening test results and ultrasound information. Specific algorithms for statistical risk assessment use some combination of biochemical markers and maternal age, the exact gestational age and nuchal translucency (in the first trimester) to estimate the risk.

Risk assessment by PRISCA 4.0 software relies upon maternal age, the exact gestational age, relevant medical history, and the measurement of multiple biochemical markers used for maternal screening. Median values throughout weeks 10 up to 13 (first trimester) and from weeks 14 up to 21 (second trimester) need to be well documented for each of the markers.

PRISCA 4.0 compares the subject's serum result with the median for that gestational age to express the result as a multiple of the median (MoM), for any of the following biochemical tests: PAPP-A and free  $\beta$ -HCG, during the first trimester; AFP, HCG (or free  $\beta$ -HCG) and unconjugated estriol, during the second trimester. Optionally, a MoM can also be calculated for nuchal translucency — an echographic measurement performed in the first trimester — and combined with the biochemical markers.

PRISCA 4.0 is able to provide MoMs corrected for variable factors such as gestational age, maternal weight, ethnicity, smoking, twin pregnancy, diabetes and IVF procedures. Once the corrected MoM

value is obtained, the likelihood ratio is calculated for each of these values, and the combination of all likelihood ratios with maternal age related risk (a priori risk) yields the final risk assessment.

PRISCA 4.0 is a user-friendly program operating under Windows, either on a stand-alone basis or in a network. It features a statistical module to allow updating of medians and the creation of a quality control program for monitoring any of the relevant maternal screening components — for example: MoMs, detection rates, false-positive rates, population parameters, and so on.

### Specimen Collection

For specimen collection requirements and sample storage instructions, please refer to the IMMULITE/IMMULITE 1000 AFP, IMMULITE 2000 AFP or IMMULITE 2500 AFP package insert.

### Warnings and Precautions

For *in vitro* diagnostic use.

### Assay Procedure

See an appropriate Operator's Manual for: preparation, setup, adjustment, assay and quality control procedures.

For the **Double Test in the second trimester**, AFP test (LKAP, L2KAP, L5KAP) should be used in conjunction with HCG test (LKCG, L2KCG, L5KCG).

For the **Triple Test in the second trimester**, AFP test (LKAP, L2KAP, L5KAP) should be used in conjunction with HCG test (LKCG, L2KCG, L5KCG) and with Unconjugated Estriol test (LKEF, L2KEF, L5KEF).

### Expected Values

For AFP median values in maternal serum for IMMULITE/IMMULITE 1000, IMMULITE 2000 and IMMULITE 2500 AFP procedures, please refer to PRISCA Statistical Guidelines for Prenatal Risk Calculation, catalog number ZS1295.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

### Limitations

**Prenatal Testing:** A reliable AFP evaluation for prenatal testing requires precise determination of the gestational age. Underestimation of the gestational age may lead to a false positive determination, while overestimation of gestational age may result in a false negative interpretation. When gestational age is uncertain, confirmation with ultrasonography is indicated.

### Performance Data

See the IMMULITE/ IMMULITE 1000 AFP, IMMULITE 2000 AFP or IMMULITE 2500 AFP package insert for information on conversion factor, calibration range, analytical sensitivity, high dose hook effect, precision, linearity, recovery, specificity, bilirubin and hemolysis.

### References

#### Down's Syndrome (Trisomy 21)

- 1) Wald, et al. Antenatal screening for Down's syndrome. *Health Technol Assess* 1998;2(1).
- 2) Wallace, et al. Calculating correct Down's syndrome risks. *Br J Obstet Gynaecol* 1999;106:371–2.
- 3) Palomaki, et al. Cigarette smoking and levels of maternal serum alpha-fetoprotein, unconjugated estriol, and hCG: impact on down syndrome screening. *Obstet Gynecol* 1993;81(5 Pt 1).
- 4) Henriques, et al. Decreased alpha-fetoprotein in amniotic fluid and maternal serum in diabetic pregnancy. *Obstet Gynecol* 1993;82(6):960–4.
- 5) Spencer, et al. Down's syndrome screening in multiple pregnancies using alpha-fetoprotein and free beta hCG. *Prenatal Diagnosis* 1994;14:537–42.
- 6) de Graaf, et al. Early pregnancy screening for fetal aneuploidy with serum markers and nuchal translucency. *Prenatal Diagnosis* 1999;19(5):458–62.
- 7) Cuckle, et al. Estimating a woman's risk of having a pregnancy associated with Down's syndrome using her age and serum alpha-fetoprotein level. *Br J Obstet Gynaecol* 1987;94:387–402.
- 8) Cuckle, et al. Improved parameters for risk estimation in Down's syndrome screening. *Prenatal Diagnosis* 1995;15:1057–65.
- 9) Snijders, et al. Maternal age- and gestation-specific risk for trisomy 21. *Ultrasound Obstet Gynecol* 1999;13:167–70.
- 10) Norgaard-Pedersen, et al. Maternal serum markers in screening for Down syndrome. *Clin Genet* 1990;37:35–43.
- 11) Wald, et al. Maternal serum screening for Down's syndrome in early pregnancy. *BMJ* 1988;297:883–7.
- 12) Wald, et al. Maternal serum unconjugated oestriol and human chorionic gonadotrophin levels in twin pregnancies: implications for screening for Down's syndrome. *Br J Obstet Gynaecol* 1991;98:905–908.
- 13) Coombes, et al. Prenatal

screening for trisomy 18 with free beta human chorionic gonadotrophin as a marker. *BMJ* 1993;307:1455–8. 14) Neveux, et al. Refinements in managing maternal weight adjustment for interpreting prenatal screening results. *Prenatal Diagnosis* 1996;16:1115–9. 15) Palomaki, et al. Risk-based prenatal screening for trisomy 18 using alpha-fetoprotein, unconjugated oestriol and human chorionic gonadotropin. *Prenatal Diagnosis* 1995;15:713–23. 16) Wald, et al. Serum screening for Down's syndrome between 8 and 14 weeks of pregnancy. *Br J Obstet Gynaecol* 1996;103:407–12. 17) Frishman, et al. Serum-triple-marker screening in in vitro fertilization and naturally conceived pregnancies. *Obstet Gynecol* 1997;90(1):98–100. 18) Spencer, et al. The influence of ethnic origin on first trimester. *Prenatal Diagnosis* 2000;20:491–4. 19) Spencer, et al. The influence of smoking on maternal serum PAPP-A and free beta hCG levels in the first trimester of pregnancy. *Prenatal Diagnosis* 1999;19(11):1065–6. 20) Cuckle, et al. Calculating correct Down's syndrome risks. *British J Obstet Gynecol* 1999;106:371–2. 21) Nicolaides, et al. Correct estimation of parameters for ultrasound nuchal translucency screening. *Prenatal Diagnosis*, 1998;18:511–23.

### Technical Assistance

Available outside the United States only. For technical assistance, contact your National Distributor.

Manufactured by EURO/DPC Ltd. under a Quality System registered to ISO 13485:2003.

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
**EURO/DPC** LTD
Diagnostic Products Corporation
EC REP DPC Biermann GmbH