Free Beta HCG Trisomy 21 Application Sheet

DPC®

Free Beta HCG Trisomy 21 Application Sheet

English

Intended Use:

IMMULITE/IMMULITE 1000 Free Beta HCG (LKFB), IMMULITE 2000 Free Beta HCG (L2KFB) and IMMULITE 2500 Free Beta HCG (L5KFB) may be used for the risk assessment for Down's syndrome (trisomy 21) and Edwards syndrome (trisomy 18) in the first and second trimesters of pregnancy. A pregnant woman's risk, as calculated by the PRISCA 4.0 software, is not confirmatory for these chromosomal abnormalities; rather, it is intended for *in vitro* diagnostic use as additional support for her decision whether or not to undergo a diagnostic procedure.

To calculate the risk assessment during the first and second trimesters of pregnancy, Free Beta HCG (LKFB, L2KFB, L5KFB) should be used in conjunction with PAPP-A (LKPC, L2KPC, L5KPC) – **Double Test for the first trimester**, or in conjunction with AFP (LKAP, L2KAP, L5KAP) – **Double Test for the second trimester**, or in conjunction with AFP (LKAP, L2KAP, L5KAP) and Unconjugated Estriol (LKEF, L2KEF, L5KEF) – **Triple Test for the second trimester**.

The concentration of free β -HCG in a given specimen determined with assays from different manufacturers can vary due to differences in assay methods and reagent specificity. The results reported by the laboratory to the physician must include the identity of the assay used. Values obtained with different Free Beta HCG assays cannot be used interchangeably. Before changing assays, the laboratory must establish a range of median values from non-affected pregnancies for the new assay based on normal sera from pregnant women with confirmed gestational age.

WARNING: This assay has not been FDA approved for the Down's syndrome (trisomy 21) or for the Edwards syndrome (trisomy 18) risk assessment in the US.

Summary and Explanation

Down's Syndrome (Trisomy 21)

Down's syndrome (trisomy 21) is the most common chromosomal abnormality in newborns, with an incidence of approximately 1:800. The only prenatal test diagnostic for chromosomal abnormalities is cytogenetic examination (karyotyping) of fetal cells using invasive procedures such as amniocentesis or chorion villus sampling (CVS).

Bases for the selection of pregnant women to undergo diagnostic procedures of this kind include: advanced maternal age, the subject's medical history, or a statistical risk assessment derived from maternal age plus some combination of biochemical screening test results and ultrasound information. Specific algorithms for statistical risk assessment use some combination of biochemical markers and maternal age, the exact gestational age and nuchal translucency (in the first trimester) to estimate the risk.

Risk assessment by PRISCA 4.0 software relies upon maternal age, the exact gestational age, relevant medical history, and the measurement of multiple biochemical markers used for maternal screening. Median values throughout weeks 10 up to 13 (first trimester) and from weeks 14 up to 21 (second trimester) need to be well documented for each of the markers.

PRISCA 4.0 compares the subject's serum result with the median for that gestational age to express the result as a multiple of the median (MoM), for any of the following biochemical tests: PAPP-A and free β -HCG, during the first trimester; AFP, HCG (or free β -HCG) and unconjugated estriol, during the second trimester. Optionally, a MoM can also be

calculated for nuchal translucency — an echographic measurement performed in the first trimester — and combined with the biochemical markers.

PRISCA 4.0 is able to provide MoMs corrected for variable factors such as gestational age, maternal weight, ethnicity, smoking, twin pregnancy, diabetes and IVF procedures. Once the corrected MoM value is obtained, the likelihood ratio is calculated for each of these values, and the combination of all likelihood ratios with maternal age related risk (a priori risk) yields the final risk assessment.

PRISCA 4.0 is a user-friendly program operating under Windows, either on a stand-alone basis or in a network. It features a statistical module to allow updating of medians and the creation of a quality control program for monitoring any of the relevant maternal screening components — for example: MoMs, detection rates, false-positive rates, population parameters, and so on.

Specimen Collection

For specimen collection requirements and sample storage instructions, please refer to the IMMULITE/IMMULITE 1000 Free Beta HCG, IMMULITE 2000 Free Beta HCG or IMMULITE 2500 Free Beta HCG package insert.

Warnings and Precautions

For *in vitro* diagnostic use.

Assay Procedure

See an appropriate Operator's Manual for: preparation, setup, adjustment, assay and quality control procedures.

For the **Double Test in the first trimester**, Free Beta HCG test (LKFB, L2KFB, L5KFB) should be used in conjunction with PAPP-A test (LKPC, L2KPC, L5KPC).

For the **Double Test in the second trimester**, Free Beta HCG test (LKFB, L2KFB, L5KFB) should be used in conjunction with AFP test (LKAP, L2KAP, L5KAP).

For the **Triple Test in the second trimester**, Free Beta HCG test (LKFB, L2KFB, L5KFB) should be used in conjunction with AFP test (LKAP, L2KAP,

L5KAP) and with Unconjugated Estriol test (LKEF, L2KEF, L5KEF).

Expected Values

For β -HCG median values in maternal serum for IMMULITE/IMMULITE 1000, IMMULITE 2000 and IMMULITE 2500 Free Beta HCG procedures, please refer to PRISCA Statistical Guidelines for Prenatal Risk Calculation, catalog number ZS1295.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Performance Data

See the IMMULITE/IMMULITE 1000 Free Beta HCG, IMMULITE 2000 Free Beta HCG or IMMULITE 2500 package insert for information on conversion factor, calibration range, analytical sensitivity, high dose hook effect, precision, linearity, recovery, specificity, bilirubin and hemolysis.

References

Down's Syndrome (Trisomy 21)

1) Wald, et al. Antenatal screening for Down's syndrome. Health Technol Assess 1998;2(1). 2) Wallace, et al. Calculating correct Down's syndrome risks. Br J Obstet Gynaecol 1999;106:371–2. 3) Palomaki, et al. Cigarette smoking and levels of maternal serum alpha-fetoprotein, unconjugated estriol, and hCG: impact on down syndrome screening. Obstet Gynecol 1993;81(5 Pt 1). 4) Henriques, et al. Decreased alpha-fetoprotein in amniotic fluid and maternal serum in diabetic pregnancy. Obstet Gynecol 1993;82(6):960–4. 5) Spencer, et al. Down's syndrome screening in multiple pregnancies using alpha-fetoprotein and free beta hCG. Prenatal Diagnosis 1994;14:537–42. 6) de Graaf, et al. Early pregnancy screening for fetal aneuploidy with serum markers and nuchal translucency. Prenatal Diagnosis 1999;19(5):458–62. 7) Cuckle, et al. Estimating a woman's risk of having a pregnancy associated with Down's syndrome using her age and serum alpha-fetoprotein level. Br J Obstet Gynaecol 1987;94:387–402. 8) Cuckle, et al. Improved parameters for risk estimation in Down's syndrome screening. Prenatal Diagnosis 1995;15:1057–65. 9) Snijders, et al. Maternal age- and gestation-specific risk for trisomy 21. Ultrasound Obstet Gynecol 1999;13:167–70. 10) Norgaard-Pedersen, et al. Maternal serum markers in screening for Down syndrome. Clin Genet 1990;37:35–43. 11) Wald, et al. Maternal serum screening for Down's syndrome in early

pregnancy. BMJ 1988;297:883–7. 12) Wald, et al. Maternal serum unconjugated oestriol and human chorionic gonadotrophin levels in twin pregnancies: implications for screening for Down's syndrome. Br J Obstet Gynaecol 1991;98:905–908 13) Coombes, et al. Prenatal screening for trisomy 18 with free beta human chorionic gonadotrophin as a marker. BMJ 1993;307:1455–8. 14) Neveux, et al. Refinements in managing maternal weight adjustment for interpreting prenatal screening results. Prenatal Diagnosis 1996;16:1115–9. 15) Palomaki, et al. Risk-based prenatal screening for trisomy 18 using alpha-fetoprotein, unconjugated oestriol and human chorionic gonadotropin. Prenatal Diagnosis 1995;15:713–23. 16) Wald, et al. Serum screening for Down's syndrome between 8 and 14 weeks of pregnancy. Br J Obstet Gynaecol 1996;103:407–12. 17) Frishman, et al. Serum-triple-marker screening in in vitro fertilization and naturally conceived pregnancies. Obstet Gynecol 1997;90(1):98–100. 18) Spencer, et al. The influence of ethnic origin on first trimester. Prenatal Diagnosis 2000;20:491–4. 19) Spencer, et al. The influence of smoking on maternal serum PAPP-A and free beta hCG levels in the first trimester of pregnancy. Prenatal Diagnosis 1999;19(11):1065–6. 20) Cuckle, et al. Calculating correct Down's syndrome risks. British J Obstet Gynecol 1999;106:371–2. 21) Nicolaides, et al. Correct estimation of parameters for ultrasound nuchal translucency screening. Prenatal Diagnosis, 1998;18:511–23.

Technical Assistance

Available outside the United States only.
For technical assistance, contact your
National Distributor.

Manufactured by EURO/DPC Ltd. under a
Quality System registered to ISO 13485:2003.

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
EURO/DPC LTD
Diagnostic Products Corporation
EC REP DPC Biermann GmbH