IMMULITE®

AlaTOP® Allergy Screen

**For use on the IMMULITE®
and IMMULITE® 1000 systems**

DPC®

IMMULITE®/IMMULITE® 1000 AlaTOP® Allergy Screen

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE and IMMULITE 1000 Analyzers — for the qualitative detection of IgE antibodies specific to inhalant allergens in serum, as an aid in the differential diagnosis of atopic allergy.

Catalog Number: **LKAT1** (100 tests)

Test Code: **TOP** Color: **Dark Gray**

CDC Analyte Identifier Code: 2805

CDC Test System Identifier Code: 10159

CLIA Complexity Category: Moderate

Summary and Explanation

In patients with extrinsic asthma, hay fever or atopic eczema, symptoms may appear immediately after exposure to specific allergens. Allergic response of this type (atopic or anaphylactic) is a function of serum reagents, principally immunoglobulin E (IgE). IgE molecules act as points of contact between the allergen and specialized cells that release histamine and other agents upon exposure to the allergen; this initiates the events which we recognize as allergic reactions.^{1,2} When evaluated in the light of other clinical and laboratory findings, *in vitro* allergy screening tests followed by allergen-specific IgE tests can help the physician identify the allergen (or allergens) to which an individual is sensitive.^{5,6,7}

Principle of the Procedure

IMMULITE/IMMULITE 1000 AlaTOP Allergy Screen is a solid-phase, sequential chemiluminescent immunoassay.

Incubation Cycles: 2 × 30 minutes.

IMMULITE/IMMULITE 1000 AlaTOP Allergy Screen is a chemiluminescent enzyme-labeled immunoassay, based on liquid ligand-labeled allergens, monoclonal antibodies, and separation by anti-ligand-coated solid phase. AlaTOP assays use a patented technology (U.S. Patent No. 4,778,751) exploiting liquid-phase kinetics in a bead format.^{3,4} It represents a

significant advance over conventional methods relying on allergens attached to a solid-phase support, such as a paper disk.

AlaTOP is unique in its use of allergens in a liquid format. The allergens are covalently bound to a soluble polymer/co-polymer matrix, which in turn is labeled with a ligand; anti-ligand is coated on the polystyrene bead to capture the ligand-labeled allergens.

Specimen Collection

Samples drawn into certain brands of gel barrier tubes have been shown to yield responses that are altered relative to those obtained with samples drawn into plain serum tubes, in the IMMULITE/IMMULITE 1000 AlaTOP Allergy Screen assay. In addition, these altered responses may be exaggerated with incomplete filling of the gel barrier tube.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Icteric or grossly contaminated samples may give erroneous results.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants.

IMMULITE/IMMULITE 1000 AlaTOP Allergy Screen has not been tested with all possible variations of tube types.

Volume Required: 50 µL serum. (Sample cup must contain at least 100 µL more than the total volume required.)

Storage: 3 days at 2–8°C^{8,11} or 6 months at –20°C.^{9,11}

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

AlaTOP Test Units (LAT1)

Each barcode-labeled unit contains one bead coated with anti-ligand. Stable at 2–8°C until expiration date.

LKAT1: 100 units.

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

AlaTOP Reagent Wedges (LATA, LATB)

With barcodes. **LATA:** 7.5 mL allergen screen. **LATB:** 7.5 mL alkaline phosphatase (bovine calf intestine) conjugated to monoclonal murine anti-IgE antibodies in buffer. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored

as indicated.

LKAT1: 1 set.

AlaTOP Adjustor (LATL)

2 mL IgE-free human serum, with preservative. Stable at 2–8°C for 30 days after opening. *Do not freeze.*

LKAT1: 1 vial.

AlaTOP Allergy Screen Positive Control (LATC2)

One vial containing 2 mL human serum with IgE specific for the IMMULITE AlaTOP Allergy Screen assay, with preservative. Refrigerate upon arrival. Stable at 2–8°C for 4 months after opening. *Do not freeze. Do not use if signs of microbiological contamination, such as a cloudy appearance, have occurred.*

LKAT1: 1 vial.

IMMULITE Multi-Allergen Negative Control (LAPNC)

One vial containing 2 mL IgE-free human serum, with preservative. Refrigerate upon arrival. Stable at 2–8°C for 4 months after opening. *Do not freeze. Do not use if signs of microbiological contamination, such as a cloudy appearance, have occurred.*

LKAT1: 1 vial.

Kit Components Supplied Separately

LSUBX: Chemiluminescent Substrate

LPWS2: Probe Wash Module

LKPM: Probe Cleaning Kit

LCHx-y: Sample Cup Holders (barcoded)

LSCP: Sample Cups (disposable)

LSCC: Sample Cup Caps (optional)

Also Required

Sample transfer pipets, distilled or deionized water.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

Note that both Reagent Wedges, A and B, must be loaded on the carousel to run this assay.

Recommended Adjustment Interval:
4 weeks.

Quality Control Samples: The control(s) supplied with the kit should be used as quality control material to monitor assay performance and the continued applicability of the stored Reactive/Nonreactive Cutoff.

Interpretation of Results

Reactive: A reactive result indicates that antibodies to one or more of the component allergens in the panel are present in the patient sample.

Nonreactive: A nonreactive result indicates nondetectable antibodies to the component allergens.

Limitations

The IMMULITE/IMMULITE 1000 AlaTOP procedure provides qualitative (reactive/nonreactive) results only. It does not provide information on the amounts of IgE specific for any component allergens.

As with all diagnostic screening tests, a definitive clinical diagnosis should not be based on the result of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1998;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in

combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in kCPS. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Specificity: The antibody is highly specific for human IgE with no crossreactivity to other human immunoglobulins.

Bilirubin: Interferes with the assay, causing a depression of values at levels of 100 mg/L and above. (See "Bilirubin" table.)

Hemolysis: Presence of packed red blood cells in concentrations up to 30 μ L/mL has no effect on results, within the precision of the assay.

Clinical Performance

The IMMULITE AlaTOP Allergy Screen procedure was used to study 205 patient specimens. Results were evaluated by clinical assessment of the patients' atopic/nonatopic status as defined by an adaptation of the criteria of Merrett & Merrett⁷, which combine clinical history, total IgE levels, and multiple individual allergen-specific IgE tests for commonly inhaled allergens. The IMMULITE AlaTOP results agreed with the clinical diagnosis of atopy status in 97.6% of the patients, as shown in the table below. The IMMULITE AlaTOP clinical sensitivity and specificity were 99.1% and 95.7%, respectively. Based on an estimated prevalence of 20%⁷ and sensitivity of 99.1% and specificity of 95.7%, the reactive and nonreactive predictive values for this data were 85.2% and 99.8%, respectively.

IMMULITE AlaTOP Allergy Screen		
	Nonreactive	Reactive
Nonatopic	88	4
Atopic	1	112

References

1) Halpern GM. Markers of human allergic disease. J Clin Immunoassay 1983;6:131-8. 2) Wide L, Bennich H, Johansson SGO. Diagnosis of allergy by an *in vitro* test for allergen antibodies. Lancet 1967;2:1105-7. 3) El Shami AS, Alaba O. Liquid-phase *in vitro* allergen-specific IgE assay with *in situ* immobilization. Adv Biosci 1989;74:191-201. 4) Alaba O, El Shami AS. Evaluation of non-specific IgE binding: comparison of two *in vitro* allergen assays. Adv Biosci 1989;74:203-14. 5) Lowenstein H. Cross reactions among pollen antigens. Allergy 1980;35:198-200. 6) Weber RW, Nelson HS. Pollen allergens and their interrelationships. Clin Rev Allergy 1985;3:291-318. 7) Merrett J, Merrett TG. Screening for IgE-mediated allergy. J Clin Immunoassay 1993;16:164-72. 8) Bazaral M, Hamburger RN. Standardization and stability of immunoglobulin E (IgE). J Allergy 1972;49:189-91. 9) Data on file.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782 or 973.927.2828
Fax: 973.927.4101. Outside the United States, contact your National Distributor.
The Quality System of Diagnostic Products Corporation is registered to ISO 13485:2003.

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O Sistema de Qualidade da Diagnostic Products Corporation está registado sob ISO 13485:2003.
EC REP DPC Biemann GmbH