IMMULITE®

Total IgE

For use on the IMMULITE®
and IMMULITE® 1000 systems

DPC®

IMMULITE®/IMMULITE® 1000 Total IgE

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE and IMMULITE 1000 Analyzers — for the quantitative measurement of total immunoglobulin type E (IgE) in serum.

Catalog Number: **LKIE1** (100 tests),
LKIE5 (500 tests)

Test Code: **TIE** Color: **Violet**

CDC Analyte Identifier Code: 2805
CDC Test System Identifier Code: 10159
CLIA Complexity Category: Moderate

Summary and Explanation

Modern laboratory methods for allergy testing are based on discoveries regarding the biology of allergy.^{4,5,7,9,10} Many allergies are mediated by IgE immunoglobulins, which act as points of contact between the allergen and specialized cells. The IgE molecules, which have a molecular mass of approximately 200,000, bind to the surface of mast cells and basophilic granulocytes.^{2,3} Subsequent binding of allergens to cell-bound IgE causes these cells to release histamines and other vasoactive substances, thereby initiating the events which we recognize as an allergic reaction.¹⁰

In deciding on a course of therapy, it is important to distinguish between IgE-mediated and non-IgE-mediated reactions.^{4,7} Measurement of the total circulating IgE level, in conjunction with other supporting diagnostic information, can aid in making this diagnosis. The supporting information should include appropriate tests for allergen-specific IgE. Measurement of the total circulating IgE level may also be of value in the early detection of allergy in infants, and as a means for predicting future atopic manifestations.^{4,6,7}

IgE levels normally show a slow increase during childhood, reaching adult levels in the second decade of life.^{1,8} In general, the total IgE level increases with the number of allergies which a person has and with the amount of exposure to relevant allergens.^{5,7}

Significant elevations may be encountered not only in sensitized individuals, but also in cases of IgE myeloma, pulmonary aspergillosis, and during the active stage of parasitic infestations.^{4,5,7,9}

Principle of the Procedure

IMMULITE/IMMULITE 1000 Total IgE is a solid-phase, chemiluminescent immunometric assay.

Incubation Cycles: 1 × 30 minutes.

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE/IMMULITE 1000 Total IgE has not been tested with all possible variations of tube types.

Volume Required: 5 µL serum. (Sample cup must contain at least 100 µL more than the total volume required.)

Storage: 3 days at 2–8°C^{13,15} or 6 months at –20°C.¹⁵

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source

materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

Total IgE Test Units (LIE1)

Each barcode-labeled unit contains one bead coated with monoclonal mouse anti-IgE antibodies. Stable at 2–8°C until expiration date.

LKIE1: 100 units. **LKIE5:** 500 units.

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

Total IgE Reagent Wedge (LIE2)

With barcode. 7.5 mL alkaline phosphatase (bovine calf intestine) conjugated to polyclonal goat anti-IgE antibody in buffer, with preservative. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

LKIE1: 1 wedge. **LKIE5:** 5 wedges.

Total IgE Adjustors (LIEL, LIEH)

Two vials (Low and High), 2 mL each, of IgE in a nonhuman serum/buffer matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

LKIE1: 1 set. **LKIE5:** 2 sets.

Kit Components Supplied Separately

Total IgE Sample Diluent (LIEZ)

For the manual dilution of patient samples.

One vial 25 mL of an IgE-free nonhuman serum/buffer matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

LSUBX: Chemiluminescent Substrate

LPWS2: Probe Wash Module

LKPM : Probe Cleaning Kit

LCHx-y: Sample Cup Holders (barcoded)

LSCP: Sample Cups (disposable)

LSCC: Sample Cup Caps (optional)

IECM: Total IgE Control Module.

Also Required
Sample transfer pipets, distilled or deionized water, controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

Recommended Adjustment Interval:
2 weeks.

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of Total IgE.

Expected Values

A reference range study was performed using the IMMULITE Total IgE assay on serum samples from nonatopic individuals. The results, tabulated below by age, are in accord with the literature.¹⁴

Age (years)	Median (IU/mL)	95%ile (IU/mL)	<i>n</i>
0–1	6.6	29	15
1–2	10.1	49	34
2–3	12.9	45	29
3–9	14.4	52	106
Adult	20.4	87	48

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

Serum total IgE levels for the majority of individuals with IgE-mediated disease can be expected to be elevated compared to the reference range for healthy adults. However, not all allergic patients exhibit elevated serum total IgE levels.¹¹

Since not all atopic reactions are IgE-mediated, a total IgE result in the reference range should always be interpreted in light of other clinical observations.

Tests for allergen-specific IgE may clarify the finding of an elevated total IgE result.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in IU/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Conversion Factor:

IU/mL × 2.4 → ng/mL

Calibration Range: up to 2,000 IU/mL (2nd IRP 75/502)

Analytical Sensitivity: 1.0 IU/mL

High-Dose Hook Effect: none up to 13,000 IU/mL.

Intraassay Precision (Within-Run):

Statistics were calculated for samples from the results of 20 replicates in a single run. (See "Intraassay Precision" table.)

Interassay Precision (Run-to-Run):

Statistics were calculated for samples assayed in 16 different runs. (See "Interassay Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.) The results show that the IMMULITE Total IgE kit maintains good linearity throughout its calibration range; but due to the heterogeneity of circulating IgE with different affinities, there may be a loss of dilutional parallelism in certain patients.

Recovery: Samples spiked 1 to 19 with three Total IgE solutions (780, 1,940 and 3,950 IU/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for human IgE, exhibiting no crossreactivity to other human immunoglobulin classes.

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of packed red blood cells in concentrations up to 30 µL/mL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Method Comparison: The assay was compared to DPC's Coat-A-Count Total IgE IRMA on 196 samples. (Concentration range: approximately 4 to 775 IU/mL. See graph.) By linear regression:

(IML) = 0.96 (Coat-A-Count IRMA) – 4.8 IU/mL
r = 0.981

Means:

133 IU/mL (IMMULITE)

144 IU/mL (Coat-A-Count IRMA)

References

- 1) Barbee RA, et al. Distribution of IgE in a community population sample: correlation with age, sex and allergen skin test reactivity. J Allergy Clin Immunol 1981;68:106–11.
- 2) Geha RS. Human IgE. J Allergy Clin Immunol 1984;74:109–20.
- 3) Gordon RR, et al. Immunoglobulin E and the eczema-asthma syndrome in early childhood. Lancet 1982;1:72–74.
- 4) Halpern GM. Markers of human allergic disease. J Clin Immunoassay 1983

Jun;6(2):131–39. 5) Hamilton RG, Adkinson NF. Clinical laboratory methods in allergic disease. Lab Management 1983 Dec;21(12):37–50. 6) Hamilton RG, Adkinson NF. Quantitation of allergen-specific IgE in serum using the radioallergosorbent test. J Clin Immunoassay 1983 Jun;6(12):147–54. 7) Homberger HA, Yunginger JW. Laboratory testing in the diagnosis and management of allergic diseases. Clin Lab Annual 1983;2:351–88. 8) Kjellman N-IM, Johansson SGO, Roth A. Serum IgE levels in healthy children. Clin Allergy 1976;6:51–59. 9) Mandy FF, Perelmutter L. Laboratory measurement of total human serum IgE. J Clin Immunoassay 1983 Jun;6(2):140–46. 10) Norman PS. The clinical significance of IgE. Hosp Pract 1975 Aug;10(8):41–9. 11) Yunginger JW. Clinical significance of IgE. In: Middleton E, et al, editors. Allergy principles and practice. Vol 2. 3rd ed. St. Louis: C.V. Mosby, 1988:849-60. 12) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 13) Bazaraal M, Hamburger RN. Standardization and stability of immunoglobulin E (IgE). J Allerg Clin Immunol 1972;49:189-91. 14) Nye L, Merrett TG, Landon J, White RJ. A detailed investigation of circulating IgE levels in a normal population. Clin Allergy 1975;1:13-24. 15) Tietz NW, editor. Clinical guide to laboratory tests. 3rd ed. Philadelphia: WB Saunders, 1995:358.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782 or 973.927.2828
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

The Quality System of Diagnostic Products Corporation is registered to ISO 13485:2003.

Interassay Precision (IU/mL)

	Mean ¹	SD ²	CV ³
1	46	3.5	7.6%
2	126	7.5	6.0%
3	327	20	6.1%
4	712	42	5.9%
5	1,081	130	12.0%
6	1,463	164	11.2%

Linearity (IU/mL)


	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	203	—	—
	4 in 8	112	102	110%
	2 in 8	62	51	122%
	1 in 8	31	25	124%
2	8 in 8	410	—	—
	4 in 8	197	205	96%
	2 in 8	112	103	109%
	1 in 8	62	51	122%
3	8 in 8	1,324	—	—
	4 in 8	648	662	98%
	2 in 8	344	331	104%
	1 in 8	160	166	96%
4	8 in 8	1,644	—	—
	4 in 8	850	822	103%
	2 in 8	377	411	92%
	1 in 8	222	206	108%

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O Sistema de Qualidade da Diagnostic Products Corporation está registado sob ISO 13485:2003.
EC REP DPC Biermann GmbH