IMMULITE®

Folic Acid

For use on the IMMULITE®
and IMMULITE® 1000 systems

DPC®

IMMULITE®/IMMULITE® 1000 Folic Acid

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE and IMMULITE 1000 Analyzers — for the quantitative measurement of folic acid in serum, heparinized plasma or ascorbic acid-treated whole blood, as an aid in clinical diagnosis and treatment of anemia.

Catalog Number: **LKFO1** (100 tests), **LKFO5** (500 tests)

Test Code: **FOL** Color: **Orange**

CDC Analyte Identifier Code: 1907
CDC Test System Identifier Code: 10159
CLIA Complexity Category: High

Summary and Explanation

Folic acid (folate) and vitamin B₁₂ are nutrients essential to hematopoiesis.¹ Megaloblastic anemia is almost always due to lack of one of these two vitamins.¹ Circulating folate levels are usually normal or elevated in vitamin B₁₂ deficiency, but red cell folate levels are frequently low in this condition.²

Folate deficiency is commonly encountered as a result of dietary deficiency (as in alcoholism) or increased demand for this vitamin (as in pregnancy).^{1,3} Unlike vitamin B₁₂, folate is a heat-labile vitamin susceptible to loss by prolonged cooking. Accordingly, the prevalence of folate deficiency exhibits major demographic variations, apparently reflecting differences in dietary and culinary habits.⁴

Circulating folate levels, being strongly influenced by recent intake, are unreliable as an index to tissue stores.^{2,5,6} Thus, folate levels measured in serum or plasma may be normal in the face of folate deficiency. Conversely, circulating levels may be low long before tissue stores have been exhausted.^{2,6} Accordingly, it is important to measure red cell folate levels whenever serum or plasma levels are measured.^{5,7,8}

Principle of the Procedure

Competitive Immunoassay.

IMMULITE Folic Acid is a boil, competitive, liquid-phase, ligand-labeled, protein binding chemiluminescent assay with *in situ* immobilization, and with an anti-ligand detection system. The solid phase, a polystyrene bead enclosed within an IMMULITE Test Unit, is coated with a murine monoclonal antibody specific for folic acid binding protein.

After the sample preparation procedure, the patient sample, ligand-labeled folic acid analog and folic acid binding protein are simultaneously introduced into the Test Unit, and incubated for approximately 30 minutes at 37°C with intermittent agitation. During this time, folic acid in the sample competes with the ligand-labeled folic acid analog for a limited amount of folic acid binding protein, and the folic acid binding protein is captured by the antibody on the bead. (Unbound analog is then removed by a centrifugal wash.)

Incubation Cycles: 2 × 30 minutes.

Specimen Collection

Patient must be in a fasting state.¹⁰

For whole blood and red cell folic acid determinations, use fresh heparinized or EDTA whole blood. The patient's hematocrit must be known in order to translate whole blood into red cell folic acid results, which are expressed as nanograms per milliliter (ng/mL) of packed red blood cells. Moreover, to calculate the red cell folic acid result exactly, the patient's serum folic acid level must also be known.

EDTA plasma should not be used as a sample type.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed serum samples are inappropriate for analysis.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving

anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE/IMMULITE 1000 Folic Acid has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 200 µL of sample — serum, plasma or hemolysate — is required for the sample pretreatment step. A single determination uses 100 µL of the *treated* sample: the sample cup should contain at least 250 µL more of the treated sample than the total volume required for the number of folic acid determinations to be performed.

Storage and Stability:

Serum or Heparinized Plasma: Storage under suitable conditions is critical for reliable folic acid results. If the sample will *not* be assayed within 8 hours, aliquot and freeze at -20°C: stable for 6–8 weeks.¹³ Thaw only once (at ambient temperature). Avoid excessive exposure of samples to direct light.

Warnings and Precautions

The Borate-KCN Buffer Solution contains cyanide. Extreme care must be taken to avoid all bodily contact with this reagent.

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chloramphenicol, at concentrations less than 0.1 g/dL has been added as a

preservative. Chloramphenicol is known to cause cancer; this disclosure is required by the state of California.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

Folic Acid Test Units (LFO1)

Each barcode-labeled unit contains one bead coated with murine monoclonal anti-folic acid binding protein antibody. Stable at 2–8°C until expiration date.

LKFO1: 100 units. **LKFO5:** 500 units.

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

Folic Acid Reagent Wedges (LFOA, LFOB)

With barcodes. **LFOA:** one wedge (7.5 mL) containing folic acid binding protein, with preservative. **LFOB:** one wedge (7.5 mL) containing alkaline phosphatase (bovine calf intestine) conjugated to anti-ligand in buffer, with preservative. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

LKFO1: 1 set. **LKFO5:** 5 sets.

Folic Acid Adjustors (LFOL, LFOH)

Two vials (Low and High) of lyophilized folic acid in a human protein-based matrix, with preservative. Reconstitute each vial with **3.0 mL** distilled or deionized water. Let stand for 30 minutes. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. Stable at 2–8°C for 30 days after reconstitution, or for 6 months (aliquotted) at -20°C.

LKFO1: 1 set. **LKFO5:** 2 sets.

Note: The Adjustors need to be boiled like patient samples (see Sample Pretreatment section).

Ligand-Labeled Folate (LLLF)
5 mL of ligand-labeled folate in buffered human protein-based solution, with preservative. Store refrigerated: stable at 2–8°C for 30 days after opening.
LKFO1: 1 vial **LKFO5:** 5 vials

Borate-KCN Buffer Solution (LBCN)
125 mL of Borate-KCN Buffer Solution, with preservative. Store refrigerated: stable at 2–8°C for 30 days after opening.

Caution! Contains cyanide. *Extreme care must be taken to avoid all bodily contact with this reagent.*
LKFO1: 1 vial **LKFO5:** 5 vials

Dithiothreitol Solution (LDTT)
3 mL of dithiothreitol solution. Store refrigerated: stable at 2–8°C for 30 days after opening.
LKFO1: 1 vial **LKFO5:** 5 vials

Kit Components Supplied Separately

Folic Acid Sample Diluent (LFOZ)
For the manual dilution of patient samples. 25 mL of human protein-based matrix. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

Ascorbic Acid
Required only for whole blood and red cell folic acid determinations. A 1% solution *must be prepared fresh*, on a daily basis. The ascorbic acid crystals should be protected from moisture so that they remain free-flowing, instead of becoming caked from excessive humidity. (Sodium ascorbate is *not* a suitable alternative.)

LSUBX: Chemiluminescent Substrate
LPWS2: Probe Wash Module
LKPM : Probe Cleaning Kit
LCHx-y: Sample Cup Holders (barcoded)
LSCP: Sample Cups (disposable)
LSCC: Sample Cup Caps (optional)

CON6: Tri-level, multi-constituent control

Also Required

Plain 12×75 mm polypropylene tubes, with loose-fitting caps — available from DPC.

Micropipets: 200 and 1,000 µL. For the 1 mL reagent addition, a reliable repeating dispenser of the Nichiryo type is suitable.

Plastic container with lid (for preparing the Working Solution).

Foil — or other means for protecting the hemolysates from light.

Covered boiling waterbath (100°C).

Sample transfer pipets, distilled or deionized water, controls.

Preparation of Hemolysates

This is required only for folic acid determinations on whole blood. Collect the sample; determine and record the hematocrit. Lyse the cells in freshly obtained EDTA or heparinized whole blood by preparing a 1-in-21 dilution of the sample (fully suspended) in freshly prepared 1% ascorbic acid solution; for example, mix 100 µL of sample with 2 mL of this solution. Vortex and let stand for 90 minutes at room temperature (15–28°C) in the dark.¹⁰ Hemolysates may be stored at 2–8°C for up to 3 hours, or frozen at –20°C for 2 weeks.^{10,11} It is important to use ascorbic acid, rather than sodium ascorbate, and to prepare the 1% solution fresh, on a daily basis.⁹

Preparation of Working Solution

Note: *Samples being tested for both vitamin B12 and folic acid should be pretreated with the Working Solution described below.*

The amount of each component depends on the number of tests to be performed. The volumes required, in microliters per test, are tabulated below. Be sure to multiply these volumes by a number slightly greater than the number of tests to be run. (The components are supplied in volumes sufficient for making up the Working Solution in approximately a 20% excess.)

	µL/test
Borate-KCN Buffer Solution	1,000
Ligand-Labeled Folate	20
Dithiothreitol Solution	20

Important: The Working Solution should be prepared on a daily basis. If not used immediately, it should be refrigerated at 2–8°C for a period of not more than 24 hours.

Sample Pretreatment

1. Pipet **200 µL** of each Adjustor, control, or patient sample — serum, plasma or hemolysate — into the tubes prepared.

For patient samples expected to have folic acid levels greater than the upper limit of the assay's calibration range, dilute the serum, plasma, or hemolysate with Folic Acid Sample Diluent.
2. Add **1,000 µL** of the Working Solution to all tubes. Vortex.
3. Loosely cap all tubes and place them in a covered, boiling waterbath (100°C) for **15 – 20 minutes**.
4. Remove the tubes from the boiling waterbath, and cool them in an ambient waterbath for **5 minutes**.
5. Pipet at least **350 µL** of the treated sample to an IMMULITE Sample Cup.

Treated samples, both serum and whole blood, are stable at room temperature (15–28°C) or refrigerated at 2–8°C for 1 hour prior to assay. Note that *commercial controls* may show *variable* stability after treatment.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

Note that both Reagent Wedges A and B must be loaded on the carousel to run this assay.

After Sample Pretreatment, process the *treated* sample according to the usual assay procedure. (Sample cup must contain at least 250 µL more than the total volume required.)

Each Sample cup holder containing pretreated sample can be followed by up to *four* test units.

Recommended Adjustment Interval:
2 weeks.

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of folic Acid.

Calculations

Whole Blood and Red Cell Results

First calculate the result *R* in nanograms per milliliter of hemolysate. Then multiply by the dilution factor 21 to obtain the whole blood folic acid concentration, in ng/mL. For an approximate measure of the packed red cell concentration, again in ng/mL, multiply the whole blood folic acid concentration by $100/H$, where *H* is the hematocrit *in percent*:

$$\text{Red Cell Folic Acid} \approx 21R \times (100 / H)$$

Strictly speaking, the serum folic acid contribution should be subtracted from the whole blood folic acid concentration before multiplying by $100/H$. Using the patient's serum folic acid level *S*, the exact equation is:

$$\text{Red Cell Folic Acid} = 21R - [S \times (100 - H) / 100] \times (100 / H)$$

The term in the square brackets is, however, in most cases so small compared to the term $21R$ that it may be justifiably neglected.

Expected Values

Because folic acid levels are strongly influenced by diet and dietary supplementation, population-based reference limits can show marked demographic differences. In the United States, for example, fortification of enriched grain products, as required by the Food and Drug Administration since 1996, has led to an estimated doubling of the mean plasma folate level among subjects not using vitamin supplements,¹⁴ and a decrease in the prevalence of low folate levels, i.e. levels below 3 ng/mL (7 nmol/L).¹⁴ Fortification programs can be expected to have a similar impact on whole blood and red cell folate levels, which are better measures of tissue stores.

Based on its relationship to DPC's Solid Phase No Boil Dualcount (see Method Comparison), the assay can be expected to have the following reference limits. Note: the Dualcount normal range study was performed in the US prior to the FDA fortification requirements.

Serum Folic Acid	3 – 17 ng/mL (6 – 39 nmol/L)
Whole Blood Folic Acid	43 – 295 ng/mL (97 – 668 nmol/L)
Red Cell Folic Acid	93 – 641 ng/mL (212 – 1,453 nmol/L)

Studies using DPC's Solid Phase No Boil assay were conducted on 100 healthy individuals for serum folic acid and on 88 healthy individuals for red cell folic acid. The groups consisted of primarily laboratory personnel, both men and women, in apparent good health. The study was conducted in the United States.

Each laboratory should establish its own reference ranges.

Limitations

The folate binding protein employed in the Reagent Wedge crossreacts to some extent with the "anti-folate" drug methotrexate (MTX), a compound used in the chemotherapy of cancer. (See Specificity section.) Since high-dose infusions of methotrexate may result in circulating levels of 45,000 ng/mL initially and levels in the order of 2,000 ng/mL 48 hours later, the kit should not be used for folic acid measurements on patients currently receiving the drug. Moreover, great caution must be exercised when interpreting folic acid levels on patients who have recently undergone therapy with methotrexate or any other drug with a similar structure.

Folic acid results near the lower limit of normal should be interpreted with caution.¹⁰

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum

products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/mL. To convert from ng/mL to S.I. units, nanomoles per liter (nmol/L), multiply by 2.266:

Conversion Factor:
ng/mL × 2.266 → nmol/L

Calibration Range: 0.5 – 24 ng/mL
(1.1 – 54 nmol/L).

Analytical Sensitivity: 0.3 ng/mL
(0.7 nmol/L).

Intraassay Precision (Within-Run):
Statistics were calculated for samples from the results of 20 replicates in a single run. (See "Intraassay Precision" table.)

Interassay Precision (Run-to-Run):
Statistics were calculated for samples assayed in 20 different runs. (See "Interassay Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three folic acid solutions (42, 137 and 211 ng/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The binding protein used in the IMMULITE Folic Acid procedure is highly specific for folic acid, with low crossreactivity to substances, other than methotrexate, which may be present in patient samples.

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay

Alternate Sample Type: Samples collected from 19 individuals (into plain and heparinized vacutainer tubes) were

assayed by the IMMULITE Folic Acid procedure.

(Heparin) = 0.98 (Serum) + 0.26 ng/mL
 $r = 0.99$

Means:
6.8 ng/mL (Serum)
7.0 ng/mL (Heparin)

Because the anticoagulant EDTA has a significant effect on results, EDTA plasma should not be used in the assay as an alternative to serum or heparinized plasma.

EDTA whole blood may be used for the determination of whole blood and red cell folic acid.

Method Comparison (Serum Samples):

The IMMULITE Folic Acid procedure was compared to DPC's Dualcount[®] Solid Phase No Boil (SPNB) assay on 71 patient serum samples. (Concentration range: approximately 1.4 to 20.3 ng/mL. See graph.) By linear regression:

(IML) = 1.05 (SPNB) – 0.57 ng/mL
 $r = 0.985$

Means:
5.9 ng/mL (IMMULITE)
6.2 ng/mL (SPNB)

Method Comparison (Hemolysates):

The IMMULITE Folic Acid procedure was compared to DPC's Dualcount[®] Solid Phase No Boil assay on 45 patient hemolysates. (Concentration range: approximately 156 to 973 ng/mL. See graph.) By linear regression:

(IML) = 1.05 (SPNB) – 41 ng/mL
 $r = 0.985$

Means:
345 ng/mL (IMMULITE)
369 ng/mL (SPNB)

References

1) Herbert V. The nutritional anemias. *Hosp Pract* 1980 Mar;15(3):65-83,87-9. 2) Lindenbaum J. Status of laboratory testing in the diagnosis of megaloblastic anemia. *Blood* 1983;61:624-7. 3) Colman N. The radioisotopic investigation of anemia. *Ligand Quarterly* 1981 Fall;4(3):24-30. 4) Grasbeck R. Biochemistry and clinical chemistry of vitamin B12 transport and the related diseases. *Clin Biochem* 1984;17:99-107. 5) Cooper BA. Folic Acid: its metabolism and utilization. *Clin Biochem* 1984;17:95-8. 6) McNeely MDD. Folic Acid assay. In: Kaplan LA, Pesce AJ, editors. *Clinical Chemistry*. St. Louis: CV Mosby, 1984: 1402-6. 7) Allen RH. Clinical role and current status of serum cobalamin (vitamin B12) assays. *Ligand*

Quarterly 1981 Fall;4(3):37-44,67. 8) Chen I-W, et al. B12. In: Kaplan LA, Pesce AJ, editors. *Clinical Chemistry*. St. Louis: CV Mosby, 1984: 1396-1400. 9) Netteland B, Bakke OM. Inadequate sample preparation as a source of error in determination of erythrocyte folate by competitive binding radioassay. *Clin Chem* 1977;23:1505-6. 10) Burtis CA, Ashwood ER, editors. *Tietz textbook of clinical chemistry*. 2nd ed. Philadelphia: WB Saunders, 1994:2056. 11) Kubasik NP, Graham M, Sine HE. Storage and stability of folate and vitamin B-12 in plasma and blood samples. *Clinica Chimica Acta* 1979;95:147-149. 12) National Committee for Clinical Laboratory Standards, Procedures for the collection of diagnostic blood specimens by venipuncture, 4th ed, NCCLS Document H3-A4, 1998. 13) Tietz NW, editor. *Clinical guide to laboratory tests*. 3rd ed. Philadelphia: WB Saunders, 1995:246. 14) Jacques PF, Selhub J, Bostom AG, et al. The effect of folic acid fortification on plasma folate and total homocysteine concentrations. *New Engl J Med* 1999;340:1449-54.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782 or 973.927.2828
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

The Quality System of Diagnostic Products Corporation is registered to ISO 13485:2003.

Tables and Graphs

Intraassay Precision (ng/mL)

	Mean ¹	SD ²	CV ³
1	2.1	0.19	9.0%
2	5.2	0.22	4.2%
3	13	0.67	5.2%

Interassay Precision (ng/mL)

	Mean ¹	SD ²	CV ³
1	1.8	0.16	8.9%
2	5.3	0.21	4.0%
3	13	0.96	7.4%

Linearity (ng/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	12	—	—
	4 in 8	5.8	6.2	94%
	2 in 8	2.9	3.1	94%
	1 in 8	1.5	1.6	94%
2	8 in 8	17	—	—
	4 in 8	8.0	8.6	93%
	2 in 8	3.9	4.3	91%
	1 in 8	2.1	2.1	100%
3	8 in 8	14	—	—
	4 in 8	7.5	7.1	106%
	2 in 8	3.8	3.5	109%
	1 in 8	2.0	1.8	111%

Recovery (ng/mL)


	Spiking Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	3.8	—	—
	A	6.2	5.7	109%
	B	9.6	10	96%
	C	16	14	114%
2	—	6.3	—	—
	A	9.3	8.1	115%
	B	13	13	100%
	C	17	17	100%
3	—	10	—	—
	A	12	12	100%
	B	19	16	119%
	C	23	20	115%

Specificity


Compound ¹	ng/mL Added ²	Apparent ng/mL ³	% Cross-reactivity ⁴
Methotrexate	100	3.4	3.4%

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O Sistema de Qualidade da Diagnostic Products Corporation está registado sob ISO 13485:2003.
EC REP DPC Biermann GmbH