IMMULITE[®]

Homocysteine

For use on the IMMULITE[®]
and IMMULITE[®] 1000 systems

DPC[®]

IMMULITE®/IMMULITE® 1000 Homocysteine

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE and IMMULITE 1000 Analyzers — for the quantitative determination of L-homocysteine in human plasma or serum. This device can assist in the diagnosis and treatment of patients suspected of having hyperhomocysteinemia or homocystinuria.

Catalog Number: **LKHO1** (100 tests)

Test Code: **HCY** Color: **Dark Gray**

Warning: Specimens from patients who are on drug therapy involving S-adenosyl-methionine may show falsely elevated levels of homocysteine.

Results on specimens obtained from patients taking methotrexate, carbamazepine, phenytoin, nitrous oxide, anti-convulsants, or 6-azauridine triacetate should be interpreted with caution as these substances interfere with homocysteine determination.

Summary and Explanation

Total homocysteine (tHcy) has emerged as an important risk factor in the assessment of cardiovascular disease.^{3-7,12} Hcy, a thiol-containing amino acid, is produced by the intracellular demethylation of methionine. Hcy thereby serves as a pool that can be later scavenged for use in the remanufacture of either methionine through the action of the folate-dependent enzyme methionine synthase or cysteine using the B6 dependent transsulphuration pathway.^{1,2} Hcy in plasma is found primarily in a protein bound form but free, oxidized and disulfide forms are also present.

Highly elevated levels of tHcy are found in patients with homocystinuria, a rare genetic disorder of the enzymes involved in Hcy metabolism^{1,2,3} Homocystinuria patients exhibit arterial thromboembolism, mental retardation and early

arteriosclerosis.¹ Less severe genetic defects are also associated with moderate levels of Hcy.^{3,4,5}

Homocysteine has been identified as a indicator of cardiovascular disease. A meta-analysis of 27 epidemiological studies has suggested that a 5 µmol/L increase in tHcy could be associated with an odds ratio for coronary artery disease (CAD) of 1.6 for men and 1.8 for women, the same increase in risk as a 0.5 mmol/L increase in cholesterol.⁷ In addition, patients with chronic renal disease complicated by arteriosclerotic cardiovascular disease show elevated tHcy due to the inability of the kidney to remove Hcy from the blood.^{1,2,3}

Principle of the Procedure

Competitive Immunoassay.

IMMULITE/IMMULITE 1000

Homocysteine involves a preliminary manual sample pretreatment step. Homocysteine in the patient plasma or serum sample is released from its binding proteins and converted to S-adenosyl-homocysteine (SAH) by an off-line 30-minute incubation at 37°C in the presence of S-adenosyl-L-homocysteine hydrolase and dithiothreitol (DTT).

The treated sample and alkaline phosphatase-labeled anti-SAH antibody are simultaneously introduced into a test unit containing an SAH-coated polystyrene bead. During a 30-minute incubation, the converted SAH from the patient sample competes with the immobilized SAH for binding the alkaline phosphatase labeled-anti-SAH antibody conjugate. Unbound enzyme conjugate is removed by centrifugal wash. Substrate is added, and the procedure continues as described for typical immunoassays in the Operator's Manual.

Incubation Cycle: 1 × 30 minutes.

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt

by the laboratory; hence the results should be interpreted with caution.

Heparinized and EDTA plasma are the samples of choice, but serum is also suitable for use. It is important to separate the plasma or serum from the cells as soon as possible after collection, as synthesis of HCY will take place in red blood cells after sampling. **Samples should be stored on ice between the time of sampling and centrifugation.** Note that storage on ice makes the use of serum samples particularly difficult.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants.

IMMULITE/IMMULITE 1000 Homocysteine has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required – untreated sample:
15 µL untreated plasma or serum.

Volume Required – treated sample:
75 µL treated sample. (Sample cup must contain at least 100 µL more than the total volume required.)

Storage: 14 days at 2–8°C, or 6 months at –20°C.

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

Homocysteine Test Units (LHO1)

Each barcode-labeled unit contains one bead coated with S-adenosyl-L-homocysteine (SAH). Stable at 2–8°C until expiration date.

LKHO1: 100 units.

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

Homocysteine Reagent Wedge (LHO2)

With barcode. 7.5 mL alkaline phosphatase (bovine calf intestine) conjugated to murine monoclonal anti-SAH in buffer. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

LKHO1: 1 wedge.

Pretreatment A Solution (LHOA)

20 mL bovine S-adenosyl-L-homocysteine hydrolase in buffer, with preservative. Stable at 2–8°C until expiration date.

Recommended usage is within 30 days after opening when stored as indicated.

LKHO1: 1 vial.

Pretreatment B Solution (LHOB)

2 mL dithiothreitol (DTT) in buffer. Stable at 2–8°C until expiration date.

Recommended usage is within 30 days after opening when stored as indicated.

LKHO1: 1 vial.

Pretreatment B Solution must be diluted 1-in-10 with distilled or deionized water before combining it with Pretreatment A

Solution. See section on Preparation of Sample Pretreatment Working Solution.

Homocysteine Adjustors (LHOL, LHOH)

Two amber vials (Low and High), 2 mL each, of synthetically-derived S-adenosyl-L-homocysteine (SAH) in a protein/buffer matrix. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

LKHO1: 1 set.

Kit Components Supplied Separately

Homocysteine Sample Diluent (LHOZ)

For the manual dilution of patient samples. 25 mL of homocysteine-free protein/buffer matrix. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

LSUBX: Chemiluminescent Substrate

LPWS2: Probe Wash Module

LKPM : Probe Cleaning Kit

LCHx-y: Sample Cup Holders (barcoded)

LSCP: Sample Cups (disposable)

LSCC: Sample Cup Caps (optional)

CCCM: A bi-level, nonhuman serum-based Cardiac Marker Control Module, containing homocysteine as one of four different constituents.

Also Required

Sample transfer pipets, distilled or deionized water, controls;

Variable micropipets: 10–500 µL.

Preparation of Sample Pretreatment Working Solution

1. Dilute sufficient volume of Pretreatment B Solution 1-in-10 with distilled or deionized water before combining it with Pretreatment A Solution in Step 2.
2. Combine equal volumes of Pretreatment A Solution and diluted Pretreatment B Solution in sufficient quantities for the number of samples to be processed.

For example, see chart below:

No. of Samples	Wkg Sol req. mL	H ₂ O mL	Pre-B µL	Pre-A mL	Total mL
5	1.5	0.90	100	1	2
10	3	1.8	200	2	4
15	4.5	2.25	250	2.5	5
20	6.0	3.15	350	3.5	7
25	7.5	3.6	400	4	8

Use Sample Pretreatment Working Solution within two hours of preparation.

Sample Pretreatment

Adjustors, Controls and patient samples must go through the Sample Pretreatment Procedure prior to analysis.

For Controls and Patient Samples

1. Label one test tube for each sample to be processed.
2. Add 300 µL of Sample Pretreatment Working Solution into the tubes prepared.
3. Add 15 µL of patient untreated plasma or serum to each tube.
4. Cover the tubes and mix well.
5. Incubate for 30 minutes at 37°C in a water bath or oven.
6. Transfer all of the treated sample into an IMMULITE Sample Cup.

For Adjustors

Because Adjustors are processed in replicates of 4, additional volume is required.

At Step 2, use 400 µL of Sample Pretreatment Working Solution.

At Step 3, use 20 µL of each of the Adjustor levels. Proceed with remaining pretreatment steps.

Treated samples, both serum and plasma, are stable at room temperature (15–28°C) or refrigerated at 2–8°C for 1 hour prior to assay.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the

IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

Recommended Adjustment Interval:
4 weeks.

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of homocysteine.

Expected Values

Homocysteine levels can vary with age, gender, geographical area and genetic factors, therefore it is important for laboratories to establish their own reference ranges based on their local populations. Literature suggests a reference range of 5–15 $\mu\text{mol/L}$ for adult males and females,¹¹ but also notes that men tend to have higher levels than women, and that postmenopausal women tend to have higher levels than premenopausal.¹²

Based on its relationship to IMMULITE 2000 Homocysteine procedure (see Method Comparison), the assay can be expected to have essentially the same reference ranges.

One hundred and twenty samples from apparently healthy adult male and female volunteers age 22–66 were analyzed using the IMMULITE 2000 Homocysteine procedure. The samples were drawn in heparinized plasma tubes and kept on ice prior to separation of the plasma from the cells. The median value was 7.7 $\mu\text{mol/L}$, with a central 95% range of

5.0 – 12 $\mu\text{mol/L}$.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

Specimens from patients who are on drug therapy involving S-adenosyl-methionine may show falsely elevated levels of homocysteine.

Although analysis of the parent compounds of carbamazepine, phenytoin, 6-azauridine, and anthopterin indicate no cross-reactivity, specimens obtained from patients treated with these drugs, as well as with methotrexate, nitrous oxide, and other anticonvulsants, should be interpreted with caution as these substances have been shown to interfere in some homocysteine assays.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in $\mu\text{mol/L}$. (Unless otherwise noted, all were generated on heparinized plasma samples.)

Calibration Range: 2–50 $\mu\text{mol/L}$

Analytical Sensitivity: 0.5 $\mu\text{mol/L}$

Precision: Samples were assayed in quadruplicate over the course of several days, for a total of 20 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three homocysteine solutions (105, 207 and 411 $\mu\text{mol/L}$) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for homocysteine. (See "Specificity" table.)

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations

up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 512 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To determine whether serum can be used in the IMMULITE Homocysteine procedure, blood was collected on ice from 20 laboratory volunteers into plain, heparinized and EDTA vacutainer tubes. Samples were separated from the cells, and matched samples were spiked with homocysteine and then assayed by the IMMULITE Homocysteine procedure, with the following results.

(EDTA) = 1.08 (Heparin) – 1.10 $\mu\text{mol/L}$
 $r = 0.983$

(Serum) = 1.11 (Heparin) – 1.79 $\mu\text{mol/L}$
 $r = 0.992$

Means:

9.18 $\mu\text{mol/L}$ (Heparin)
8.82 $\mu\text{mol/L}$ (EDTA)
8.38 $\mu\text{mol/L}$ (Serum)

To assess the effect of storage temperature, heparinized, EDTA and plain vacutainer tubes were collected from five volunteers for each tube type. Some tubes were stored at room temperature, while other tubes were kept on ice for various time periods *prior to separation*. The graphs below show the effect of storage time and temperature for heparin, EDTA and serum. (See graphs 1-3).

Method Comparison: The assay was compared to DPC's IMMULITE 2000 Homocysteine on 257 plasma samples (Concentration range: approximately 4 to 50 $\mu\text{mol/L}$, as measured by IMMULITE 2000. See graph.) By linear regression:

(IML) = 0.97 (IML 2000) + 0.86 $\mu\text{mol/L}$
 $r = 0.975$

Means:

13.5 $\mu\text{mol/L}$ (IMMULITE)
13.1 $\mu\text{mol/L}$ (IMMULITE 2000)

References

1) Selhub J, Miller JW. The pathogenesis of homocysteinemia: interruption of the coordinate regulation by s-adenosylmethionine of the

remethylation and transsulfuration of homocysteine. Am J Clin Nutr. 1992;55:131-8. 2) Jacques PF, Bostom AG, Williams RR, et al. Relation between folate status, a common mutation in methylenetetrahydrofolate reductase and plasma homocysteine concentrations. Circulation 1996;93:7-9. 3) Malinow MR. Plasma homocysteine and arterial occlusive diseases: A mini review. Clin Chem. 1995;40:173-76. 4) Clarke R, Daly L, Robinson K, Naughten E, et al. Hyperhomocysteinemia: An independent risk factor for vascular disease. New Eng J Med 1991;324:1149-55. 5) Deloughery TG, Evans A, Sadeghi A, et al. Common mutation in methylenetetrahydrofolate reductase: Correlation with homocysteine metabolism and late-onset vascular disease. Circulation 1996;94:3074-78. 6) Bostom AG, Selhub J. Homocysteine and arteriosclerosis: Subclinical and clinical disease associations. Circulation 1999;99:2361-63. 7) Boushey CJ, Beresford SAA, et al. A quantitative assessment of plasma homocysteine as a risk factor for vascular disease. JAMA 1995;274(13):1049-57. 8) Guttormsen AB, Svarstad E, et al. Elimination of homocysteine in subjects with end-stage renal failure. Irish J Med Sci 1995;164:8. 9) Bostom AG, Lathrop, L. Hyper-homocysteinemia in end-stage renal disease (ESDR): Prevalence, etiology and potential relationship to arteriosclerotic outcomes. Kidney Int 1997;52:10-20. 10) Refsum H, Ueland PM. Clinical significance of pharmacological modulation of homocysteine metabolism. TIPS 1990;11:411-16. 11) Ueland PM, Refsum H, Stabler SP, et al. Total homocysteine in plasma or serum: Methods and clinical applications. Clin Chem 1993;39:1764-79. 12) Nehler MR, Taylor LM, Porter JM. Homocysteine as a risk factor for atherosclerosis: A Review. Cardiovasc Pathol 1997;6:1-9.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782 or 973.927.2828
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

The Quality System of Diagnostic Products Corporation is registered to ISO 13485:2003.

Tables and Graphs

Precision (µmol/L)

	Mean	Within-Run		Total	
		SD	CV	SD	CV
1	11.05	1.01	9.1%	1.09	9.9%
2	19.06	1.37	7.2%	1.82	9.5%
3	25.33	1.79	7.1%	1.93	7.6%
4	27.49	1.85	6.7%	3.00	10.9%
5	35.99	2.25	6.3%	3.48	9.7%

Specificity

Compound ¹	µmol/L Added ²	% Cross reactivity ³
Adenosine	5,000	ND
S-adenosyl-L-methionine	500	0.6%
Cystathionine	500	5.7%
L-Cysteine	100,000	ND
Gluthathione	100,000	ND

ND: Not detectable.⁴

Linearity (µmol/L)


	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8	7.28	—	—
	4 in 8	3.92	3.64	108%
	2 in 8	1.63	1.82	90%
	1 in 8	0.73	0.91	80%
2	8 in 8	8.68	—	—
	4 in 8	4.13	4.34	95%
	2 in 8	1.87	2.17	86%
	1 in 8	0.94	1.09	86%
3	8 in 8	28.33	—	—
	4 in 8	12.68	14.17	89%
	2 in 8	7.03	7.08	99%
	1 in 8	2.81	3.54	79%
4	8 in 8	32.39	—	—
	4 in 8	16.88	16.20	104%
	2 in 8	8.17	8.10	101%
	1 in 8	3.94	4.05	97%

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
5	8 in 8	36.91	—	—
	4 in 8	17.50	18.46	95%
	2 in 8	8.62	9.23	93%
	1 in 8	4.18	4.61	91%

Recovery (µmol/L)

	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	6.09	—	—
	A	12.65	11.04	115%
	B	17.39	16.14	108%
	C	30.02	26.34	114%
2	—	9.38	—	—
	A	15.75	14.16	111%
	B	20.88	19.26	108%
	C	32.87	29.46	112%
3	—	20.88	—	—
	A	26.89	25.09	107%
	B	30.09	30.19	100%
	C	40.77	40.39	101%
4	—	30.25	—	—
	A	31.51	33.99	93%
	B	40.02	39.09	102%
	C	49.77	44.29	101%
5	—	34.60	—	—
	A	38.32	38.12	101%
	B	42.65	43.22	99%
	C	57.27	53.42	107%

Method Comparison


$$(IML) = 0.97 (IML 2000) + 0.86 \mu\text{mol/L}$$


$$r = 0.975$$

Effect of Storage Temperature on Different Sample Types¹

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
Homocysteine: Homoisteina. **Effect of Storage
Temperature on Different Sample Types:**

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O Sistema de Qualidade da Diagnostic Products Corporation está registado sob ISO 13485:2003.
EC REP DPC Biermann GmbH