# IMMULITE®

## Albumin

For use on the IMMULITE®  
and IMMULITE® 1000 systems

**DPC®**

# IMMULITE/IMMULITE 1000<sup>®</sup> Albumin

---

## English

---

**Intended Use:** For *in vitro* diagnostic use with the IMMULITE and IMMULITE 1000 Analyzers — for the detection and quantitation of microalbuminuria, that is, concentrations of albumin in human urine not detectable by conventional dipstick methods. It is intended for use as an aid in monitoring the development of incipient diabetic nephropathy.

Catalog Number: **LKHA1** (100 tests),  
**LKHA5** (500 tests)

Test Code: **ALB** Color: **Brown**

CDC Test System Identifier Code: 10159  
CLIA Complexity Category: Moderate

### Summary and Explanation

Microalbuminuria is a condition characterized by increased urinary excretion of albumin in the absence of overt nephropathy.<sup>9,10,21</sup> Conventional dipstick and acid precipitation tests for detecting protein in urine lack the sensitivity required to delineate this condition. Dip-sticks may yield negative or trace results even when the albumin excretion rate is 10 or 20 times normal; and the rate must increase to 200 or 300 micrograms per minute ( $\mu\text{g}/\text{min}$ ) before nephropathy becomes clinically apparent as persistent proteinuria.<sup>14,21</sup> Published studies based on highly sensitive immunoassays for albumin have established the upper limit of normal for adults as approximately 12 or 15  $\mu\text{g}/\text{min}$ ; accordingly, investigators have typically fixed upon a range extending from 15 or 30  $\mu\text{g}/\text{min}$  to about 150  $\mu\text{g}/\text{min}$  as an operational definition of microalbuminuria.<sup>11,12,20</sup>

Measurement of urinary albumin has been used in conjunction with other assays to help localize the site of renal impairment. Thus, elevated urinary albumin levels are to be expected when the glomerular filtration process is overactive; impaired tubular resorption, on the other hand, leads to increased urinary levels of smaller molecules like  $\beta_2$ -microglobulin.<sup>16</sup>

A variety of factors other than glomerular integrity — notably exercise — may have an influence on reference ranges for urinary albumin. Some reports indicate that ranges based on overnight specimens are lower than ranges based on daytime or 24-hour specimens.<sup>2,6,17</sup> It is tempting to explain this in terms of differences in posture and activity which the collection times entail. But other studies have demonstrated that exercise, within limits, has no significant effect on the albumin excretion rate in healthy adults.<sup>1,7,15</sup>

This does not hold true for individuals with diabetes or essential hypertension.<sup>1,8,13,15,18,24</sup> Although newly diagnosed diabetics may have baseline levels within the reference range, their urinary albumin levels are likely to show an increased responsiveness to exercise, a responsiveness which can be eliminated in the short term by careful metabolic control.<sup>4,5,19,25</sup>

Interest in measuring subclinical elevations in the albumin excretion rate has focused on individuals with an already established diagnosis of diabetes or essential hypertension.<sup>3</sup> Providing proper care is taken to minimize the influence of exercise and poor metabolic control on the albumin excretion rate, the urinary albumin level has proved to be an excellent predictor of progression to overt nephropathy in both insulin-dependent and non-insulin-dependent diabetes.<sup>11,12,14,20,23</sup> Thus, according to the literature, it is very unlikely that a diabetic with a resting level of less than 15  $\mu\text{g}/\text{min}$  will progress within a decade to overt clinical nephropathy; whereas diabetics who have resting albumin levels greater than 15  $\mu\text{g}/\text{min}$  even while in good metabolic control can be expected to show yearly increases of (on average) 25  $\mu\text{g}/\text{min}$ .<sup>12</sup>

### Principle of the Procedure

IMMULITE/IMMULITE 1000 Albumin is a solid-phase, competitive chemiluminescent enzyme immunoassay.

**Incubation Cycles:** 1 × 30 minutes.

## Specimen Collection

Collect a timed 3-hour, overnight or 24-hour urine sample — without preservative — recording the time, duration and total volume of the collection. Mix well and retain an aliquot for analysis. Samples should not be acidified.

The urine should be cleared by filtration or centrifugation before use. (Bloody specimens are unsuitable for use.) Remove the supernatant urine from the sediment and mix by *gentle* swirling or inversion. Use disposable-tip pipets.

**Volume Required:** 50  $\mu$ L urine. (Sample cup must contain at least 100  $\mu$ L more than the total volume required.)

**Storage:** 14 days at 2–8°C<sup>28</sup>  
Do *not* freeze.<sup>27,28</sup>

Before assay, allow the samples to come to room temperature (15–28°C). Samples which are cloudy or have particulate material should be clarified by low-speed centrifugation.

Consult the Summary and Explanation section and the References for information regarding the effect of exercise and other factors on urinary albumin determinations.

## Warnings and Precautions

For *in vitro* diagnostic use.

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight. (See insert.)

**Water:** Use distilled or deionized water.

## Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

### Albumin Test Units (LHA1)

Each barcode-labeled unit contains one bead coated with monoclonal murine anti-HSA (human serum albumin). Stable at 2–8°C until expiration date.

**LKHA1:** 100 units. **LKHA5:** 500 units.

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

### Albumin Reagent Wedge (LHA2)

With barcode. 7.5 mL containing alkaline phosphatase (bovine calf intestine) conjugated to HSA (human serum albumin) in buffer, with preservative. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

**LKHA1:** 1 wedge. **LKHA5:** 5 wedges.

### Albumin Adjustors (LHAL, LHAH)

Two vials (Low and High) of lyophilized albumin in a protein/buffer matrix, with preservative. At least 30 minutes before use: Reconstitute each vial with **2.0 mL** distilled or deionized water. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. Stable at 2–8°C for 30 days after reconstitution, or for 6 months (aliquotted) at –20°C.

**LKHA1:** 1 set. **LKHA5:** 2 sets.

## Kit Components Supplied Separately

### Albumin Sample Diluent (LHAZ)

For the manual dilution of patient samples. One vial containing 25 mL of a protein-based matrix, free of human albumin, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

**LSUBX:** Chemiluminescent Substrate

**LPWS2:** Probe Wash Module

**LKPM :** Probe Cleaning Kit

**LCHx-y:** Sample Cup Holders (barcoded)

**LSCP:** Sample Cups (disposable)

**LSCC:** Sample Cup Caps (optional)

**LHACM:** Bi-level, protein-based albumin control.

Also Required  
Sample transfer pipets, distilled or deionized water, controls.

### Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

**Recommended Adjustment Interval:**  
2 weeks.

**Quality Control Samples:** Use controls or urine pools with at least two levels (low and high) of albumin.

### Expected Values

Results for albumin are usually reported as excretion rates, in micrograms of albumin per minute ( $\mu\text{g}/\text{min}$ ). This is obtained by multiplying the concentration in  $\mu\text{g}/\text{mL}$  by the total volume of the urine collection (in milliliters) and then dividing by the duration of the collection (in minutes).

Based on its relationship to DPC's Double Antibody Albumin kit (see Method Comparison), the IMMULITE Albumin assay can be expected to have approximately the same reference range limits for urine samples, expressed as micrograms of albumin excreted per minute ( $\mu\text{g}/\text{min}$ ):

nondetectable to 18  $\mu\text{g}/\text{min}$ ,

with a median of approximately 6  $\mu\text{g}/\text{min}$ . A total of 55 timed, 3-hour morning urine samples were collected from healthy laboratory volunteers at two separate laboratories and assayed by the Double Antibody Albumin kit. 95% had excretion rates of less than 13  $\mu\text{g}/\text{min}$ . 24-hour as well as 3-hour specimens were collected from the 20 volunteers at one of the two laboratories. 90% of the results fell below

11.5 and 12.5  $\mu\text{g}/\text{min}$  for the 24-hour and 3-hour specimens, respectively.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

### Limitations

There are other potential causes of microalbuminuria besides incipient diabetic nephropathy. Consideration must be given to exercise, poor diabetic control and nondiabetic renal or systemic diseases, including hypertension. Urinary tract infections and congestive heart disease are also possible causes of subclinical elevations in the albumin excretion rate.

Samples with a pH less than 4 or greater than 8 may yield results which are respectively too high or too low. Acidified samples are accordingly unsuitable for use.

The test should not be performed if the sample exhibits significant bacterial growth or if the patient shows signs of a urinary tract infection.

Bloody specimens are unsuitable for use, even if cleared by centrifugation, since blood is a likely sign of contamination — e.g. from menstrual flow — and albumin circulates in plasma at levels approximately 4,000 times those normally encountered in urine.

### Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in  $\mu\text{g}/\text{mL}$ . (Unless otherwise noted, all results were generated on urine samples.)

**Calibration Range:** 2.5 – 60  $\mu\text{g}/\text{mL}$

**Analytical Sensitivity:** 1.0  $\mu\text{g}/\text{mL}$

**Precision:** Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

**Linearity:** Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

**Recovery:** Samples spiked 1 to 19 with three albumin solutions (200, 300 and 400  $\mu\text{g}/\text{mL}$ ) were assayed. (See "Recovery" table for representative data.)

**Specificity:** The antibody is highly specific for human albumin. (See "Specificity" table.)

**Method Comparison:** The assay was compared to DPC's Double Antibody Albumin on 124 urine samples (Concentration range: approximately 5 to 60 µg/mL. See graph.) By linear regression:

(IML) = 0.92 (DAb) + 1.9 µg/mL  
r = 0.972

Means:

24.7 µg/mL (IMMULITE)  
24.9 µg/mL (Double Antibody)

## References

1) Dahlquist G, et al. Effect of metabolic control and duration on exercise-induced albuminuria in diabetic teenagers. *Acta Paediatr Scand* 1983;72:895-902. 2) Davies AG, et al. Urinary albumin excretion in school children. *Arch Dis Child* 1984;59:625-30. 3) Hostetter TH. Diabetic nephropathy. *New Engl J Med* 1985;312:642-4. 4) Kroc Collaborative Study Group. Blood glucose control and the evolution of diabetic retinopathy and albuminuria; a preliminary multicenter trial. *New Engl J Med* 1984;311:365-72. 5) Microalbuminuria and renal hyperfiltration. In: Gambino R, editor. *Lab Report for Physicians* 1984; 6(8): 57-61. 6) Mogensen CE. Urinary albumin excretion in early and long-term juvenile diabetes. *Scand J Clin Lab Invest* 1971;28:183-93. 7) Mogensen CE, Vittinghus E, Sølling K. Abnormal albumin excretion after two provocative renal tests in diabetes: physical exercise and lysine injection. *Kidney Int* 1979;16:385-93. 8) Mogensen CE. Antihypertensive treatment inhibiting the progression of diabetic nephropathy. *Acta Endocrinol* 1980;94 (suppl. 238):103-11. 9) Mogensen CE, et al. Early changes in kidney function, blood pressure and the stages in diabetic nephropathy. In: Keen H, Legrain M, editors. *Prevention and treatment of diabetic nephropathy*. Boston: MTP Press, 1983: 57-83, 99-105. 10) Mogensen CE, Christensen CK, Vittinghus E. The stages in diabetic renal disease, with emphasis on the stage of incipient diabetic nephropathy. *Diabetes* 1983;32 (Suppl. 2):64-78. 11) Mogensen CE. Microalbuminuria predicts clinical proteinuria and early mortality in maturity-onset diabetes. *New Engl J Med* 1984;310:356-60. 12) Mogensen CE, Christensen CK. Predicting diabetic nephropathy in insulin-dependent patients. *New Engl J Med* 1984;311:89-93. 13) Mohamed A, et al. Response of urinary albumin to submaximal exercise in newly diagnosed non-insulin dependent diabetes. *Br Med J* 1984;288:1342-3. 14) Parving H-H, et al. Early detection of patients at risk of developing diabetic nephropathy; a longitudinal study of urinary albumin excretion. *Acta Endocrinol*

1982;100:550-5. 15) Pedersen EB, Mogensen CE, Larsen JS. Effects of exercise on urinary excretion of albumin and β<sub>2</sub>-microglobulin in young patients with mild essential hypertension without treatment and during long-term propranolol treatment. *Scand J Clin Lab Invest* 1981;41:493-8. 16) Ratge D, Wisser H. Urinary protein profiling by high performance gel permeation chromatography. *J Chromatogr* 1982;230:47-56. 17) Rowe DJF, et al. Effect of glycaemic control and duration of disease on overnight albumin excretion in diabetic children. *Br Med J* 1984;289:957-9. 18) Viberti GC, et al. Increased glomerular permeability to albumin induced by exercise in diabetic subjects. *Diabetologia* 1978;14:293-300. 19) Viberti GC, et al. Effect of control of blood glucose on urinary excretion of albumin and β<sub>2</sub>-microglobulin in insulin-dependent diabetes. *New Engl J Med* 1979;300:638-41. 20) Viberti GC, et al. Microalbuminuria as a predictor of clinical nephropathy in insulin-dependent diabetes mellitus. *Lancet* 1982;ii(26 Jun):1430-2. See also *ibid* 1982;ii(7 Aug):330, 1982;ii(30 Oct):985, 1982;ii(18 Dec):1407 and 1983;ii(12 Feb):352. 21) Viberti GC, et al. Proteinuria in diabetes mellitus: role of spontaneous and experimental variation of glycemia. *Kidney Int* 1982;21:714-20. 22) Viberti GC, Mogensen CE, et al. Urinary excretion of albumin in normal man: the effect on water loading. *Scand J Clin Lab Invest* 1982;42:147-51. 23) Viberti GC, et al. Long term correction of hyperglycemia and progression of renal failure in insulin dependent diabetes. *Br Med J* 1983;286:598-602. 24) Vittinghus E, Mogensen CE. Graded exercise and protein excretion in diabetic man and the effect of insulin treatment. *Kidney Int* 1982;21:725-9. 25) Wiseman MJ, et al. Effect of blood glucose control on increased glomerular filtration rate and kidney size in insulin-dependent diabetes. *New Engl J Med* 1985;312:617-21. 26) Rosenstock J, Raskin P. Early diabetic nephropathy: assessment and potential therapeutic interventions. *Diabetes Care* 1986;9:529-45. 27) Erving LD, et al. Screening of microalbuminuria in patients with diabetes mellitus: frozen storage of urine samples decreases their albumin content. *Clin Chem* 1989;35(2):308-10. 28) Osberg I, et al. Effects of storage time and temperature on measurement of small concentrations of albumin in urine. *Clin Chem* 1990;36(8):1428-30.

## Technical Assistance

In the United States, contact DPC's Technical Services department.  
Tel: 800.372.1782 or 973.927.2828  
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

Manufactured by EURO/DPC Ltd. under a Quality System registered to ISO 13485:2003.

## Tables and Graphs

### Precision (µg/mL)

	Mean <sup>3</sup>	Within-Run <sup>1</sup>		Total <sup>2</sup>	
		SD <sup>4</sup>	CV <sup>5</sup>	SD	CV
1	4.2	0.16	3.8%	0.24	5.7%
2	9.1	0.31	3.4%	0.44	4.8%
3	26.9	0.91	3.4%	1.10	4.1%
4	35.5	1.31	3.7%	2.39	6.7%
5	51.4	2.28	4.4%	3.20	6.2%

### Linearity (µg/mL)

	Dilution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	4 in 4 <sup>5</sup>	9.5	—	—
	2 in 4	5.0	4.8	104%
	1 in 4	2.5	2.4	104%
2	4 in 4	15	—	—
	2 in 4	6.9	7.5	92%
	1 in 4	3.8	3.8	100%
3	8 in 8	22	—	—
	4 in 8	11	11	100%
	2 in 8	5.6	5.5	102%
	1 in 8	2.7	2.8	96%
4	8 in 8	36	—	—
	4 in 8	19	18	106%
	2 in 8	8.7	9.0	97%
	1 in 8	4.0	4.5	89%
5	8 in 8	41	—	—
	4 in 8	20	21	95%
	2 in 8	8.5	10	85%
	1 in 8	4.7	5.1	92%

### Recovery (µg/mL)


	Spiking Solution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	—	5.6	—	—
	A	15.4	15.6	99%
	B	20.2	20.6	98%
	C	23.2	25.6	91%
2	—	8.9	—	—
	A	17.9	18.9	95%
	B	21.4	23.9	90%
	C	27.0	28.9	93%
3	—	9.5	—	—
	A	19.0	19.5	97%
	B	24.1	24.5	98%
	C	28.7	29.5	97%
4	—	18.9	—	—
	A	27.9	28.9	97%
	B	30.6	33.9	90%
	C	36.7	38.9	94%
5	—	23.7	—	—
	A	34.4	33.7	102%
	B	39.1	38.7	101%
	C	39.3	43.7	90%
6	—	25.6	—	—
	A	32.3	35.6	91%
	B	37.7	40.6	93%
	C	40.8	45.6	90%

### Specificity

Compound <sup>1</sup>	µg/mL Added <sup>2</sup>	% Cross-reactivity <sup>3</sup>
Alpha-fetoprotein	20	ND
Urea	625	ND
Human IgG	500	ND
Creatinine	2,500	ND
Uroporphyrin ester	0.10	ND
Bovine serum albumin	10,000	ND
Transferrin	100	ND
Bilirubin	400	ND
Gentamycin sulfate	120	ND
Glucose	5,000	ND

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
**EURO/DPC** LTD
Diagnostic Products Corporation
EC REP DPC Biermann GmbH