IMMULITE®

Insulin

For use on the IMMULITE®
and IMMULITE® 1000 systems

DPC®

IMMULITE®/IMMULITE 1000 Insulin

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE and IMMULITE 1000 Analyzers — for the quantitative measurement of insulin in, serum or heparinized plasma, for the management of diabetes.

Catalog Number: **LKIN1** (100 tests), **LKIN5** (500 tests)

Test Code: **INS** Color: **Orange**

Summary and Explanation

Human insulin is a polypeptide hormone originating in the beta cells of the pancreas and serving as a principal regulator for the storage and production of carbohydrates. Its secretion is normally stimulated by increases in the amount of glucose in circulation. This leads to higher insulin levels and more rapid tissue-assimilation of glucose — followed by a decline in the insulin level as the glucose level subsides.

In a number of conditions, notably insulinoma and diabetes, this relationship is impaired. Insulin tends to circulate at inappropriately high levels in patients with insulin-secreting pancreatic tumors; such tumors can thus be a cause of hypoglycemia. Accordingly, insulin immunoassays — used sometimes in connection with provocative doses of tolbutamide or calcium — play an essential role in the identification (and localization) of insulinomas. The finding of fasting hypoglycemia in association with an *inappropriately high* serum insulin concentration is considered diagnostic.

Insulin levels do not figure in the subclassification of diabetes worked out by the National Diabetes Data Group. Nevertheless, when obtained in the course of a glucose tolerance test, they appear to be of some prognostic value in predicting the benefits of insulin therapy and the likelihood of progression to insulin-dependence and the complications (such as retinopathy) characteristic of diabetes.

The application of insulin immunoassay to patients already undergoing insulin therapy is complicated by the fact that such therapy typically leads to the formation of anti-insulin antibodies capable of interfering with the assay. Some investigators have sought therefore to measure insulin in urine, or in serum samples subjected to column chromatography or PEG precipitation. But the measurement of "free" insulin remains of limited interest as a technique for monitoring insulin therapy in the absence of statistics establishing therapeutic or toxic ranges. So far, it appears that glucose control in diabetics cannot in general be achieved by normalizing the insulin profile. Nor is it known at what point abnormally high insulin levels become dangerous.

Principle of the Procedure

IMMULITE/IMMULITE 1000 Insulin is a solid-phase, two-site chemiluminescent immunometric assay.

Incubation Cycles: 1 × 60 minutes.

Specimen Collection

EDTA tubes should *not* be used in the IMMULITE/IMMULITE 1000 Insulin procedure.

In collecting samples for insulin determinations, it is important to avoid hemolysis, which can lead to spuriously low values.^{22,23}

Icteric, lipemic, or grossly contaminated samples may give erroneous results.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants.

IMMULITE/IMMULITE 1000 Insulin has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 100 μ L serum or heparinized plasma. (Sample cup must contain at least 250 μ L more than the total volume required.)

Storage: 7 days at 2–8°C, or 3 months at –20°C.²¹

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

Insulin Test Units (LIN1)

Each barcode-labeled unit contains one bead coated with monoclonal murine anti-insulin. Stable at 2–8°C until expiration date.

LKIN1: 100 units. **LKIN5:** 500 units.

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

Insulin Reagent Wedge (LIN2)

With barcode. 7.5 mL alkaline phosphatase (bovine calf intestine) conjugated to polyclonal sheep anti-insulin antibody and alkaline phosphatase (bovine calf intestine) conjugated to monoclonal murine anti-insulin antibody in buffer, with preservative. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

LKIN1: 2 wedges. **LKIN5:** 10 wedges.

Insulin Adjustors (LINL, LINH)

Two vials (Low and High) of lyophilized insulin in a nonhuman serum matrix, with preservative. At least 30 minutes before use, reconstitute each vial by adding **4.0 mL** distilled or deionized water. Mix by *gentle* swirling or inversion. Aliquot and freeze: stable at –20°C for 60 days after reconstitution.

LKIN1: 1 set. **LKIN5:** 2 sets.

Insulin Controls (LINC1, LINC2)

Two vials of lyophilized insulin in a nonhuman serum matrix, with preservative. Reconstitute each vial by adding **4.0 mL** distilled or deionized water. Mix by *gentle* swirling or inversion. Aliquot and freeze: stable at –20°C for 60 days after reconstitution. Refer to Insulin Control insert for the respective concentrations in μ IU/mL.

LKIN1: 1 set. **LKIN5:** 2 sets.

Kit Components Supplied Separately

Insulin Sample Diluent (LINZ)

For the manual dilution of samples. One vial with 25 mL of ready-to-use, insulin-free nonhuman serum matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

LSUBX: Chemiluminescent Substrate

LPWS2: Probe Wash Module

LKPM : Probe Cleaning Kit

LCHx-y: Sample Cup Holders (barcoded)

LSCP: Sample Cups (disposable)

LSCC: Sample Cup Caps (optional)

Also Required

Sample transfer pipets, distilled or deionized water.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for: preparation, setup, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

Sample cup must contain at least 250 μL more than the total volume required. Each Sample cup holder can be followed by up to four test units.

Recommended Adjustment Interval:
4 weeks.

Quality Control Samples: Use Insulin Controls supplied with the kit.

Expected Values

A study performed on 83 fasting, apparently healthy laboratory volunteers yielded a median of 8.9 $\mu\text{IU/mL}$ and a lower 95% range of up to 28.4 $\mu\text{IU/mL}$.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

Circulating anti-insulin antibodies are often found in patients who have been treated with nonhuman forms of insulin. If present, these antibodies may interfere with the assay.

For individuals who are significantly overweight, fasting insulin levels are typically somewhat higher than for adults of normal weight.

EDTA plasma has an effect on the measurement of insulin in the IMMULITE/IMMULITE 1000 Insulin procedure.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum

products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in $\mu\text{IU/mL}$. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Calibration Range: Up to 300 $\mu\text{IU/mL}$ (WHO NIBSC 1st IRP 66/304).

Analytical Sensitivity: 2 $\mu\text{IU/mL}$.

High-dose Hook Effect:
None up to 70,000 $\mu\text{IU/mL}$.

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three insulin solutions (280, 960 and 2,000 $\mu\text{IU/mL}$) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is specific for insulin. (See "Specificity" table.)

Bilirubin: Severe icterus (bilirubin up to 100 or 200 mg/L) may cause a depression of values. (See "Bilirubin" table.)

Hemolysis: May cause a depression of values. (See "Hemolysis" table.)

Lipemia: Presence of triglycerides may cause a depression of values. (See "Lipemia" table.)

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 18 volunteers into plain glass and plastic serum tubes, plastic heparinized and EDTA tubes, sodium fluoride and plastic gel barrier tubes (SST[®]). All tubes were from Becton

Dickinson. Eight samples were spiked with insulin to provide values throughout the calibration range of the assay. All samples were assayed by the IMMULITE Insulin procedure, with the following results. By linear regression:

(Serum Glass) = 1.06 (Serum Plastic) – 0.6 µIU/mL
r = 0.997

(Na Heparin Plastic) = 1.06 (Serum Plastic) + 0.98 µIU/mL
r = 0.995

(SST Plastic) = 0.996 (Serum Plastic) + 1.2 µIU/mL
r = 0.996

(NaFI Plastic) = 0.892 (Serum Plastic) – 8.0 µIU/mL
r = 0.994

(EDTA Plastic) = 0.31 (Serum Plastic) – 1.4 µIU/mL
r = 0.997

Means:

87.5 µIU/mL (Serum Plastic)
92.4 µIU/mL (Serum Glass)
93.5 µIU/mL (Na Heparin Plastic)
88.3 µIU/mL (SST Plastic)
70.0 µIU/mL (NaFI Plastic)
25.7 µIU/mL (EDTA Plastic)

EDTA plasma is not suitable for use.
Sodium fluoride causes a depression in values.

Method Comparison: The assay was compared to DPC's IMMULITE 2000 Insulin on 270 samples. (Concentration range: approximately 2 to 291 µIU/mL. See graph.) By linear regression:

(IML) = 1.04 (IML 2000) + 0.2 µIU/mL
r = 0.989

Means:

23.6 µIU/mL (IMMULITE)
22.5 µIU/mL (IMMULITE 2000)

References

1) Bates HM. Insulinoma and pheochromocytoma. *Lab Management* 1983 (May);21:11–12, 15. 2) Bennett PH. The diagnosis of diabetes; new international classification and diagnostic criteria. *Ann Rev Med* 1983;34:295–309. 3) Berntorp K et al. Relation between plasma insulin and blood glucose in a cross-sectional population study of the oral glucose tolerance test. *Acta Endocrinol* 1983;102:549–56. 4) Boehm TM, Lebovitz HE. Statistical analysis of glucose and insulin responses to intravenous tolbutamide: evaluation of hypoglycemic and hyperinsulinemic states. *Diabetes Care* 1979;2:479–90. 5) Fajans SS, Floyd JC. Diagnosis and medical management of insulinomas. *Ann Rev Med* 1979;30:313–29. 6) Hirata Y et al. Autoimmunity in diabetes mellitus. *Gunma Symposia Endocrinol* 1980;17:81–91. 7) Kahn CR, Rosenthal AS. Immunologic reactions

to insulin: insulin allergy, insulin resistance and the autoimmune insulin syndrome. *Diabetes Care* 1979;2:283–95. 8) Kaplan EL, Rubinstein AH et al. Calcium infusion: a new provocative test for insulinomas. *Annals Surgery* 1979;190:501–7. 9) Malone JI, Root AW. Plasma free insulin concentrations: keystone to effective management of diabetes mellitus in children. *J Pediatr* 1981;99:862–7. See also *ibid* 1983;102:800–2. 10) Marschner I et al. Group experiments on the radioimmunological insulin determination. *Horm Metabol Res* 1974;6:293–6. 11) National Diabetes Data Group. Classification and diagnosis of diabetes mellitus and other categories of glucose intolerance. *Diabetes* 1979;28:1039–57. 12) Pfeifer MA et al. Insulin secretion in diabetes mellitus. *Am J Med* 1981;70:579–88. 13) Reeves WG. Insulin antibody determination: theoretical and practical considerations. *Diabetologia* 1983;24:399–403. 14) Sachse G, Willms B. Klinische Bedeutung der Seruminsulinbestimmung: Diagnose und Prognose bei drohendem Sekundärversagen der Sulfonylharnstofftherapie. *Medizin Klinik* 1979;74:1635–40. 15) Schade DS, Eaton RP, Spencer W. Normalization of plasma insulin profiles in diabetic subjects with programmed insulin delivery. *Diabetes Care* 1980;3:9–14. 16) Service FJ et al. Insulinoma: clinical and diagnostic features of 60 consecutive cases. *Mayo Clinic Proceed* 1976;51:417–29. 17) Turkington RW, Estkowski A, Link M. Secretion of insulin or connecting peptide; a predictor of insulin dependence of obese diabetics. *Arch Intern Med* 1982;142:1102–5. 18) Turkington RW, Weindling HK. Insulin secretion in the diagnosis of adult-onset diabetes mellitus. *JAMA* 1978;240:833–36. See also *ibid* 1979;241:1462–3. 19) Yalow R, Bauman WA. Plasma insulin in health and disease. In: Ellenberg M, Rifkin H, editors. *Diabetes mellitus; theory and practice*. New York: Excerpta Medica, 1983: 119–50. 20) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 21) Burtis CA, Ashwood ER, editors. *Tietz textbook of clinical chemistry*. 2nd ed. Philadelphia: WB Saunders, 1994:943–4. 22) Sapin R, Ongagna JC, Gasser F, Grucker D. Insulin measurements in haemolysed serum; influence of insulinase inhibitors. *Clin Chim Acta* 1998;274:111–7. 23) Chevenne D, Letailleur A, Trivin F, Porquet D. Effect of hemolysis on the concentration of insulin in serum determined by RIA and IRMA. *Clin Chem* 1998;44:354–6.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782 or 973.927.2828
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

Tables and Graphs

Precision (µIU/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	7.39	0.47	6.4%	0.59	8.0%
2	12.3	0.65	5.3%	0.74	6.0%
3	17.8	1.08	6.1%	1.26	7.1%
4	25.5	1.46	5.7%	1.50	5.9%
5	102	5.26	5.2%	6.20	6.1%
6	300	15.8	5.3%	21.0	7.0%

Linearity (µIU/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8	118	—	—
	4 in 8	51.7	59.0	88%
	2 in 8	27.8	29.5	94%
	1 in 8	14.8	14.8	100%
2	8 in 8	119	—	—
	4 in 8	61.5	59.5	103%
	2 in 8	34.7	29.8	116%
	1 in 8	17.6	14.9	118%
3	8 in 8	120	—	—
	4 in 8	63.8	60.0	106%
	2 in 8	33.5	30.0	112%
	1 in 8	17.3	15.0	115%
4	8 in 8	125	—	—
	4 in 8	57.0	62.5	91%
	2 in 8	30.8	31.3	98%
	1 in 8	15.6	15.6	100%
5	8 in 8	128	—	—
	4 in 8	54.2	64.0	85%
	2 in 8	28.7	32.0	90%
	1 in 8	14.7	16.0	92%

Recovery (µIU/mL)

	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	13.2	—	—
	A	28.5	26.5	108%
	B	59.7	60.5	99%
	C	96.0	113	85%
2	—	16.9	—	—
	A	29.5	30.1	98%
	B	58.6	64.1	91%
	C	98.0	116	84%
3	—	20.0	—	—
	A	33.1	33.0	100%
	B	63.7	67.0	95%
	C	108	119	91%
4	—	26.2	—	—
	A	38.0	38.9	98%
	B	66.6	72.9	91%
	C	108	125	86%
5	—	32.1	—	—
	A	43.6	44.5	98%
	B	79.4	78.5	101%
	C	121	130	93%
6	—	43.6	—	—
	A	51.7	55.4	93%
	B	84.8	89.4	95%
	C	125	141	89%

Specificity

Compound ¹	ng/mL Added ²	%Cross-reactivity ³
C-Peptide	1,000	ND
Glucagon	1,000	ND
Proinsulin	10	8.5%

ND: not detectable.⁴

Bilirubin

	Unspiked ¹	Conjugated ²		Unconjugated ³	
		100 mg/L	200 mg/L	100 mg/L	200 mg/L
1	6.90	5.25	5.30	5.35	4.55
2	15.5	12.2	12.7	11.5	11.0
3	20.6	16.1	15.3	15.1	14.4
4	70.8	58.6	57.4	54.0	50.7
5	285	189	208	202	193


Hemolysis

	Unspiked	Hemoglobin			
		157 mg/dL	314 mg/dL	550 mg/dL	550 mg/dL
1	7.15	6.30	6.25	5.50	
2	15.1	14.9	14.2	12.7	
3	20.4	19.0	17.7	14.4	
4	71.1	66.0	61.3	52.9	
5	279	264	240	201	

Lipemia

	Triglycerides Added mg/dL	Observed	Expected	%O/E
1	1,000	5.0	5.7	88%
	2,000	4.4	5.4	81%
	3,000	4.0	5.1	78%
2	—	17.0	—	—
	1,000	15.1	16.2	93%
	2,000	12.8	15.3	84%
3	—	24.2	—	—
	1,000	22.8	23.0	99%
	2,000	21.6	21.8	99%
4	—	96.0	—	—
	1,000	92.7	91.2	102%
	2,000	89.6	86.4	104%
5	—	287	—	—
	1,000	256	273	94%
	2,000	259	258	100%
6	—	287	—	—
	1,000	256	273	94%
	2,000	259	258	100%
7	—	287	—	—
	1,000	256	273	94%
	2,000	259	258	100%

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
EURO/DPC LTD
Diagnostic Products Corporation
EC REP DPC Biermann GmbH