IMMULITE®

C-Peptide

**For use on the IMMULITE®
and IMMULITE® 1000 systems**

DPC®

IMMULITE®/IMMULITE® 1000 C-Peptide

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE and IMMULITE 1000 Analyzers — for the quantitative measurement of C-peptide in serum, heparinized plasma, or urine, as an aid in the diagnosis and treatment of patients with abnormal insulin secretion.

Catalog Number: **LKPE1** (100 tests),
LKPE5 (500 tests)

Test Code: **CPE** Color: **Dark Gray**

CDC Analyte Identifier Code: 1040
CDC Test System Identifier Code: 10159
CLIA Complexity Category: Moderate

Summary and Explanation

Human C-peptide is a 31 amino acid chain with a molecular mass of approximately 3,020 daltons. Metabolically inert, it originates in the pancreatic β -cells as a by-product of the enzymatic cleavage of proinsulin to insulin.^{1,2,5} In this process, insulin and C-peptide are split from the prohormone and secreted into the portal circulation in equimolar concentrations.^{4,5,7} It is this fact which underlies the clinical interest in plasma determinations of C-peptide.

Within limits, C-peptide levels can serve as a valuable index to insulin secretion. Thus, low C-peptide levels are to be expected where insulin secretion is diminished, as in insulin-dependent diabetes, or suppressed, as a normal response to exogenous insulin; whereas elevated C-peptide levels may result from the increased β -cell activity observed in insulinomas.^{3,4,6,9}

Accordingly, in the differential diagnosis of hypoglycemia, C-peptide determinations can be used to supplement insulin measurements as an index to pancreatic activity in the classic 72-hour fasting test, and as the sole indicator of pancreatic activity where insulin itself is administered to check for suppressibility.^{1,8} In addition, covert self-administration of insulin can be virtually ruled out as the cause of hyperinsulinemia by the finding of an elevated C-peptide level.^{2,3,8,9}

Circulating anti-insulin antibodies are commonly encountered in patients who have undergone insulin therapy. These would typically interfere with immunoassays for insulin, making it impossible to use insulin measurements in this context to check on residual β -cell activity, even if treatment were temporarily suspended. C-peptide measurements have therefore been used as an alternative in this context, to yield information on the natural history of insulin-dependent diabetes, to indirectly monitor insulin secretion in the presence of anti-insulin antibodies, and to help settle on an appropriate course of treatment.^{3,6,7,10}

C-peptide has also been measured as an additional means for evaluating glucose tolerance and glibenclamide-glucose tests.^{2,3,10}

Principle of the Procedure

IMMULITE/IMMULITE 1000 C-Peptide is a solid-phase, competitive chemiluminescent enzyme immunoassay.

Incubation Cycles: 1 × 60 minutes.

Specimen Collection

Serum and Heparinized Plasma:

The patient should be fasting. Collect blood by venipuncture,¹² avoiding hemolysis, into plain tubes (without anticoagulant) or heparinized tubes, noting the time of collection, and separate the serum or plasma from the cells.

EDTA plasma and sodium fluoride plasma are unsuitable for use.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Icteric or grossly contaminated samples may give erroneous results.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin,

ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE/IMMULITE 1000 C-Peptide has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Storage: Assay within 2–3 hours, or store frozen at –20°C for 1 week.¹¹

Urine – Collection and Storage:

Collect a 24-hour urine, without preservative, keeping the specimen refrigerated at 2–8°C during collection. Record the total volume of the collection and retain a well-mixed aliquot for analysis. Before assay, clear the sample by centrifugation or filtration.

Storage: For longer storage, aliquot and freeze: stable at –20°C for 30 days.

Predilution of Urine Samples: All urine samples should be diluted *at least* 1-in-4 with C-Peptide Sample Diluent before assay. (Urine samples *should not* be assayed without prior dilution.) A 1-in-20 dilution will bring normal urine samples within the assay's calibration range.

Volume Required:

50 µL serum, plasma, or urine (prediluted). (Sample cup must contain at least 100 µL more than the total volume required.)

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

C-Peptide Test Units (LPE1)

Each barcode-labeled unit contains one bead coated with polyclonal rabbit anti-C-Peptide. Stable at 2–8°C until expiration date.

LKPE1: 100 units. **LKPE5:** 500 units.

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

C-Peptide Reagent Wedge (LPE2)

With barcode. 7.5 mL alkaline phosphatase (bovine calf intestine) conjugated to C-peptide in buffer, with preservative. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

LKPE1: 1 wedge. **LKPE5:** 5 wedges.

C-Peptide Adjustors (LPEL, LPEH)

Two vials (Low and High) of lyophilized C-peptide in buffered human albumin, with preservative. Reconstitute each vial with **4.0 mL** distilled or deionized water. Let stand for 30 minutes. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. After reconstitution, aliquot and freeze. Stable at –20°C for 6 months. Discard aliquots after use.

LKPE1: 1 set. **LKPE5:** 2 sets.

Kit Components Supplied Separately

C-Peptide Sample Diluent (LPEZ)

For the manual dilution of urine samples and high serum/plasma samples. One vial containing 25 mL of processed, C-peptide-free buffered human albumin, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

LSUBX: Chemiluminescent Substrate

LPWS2: Probe Wash Module

LKPM: Probe Cleaning Kit

LCHx-y: Sample Cup Holders (barcoded)

LSCP: Sample Cups (disposable)

LSCC: Sample Cup Caps (optional)

PECM: Tri-level C-Peptide Control Module

Also Required

Sample transfer pipets, distilled or deionized water, controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

Predilution of Urine Samples: All urine samples should be diluted *at least* 1-in-4 with C-Peptide Sample Diluent before assay.

Recommended Adjustment Interval: 2 weeks.

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of C-Peptide.

Calculation of Results

The results for *urine* samples, as returned by the IMMULITE Analyzer, must be corrected for dilution by multiplying by the dilution factor. Results are usually reported as micrograms of C-peptide per 24 hours ($\mu\text{g/day}$). This can be obtained by multiplying the result in ng/mL ($= \mu\text{g/L}$)

by the total volume of the 24-hour collection *in liters*.

Expected Values

Serum and Heparinized Plasma:

Based upon the agreement with DPC's IMMULITE 2000 C-Peptide assay, the IMMULITE C-Peptide procedure can be expected to yield a median of 2.1 ng/mL (695 pmol/L) and a nonparametric central 95% range of:

1.1 – 5 ng/mL (364 – 1,655 pmol/L)

for fasting individuals. The reference range study was performed on serum samples from 71 fasting laboratory volunteers.

Urine:

Based upon the agreement with DPC's Double Antibody C-Peptide assay, the IMMULITE C-Peptide procedure can be expected to yield a mean \pm SD of $79 \pm 56 \mu\text{g/day}$, with a range of 2 to 260 $\mu\text{g/day}$, representing the central 95% of the observations.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

Because metabolism of C-peptide differs from that of insulin, C-peptide levels are at best a semi-quantitative index of insulin secretion. The half-life of C-peptide in plasma has been estimated as approximately 30 minutes, compared to approximately 5 minutes for insulin. Because of the difference in half-life, C-peptide circulates in plasma at a level roughly five times that of insulin, even though the two molecules are secreted in an equimolar ratio. Again, the liver plays a major role in clearing insulin, whereas C-peptide is removed by degradation and elimination mainly through the kidneys. Hepatic and renal complications will therefore affect the circulating C-peptide/insulin ratio.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all

immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/mL. (Unless otherwise specified, all results were generated on serum samples collected in tubes without anticoagulants, gel barriers, or clot-promoting additives.)

Conversion Factor:
ng/mL × 331 → pmol/L

Calibration Range: 0.5 – 7 ng/mL
(166 to 2,317 pmol/L)
WHO 1st IRP 84/510.

Analytical Sensitivity: 0.3 ng/mL
(99 pmol/L) for serum and plasma samples.

Precision: Samples were repeatedly assayed in quadruplicate over the course of several days, for a total of 20 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with two C-peptide solutions (18 and 72 ng/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for C-peptide. (See "Specificity" table.)

Bilirubin: Interferes with the assay at levels of 100 mg/L and 200 mg/L, causing elevation of values. (See "Bilirubin" table)

Hemolysis: Presence of packed red blood cells in concentrations up to 30 µL/mL has no effect on results, within the precision of the assay.

Alternate Sample Types: No significant effect for heparinized plasma, Becton Dickinson SST[®] vacutainer tubes. EDTA

plasma and sodium fluoride plasma are unsuitable for use. By linear regression:

(Heparin) = 1.03 (Serum) – 0.04 ng/mL
r = 0.99

(SST) = 1.16 (Plain Tubes) + 0.81 ng/mL
r = 0.98

Means:
4.4 ng/mL (Serum)
4.5 ng/mL (Heparin)
5.9 ng/mL (SST Tubes)

Method Comparison–Serum: The assay was compared to DPC's IMMULITE 2000 C-Peptide on 66 serum samples. (Concentration range: approximately 0.5 to 6.5 ng/mL. See graph.) By linear regression:

(IML) = 1.01 (IML 2000) – 0.1 ng/mL
r = 0.990

Means:
2.8 ng/mL (IMMULITE)
2.8 ng/mL (IMMULITE 2000)

Method Comparison–Urine: The assay was compared to DPC's Double Antibody C-Peptide on 24 urine samples. (Concentration range: up to approximately 6.7 ng/mL. See "Method Comparison: Urine" graph.) By linear regression:

(IMMULITE) = 0.89 (DA) – 0.01 ng/mL
r = 0.983

Means:
2.51 ng/mL (IMMULITE)
2.81 ng/mL (Double Antibody)

References

- 1) Beischer W. Proinsulin and C-Peptide in humans. In: Fotherby K, Pal S, editors. Hormones in Normal and Abnormal Human Tissues. Vol 3. Berlin: Walter DeGruyter, 1983: 1-43. 2) Beyer J, Krause U, Cordes U. C-Peptide: its biogenesis, structure, determination and clinical significance. Giornale Ital Chem Clin 1979;4 Suppl 1:9-22. 3) Bonser A, Garcia-Webb P. C-Peptide measurement: methods and clinical utility. CRC Crit Rev Clin Lab Sci 1984;19:297-352. 4) Gonen G, Rubenstein AH, Horwitz DL, Blix PM. Clinical significance of C-Peptide. In: Baba S, Kaneko T, Yanaihara N, editors. Proinsulin, Insulin, C-Peptide. Amsterdam: Excerpta Medica, 1979: 246-53. 5) Horwitz D, et al. Proinsulin, insulin and C-Peptide concentrations in human portal and peripheral blood. J Clin Invest 1975;55:1278-83. 6) Horwitz D, Kuzuya H, Rubenstein AH. Circulating serum C-Peptide. N Engl J Med 1976;295:207-9. 7) Polonsky KS, Rubenstein AH. C-Peptide as a measure of the secretion and hepatic extraction of insulin; pitfalls and limitations. Diabetes 1984;33:486-93. 8) Rendell M. Expanding clinical use of

C-Peptide radioimmunoassay. Acta Diabetol Lat 1983;20:105-13. 9) Rubenstein AH, Kuzuya H, Horwitz DL. Clinical significance of circulating C-Peptide in diabetes mellitus and hypoglycemic disorders. Arch Intern Med 1977;137:625-32. 10) Turkington RW, Estkowski A, Link M. Secretion of insulin or connecting peptide; a predictor of insulin dependence of obese diabetics. Arch Intern Med 1982;142:1102-5. 11) Data on file. 12) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782 or 973.927.2828
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

Manufactured by EURO/DPC Ltd. under a Quality System registered to ISO 13485:2003.

Tables and Graphs

Precision (ng/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	0.97	0.1	10.3%	0.14	14.4%
2	2.1	0.16	7.6%	0.22	10.5%
3	3.5	0.25	7.1%	0.29	8.3%
4	4.6	0.3	6.5%	0.36	7.8%
5	5.7	0.39	6.8%	0.43	7.5%

Linearity (ng/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	4 in 4 ⁵	3.4	—	—
	2 in 4	1.5	1.7	88%
	1 in 4	0.77	0.85	91%
2	8 in 8	3.8	—	—
	4 in 8	2.1	1.9	111%
	2 in 8	1.0	0.95	105%
	1 in 8	0.52	0.48	108%

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
3	8 in 8	5.3	—	—
	4 in 8	2.4	2.7	89%
	2 in 8	1.4	1.3	108%
	1 in 8	0.57	0.66	86%

Recovery (ng/mL)

	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	0.27	—	—
	A	1.1	1.2	92%
	B	3.7	3.9	95%
2	—	0.60	—	—
	A	1.6	1.5	107%
	B	4.6	4.2	110%
3	—	0.83	—	—
	A	1.6	1.7	94%
	B	4.4	4.4	100%

Specificity


Compound ¹	Amount Added ²	Apparent ng/mL ³	% Cross-reactivity ⁴
Insulin	200 µIU/mL	ND	ND
Glucagon	15,000 ng/mL	ND	ND
Secretin	15,000 ng/mL	ND	ND
Proinsulin	10.0 ng/mL	1.25	13%

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
Method Comparison: C-Peptide: C-Peptide

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
EURO/DPC LTD
Diagnostic Products Corporation
EC REP DPC Biermann GmbH