IMMULITE®

Growth Hormone (hGH)

For use on the IMMULITE®
and IMMULITE® 1000 systems

DPC®

IMMULITE®/IMMULITE 1000® Growth Hormone (hGH)

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE and IMMULITE 1000 Analyzers — for the quantitative measurement of human growth hormone (hGH) in serum, as an aid in patient management.

Catalog Number: **LKGH1** (100 tests), **LKGH5** (500 tests)

Test Code: **GH** Color: **Light Gray**

CDC Analyte Identifier Code: **2547**

CDC Test System Identifier Code: **10159**

CLIA Complexity Category: **Moderate**

Summary and Explanation

Human growth hormone (hGH, somatotropin) is a polypeptide originating in the anterior pituitary. It is 191 amino acids in length and has a molecular mass of approximately 22,000 daltons. Its metabolic effects are primarily anabolic. It promotes protein conservation and engages a wide range of mechanisms for protein synthesis. It also enhances glucose transport and facilitates the buildup of glycogen stores. Another family of peptide hormones, the somatomedins, mediates its cascade of growth-promoting actions.¹

Measurement of hGH is primarily of interest in the diagnosis and treatment of various forms of inappropriate growth hormone secretion. Clinical disorders of hyposecretion include dwarfism and unattained growth potential. Hypersecretion is associated with gigantism and acromegaly.

Caution must be exercised in the clinical interpretation of growth hormone levels. These vary throughout the day, making it difficult to define a reference range or to judge an individual's status based on single determinations. Many factors are known to influence the rate of growth hormone secretion, including periods of sleep and wakefulness, exercise, stress, hypoglycemia, estrogens, corticosteroids, L-dopa, and others.¹ Because of its similarity to prolactin and placental lactogen, earlier growth hormone

immunoassays were often plagued with falsely high values in pregnant and lactating women.

Because not all acromegalic individuals have elevated baseline levels, suppression tests based on glucose loading are of value in this context. In spite of the induced hyperglycemia, there is rarely a decrease from baseline levels in acromegaly.²

Growth hormone-deficient individuals have fasting/resting levels similar to those found in healthy individuals. Various challenge tests have therefore been devised to differentiate these groups. Thus, with the onset of deep sleep or after 15 to 20 minutes of vigorous exercise, growth hormone levels normally show a rise. Other tests of growth hormone responsiveness are based on the administration of L-dopa, arginine and insulin. Propranolol or estrogen are sometimes given in conjunction with the primary stimulus to accentuate the response.^{4,5}

A small number of cases of dwarfism have been documented in which both the basal level and the response to challenge testing were normal. Such cases may involve tissue insensitivity to either growth hormone or the somatomedins, or the presence of antibodies or immunoreactive but biologically inactive growth hormone.⁴

Principle of the Procedure

IMMULITE/IMMULITE 1000 Growth Hormone is a solid-phase, two-site chemiluminescent immunometric assay.

The solid phase (bead) is coated with murine monoclonal anti-hGH antibody. The reagent contains alkaline phosphatase (bovine calf intestine) conjugated to a rabbit anti-hGH polyclonal antibody. The reagent and hGH in the sample are incubated together with a bead coated with a murine anti-hGH monoclonal antibody to form an antibody-sandwich complex. Unbound enzyme conjugate is then removed by a centrifugal wash. Finally, chemiluminescent substrate is added to the bead and signal is generated in proportion to the bound enzyme.

Incubation Cycles: 1 × 30 minutes.

Specimen Collection

Patient must be fasting and at complete rest 30 minutes before blood collection.^{6,7}

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Icteric or grossly contaminated samples may give erroneous results.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE/IMMULITE 1000 Growth Hormone has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 50 µL serum. (Sample cup must contain at least 100 µL more than the total volume required.)

Storage: Refrigerate immediately. Stable at 2–8°C for 8 hours,^{6,7} or for 2 months (aliquotted) at –20°C.

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

GH Test Units (LGH1)

Each barcode-labeled unit contains one bead coated with murine monoclonal anti-hGH antibody. Stable at 2–8°C until expiration date.

LKGH1: 100 units. **LKGH5:** 500 units.

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

GH Reagent Wedge (LGH2)

With barcode. 7.5 mL alkaline phosphatase (bovine calf intestine) conjugated to rabbit polyclonal anti-hGH antibody in buffer, with preservative. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

LKGH1: 1 wedge. **LKGH5:** 5 wedges.

GH Adjustors (LGHL, LGHH)

Two vials (Low and High) containing lyophilized hGH in nonhuman serum, with preservative. Reconstitute each vial with **3.0 mL** distilled or deionized water. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. Stable at 2–8°C for 30 days after reconstitution, or for 6 months (aliquotted) at –20°C.

LKGH1: 1 set. **LKGH5:** 2 sets.

Kit Components Supplied Separately

GH Sample Diluent (LGHZ)

For the manual dilution of patient samples. 25 mL of processed nonhuman serum,

with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

LSUBX: Chemiluminescent Substrate
LPWS2: Probe Wash Module
LKPM : Probe Cleaning Kit
LCHx-y: Sample Cup Holders (barcoded)
LSCP: Sample Cups (disposable)
LSCC: Sample Cup Caps (optional)

CON6: Tri-level, multi-constituent control

Also Required

Sample transfer pipets, distilled or deionized water, controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for: preparation, setup, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

Recommended Adjustment Interval:
2 weeks.

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of growth hormone.

Expected Values

A study performed at a large reference laboratory in the western United States on 59 healthy adult males and 60 healthy adult females yielded the following reference intervals:

Males: Up to 1 ng/mL

Females: Up to 10 ng/mL

Secretion of growth hormone by the pituitary gland is episodic and pulsatile, and transient levels up to 40 ng/mL have been observed in healthy subjects. Because serum levels are rather low between pulses in healthy individuals, immunoassays with sufficient sensitivity are necessary to distinguish patients with abnormally low values from healthy subjects who have values that happen to fall in the low-normal range. Thus, basal or random levels of hGH provide little

diagnostic information. In some individuals, spontaneous growth hormone secretion is better monitored by using a continuous withdrawal pump, or by drawing specimens for hGH assay every 20 to 30 minutes over a 12-to-24-hour period. A number of provocative tests designed to stimulate or suppress release are also used to assess growth hormone secretion; the definition of abnormal responses varies with the laboratory.¹

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

Patients on hGH therapy may develop human antibodies to hGH that will interfere in the assay and cause falsely low values.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Please note that beginning with **software versions 4.13 (IMMULITE) and 5.7b (IMMULITE 1000)**, the alternate units are in terms of WHO NIBSC 2nd International Standard 98/574 for r-DNA Growth Hormone, with a conversion factor of 2.4.

For previous software versions, the alternate units are in terms of WHO

NIBSC 1st International Standard 80/505, with a conversion factor of 2.6.

If customers reprint results generated on previous software versions, results will be expressed relative to WHO NIBSC 2nd International Standard 98/574 for r-DNA Growth Hormone.

Conversion Factor:

ng/mL × 2.6 → mIU/L WHO 1st IS 80/505
ng/mL × 2.4 → mIU/L WHO 2nd IS 98/574

Calibration Range: Up to 40 ng/mL
(Up to 104 mIU/L, WHO 1st IS 80/505)
(Up to 96 mIU/L, WHO 2nd IS 98/574)

Analytical Sensitivity: Up to 0.01 ng/mL
(0.026 mIU/L).

High-dose Hook Effect:

None up to 17,000 ng/mL.

Intraassay Precision (Within-Run):

Statistics were calculated for samples from the results of 20 replicates in a single run. (See "Intraassay Precision" table.)

Interassay Precision (Run-to-Run):

Statistics were calculated for samples assayed in 11 different runs. (See "Interassay Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three growth hormone solutions (115, 230 and 476 ng/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody used in the IMMULITE Growth Hormone procedure is highly specific for hGH, with low crossreactivity to other naturally occurring hormones. (See "Specificity" table.)

A second, smaller form of growth hormone (approximately 20,000 daltons) accounts for 15% of the total growth hormone content of the human pituitary. This form has the normal growth promoting activity, but lacks the insulin-like activity associated with the 22,000 dalton form. The 20 kDa form was tested for crossreactivity in the IMMULITE Growth Hormone assay and was found to crossreact to the extent of 63%.

Bilirubin: Presence of bilirubin in concentrations up to 100 mg/L has no effect on results, within the precision of the assay. Levels of 200 mg/L interfere with

the assay, causing a slight depression of values. (See "Bilirubin" table.)

Hemolysis: Presence of packed red blood cells in concentrations up to 30 µL/mL has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 15 volunteers into plain, heparinized, and Becton Dickinson SST[®] vacutainer tubes. Equal volumes of the matched samples were spiked with various concentrations of hGH, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE Growth Hormone procedure.

(Heparin) = 0.99 (Serum) – 0.55 ng/mL
r = 0.982

(SST) = 0.92 (Plain Tubes) + 2.3 ng/mL
r = 0.975

Means:

26.1 ng/mL (Heparin)
25.8 ng/mL (Serum)
26.0 ng/mL (SST Tubes)

Method Comparison: The IMMULITE Growth Hormone procedure was compared to an immunoradiometric assay (Kit A) on 84 patient samples. (Concentration range: approximately 0.04 to 27 ng/mL. See graph.) By linear regression:

(IML) = 0.92 (Kit A) + 0.47 ng/mL
r = 0.937

Means:

6.0 ng/mL (IMMULITE)
6.0 ng/mL (Kit A)

References

- 1) Whitley RJ, Meikle AW, Watts NB. Endocrinology. Part 2: Protein hormones. In: Burtis CA, Ashwood ER, editors. Tietz textbook of clinical chemistry. 2nd ed. Philadelphia: Saunders, 1994: 1665-70. 2) Daughaday WM, Cryer PE. Growth hormone hypersecretion and acromegaly. Hospital Practice. no. 8 August 1978;13:75-80. 3) Raite S. The standards for human growth hormone assays. In: Laron Z, Butenandt O, editors. Evaluation of Growth Hormone Secretion. Basel: Karger, 1983: 162-9. 4) Van Wyk JJ. Diagnosis of growth hormone deficiency and other causes of growth failure. The Endocrine Society: Syllabus for the 35th Annual Postgraduate Assembly October 1983: 90-106. 5) Van Wyk JJ, Underwood L. Growth hormone, somatomedins and growth failure. Hospital Practice no. 8 August 1978;13:57-67. 6) Tietz NW, editor. Clinical guide to laboratory tests. 3rd ed. Philadelphia: W.B. Saunders,

1995: 300. 7) Burtis CA, Ashwood ER, editors. Tietz textbook of clinical chemistry. 2nd ed. Philadelphia: W.B. Saunders, 1994: 1670. 8) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782 or 973.927.2828
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

Manufactured by EURO/DPC Ltd. under a Quality System registered to ISO 13485:2003.

Tables and Graphs

Intraassay Precision (ng/mL)

	Mean ¹	SD ²	CV ³
1	1.7	0.09	5.3%
2	3.1	0.18	5.8%
3	7.8	0.47	6.0%
4	15	0.93	6.2%
5	31	2.0	6.5%

Interassay Precision (ng/mL)

	Mean ¹	SD ²	CV ³
1	3.0	0.17	5.7%
2	7.1	0.39	5.5%
3	9.3	0.58	6.2%
4	16	0.90	5.6%
5	18	1.1	6.1%

Linearity (ng/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	12	—	—
	4 in 8	5.1	6.0	85%
	2 in 8	2.6	3.0	87%
	1 in 8	1.5	1.5	100%
2	8 in 8	23	—	—
	4 in 8	11	12	92%
	2 in 8	5.1	5.8	88%
	1 in 8	2.7	2.9	93%
3	8 in 8	29	—	—
	4 in 8	14	15	93%
	2 in 8	6.3	7.3	86%
	1 in 8	3.2	3.6	89%
4	8 in 8	40	—	—
	4 in 8	21	20	105%
	2 in 8	9.9	10	99%
	1 in 8	4.5	5.0	90%

Recovery (ng/mL)

Sample	Spiking		Expected ³	%O/E ⁴
	Solution ¹	Observed ²		
1	—	7.5	—	—
	A	15	13	115%
	B	21	19	111%
	C	33	31	106%
2	—	8.2	—	—
	A	14	14	100%
	B	23	20	115%
	C	35	32	109%
3	—	23	—	—
	A	31	28	111%
	B	36	34	106%
	C	>40	46	—
4	—	26	—	—
	A	31	31	100%
	B	38	37	103%
	C	>40	49	—

Specificity


Compound ¹	ng/mL Added ²	Apparent Conc. ng/mL ³	% Cross-reactivity ⁴
FSH	10,000	ND	ND
β-subunit HCG	10,000	ND	ND
HPL	10,000	ND	ND
LH	3,663	ND	ND
Prolactin	40,000	ND	ND
	10,000	ND	ND
	1,000	ND	ND
	200	ND	ND
TSH	11,125	0.941	0.0085%

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
EURO/DPC LTD
Diagnostic Products Corporation
EC REP DPC Biermann GmbH