IMMULITE®

Carbamazepine

For use on the IMMULITE®
and IMMULITE® 1000 systems

DPC®

IMMULITE®/IMMULITE® 1000 Carbamazepine

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE and IMMULITE 1000 Analyzers — for the quantitative measurement of carbamazepine in serum or heparinized plasma, as an aid in monitoring the therapeutic administration of this drug.

Catalog Number: **LKCB1** (100 tests), **LKCB5** (500 tests)

Test Code: **CMP** Color: **Light Gray**

Summary and Explanation

Carbamazepine is an iminostilbene derivative used for the treatment of epilepsy; trigeminal neuralgia; and simple and complex, partial and generalized convulsive seizures.¹⁻³ It is administered either alone or in combination with other antiepileptic drugs.

Carbamazepine is rapidly absorbed into the bloodstream, where it is highly protein bound (60 to 80 percent).³ The drug is metabolized by hepatic oxidation to the 10,11-epoxide form, which is as active as the parent drug. The epoxide form is further metabolized to the 10,11-dihydroxide, which is eliminated in the urine.^{1,4} The presence of other drugs can significantly alter the rate of epoxide formation: phenytoin, phenobarbital and felbamate induce the metabolic enzyme activity, whereas erythromycin and propoxyphene inhibit this activity.¹⁻⁴

The therapeutic range for carbamazepine is between 4 and 12 µg/mL (between 17 and 51 µmol/L), with toxic levels above 15 µg/mL.³ However, individual differences in absorption and metabolism can result in highly variable responses to carbamazepine. Monitoring levels assists the physician in optimizing dosage and minimizing toxic side effects for each patient.

Principle of the Procedure

IMMULITE/IMMULITE 1000 Carbamazepine is a solid-phase, competitive chemiluminescent enzyme immunoassay.

Incubation Cycles: 1 × 30 minutes.

Specimen Collection

The time of collection will depend on whether the desired measurement is for peak or trough values.⁶

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants.

IMMULITE/IMMULITE 1000 Carbamazepine has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 20 µL serum or heparinized plasma. (Sample cup must contain at least 100 µL more than the total volume required.)

Storage: 2 days at 2–8°C or 2 months at –20°C.

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis

B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

Carbamazepine Test Units (LCB1)

Each barcode-labeled unit contains one bead coated with monoclonal murine anti-carbamazepine antibody. Stable at 2–8°C until expiration date.

LKCB1: 100 units. **LKCB5:** 500 units.

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

Carbamazepine Reagent Wedge (LCB2)

With barcode. 7.5 mL alkaline phosphatase (bovine calf intestine) conjugated to carbamazepine in buffer. Store capped and refrigerated: stable at 2–8°C until expiration date.

Recommended usage is within 30 days after opening when stored as indicated.

LKCB1: 1 wedge. **LKCB5:** 5 wedges.

Carbamazepine Adjustors (LCBL, LCBH)

Two vials (Low and High), 2 mL each, of carbamazepine in processed human serum, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

LKCB1: 1 set. **LKCB5:** 2 sets.

Kit Components Supplied Separately

LSUBX: Chemiluminescent Substrate

LPWS2: Probe Wash Module

LKPM: Probe Cleaning Kit

LCHx-y: Sample Cup Holders (barcoded)

LSCP: Sample Cups (disposable)

LSCC: Sample Cup Caps (optional)

SDCM: Bi-level, human serum-based control containing carbamazepine.

Also Required

Sample transfer pipets, distilled or deionized water, controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

Recommended Adjustment Interval:

2 weeks.

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of carbamazepine.

Therapeutic Range

The following results are reported in the literature, with the therapeutic range for carbamazepine between 4 and 12 µg/mL (17 and 51 µmol/L) and toxic levels above 15 µg/mL (63 µmol/L) for trough level in adults and children.³

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

A complete understanding of the information contained in this package insert is necessary for proper use of this product.

The specificity table should be consulted for a list of substances which have measurable crossreactivity in this assay. Other medications and/or food not found on this list may nevertheless crossreact in the assay.

The following should be considered in the interpretation of the carbamazepine concentration:

- All concurrent medications taken by the patient
- The time of sampling in relation to the dosing history
- Existing diseases or conditions

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1998;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in µg/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Conversion Factor:
µg/mL × 4.23 → µmol/L

Calibration Range: 1.25 – 20 µg/mL
(5.29 – 85 µmol/L).

Analytical Sensitivity: 0.2 µg/mL
(0.9 µmol/L).

Intraassay Precision (Within-Run):
Statistics were calculated for samples from the results of 20 replicates in a single run. (See "Intraassay Precision" table.)

Interassay Precision (Run-to-Run):
Statistics were calculated for samples assayed in 10 different runs. (See "Interassay Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with two carbamazepine solutions (96 and

202 µg/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for carbamazepine, with less than 20% crossreactivity to the major metabolite (carbamazepine-10,11-epoxide) and no detectable crossreactivity to other compounds and therapeutic drugs that may be present in patient samples. (See "Specificity" table.)

Effect of Epoxide: Various amounts of carbamazepine-10,11-epoxide were spiked into samples containing various concentrations of carbamazepine. (See "Effect of Epoxide" table.)

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 508 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 5,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 20 volunteers into plain and EDTA tubes. Equal volumes of the matched samples were spiked with various concentrations of carbamazepine, to obtain values throughout the calibration range of the assay, and then assayed by the Carbamazepine I procedure.

(EDTA) = 0.94 (Serum) + 0.02 µg/mL
r = 0.999

Means:
6.7 µg/mL (Serum)
6.3 µg/mL (EDTA)

Method Comparison: The assay was compared to Kit A on 38 patient samples. (Concentration range: approximately 0.7 to 14.7 µg/mL. See first graph.) By linear regression:

(IMMULITE) = 1.1 (Kit A) + 0.12 µg/mL
r=0.972

Means:
7.2 µg/mL (IMMULITE)
6.4 µg/mL (Kit A)

In a second study, the assay was compared to Kit B on 47 patient samples (Concentration range: approximately 1.6

to 11 µg/mL. See second graph.) By linear regression:

(IMMULITE) = 1.1 (Kit B) – 0.3 µg/mL
r = 0.973

Means:
7.3 µg/mL (IMMULITE)
7.1 µg/mL (Kit B)

References

- 1) Drug monitoring data pocket guide II. 2nd ed. Washington, DC: AACCC Press, 1994: 49-51.
- 2) McNamara JO. Drugs effective in the treatment of the epilepsies. In: GL Hardman, LL Limbird, PB Molinoff, RW Ruddon, editors. Goodman & Gilman's the pharmacological basis of therapeutics. 9th ed. New York: McGraw-Hill, 1996: 461-86.
- 3) Moyer TP, Pippenger CE. Therapeutic drug monitoring. In: Burtis CA, Ashwood ER, editors. Tietz textbook of clinical chemistry. 2nd ed. Philadelphia: WB Sanders, 1994: 1094-154.
- 4) Kutt H. Carbamazepine: chemistry and methods of determination. Complex partial seizures and their treatment. Adv Neurol 1975;11:249-61.
- 5) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.
- 6) Jacobs DS, Kaster BL Jr, Demott WR, Wolfson WL. Laboratory test handbook. Stowe, OH: Lexi-Compl. Mosby, 1990: 782.

Technical Assistance

In the United States, contact DPC's Technical Services department.
Tel: 800.372.1782 or 973.927.2828
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

Manufactured by EURO/DPC Ltd. under a Quality System registered to ISO 13485:2003.

Tables and Graphs

Intraassay Precision (µg/mL)

	Mean ¹	SD ²	CV ³
1	1.4	0.07	5.0%
2	3.0	0.10	3.3%
3	5.0	0.16	3.2%
4	8.3	0.23	2.8%
5	9.8	0.34	3.5%
6	15.5	0.40	2.6%
7	19.8	0.53	2.7%

Interassay Precision (µg/mL)

	Mean ¹	SD ²	CV ³
1	1.3	0.11	8.5%
2	2.2	0.18	8.2%
3	4.8	0.22	4.6%
4	7.9	0.25	3.2%
5	9.5	0.43	4.5%
6	14.8	0.49	3.3%
7	18.6	1.01	5.4%

Linearity (µg/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	7.8	—	—
	4 in 8	3.7	3.9	95%
	2 in 8	1.7	2.0	85%
	1 in 8	0.95	0.98	97%
2	8 in 8	11.6	—	—
	4 in 8	5.9	6.0	98%
	2 in 8	3.3	3.0	110%
	1 in 8	1.5	1.5	100%
3	8 in 8	15.4	—	—
	4 in 8	7.7	7.7	100%
	2 in 8	3.5	3.9	90%
	1 in 8	1.8	1.9	95%
4	8 in 8	16.4	—	—
	4 in 8	8.3	8.2	101%
	2 in 8	4.1	4.1	100%
	1 in 8	1.9	2.1	90%
5	8 in 8	18.5	—	—
	4 in 8	8.9	9.3	96%
	2 in 8	4.7	4.6	102%
	1 in 8	2.1	2.3	91%

Recovery (µg/mL)

	Solution ¹	Observed ²	Expected ³	% O/E ⁴
1	—	1.8	—	—
	A	6.7	6.5	103%
	B	11.5	11.8	97%
2	—	2.0	—	—
	A	6.9	6.7	103%
	B	12.2	12.0	102%
3	—	4.9	—	—
	A	10.2	9.5	107%
	B	15.2	14.8	103%
4	—	5.3	—	—
	A	9.9	9.8	101%
	B	15.6	15.1	103%
5	—	10.1	—	—
	A	14.6	14.4	101%
	B	19.6	19.7	99%

Specificity


Compound ¹	µg/mL Added ²	Apparent Conc. ³	% Cross-reactivity ⁴
Carbamazepine-10,11-epoxide	20	3.7	19%
Clonazepam	0.12	ND	ND
Diazepam	2.5	ND	ND
Ethosuximide	100	ND	ND
Imipramine	1.0	ND	ND
Phenobarbital	40	ND	ND
Phenytoin	20	ND	ND
Phenylbutazone	45	ND	ND
Primidone	12	ND	ND
Secobarbital	15	ND	ND
Theophylline	20	ND	ND
Valproic acid	100	ND	ND

ND: not detectable.⁵

Effect of Epoxide (µg/mL)


	Epoxide Added ¹	Observed ²
1	Unspiked	5.1
	+ 3.0	5.0
	+ 6.0	6.0
	+12.0	5.9
	+ 24.0	7.8
2	Unspiked	7.8
	+ 3.0	7.7
	+ 6.0	7.1
	+12.0	8.1
	+ 24.0	9.5
3	Unspiked	10.3
	+ 3.0	10.6
	+ 6.0	10.2
	+12.0	11.5
	+ 24.0	11.5
4	Unspiked	12.5
	+ 3.0	11.7
	+ 6.0	13.2
	+12.0	12.6
	+ 24.0	15.0

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
EURO/DPC LTD
Diagnostic Products Corporation
EC REP DPC Biermann GmbH