# IMMULITE®

## Digitoxin

For use on the IMMULITE®  
and IMMULITE® 1000 systems

**DPC®**

# IMMULITE®/IMMULITE® 1000 Digitoxin

## English

**Intended Use:** For *in vitro* diagnostic use with the IMMULITE and IMMULITE 1000 Analyzers — for the quantitative measurement of digitoxin in serum, as an aid in the therapeutic administration of this cardioglycoside.

Catalog Number: **LKDG1** (100 tests), **LKDG5** (500 tests)

Test Code: **DGT** Color: **Red**

CDC Analyte Identifier Code: 1303

CDC Test System Identifier Code: 10159

CLIA Complexity Category: Moderate

### Summary and Explanation

Digitoxin is a long-acting cardiac glycoside used in the treatment of congestive heart failure and certain arrhythmias. It is structurally identical to digoxin except for the absence of a single hydroxyl group. Being less polar, digitoxin is absorbed more completely, and circulates to a much greater extent protein-bound. In the liver, it is extensively metabolized to both active and inactive forms. Digitoxin has a comparatively long half-life of 5 to 7 days and, except in the nephrotic syndrome, is largely unaffected by renal dysfunction.

A transient peak is associated with the distribution phase; thereafter the circulating level exhibits a steady and very gradual decline. Accordingly, blood samples for digitoxin analysis should be collected just before the next dose — or in any case not less than 6 hours after the previous dose — to allow time for plasma and tissue levels to reach equilibrium.

Serum determinations can assist in verifying compliance and in adjusting the dose to achieve a desired plasma concentration. Monitoring circulating digitoxin levels is especially valuable when there are enterohepatic or renal complications, and during concomitant therapy with drugs like quinidine, which can affect the pharmacokinetics of digitoxin.

It must be remembered, however, that the relation between circulating level and

clinical response varies greatly from one individual to another, due to the many factors which are at work. This is reflected in the considerable overlap between therapeutic and toxic ranges summarized in the literature.

Thus, circulating levels of approximately 10 to 25 ng/mL are observed in the majority of patients exhibiting a satisfactory therapeutic response to digitoxin, while levels greater than 45 ng/mL are typical of patients with symptoms of digitoxin toxicity.<sup>2,4,16</sup> But concentrations higher than the conventional therapeutic range may be required to control certain arrhythmias. Conversely, circulating levels well below 45 ng/mL may give rise to symptoms of toxicity in patients rendered more sensitive to the drug's action by, for example, a lowered potassium level. Even so, the higher the circulating level, the more probable it is that clinical symptoms of toxicity can be attributed to the drug — and the lower the level, the less likely.

Serum determinations should supplement clinical judgment, not substitute for it. An isolated value should not be used as the sole criterion for the assessment of efficacy or toxicity; rather it should always be interpreted in the total clinical context.

### Principle of the Procedure

IMMULITE/IMMULITE 1000 Digitoxin is a solid-phase, competitive chemiluminescent enzyme-labeled immunoassay.

**Incubation Cycles:** 1 × 30 minutes.

### Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Collect blood into plain (without anticoagulant) or heparinized tubes, and separate the serum from the cells. The time of collection should be noted.

To avoid misleadingly high digitoxin results, patient samples should be drawn

just before the next dose — or in any case no less than six hours after the administration of this cardiac glycoside.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE/IMMULITE 1000 Digitoxin has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

**Volume Required:** 20  $\mu$ L serum. (Sample cup must contain at least 100  $\mu$ L more than the total volume required.)

**Storage:** 7 days at 2–8°C or 2 months at –20°C.

## Warnings and Precautions

For *in vitro* diagnostic use.

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight. (See insert.)

**Water:** Use distilled or deionized water.

## Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

### Digitoxin Test Units (LDG1)

Each barcode-labeled unit contains one bead coated with polyclonal rabbit digitoxin antibody. Stable at 2–8°C until expiration date.

**LKDG1:** 100 units. **LKDG5:** 500 units.

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

### Digitoxin Reagent Wedge (LDG2)

With barcode. 7.5 mL alkaline phosphatase (bovine calf intestine) conjugated to digitoxin in buffer, with preservative. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

**LKDG1:** 1 wedge. **LKDG5:** 5 wedges.

### Digitoxin Adjustors (LDGL, LDGH)

Two vials (Low and High), 4 mL each, of digitoxin in processed human serum, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

**LKDG1:** 1 set. **LKDG5:** 2 sets.

## Kit Components Supplied Separately

**LSUBX:** Chemiluminescent Substrate

**LPWS2:** Probe Wash Module

**LKPM :** Probe Cleaning Kit

**LCHx-y:** Sample Cup Holders (barcoded)

**LSCP:** Sample Cups (disposable)

**LSCC:** Sample Cup Caps (optional)

**SDCM:** Bi-level, multi-constituent control

Also Required

Sample transfer pipets, distilled or deionized water, controls.

## Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

**Recommended Adjustment Interval:**  
2 weeks.

**Quality Control Samples:** Use controls or sample pools with at least two levels (low and high) of digitoxin.

### Therapeutic Range

The following results are reported in the literature for therapeutic and toxic circulating levels of digitoxin.<sup>16</sup>

	ng/mL	nmol/L
Therapeutic Range	10 – 25	13 – 33
Toxic Level	> 45	> 59

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

### Limitations

Many physiological, pharmacological, pathological and genetic factors may affect the interpretation of results (see the Summary and Explanation section).

In cases of overdosage, specimens should be obtained *before* administering digoxin immune Fab antidote preparations (e.g. Digibind<sup>®</sup>) since these can be expected to interfere with any digitoxin immunoassay procedure. Once such a preparation has been administered, accurate determinations of digitoxin can be made only by specialized techniques (not routinely available in most clinical settings) until the Fab fragments have been eliminated from the body: this takes several days for patients with intact renal function, and a week or longer for patients with renal impairment. In addition to *in vitro* effects, these preparations, administered intravenously, cause a rapid rise in pharmacologically inactive, Fab-bound, total circulating digitoxin that may be detected immunologically.<sup>15</sup>

The presence of Uzara<sup>®</sup> in the serum of patients can result in falsely elevated digitoxin results.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

### Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

**Conversion Factor:**  
ng/mL × 1.307 → nmol/L

**Calibration Range:** 5 – 60 ng/mL  
(6.5 – 78 nmol/L).

**Analytical Sensitivity:** 0.65 ng/mL  
(0.85 nmol/L).

**Intraassay Precision (Within-Run):**  
Statistics were calculated for samples from the results of 20 replicates in a single run. (See "Intraassay Precision" table.)

**Interassay Precision (Run-to-Run):**  
Statistics were calculated for samples assayed in 25 different runs. (See "Interassay Precision" table.)

**Specificity:** The antibody is highly specific for digitoxin, with low crossreactivity to other compounds and therapeutic drugs that may be present in patient samples. (See "Specificity" table.)

**Linearity:** Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

**Recovery:** Samples spiked 1 to 19 with three digitoxin solutions (154, 594 and 1,218 ng/mL) were assayed. (See "Recovery" table for representative data.)

**Bilirubin:** Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

**Hemolysis:** Presence of packed red blood cells in concentrations up to 30  $\mu$ L/mL has no effect on results, within the precision of the assay.

**Alternate Sample Type:** To assess the effect of alternate sample types, blood was collected from 13 volunteers into plain, heparinized and Becton Dickinson SST<sup>®</sup> vacutainer tubes. Equal volumes of the matched samples were spiked with various concentrations of digitoxin, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE Digitoxin procedure.

(Heparin) = 0.89 (Serum) + 0.96 ng/mL  
 $r = 0.982$

(SST) = 0.99 (Plain Tubes) + 0.22 ng/mL  
 $r = 0.954$

Means:

14.8 ng/mL (Serum)  
14.2 ng/mL (Heparin)  
14.9 ng/mL (SST)

**Method Comparison:** The assay was compared to DPC's Double Antibody Digitoxin on 139 patient samples. (Concentration range: approximately 5.8 to 68 ng/mL. See graph.) By linear regression:

(IML) = 0.96 (DAb) – 0.55 ng/mL  
 $r = 0.961$

Means:

26 ng/mL (IMMULITE)  
28 ng/mL (DAb)

## References

- 1) Aronson JK. Clinical pharmacokinetics of cardiac glycosides in patients with renal dysfunction. Clin Pharmacokinet 1983;8:155-78.
- 2) Baer DM. Key points for interpretation of drug concentrations: digoxin and digitoxin. In: Baer DM, Dito WR, editors. Interpretations in therapeutic drug monitoring. Chicago: American Society of Clinical Pathologists, 1982: 351-4.
- 3) Bresnahan JF, Vlietstra RE. Digitalis glycosides. Mayo Clin Proc 1979;54:675-84.
- 4) Dito WR. Cardiac glycosides and therapeutic drug monitoring. In: Baer DM, Dito WR, editors. Interpretations in therapeutic drug monitoring. Chicago: American Society of Clinical Pathologists, 1981: 61-73.
- 5) Fenster PE, et al. Digitoxin-quinidine interaction: pharmacokinetic evaluation. Ann Int Med 1980;93:698-701.
- 6) Garty M. Digitoxin elimination reduced during quinidine therapy. Ann Int Med 1981;94:35-7.
- 7)

- 8) Kuemmel HC, et al. Studies to determine individual toxic doses of digitoxin. Munchener Med Wochenschr 1982;124:545-9.
- 9) Kuhlmann J. Sources of digitalis intoxication. Med Klinik 1982;77:39-44.
- 10) Larsen A, Storstein L. Digitoxin kinetics and renal excretion in children. Clin Pharmacol Ther 1983;33:717-27.
- 11) Maier C, Kolenda KD. Digitoxin in the presence of hepatorenal failure. Deutsche Med Wochenschr 1983;108:1475-80.
- 12) Rietbrock N, Alken RG. The therapy of heart failure with digitalis. Deutsche Med Wochenschr 1980;105:1622-8.
- 13) Smith TW. Digitalis in the management of heart failure. Hosp Pract 1984 Mar;19(3):67-92.
- 14) Smith TW, et al. Digitalis glycosides: mechanisms and manifestations of toxicity. Prog Cardiovasc Dis 1984;26:413-58 and 495-540.
- 15) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.
- 16) Digoxin immune Fab. In: McEvoy GK, editor. AHFS Drug Information 92. American Society of Hospital Pharmacists, 1992: 1977-81.
- 17) National Academy of Clinical Biochemistry. Guidelines for therapeutic drug monitoring services. Printed by Durik Advertising, Inc., 1999: 82-83.

## Technical Assistance

In the United States, contact DPC's Technical Services department.  
Tel: 800.372.1782 or 973.927.2828  
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

Manufactured by EURO/DPC Ltd. under a Quality System registered to ISO 13485:2003.

## Tables and Graphs

### Intraassay Precision (ng/mL)

	Mean <sup>1</sup>	SD <sup>2</sup>	CV <sup>3</sup>
1	8.6	0.43	5.0%
2	22	0.86	3.9%
3	46	2.1	4.6%

### Interassay Precision (ng/mL)

	Mean <sup>1</sup>	SD <sup>2</sup>	CV <sup>3</sup>
1	8.9	0.74	8.3%
2	24	1.5	6.3%
3	48	2.8	5.8%

### Specificity


Compound <sup>1</sup>	ng/mL Added <sup>2</sup>	% Cross-reactivity <sup>3</sup>
Digitoxin	—	100%
Digoxin	250	6.9%
Aldosterone	100,000	0.011%
Amikacin	100,000	0.003%
Corticosterone	100,000	0.012%
Cortisol	100,000	0.014%
Cortisone	100,000	0.008%
Dexamethasone	100,000	0.003%
11-Deoxycorticosterone	100,000	0.026%
Hydrocortisone	100,000	0.007%
Estradiol	100,000	0.004%
Estriol	100,000	ND
Methotrexate	100,000	ND
Prednisolone	100,000	0.005%
Prednisone	100,000	0.003%
Pregnenolone	100,000	0.007%
Progesterone	100,000	0.026%
Spirolactone	100,000	0.008%
Testosterone	100,000	0.027%
Tobramycin	100,000	0.005%
Uzara	1,397	2.11%
	2,794	1.96%

ND: not detectable.<sup>4</sup>

### Linearity (ng/mL)

	Dilution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	5 in 5	27	—	—
	4 in 5	19	22	86%
	3 in 5	15	16	94%
	2 in 5	9.9	11	90%
	1 in 5	5.3	5.4	98%
2	5 in 5	27	—	—
	4 in 5	21	22	95%
	3 in 5	18	16	113%
	2 in 5	11	11	100%
3	5 in 5	28	—	—
	4 in 5	23	22	105%
	3 in 5	16	17	94%
	2 in 5	12	11	109%
4	5 in 5	31	—	—
	4 in 5	22	25	88%
	3 in 5	17	19	89%
	2 in 5	12	12	100%
	5 in 5	6.8	6.2	110%

### Method Comparison

$$(IML) = 0.96 (DAb) - 0.55 \text{ ng/mL}$$

$$r = 0.961$$

## Recovery (ng/mL)

	Solution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	—	13	—	—
	A	18	20	90%
	B	34	42	81%
	C	> 60	73	N/A
2	—	15	—	—
	A	22	22	100%
	B	40	44	91%
	C	> 60	75	N/A
3	—	15	—	—
	A	22	22	100%
	B	44	44	100%
9	C	> 60	75	N/A
4	—	18	—	—
	A	22	25	88%
	B	42	47	89%
	C	> 60	78	N/A
5	—	24	—	—
	A	30	31	97%
	B	48	53	91%
	C	> 60	84	N/A

N/A: not applicable, i.e. observed value exceeds 60 ng/mL.<sup>5</sup>

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applicable, i.e. observed value exceeds 60 ng/mL. **Method Comparison:** Digitoxin: Digitoxine.

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
**EURO/DPC** LTD
Diagnostic Products Corporation
EC REP DPC Biermann GmbH