# IMMULITE®

## Digoxin

For use on the IMMULITE®  
and IMMULITE® 1000 systems

**DPC®**

# IMMULITE®/IMMULITE® 1000 Digoxin

## English

**Intended Use:** For *in vitro* diagnostic use with the IMMULITE and IMMULITE 1000 Analyzers — for the quantitative measurement of digoxin in serum or heparinized plasma, as an aid in monitoring the therapeutic administration of this cardioglycoside, while avoiding toxicity.

Catalog Number: **LKDI1** (100 tests), **LKDI5** (500 tests)

Test Code: **DGX** Color: **Light Green**

### Summary and Explanation

Immunoassay has proved particularly useful in determining drug overdoses in patients treated with cardiac glycosides. The technique is also useful (a) in clarifying situations where a patient's symptoms might be due either to intrinsic heart disease or to digitalis intoxication; (b) where there is doubt concerning the type of digitalis preparation the patient is taking — in this case, a digitoxin immunoassay is also necessary; (c) for measuring the digoxin ingestion of patients with an inadequate history of previous dosage; (d) in documenting cases of underdigitalization as well as digitalis (digoxin) excess; (e) in monitoring the toxic response in patients with myocardial disease associated with hypokalemia, hypomagnesemia, hypercalcemia, hypoxia and alkalosis, which are particularly sensitive to digitalis; and (f) in preventing overdigitalization, particularly in patients whose renal function is deteriorating or for whom an increased digoxin dosage is contemplated.

The high sensitivity of digoxin immunoassay is especially necessary in view of the small differences and occasional overlap that exist between therapeutic and toxic levels of circulating digoxin.<sup>13</sup> (Intoxication is defined in terms of arrhythmias and disturbances of cardiac conduction due to the drug's presence.) Smith et al reported serum digoxin concentrations of 0.8–2.4 ng/mL (1.0–3.1 nmol/L) in nontoxic patients, and 2.1–

8.7 ng/mL (2.7–11.1 nmol/L) in toxic patients, based on data taken six hours post-dose in patients with normal renal function.<sup>8</sup> More recent clinical investigations have confirmed the association of toxicity with serum levels above 2 ng/mL.

### Principle of the Procedure

IMMULITE/IMMULITE 1000 Digoxin is a solid-phase, competitive chemiluminescent enzyme immunoassay.

**Incubation Cycles:** 1 × 30 minutes.

### Specimen Collection

Collect blood in plain tubes without gel barriers, noting the time of collection, just before the next dose — or, in any case, no less than six hours after the administration of this cardioglycoside, to avoid misleadingly high digoxin results.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants.

IMMULITE/IMMULITE 1000 Digoxin has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

**Volume Required:** 50 µL serum. (Sample cup must contain at least 100 µL more than the total volume required.)

**Storage:** 7 days at 2–8°C or 2 months at –20°C.

**SST Tubes:** Do not use serum separator tubes. Manufacturer's studies have shown that Becton Dickinson SST<sup>®</sup> vacutainer tubes may cause elevation of values of approximately 30%.

## Warnings and Precautions

For *in vitro* diagnostic use.

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight. (See insert.)

**Water:** Use distilled or deionized water.

## Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

### Digoxin Test Units (LDI1)

Each barcode-labeled unit contains one bead coated with monoclonal murine anti-digoxin. Stable at 2–8°C until expiration date.

**LKDI1:** 100 units. **LKDI5:** 500 units.

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

### Digoxin Reagent Wedge (LDI2)

With barcode. 7.5 mL alkaline phosphatase (bovine calf intestine) conjugated to digoxin in buffer, with preservative. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is

within 30 days after opening when stored as indicated.

**LKDI1:** 1 wedge. **LKDI5:** 5 wedges.

### Digoxin Adjustors (LDIL, LDIH)

Two vials (Low and High), 4 mL each, of digoxin in processed human serum, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

**LKDI1:** 1 set. **LKDI5:** 2 sets.

## Kit Components Supplied Separately

**LSUBX:** Chemiluminescent Substrate

**LPWS2:** Probe Wash Module

**LKPM:** Probe Cleaning Kit

**LCHx-y:** Sample Cup Holders (barcoded)

**LSCP:** Sample Cups (disposable)

**LSCC:** Sample Cup Caps (optional)

**SDCM:** Bi-level, multi-constituent control

Also Required

Sample transfer pipets, distilled or deionized water, controls.

## Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

**Recommended Adjustment Interval:** 2 weeks.

**Quality Control Samples:** Use controls or sample pools with at least two levels (low and high) of digoxin.

## Therapeutic Range

The following results are reported in the literature, based on data taken six hours post-dose in patients with normal renal function.<sup>8</sup>

Group	Digoxin, ng/mL	Digoxin, nmol/L
Nontoxic patients	0.8 – 2.4	1.0 – 3.1
Toxic patients	2.1 – 8.7	2.7 – 11.1

Consider these limits as *guidelines* only.

### Limitations

Sole reliance on digoxin concentrations for determining digitalis toxicity is not warranted, and must be supplemented with additional clinical and electrocardiographic information.

Many physiological, pharmacological, pathological and genetic factors may affect the interpretation of results (see Summary and Explanation section).

The presence of endogenous digoxin-like immunoreactive factors (DLIF) has been reported in the serum of patients with renal and hepatic failure, and in newborns and women in the third trimester of pregnancy.<sup>14</sup> These factors can result in falsely elevated digoxin results in a number of commercially available immunoassays.

In cases of overdosage, specimens should be obtained *before* administering digoxin immune Fab antidote preparations (e.g. Digibind<sup>®</sup>) since these can be expected to interfere with any digoxin immunoassay procedure. Once such a preparation has been administered, accurate determinations of digoxin can be made only by specialized techniques (not routinely available in most clinical settings) until the Fab fragments have been eliminated from the body: this takes several days for patients with intact renal function, and a week or longer for patients with renal impairment. In addition to *in vitro* effects, these preparations, administered intravenously, cause a rapid rise in pharmacologically inactive, Fab-bound, total circulating digoxin that may be detected immunologically.<sup>16</sup>

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum

products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

### Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

**Conversion Factor:**  
ng/mL × 1.281 → nmol/L

**Calibration Range:** 0.5 – 8 ng/mL  
(0.6 – 10 nmol/L).

**Analytical Sensitivity:** 0.1 ng/mL  
(0.13 nmol/L).

**Precision:** Samples were processed in quadruplicate over the course of 5 days, 4 runs per day, for a total of 20 runs and 80 replicates. (See "Precision" table.)

**Linearity:** Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

**Recovery:** Samples spiked 1 to 19 with three digoxin solutions (13, 26 and 52 ng/mL) were assayed. (See "Recovery" table for representative data.)

**Specificity:** The antibody is highly specific for digoxin, though Lanatoside C and Methyl digoxin do exhibit significant crossreactivity. (See "Specificity" table.)

**Bilirubin:** Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

**Hemolysis:** Presence of hemoglobin in concentrations up to 513 mg/dL has no effect on results, within the precision of the assay.

**Lipemia:** Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

**Alternate Sample Type:** Matched samples from 18 volunteers, collected in plain serum tubes, and EDTA and heparinized plasma tubes were supplemented with digoxin to produce concentrations ranging from approximately 0.7 to 6.5 ng/mL. Average recovery for each sample type, relative to plain serum, was as follows:

EDTA plasma = 97%  
 Heparinized plasma = 95%

**Method Comparison 1:** The assay was compared to DPC's IMMULITE 2000 Digoxin on 128 serum samples. (Concentration range: approximately 0.5 to 6.5 ng/mL. See graph 1.) By linear regression:

(IML) = 1.04 (IML 2000) – 0.06 ng/mL  
 r = 0.985

Means:  
 1.5 ng/mL (IMMULITE)  
 1.5 ng/mL (IMMULITE 2000)

**Method Comparison 2:** The assay was compared to DPC's Coat-A-Count Digoxin on 116 serum samples. (Concentration range: approximately 0.5 to 3.3 ng/mL. See graph 2.) By linear regression:

(IML) = 1.0 (Coat-A-Count) + 0.12 ng/mL  
 r = 0.945

Means:  
 1.2 ng/mL (IMMULITE)  
 1.1 ng/mL (Coat-A-Count)

## References

1) Butler VP, Jr., Lindenbaum J. Serum digitalis measurements in the assessment of digitalis resistance and sensitivity. *Am J Med* 1975;58:460-9. 2) Doherty JE. How and when to use the digitalis serum levels. *J Am Med Assoc* 1978;239:2594-6. 3) Gault MH, et al. Studies of digoxin dosage, kinetics and serum concentrations in renal failure and review of the literature. *Nephron* 1976;17:161-87. 4) Huffman DH, et al. Association between clinical cardiac status, laboratory parameters, and digoxin usage. *Am Heart J* 1976;91:28-34. 5) Rietbrock N, Alken RG. The therapy of heart failure with digitalis. *Deutsche Med Wochenschr* 1980;105:1622-8. 6) Shapiro B, et al. Pitfalls in the application of digoxin determinations. *Sem Nucl Med* 1975;5:205-20. 7) Shapiro W. Current considerations in digoxin usage. *CRC Crit Rev Clin Lab Sci* 1978;9:321-46. 8) Smith TW, et al. Determination of therapeutic and toxic serum digoxin concentrations by radioimmunoassay. *N Engl J Med* 1969;281:1212-16. 9) Smith TW. Medical intelligence. *N Engl J Med* 1973;288:719-22 and 942-6. 10) Smith TW. Digitalis toxicity: epidemiology and clinical use of

serum concentration measurements. *Am J Med* 1975;58:470-6. 11) Weintraub M. Interpretation of the serum digoxin concentration. *Clin Pharmacokinetics* 1977;2:205-19. 12) Wettrell G. Digoxin therapy in infants. *Acta Paediatr Scand* 1976; Suppl 257:7-28. 13) Hallworth MJ. Practical therapeutic drug monitoring: digoxin. *Lab Medica* 1991 Mar/Apr;19:22. 14) Valdes R Jr. Endogenous digoxin-like immunoreactive factors: impact on digoxin measurements and potential physiological implications. *Clin Chem* 1985;9:1525-32. 15) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 16) Digoxin immune Fab. In: McEvoy GK, editor. *AHFS Drug Information* 92. American Society of Hospital Pharmacists, 1992: 1977-81.

## Technical Assistance

In the United States, contact DPC's Technical Services department.  
 Tel: 800.372.1782 or 973.927.2828  
 Fax: 973.927.4101. Outside the United States, contact your National Distributor.

Manufactured by EURO/DPC Ltd. under a Quality System registered to ISO 13485:2003.

## Tables and Graphs

### Precision (ng/mL)

	Mean <sup>3</sup>	Within-Run <sup>1</sup>			Total <sup>2</sup>	
		SD <sup>4</sup>	CV <sup>5</sup>	SD	CV	
1	0.7	0.04	5.7%	0.06	8.6%	
2	1.2	0.06	5.0%	0.08	6.7%	
3	2.1	0.09	4.3%	0.12	5.7%	
4	3.7	0.13	3.5%	0.23	6.2%	
5	6.4	0.23	3.6%	0.44	6.9%	

**Linearity (ng/mL)**



	Dilution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	8 in 8 <sup>5</sup>	4.8	—	—
	4 in 8	2.4	2.4	100%
	2 in 8	1.0	1.2	83%
	1 in 8	0.48	0.60	80%
2	8 in 8	5.7	—	—
	4 in 8	2.7	2.9	93%
	2 in 8	1.3	1.4	93%
	1 in 8	0.62	0.71	87%
3	8 in 8	7.0	—	—
	4 in 8	3.1	3.5	89%
	2 in 8	1.4	1.8	78%
	1 in 8	0.75	0.88	85%
4	8 in 8	8.0	—	—
	4 in 8	3.6	4.0	90%
	2 in 8	1.6	2.0	80%
	1 in 8	0.76	1.0	76%

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Fabricado pela EURO/DPC Ltd. de acordo com o Sistema de Qualidade registado segundo a norma ISO 13485:2003.
**EURO/DPC** LTD
Diagnostic Products Corporation
EC REP DPC Biermann GmbH