# IMMULITE®

## Phenytoin

For use on the IMMULITE®  
and IMMULITE® 1000 systems

**DPC®**

# IMMULITE®/IMMULITE® 1000 Phenytoin

## English

**Intended Use:** For *in vitro* diagnostic use with the IMMULITE and IMMULITE 1000 Analyzers — for the quantitative measurement of phenytoin in serum or heparinized plasma, as an aid in monitoring drug therapy.

Catalog Number: **LKPN1** (100 tests),  
**LKPN5** (500 tests)

Test Code: **PHN** Color: **Dark Gray**

### Summary and Explanation

Phenytoin (5,5-diphenylhydantoin; MW: 252.3 Da) is an antiepileptic drug widely used for treating several types of seizure. It is typically administered orally as phenytoin (Dilantin®) itself, or else by intravenous infusion (IV) or intramuscular injection (IM) of the prodrug fosphenytoin (Cerebyx®), a relatively water-soluble, phosphate ester of phenytoin. Depending on the route of administration, conversion of fosphenytoin to phenytoin is expected to be essentially complete within 2 to 4 hours.<sup>1-4</sup>

Most of the phenytoin in circulation is tightly bound to plasma proteins, chiefly albumin. Normally, only about 10% circulates in the pharmacologically active free form, though in some conditions, such as renal failure, this fraction may increase to 30% or more.<sup>5</sup> The amount circulating in free form is expected to be increased in hypoalbuminemia and in conditions where other compounds compete with phenytoin for protein binding sites. Examples include valproic acid and other anticonvulsants which may be administered concurrently with phenytoin; salicylates and nonsteroidal anti-inflammatory drugs such as oxaprozin (Daypro®); as well as many other drugs which bind strongly to albumin. Likewise, high free fatty acid levels may displace phenytoin from plasma proteins.

Phenytoin is cleared through the liver, where it is converted primarily to 5-(*p*-hydroxyphenyl)-5-phenylhydantoin, which is subsequently glucuronidated and excreted in the urine.<sup>2,5</sup> Though inactive,

these metabolites — HPPH and its glucuronide HPPG — are structurally similar to phenytoin and may compete with phenytoin for protein binding sites. The clearance rate declines somewhat with age, and is reduced in renal insufficiency. Drugs enhancing or inhibiting the liver's ability to metabolize phenytoin can be expected to affect circulating levels. Moreover, metabolic capacity is limited; hence the relationship between dosage and circulating levels is not uniform: at higher levels, even within the therapeutic range, small increases in administered dose may result in large increases in circulating levels.

Phenytoin is commonly monitored by immunoassay because (total) circulating levels of the drug correlate more strongly with therapeutic effect and toxicity than does the administered dose.<sup>6</sup> The generally accepted therapeutic range for total circulating phenytoin is 10 to 20 µg/mL (40 to 80 µmol/L), but these limits should be considered as guidelines only.<sup>1,7-9</sup> (There is no established range for free circulating phenytoin.) Dosage optimization must take into account the patient's clinical response to the drug. Low levels may be associated with inadequate control, as manifested in the frequency and/or severity of the seizures. High levels may be associated with symptoms of neurotoxicity — including involuntary rhythmic eye movements (nystagmus), nausea, vomiting, impaired coordination (ataxia), speech interference (dysarthria), severe lethargy, and so on — or even with a worsening of the seizures.

### Principle of the Procedure

IMMULITE/IMMULITE 1000 Phenytoin is a solid-phase, competitive chemiluminescent enzyme immunoassay.

**Incubation Cycles:** 1 × 30 minutes.

### Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Do not use gel barrier tubes.<sup>1,10-11</sup> For patients being treated with fosphenytoin (Cerebyx<sup>®</sup>), it is important not to collect samples for phenytoin analysis until at least 2 hours after (the completion of) intravenous infusion, or 4 hours after intramuscular injection, when conversion of the prodrug to phenytoin can be expected to be essentially complete.<sup>1-3</sup>

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE/IMMULITE 1000 Phenytoin has not been tested with all possible variations of tube types.

**Volume Required:** 20 µL serum or heparinized plasma. (Sample cup must contain at least 100 µL more than the total volume required.)

**Storage:** 7 days at 2–8°C or 1 month at –20°C.<sup>12</sup>

## Warnings and Precautions

For *in vitro* diagnostic use.

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight. (See insert.)

**Water:** Use distilled or deionized water.

## Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

### Phenytoin Test Units (LPN1)

Each barcode-labeled unit contains one bead coated with polyclonal rabbit anti-phenytoin antibody. Stable at 2–8°C until expiration date.

**LKPN1:** 100 units. **LKPN5:** 500 units.

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

### Phenytoin Reagent Wedge (LPN2)

With barcode. 7.5 mL alkaline phosphatase (bovine calf intestine) conjugated to phenytoin in buffer, with preservative. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

**LKPN1:** 2 wedges. **LKPN5:** 10 wedges.

### Phenytoin Adjustors (LPNL, LPNH)

Two vials (Low and High), 2 mL each, of phenytoin in processed human serum, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

**LKPN1:** 1 set. **LKPN5:** 2 sets.

## Kit Components Supplied Separately

**LSUBX:** Chemiluminescent Substrate

**LPWS2:** Probe Wash Module

**LKPM:** Probe Cleaning Kit

**LCHx-y:** Sample Cup Holders (barcoded)

**LSCP:** Sample Cups (disposable)

**LSCC:** Sample Cup Caps (optional)

**SDCM:** Bi-level, human serum-based control

Also Required  
Sample transfer pipets, distilled or deionized water, controls.

## Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

**Recommended Adjustment Interval:**  
2 weeks.

**Quality Control Samples:** Use controls or sample pools with at least two levels (low and high) of phenytoin.

## Therapeutic Range

The commonly accepted therapeutic range for phenytoin is 10 to 20 µg/mL (40 to 80 µmol/L).<sup>1,7-9</sup> Toxicity has been reported at concentrations greater than 20 µg/mL.

Consider these limits as guidelines only.

## Limitations

The following should be considered in the interpretation of the phenytoin concentration:

All concurrent medications taken by the patient

The time of sampling in relation to the dosing history

Existing diseases or conditions

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this

assay should always be used in combination with the clinical examination, patient medical history, and other findings.

## Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in µg/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

**Conversion Factor:**  
µg/mL × 3.96 → µmol/L

**Calibration Range:** 2.5 – 40 µg/mL  
(9.9 – 158.4 µmol/L).

**Analytical Sensitivity:** 0.03 µg/mL  
(0.12 µmol/L).

**Precision:** Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

**Linearity:** Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

**Recovery:** Samples spiked 1 to 19 with three phenytoin solutions (100, 200 and 400 µg/mL) were assayed. (See "Recovery" table for representative data.)

**Specificity:** The antibody is highly specific for phenytoin. There is an approximate 14% crossreactivity with the principal metabolite, HPPH, but no significant crossreactivity with the glucuronide ester secreted form, HPPG. (See "Specificity" table.)

**Bilirubin:** Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

**Hemolysis:** Presence of hemoglobin in concentrations up to 508 mg/dL has no effect on results, within the precision of the assay.

**Lipemia:** Presence of triglycerides in concentrations up to 3 000 mg/dL has no effect on results, within the precision of the assay.

**Method Comparisons:** The assay was compared to a commercially available phenytoin assay (Kit A) on 95 samples. (Concentration range: approximately 4.8 to 39 µg/mL. See "Method Comparison: Kit A" graph.) By linear regression:

(IML) = 1.04 (Kit A) + 0.87 µg/mL  
r = 0.988

Means:  
21.1 µg/mL (IMMULITE)  
18.6 µg/mL (Kit A)

The assay was compared to a commercially available phenytoin assay (Kit B) on 44 samples. (Concentration range: approximately 3 to 30 µg/mL. See "Method Comparison: Kit B" graph.) By linear regression:

(IMMULITE) = 1.1 (Kit B) + 0.39 µg/mL  
r = 0.987

Means:  
13.3 µg/mL (IMMULITE)  
11.6 µg/mL (Kit B)

## References

1) Warner A, Privitera M, Bates D. Standards of laboratory practice: antiepileptic drug monitoring. National Academy of Clinical Biochemistry. Clin Chem 1998;44:1085-95. 2) Roberts WL, De BK, Coleman JP, Annesley TM. Falsely increased immunoassay measurements of total and unbound phenytoin in critically ill uremic patients receiving fosphenytoin. Clin Chem 1999;45:829-37. 3) Kugler AR, Annesley TM, Nordblom GD, et al. Cross-reactivity of fosphenytoin in two human plasma phenytoin immunoassays. Clin Chem 1998;44:1474-80. 4) Aweeka FT, Gottwald MD, Gambertoglio JG, et al. Pharmacokinetics of fosphenytoin in patients with hepatic or renal disease. Epilepsia 1999;40:777-82. 5) Rainey PM, Rogers KE, Roberts WL. Metabolite and matrix interference in phenytoin immunoassays. Clin Chem 1996;42:1645-53. 6) Yukawa E. Optimisation of antiepileptic drug therapy. The importance of serum drug concentration monitoring. Clin Pharmacokinet 1996;31:120-30. 7) Burtis CA, Ashwood ER, editors. Tietz textbook of clinical chemistry. 3rd ed. Philadelphia: WB Saunders, 1999: 874-5. 8) Hayes G, Kootsikis ME. Reassessing the lower end of the phenytoin therapeutic range: a review of the literature. Ann Pharmacother 1993;27:1389-92. 9) Schumacher GE, Barr JT, Browne TR, Collins JF. Test performance characteristics of the serum phenytoin concentration (SPC): the relationship between SPC and patient response. The Veterans Administration Epilepsy Cooperative Study Group. Ther Drug Monit 1991;13:318-24. 10) Cai WM, Leader WG, Porter WH, Chandler MH. Influence of serum separator tubes on total and free phenytoin concentrations and dosages. Ther Drug Monit 1993;15:427-30. 11) Mauro LS, Mauro VF. Effect of serum separator tubes on free and total phenytoin and carbamazepine serum concentrations. Ther Drug Monit 1991;13:240-3. 12) Quest Diagnostics. The 1997/1998 Quest Diagnostics Nichols Institute reference manual. San Juan Capistrano: Quest

Diagnostics, 1997: 228. 13) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.

## Technical Assistance

In the United States, contact DPC's Technical Services department.  
Tel: 800.372.1782 or 973.927.2828  
Fax: 973.927.4101. Outside the United States, contact your National Distributor.

The Quality System of Diagnostic Products Corporation is registered to ISO 13485:2003.

## Tables and Graphs

### Precision (µg/mL)

|   | Mean <sup>3</sup> | Within-Run <sup>1</sup> |                 | Total <sup>2</sup> |      |
|---|-------------------|-------------------------|-----------------|--------------------|------|
|   |                   | SD <sup>4</sup>         | CV <sup>5</sup> | SD                 | CV   |
| 1 | 6.8               | 0.40                    | 5.9%            | 0.42               | 6.2% |
| 2 | 9.3               | 0.51                    | 5.5%            | 0.57               | 6.1% |
| 3 | 14.1              | 0.70                    | 5.0%            | 0.76               | 5.4% |
| 4 | 24.1              | 1.4                     | 5.8%            | 1.5                | 6.2% |
| 5 | 29.6              | 1.8                     | 6.1%            | 2.1                | 7.1% |

### Linearity (µg/mL)

|   | Dilution <sup>1</sup> | Observed <sup>2</sup> | Expected <sup>3</sup> | %O/E <sup>4</sup> |
|---|-----------------------|-----------------------|-----------------------|-------------------|
| 1 | 8 in 8 <sup>5</sup>   | 22.4                  | —                     | —                 |
|   | 4 in 8                | 11.0                  | 11.2                  | 98%               |
|   | 2 in 8                | 5.3                   | 5.6                   | 94%               |
|   | 1 in 8                | 2.5                   | 2.8                   | 88%               |
| 2 | 8 in 8                | 24.4                  | —                     | —                 |
|   | 4 in 8                | 13.1                  | 12.2                  | 108%              |
|   | 2 in 8                | 6.8                   | 6.1                   | 111%              |
|   | 1 in 8                | 2.7                   | 3.0                   | 87%               |
| 3 | 8 in 8                | 32.8                  | —                     | —                 |
|   | 4 in 8                | 17.9                  | 16.4                  | 109%              |
|   | 2 in 8                | 7.6                   | 8.2                   | 93%               |
|   | 1 in 8                | 3.6                   | 4.1                   | 88%               |
| 4 | 8 in 8                | 35.7                  | —                     | —                 |
|   | 4 in 8                | 18.3                  | 17.9                  | 103%              |
|   | 2 in 8                | 8.6                   | 8.9                   | 96%               |
|   | 1 in 8                | 4.2                   | 4.5                   | 94%               |

### Recovery (µg/mL)


|   | Solution <sup>1</sup> | Observed <sup>2</sup> | Expected <sup>3</sup> | %O/E <sup>4</sup> |
|---|-----------------------|-----------------------|-----------------------|-------------------|
| 1 | —                     | 4.1                   | —                     | —                 |
|   | A                     | 10.2                  | 8.9                   | 115%              |
|   | B                     | 15.8                  | 13.9                  | 114%              |
|   | C                     | 23.6                  | 23.9                  | 99%               |
| 2 | —                     | 4.0                   | —                     | —                 |
|   | A                     | 9.5                   | 8.8                   | 108%              |
|   | B                     | 14.1                  | 13.8                  | 102%              |
|   | C                     | 25.8                  | 23.8                  | 108%              |
| 3 | —                     | 8.6                   | —                     | —                 |
|   | A                     | 13.4                  | 13.2                  | 101%              |
|   | B                     | 17.9                  | 18.2                  | 98%               |
|   | C                     | 28.1                  | 28.2                  | 100%              |
| 4 | —                     | 9.2                   | —                     | —                 |
|   | A                     | 13.1                  | 13.7                  | 96%               |
|   | B                     | 18.1                  | 18.7                  | 97%               |
|   | C                     | 28.0                  | 28.7                  | 98%               |

### Specificity (µg/mL)


| Compound <sup>1</sup> | Amount Added <sup>2</sup> | Apparent Conc. <sup>3</sup> | % Cross-reactivity <sup>4</sup> |
|-----------------------|---------------------------|-----------------------------|---------------------------------|
| Clonazepam            | 0.12                      | ND                          | ND                              |
| Diazepam              | 2.5                       | ND                          | ND                              |
| Ethosuximide          | 100                       | ND                          | ND                              |
| Imipramine            | 1.0                       | ND                          | ND                              |
| Phenylbutazone        | 45                        | ND                          | ND                              |
| Primidone             | 12                        | ND                          | ND                              |
| Secobarbital          | 15                        | ND                          | ND                              |
| Carbamazepine         | 20                        | ND                          | ND                              |
| Phenobarbital         | 40                        | ND                          | ND                              |
| Theophylline          | 20                        | ND                          | ND                              |
| Valproic acid         | 100                       | ND                          | ND                              |
| HPPH                  | 20                        | 2.8                         | 14.2%                           |
| HPPG                  | 100                       | 0.03                        | 0.03%                           |

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O Sistema de Qualidade da Diagnostic Products Corporation está registado sob ISO 13485:2003.
EC REP DPC Biemann GmbH